

May 19, 2015

URGENT PRODUCT CORRECTION NOTIFICATION
Software Anomaly Affecting User Defined Protocol Feature on ORTHO VISION™ Analyzer for ORTHO BioVue® Cassettes

Please distribute this information to the appropriate personnel at your facility

Dear Valued Customer,

This notification is to inform you of an Urgent Product Correction Notification involving the following product:

Affected System	Affected Software Version	Product Code
ORTHO VISION™ Analyzer	1.0.4	6904579

Description of the Issue

Ortho-Clinical Diagnostics, Inc. (OCD) has identified an anomaly in ORTHO VISION™ Analyzer software that may occur when an operator modifies an assay column name when using the User Defined Protocol (UDP) feature. When a User Defined Protocol (UDP) is created, if the operator modifies the cassette selection and changes the column names in the UDP test, under specific conditions the column name for the test template may get changed until the system is restarted. After the system is restarted, the name that was modified by the operator will be cleared and the column name of the test template will return to the default configuration for the standard assay template.

Impact to Results

If a result is reported solely based on a column result of a test template affected by this issue using the Graphical User Interface (GUI) or Order Report, an erroneous test result could be reported.

To date, no customer complaints or patient injury due to this issue have been reported to OCD.

Resolution

This issue has been resolved in ORTHO VISION™ Analyzer Software Version 1.2.10, which will be released imminently. The UDP feature should not be used until Software Version 1.2.10 is installed and validated according to your local requirements.

Actions Required from You

1. Do not use the User Defined Protocol (UDP) feature of your ORTHO VISION™ Analyzer until Software Version 1.2.10 is installed and validated according to your local requirements.
2. Post this notification by each ORTHO VISION™ Analyzer or with the user documentation.
3. Complete and return the **Customer Acknowledgement Form** within 2 business days to acknowledge your reading and understanding of this notice.

Ortho Clinical Diagnostics

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

We sincerely apologize for the inconvenience this may cause your laboratory. If you have any additional questions, please contact Customer Technical Services at **1800 5646 766**.

Yours sincerely,



Cherie Yip
Regulatory Associate