



2 Science Park Drive, #07-13 Ascent
Singapore Science Park 1
Singapore 118222

URGENT
FIELD SAFETY NOTICE – VOLUNTARY RECALL (FINAL)

18 MARCH 2019

To Whom It May Concern:

RE: Voluntary Field Action of specific ACUVUE® OASYS for ASTIGMATISM Brand Contact Lenses

Lot NUMBER: B00RWF0

On March 15, 2019, Johnson & Johnson Vision Care, Inc. initiated a recall in United States of America associated with two failures found during bounding testing for a package integrity issue. We are voluntarily recalling one (1) lots of ACUVUE® OASYS for ASTIGMATISM Contact Lenses because they do not meet our quality standards. This recall **does not** impact product distributed in **Asia Pacific countries including Singapore**.


In United States of America, limited quantity of 300 lenses shipped to one distributor, all product to be returned and scheduled for destruction. **No impacted lot** was in the hands of end consumers.

The cause of the issue was identified as a build-up on the saline dose tips caused saline to be dosed sideways onto packages wetting the seal area. Countermeasures have been put in place to reduce the likelihood of this occurring again.

At Johnson & Johnson Vision Care, Inc., we continue to strive for the highest quality products for our customers. ACUVUE® brand products are safe when used as directed and can continue to be used with confidence. Our contact lenses undergo rigorous evaluation to provide assurance that they meet our own high quality standards as well as regulatory standards in all markets where they are sold. We apologize for any inconvenience this causes and thank you for your understanding.

Sincerely,

Christopher Chong
Meng Lim



Christopher Lim
Director, Business Regional Quality Asia Pacific
Johnson & Johnson Vision

CC:

Lex Tan – QA Specialist DID: +65 6918 6254 | Email: [REDACTED]

Alvin Loh – Business Regional Quality, Manager ASEAN DID (65) 6918 7295 | Email: [REDACTED]

Eric Chin – Regulatory Affairs, Senior Manager_| Email: [REDACTED]