

URGENT: FIELD SAFETY NOTICE

Trauma Sterile Tube (M size)

2 January, 2020

Dear Operating Room Supervisors, Materials Management Personnel, and Chief of Surgery:

Synthes GmbH is initiating a voluntary field safety notice for all lots of Sterile Screw implants packed in Sterile Tubes with part numbers ending in suffix TS with inner tube components manufactured by Flex Precision Plastics Solutions since June 17, 2019.

Product Info:

Sterile screw implants packed in Sterile Tubes (size M) with part numbers ending in suffix TS. Sterile Tube Packaging is a packaging system to deliver single, sterile screws to the sterile field. See attached list of affected product codes.

Reason for Field Safety Notification:

There is a possibility that

- 1) the inner cap could open when the outer cap is removed (picture 1) or
- 2) the entire inner tube remains in the outer tube (picture 2)

In both instances this could prevent the screw from being removed and used causing a surgical delay. The inner tube should remain closed and retained in the outer cap when presented to the sterile field.

Picture 1



Picture 2



Potential Patient Impact:

The possible clinical implications related to this issue may include:

- A surgical delay may occur while a replacement screw implant packaged in a sterile tube, blister or non-sterile variant is located, depending on availability.

Instructions to follow:

To reduce the possibility of the inner cap coming loose or the inner tube remaining in the outer tube while removing the outer cap follow these techniques:

- Option 1: Suggested opening technique for all affected lots

1. When opening the outer tube, ensure that the inner cap and inner tube are intact and turning freely in the outer tube



2. If the inner cap is not turning freely, push the outer cap towards the tube and keep turning counter-clockwise until you see the inner tube and the inner cap turning in the outer tube (remove the labels to increase visibility).



3. Once you see that the inner cap is turning in the outer tube then pull out the outer cap with the inner tube and inner cap intact. Hold the outer cap at an upward angle as you separate the outer cap from the outer tube.



4. Continue as usual by following steps 2-4 in the Sterile Tube User Guide (106926-190205 DSEM).

- Option 2: In case Option 1 didn't work

1. Reject the screw and use a replacement sterile tube, blister or non-sterile variant, depending on availability.
2. Please ensure a product complaint is filed and product returned for the issue according to your local procedure.

ACTION REQUIRED:

Our records show that your facility has received products that are subject to this field safety notification. We ask that you review the information contained in this notification and complete the Verification Section located on page 4 of this letter.

1. Review the instructions listed in this notification.
2. Forward this notice to anyone in your facility that needs to be informed.
3. If any of the subject product has been forwarded to another facility, contact that facility.
4. Keep a copy of this notice with the subject product.
5. Review, complete, sign and return the attached business reply form of this notice and hand it to your Depuy Synthes sales representative within (3) business days.


If you require any assistance, please contact your Depuy Synthes sales representative.

Attachments:

Attachment 1: Business Reply Form

Attachment 2: List of product part & lot number

Yours sincerely,



Lee Ching Hwee
Manager, Regulatory Affairs

cc: Chairman Medical Board
Relevant Head of Departments

ATTACHMENT 1: Business Reply Form (BRF)

Your timely response to this customer notification is requested. Please complete this form and hand it to your Deputy Synthes sales representative **within 3 business days**.

☐ I have read and understand the notification

Print Name of Person Completing Business Reply Form:	Sign & Date:
Customer Account Details (Name and Address)	

Attachment 2: List of product part & lot numbers:

Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.
02.125.124TS	5L55742	02.200.012TS	6L42649	02.200.022TS	6L30486	02.200.044TS	6L28353	02.210.942TS	6L11679
02.125.124TS	5L66887	02.200.012TS	6L42144	02.200.022TS	6L31040	02.200.044TS	6L30704	02.210.942TS	6L11682
02.127.130TS	5L47705	02.200.012TS	6L42854	02.200.022TS	6L41824	02.200.045TS	6L27690	02.211.010TS	5L57889
02.127.130TS	5L47706	02.200.014TS	5L56114	02.200.022TS	6L42857	02.200.045TS	6L29730	02.211.010TS	6L18450
02.127.130TS	5L47707	02.200.014TS	6L16422	02.200.024TS	5L57886	02.200.045TS	6L31442	02.211.010TS	6L38921
02.127.130TS	5L49910	02.200.014TS	6L20256	02.200.024TS	6L16424	02.200.046TS	6L14338	02.211.012TS	5L56127
02.127.130TS	5L52128	02.200.014TS	6L21149	02.200.024TS	6L22718	02.200.046TS	6L28354	02.211.014TS	5L49917
02.127.130TS	5L53660	02.200.014TS	6L25145	02.200.024TS	6L26644	02.200.046TS	6L30705	02.211.014TS	5L57891
02.127.130TS	6L33811	02.200.014TS	6L27688	02.200.024TS	6L30487	02.200.048TS	6L30492	02.211.014TS	6L28356
02.127.130TS	6L35523	02.200.014TS	6L26639	02.200.024TS	6L32139	02.200.048TS	6L42656	02.211.016TS	5L56128
02.127.130TS	6L35539	02.200.014TS	6L28036	02.200.024TS	6L31041	02.200.050TS	6L25149	02.211.016TS	6L25150
02.127.130TS	6L36827	02.200.014TS	6L18448	02.200.026TS	5L56119	02.200.050TS	6L28355	02.211.016TS	6L28053
02.127.130TS	6L37683	02.200.014TS	6L28349	02.200.028TS	5L49916	02.200.050TS	6L40751	02.211.018TS	5L53662
02.127.130TS	6L39067	02.200.014TS	6L30698	02.200.028TS	6L20260	02.200.050TS	6L30706	02.211.018TS	6L16426
02.127.130TS	6L39619	02.200.014TS	6L29726	02.200.028TS	6L22719	02.200.050TS	6L42654	02.211.018TS	6L25151
02.127.130TS	6L41717	02.200.014TS	6L31038	02.200.028TS	6L26645	02.200.055TS	5L56121	02.211.018TS	6L28054
02.127.130TS	6L42372	02.200.014TS	6L41823	02.200.028TS	6L30488	02.200.055TS	6L28044	02.211.022TS	6L25152
02.127.130TS	6L42718	02.200.014TS	6L42145	02.200.028TS	6L32140	02.200.055TS	6L30493	02.211.024TS	5L53663
02.127.130TS	6L42719	02.200.014TS	6L42647	02.200.028TS	6L41825	02.200.055TS	6L31443	02.211.024TS	6L25153
02.127.130TS	6L44267	02.200.014TS	6L42650	02.200.030TS	5L57887	02.210.106TS	6L31485	02.211.024TS	6L28374
02.127.130TS	6L44253	02.200.014TS	6L42856	02.200.030TS	6L20261	02.210.110TS	6L30494	02.211.026TS	6L28375
02.127.136TS	6L34623	02.200.016TS	5L57884	02.200.030TS	6L22734	02.210.112TS	6L35133	02.211.026TS	6L31043
02.130.324TS	5L49913	02.200.016TS	6L16423	02.200.030TS	6L28037	02.210.112TS	6L35429	02.211.030TS	6L28376
02.130.324TS	5L49912	02.200.016TS	6L20257	02.200.030TS	6L30489	02.210.112TS	6L42150	02.211.032TS	6L20264
02.130.324TS	5L49915	02.200.016TS	6L21151	02.200.030TS	6L29727	02.210.112TS	6L41829	02.211.034TS	5L49918
02.130.324TS	6L21133	02.200.016TS	6L25146	02.200.030TS	6L32141	02.210.114TS	5L52594	02.211.038TS	6L18451
02.130.324TS	6L21135	02.200.018TS	5L56116	02.200.030TS	6L31438	02.210.114TS	6L37348	02.211.040ST	5L53664
02.130.324TS	6L21136	02.200.018TS	6L20258	02.200.030TS	6L42651	02.210.114TS	6L38918	02.211.040ST	4L56097
02.130.324TS	6L21137	02.200.018TS	6L21152	02.200.032TS	6L22735	02.210.116TS	5L56122	02.211.042ST	5L37907
02.130.324TS	6L21138	02.200.018TS	6L22717	02.200.032TS	6L28039	02.210.116TS	6L35430	02.211.042ST	6L20617
02.130.324TS	6L21139	02.200.018TS	6L25147	02.200.032TS	6L29729	02.210.116TS	6L42156	02.211.042ST	6L20593
02.130.324TS	6L21141	02.200.018TS	6L26640	02.200.034TS	5L53661	02.210.116TS	6L41830	02.211.044TS	6L42151
02.130.324TS	6L21132	02.200.018TS	6L30103	02.200.034TS	6L22736	02.210.118TS	6L18449	02.211.048TS	5L57892
02.130.324TS	6L26636	02.200.018TS	6L30484	02.200.034TS	6L28040	02.210.118TS	6L35134	02.211.048TS	6L20265
02.130.324TS	6L26632	02.200.018TS	6L30699	02.200.034TS	6L31439	02.210.118TS	6L37349	02.211.050TS	5L60577
02.130.324TS	6L26634	02.200.018TS	6L27689	02.200.034TS	6L31042	02.210.118TS	6L39440	02.211.050TS	6L31044
02.130.324TS	6L26633	02.200.018TS	6L40749	02.200.036TS	5L56120	02.210.118TS	6L42146	02.211.052TS	6L28377
02.130.324TS	6L26635	02.200.020TS	5L57885	02.200.036TS	6L16425	02.210.118TS	6L41839	02.211.054TS	5L49920
02.200.010TS	6L14332	02.200.020TS	6L21153	02.200.036TS	6L22744	02.210.120TS	5L61650	02.211.056TS	6L28378
02.200.010TS	6L26637	02.200.020TS	6L25148	02.200.036TS	6L28041	02.210.120TS	6L37352	04.130.326TS	5L47708
02.200.010TS	6L30482	02.200.020TS	6L26641	02.200.036TS	6L30490	02.210.120TS	6L42147	04.130.326TS	5L47709
02.200.010TS	6L32135	02.200.020TS	6L28350	02.200.036TS	6L31440	02.210.120TS	6L41840	04.130.326TS	6L31049
02.200.012TS	5L61648	02.200.020TS	6L20259	02.200.036TS	6L41827	02.210.122TS	5L56123	04.130.326TS	6L31459
02.200.012TS	6L20255	02.200.020TS	6L30485	02.200.038TS	6L28042	02.210.122TS	5L65722	04.130.326TS	6L31454
02.200.012TS	6L21147	02.200.020TS	6L30700	02.200.038TS	6L28351	02.210.122TS	6L42148	04.130.326TS	6L31455
02.200.012TS	6L22716	02.200.020TS	6L29728	02.200.038TS	6L30701	02.210.122TS	6L41841	04.130.326TS	6L31456
02.200.012TS	6L25144	02.200.020TS	6L32137	02.200.040TS	5L57888	02.210.124TS	6L38919	04.130.326TS	6L31457
02.200.012TS	6L26638	02.200.020TS	6L31039	02.200.040TS	6L22746	02.210.124TS	6L42149	04.130.326TS	6L31458
02.200.012TS	6L30102	02.200.020TS	6L42648	02.200.040TS	6L28043	02.210.126TS	6L30104	04.130.326TS	6L31484
02.200.012TS	6L30483	02.200.020TS	6L42652	02.200.040TS	6L30702	02.210.128TS	6L42858	04.130.326TS	6L32159
02.200.012TS	6L31037	02.200.020TS	6L42655	02.200.040TS	6L30491	02.210.942TS	6L11674	04.130.326TS	6L32160
02.200.012TS	6L33565	02.200.022TS	5L56118	02.200.040TS	6L31441	02.210.942TS	6L11675	04.130.326TS	6L32164
02.200.012TS	6L32136	02.200.022TS	6L21154	02.200.040TS	6L41828	02.210.942TS	6L11677	04.130.326TS	6L32142
02.200.012TS	6L30697	02.200.022TS	6L26643	02.200.042TS	6L30703	02.210.942TS	6L11681	04.130.326TS	6L44066
02.200.012TS	6L29724	02.200.022TS	6L32138	02.200.042TS	6L28352	02.210.942TS	6L11685	04.200.010TS	5L54914
02.200.012TS	6L40748	02.200.022TS	6L40750	02.200.042TS	6L42653	02.210.942TS	6L11684	04.200.010TS	5L57903

Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.
04.200.010TS	5L57942	04.200.012TS	6L42721	04.200.014TS	6L39292	04.200.016TS	6L35652	04.200.018TS	6L35653
04.200.010TS	5L57948	04.200.012TS	6L42759	04.200.014TS	6L39293	04.200.016TS	6L35624	04.200.018TS	6L35625
04.200.010TS	6L16427	04.200.012TS	6L43775	04.200.014TS	6L39591	04.200.016TS	6L36609	04.200.018TS	6L36610
04.200.010TS	6L26379	04.200.014TS	5L46304	04.200.014TS	6L39592	04.200.016TS	6L37242	04.200.018TS	6L37243
04.200.010TS	6L34081	04.200.014TS	5L49941	04.200.014TS	6L41808	04.200.016TS	6L37443	04.200.018TS	6L37444
04.200.010TS	6L36630	04.200.014TS	5L52151	04.200.014TS	6L41810	04.200.016TS	6L37444	04.200.018TS	6L37447
04.200.010TS	6L36606	04.200.014TS	5L52131	04.200.014TS	6L41812	04.200.016TS	6L37445	04.200.018TS	6L37448
04.200.010TS	6L37241	04.200.014TS	5L56131	04.200.014TS	6L42036	04.200.016TS	6L37699	04.200.018TS	6L37700
04.200.010TS	6L37436	04.200.014TS	5L57949	04.200.014TS	6L43776	04.200.016TS	6L37693	04.200.018TS	6L37701
04.200.010TS	6L37437	04.200.014TS	5L56130	04.200.014TS	6L42749	04.200.016TS	6L37694	04.200.018TS	6L37702
04.200.010TS	6L39589	04.200.014TS	5L56133	04.200.014TS	6L42754	04.200.016TS	6L37695	04.200.018TS	6L38750
04.200.010TS	6L39590	04.200.014TS	5L57907	04.200.014TS	6L42722	04.200.016TS	6L37696	04.200.018TS	6L38751
04.200.010TS	6L44495	04.200.014TS	5L57944	04.200.014TS	6L43809	04.200.016TS	6L37697	04.200.018TS	6L38752
04.200.010TS	6L44496	04.200.014TS	5L60134	04.200.014TS	6L43815	04.200.016TS	6L37698	04.200.018TS	6L38753
04.200.012TS	5L47710	04.200.014TS	5L60160	04.200.014TS	6L44067	04.200.016TS	6L39294	04.200.018TS	6L40757
04.200.012TS	5L52150	04.200.014TS	5L60163	04.200.014TS	6L44078	04.200.016TS	6L40754	04.200.018TS	6L40763
04.200.012TS	5L52130	04.200.014TS	6L14155	04.200.014TS	6L44082	04.200.016TS	6L40762	04.200.018TS	6L40767
04.200.012TS	5L56129	04.200.014TS	6L16429	04.200.016TS	5L49942	04.200.016TS	6L40766	04.200.018TS	6L41720
04.200.012TS	5L57904	04.200.014TS	6L27595	04.200.016TS	5L49943	04.200.016TS	6L41719	04.200.018TS	6L41798
04.200.012TS	6L16428	04.200.014TS	6L27596	04.200.016TS	5L53665	04.200.016TS	6L41797	04.200.018TS	6L41801
04.200.012TS	6L18220	04.200.014TS	6L27597	04.200.016TS	5L54958	04.200.016TS	6L41800	04.200.018TS	6L42067
04.200.012TS	6L18221	04.200.014TS	6L27598	04.200.016TS	5L54960	04.200.016TS	6L42066	04.200.018TS	6L42070
04.200.012TS	6L18222	04.200.014TS	6L27599	04.200.016TS	5L54915	04.200.016TS	6L42069	04.200.018TS	6L42053
04.200.012TS	6L18223	04.200.014TS	6L28011	04.200.016TS	5L57908	04.200.016TS	6L42052	04.200.018TS	6L42751
04.200.012TS	6L28008	04.200.014TS	6L28012	04.200.016TS	5L56134	04.200.016TS	6L42750	04.200.018TS	6L42756
04.200.012TS	6L28009	04.200.014TS	6L28013	04.200.016TS	5L57950	04.200.016TS	6L42755	04.200.018TS	6L42724
04.200.012TS	6L28010	04.200.014TS	6L34983	04.200.016TS	5L60135	04.200.016TS	6L42723	04.200.018TS	6L43778
04.200.012TS	6L29338	04.200.014TS	6L35313	04.200.016TS	5L60161	04.200.016TS	6L43777	04.200.018TS	6L43812
04.200.012TS	6L34082	04.200.014TS	6L35307	04.200.016TS	5L60164	04.200.016TS	6L43810	04.200.018TS	6L43820
04.200.012TS	6L34083	04.200.014TS	6L35308	04.200.016TS	6L14156	04.200.016TS	6L43816	04.200.018TS	6L44069
04.200.012TS	6L35100	04.200.014TS	6L35309	04.200.016TS	6L14157	04.200.016TS	6L44068	04.200.018TS	6L44080
04.200.012TS	6L35301	04.200.014TS	6L35310	04.200.016TS	6L16431	04.200.016TS	6L44079	04.200.020TS	5L54917
04.200.012TS	6L35303	04.200.014TS	6L35311	04.200.016TS	6L16430	04.200.016TS	6L44083	04.200.020TS	5L60136
04.200.012TS	6L35302	04.200.014TS	6L35312	04.200.016TS	6L16432	04.200.018TS	5L47711	04.200.020TS	5L60162
04.200.012TS	6L35304	04.200.014TS	6L35314	04.200.016TS	6L16433	04.200.018TS	5L53686	04.200.020TS	6L12442
04.200.012TS	6L35305	04.200.014TS	6L35623	04.200.016TS	6L16434	04.200.018TS	5L53666	04.200.020TS	6L14162
04.200.012TS	6L35306	04.200.014TS	6L36608	04.200.016TS	6L16435	04.200.018TS	5L54959	04.200.020TS	6L18225
04.200.012TS	6L36607	04.200.014TS	6L36937	04.200.016TS	6L18224	04.200.018TS	5L54916	04.200.020TS	6L18226
04.200.012TS	6L36935	04.200.014TS	6L36990	04.200.016TS	6L26381	04.200.018TS	6L12441	04.200.020TS	6L20120
04.200.012TS	6L36989	04.200.014TS	6L37441	04.200.016TS	6L26382	04.200.018TS	6L14159	04.200.020TS	6L26384
04.200.012TS	6L37439	04.200.014TS	6L37442	04.200.016TS	6L27601	04.200.018TS	6L14160	04.200.020TS	4L55485
04.200.012TS	6L37440	04.200.014TS	6L37689	04.200.016TS	6L26383	04.200.018TS	6L16438	04.200.020TS	6L34972
04.200.012TS	6L37438	04.200.014TS	6L37687	04.200.016TS	6L27600	04.200.018TS	6L16439	04.200.020TS	6L34986
04.200.012TS	6L37684	04.200.014TS	6L37688	04.200.016TS	6L27602	04.200.018TS	6L16440	04.200.020TS	6L34992
04.200.012TS	6L37685	04.200.014TS	6L37690	04.200.016TS	6L26380	04.200.018TS	6L16436	04.200.020TS	6L35320
04.200.012TS	6L37686	04.200.014TS	6L37691	04.200.016TS	6L28014	04.200.018TS	6L16437	04.200.020TS	6L35321
04.200.012TS	6L38863	04.200.014TS	6L37692	04.200.016TS	6L34084	04.200.018TS	6L27603	04.200.020TS	6L35322
04.200.012TS	6L40760	04.200.014TS	6L38746	04.200.016TS	6L34851	04.200.018TS	6L27604	04.200.020TS	6L35630
04.200.012TS	6L40752	04.200.014TS	6L38747	04.200.016TS	6L34970	04.200.018TS	6L27605	04.200.020TS	6L36611
04.200.012TS	6L41718	04.200.014TS	6L38748	04.200.016TS	6L34991	04.200.018TS	6L27606	04.200.020TS	6L37244
04.200.012TS	6L41796	04.200.014TS	6L38749	04.200.016TS	6L34984	04.200.018TS	6L27607	04.200.020TS	6L37245
04.200.012TS	6L42051	04.200.014TS	6L39068	04.200.016TS	6L35317	04.200.018TS	6L27608	04.200.020TS	6L37449
04.200.012TS	6L42748	04.200.014TS	6L39075	04.200.016TS	6L35318	04.200.018TS	6L28015	04.200.020TS	6L37450
04.200.012TS	6L42753	04.200.014TS	6L39069	04.200.016TS	6L35315	04.200.018TS	6L34971	04.200.020TS	6L38864
04.200.012TS	6L42758	04.200.014TS	6L39290	04.200.016TS	6L35316	04.200.018TS	6L34985	04.200.020TS	6L40764
04.200.012TS	6L42757	04.200.014TS	6L39291	04.200.016TS	6L35316	04.200.018TS	6L35319	04.200.020TS	6L40756

Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.
04.200.020TS	6L40768	04.200.024TS	6L43818	04.200.028TS	6L43781	04.200.034TS	6L29340	04.200.042TS	6L40669
04.200.020TS	6L41721	04.200.024TS	6L44071	04.200.028TS	6L42770	04.200.034TS	6L34990	04.200.042TS	6L41703
04.200.020TS	6L41799	04.200.026TS	6L12151	04.200.028TS	6L42767	04.200.034TS	6L34979	04.200.042TS	6L40670
04.200.020TS	6L41806	04.200.026TS	6L12153	04.200.028TS	6L42764	04.200.034TS	6L38764	04.200.042TS	6L44500
04.200.020TS	6L42054	04.200.026TS	6L16442	04.200.028TS	6L42765	04.200.034TS	6L39308	04.200.044TS	5L49946
04.200.020TS	6L43779	04.200.026TS	6L16444	04.200.028TS	6L43814	04.200.034TS	6L41699	04.200.044TS	5L60138
04.200.020TS	6L42752	04.200.026TS	6L16445	04.200.028TS	6L43819	04.200.034TS	6L41700	04.200.044TS	6L18239
04.200.020TS	6L42725	04.200.026TS	6L18231	04.200.030TS	5L56138	04.200.036TS	5L53671	04.200.044TS	6L24944
04.200.020TS	6L42760	04.200.026TS	6L18232	04.200.030TS	6L14170	04.200.036TS	6L12156	04.200.044TS	6L26376
04.200.020TS	6L43811	04.200.026TS	6L18233	04.200.030TS	6L16447	04.200.036TS	6L16451	04.200.044TS	6L35327
04.200.020TS	6L43817	04.200.026TS	6L18234	04.200.030TS	6L18237	04.200.036TS	6L16450	04.200.044TS	6L36741
04.200.020TS	6L44070	04.200.026TS	6L18235	04.200.030TS	6L18238	04.200.036TS	6L24940	04.200.044TS	6L37246
04.200.022TS	5L47712	04.200.026TS	6L26386	04.200.030TS	6L24937	04.200.036TS	6L24941	04.200.044TS	6L37247
04.200.022TS	5L54918	04.200.026TS	6L26373	04.200.030TS	6L24938	04.200.036TS	6L34980	04.200.044TS	6L37467
04.200.022TS	5L61654	04.200.026TS	6L27613	04.200.030TS	6L34977	04.200.036TS	6L37459	04.200.044TS	6L39599
04.200.022TS	6L12149	04.200.026TS	6L34975	04.200.030TS	6L34989	04.200.036TS	6L38766	04.200.044TS	6L39600
04.200.022TS	6L18227	04.200.026TS	6L34988	04.200.030TS	6L35654	04.200.036TS	6L38767	04.200.044TS	6L39601
04.200.022TS	6L26385	04.200.026TS	6L35324	04.200.030TS	6L35627	04.200.036TS	6L42039	04.200.044TS	6L39602
04.200.022TS	6L27609	04.200.026TS	6L35325	04.200.030TS	6L36938	04.200.036TS	6L42040	04.200.044TS	6L42043
04.200.022TS	6L27610	04.200.026TS	6L39297	04.200.030TS	6L37456	04.200.038TS	5L52133	04.200.044TS	6L42044
04.200.022TS	6L28016	04.200.026TS	6L39298	04.200.030TS	6L38762	04.200.038TS	5L57910	04.200.045TS	6L14174
04.200.022TS	6L28017	04.200.026TS	6L39299	04.200.030TS	6L39303	04.200.038TS	6L28020	04.200.045TS	6L16452
04.200.022TS	6L28018	04.200.026TS	6L39593	04.200.030TS	6L39304	04.200.038TS	4L55488	04.200.045TS	6L29341
04.200.022TS	6L34973	04.200.026TS	6L41697	04.200.030TS	6L39305	04.200.038TS	6L35326	04.200.045TS	6L37248
04.200.022TS	6L37451	04.200.026TS	6L42068	04.200.030TS	6L39306	04.200.038TS	6L37460	04.200.045TS	6L38865
04.200.022TS	6L37452	04.200.026TS	6L42071	04.200.030TS	6L39598	04.200.038TS	6L40695	04.200.045TS	6L42056
04.200.022TS	6L37703	04.200.026TS	6L42055	04.200.030TS	6L39595	04.200.038TS	6L41701	04.200.045TS	5L54921
04.200.022TS	6L38754	04.200.026TS	6L42763	04.200.030TS	6L39596	04.200.038TS	6L42041	04.200.046TS	5L57911
04.200.022TS	6L38755	04.200.026TS	6L42373	04.200.030TS	6L39597	04.200.038TS	6L42042	04.200.046TS	5L60139
04.200.022TS	6L38756	04.200.026TS	6L42374	04.200.030TS	6L41698	04.200.040TS	5L47721	04.200.046TS	6L14175
04.200.022TS	6L39070	04.200.028TS	5L54920	04.200.030TS	6L42037	04.200.040TS	5L49945	04.200.046TS	6L26377
04.200.022TS	6L39295	04.200.028TS	5L56135	04.200.030TS	6L42038	04.200.040TS	5L52134	04.200.046TS	6L35328
04.200.024TS	5L54919	04.200.028TS	6L14168	04.200.032TS	6L42375	04.200.040TS	6L12157	04.200.046TS	6L35628
04.200.024TS	6L14165	04.200.028TS	6L12155	04.200.030TS	6L42376	04.200.040TS	6L12158	04.200.046TS	6L37249
04.200.024TS	6L14166	04.200.028TS	6L14169	04.200.030TS	6L43782	04.200.040TS	6L28021	04.200.046TS	6L37468
04.200.024TS	6L14164	04.200.028TS	6L16446	04.200.030TS	6L42768	04.200.040TS	6L28318	04.200.046TS	6L37470
04.200.024TS	6L16441	04.200.028TS	6L18236	04.200.032TS	5L53667	04.200.040TS	6L36739	04.200.046TS	6L37469
04.200.024TS	6L18228	04.200.028TS	6L21022	04.200.032TS	5L56140	04.200.040TS	6L36742	04.200.046TS	6L39309
04.200.024TS	6L18229	04.200.028TS	6L20121	04.200.032TS	5L57909	04.200.040TS	6L37461	04.200.046TS	6L41704
04.200.024TS	6L18230	04.200.028TS	6L20122	04.200.032TS	5L56139	04.200.040TS	6L38765	04.200.046TS	6L42045
04.200.024TS	6L27611	04.200.028TS	6L26374	04.200.032TS	6L16448	04.200.040TS	6L38768	04.200.046TS	6L42046
04.200.024TS	6L27612	04.200.028TS	6L34976	04.200.032TS	6L24939	04.200.040TS	6L38769	04.200.048TS	6L12159
04.200.024TS	6L29339	04.200.028TS	6L35626	04.200.032TS	6L27614	04.200.040TS	6L41702	04.200.048TS	6L27615
04.200.024TS	6L34974	04.200.028TS	6L37455	04.200.032TS	6L34978	04.200.042TS	5L60137	04.200.048TS	6L26378
04.200.024TS	6L34987	04.200.028TS	6L37454	04.200.032TS	6L36939	04.200.042TS	6L24942	04.200.048TS	6L28023
04.200.024TS	6L35323	04.200.028TS	6L37706	04.200.032TS	6L37457	04.200.042TS	6L24943	04.200.048TS	6L35329
04.200.024TS	6L37453	04.200.028TS	6L38761	04.200.032TS	6L37458	04.200.042TS	6L26375	04.200.048TS	6L35330
04.200.024TS	6L37704	04.200.028TS	6L38759	04.200.032TS	6L37707	04.200.042TS	6L28022	04.200.048TS	6L37471
04.200.024TS	6L37705	04.200.028TS	6L38760	04.200.032TS	6L38763	04.200.042TS	6L34981	04.200.048TS	6L37708
04.200.024TS	6L38757	04.200.028TS	6L39072	04.200.032TS	6L39307	04.200.042TS	6L36740	04.200.048TS	6L39603
04.200.024TS	6L38758	04.200.028TS	6L39073	04.200.032TS	6L44072	04.200.042TS	6L37463	04.200.048TS	6L39604
04.200.024TS	6L39071	04.200.028TS	6L39074	04.200.034TS	5L47720	04.200.042TS	6L37464	04.200.050TS	5L47722
04.200.024TS	6L39296	04.200.028TS	6L39301	04.200.034TS	5L61655	04.200.042TS	6L37465	04.200.050TS	5L52135
04.200.024TS	6L42762	04.200.028TS	6L39302	04.200.034TS	6L16449	04.200.042TS	6L37466	04.200.050TS	5L60140
04.200.024TS	6L43780	04.200.028TS	6L39300	04.200.034TS	6L28019	04.200.042TS	6L37462	04.200.050TS	6L14176
04.200.024TS	6L43813	04.200.028TS	6L39594	04.200.034TS	4L55486	04.200.042TS	6L38770	04.200.050TS	6L12160

Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.
04.200.050TS	6L16453	04.210.120TS	6L30109	04.211.030TS	4L56098	201.370TS	6L29351	202.210TS	5L54946
04.200.050TS	6L24945	04.210.120TS	6L30709	04.211.030TS	6L44671	201.372TS	6L29352	202.210TS	5L60583
04.200.050TS	6L24946	04.210.120TS	6L31055	04.211.030TS	6L44658	201.374TS	6L31489	202.210TS	5L61676
04.200.050TS	6L34982	04.210.120TS	6L32162	04.211.032TS	4L56099	201.375TS	6L31490	202.210TS	5L61677
04.200.050TS	6L35331	04.210.120TS	6L33567	04.211.032TS	6L44268	201.378TS	6L29353	202.210TS	5L60149
04.200.050TS	6L37472	04.210.120TS	6L40759	04.211.032TS	6L44254	201.756TS	6L20267	202.210TS	5L60150
04.200.050TS	6L37473	04.210.120TS	6L42154	04.211.032TS	6L44596	201.760TS	6L16461	202.210TS	6L14178
04.200.050TS	6L37474	04.210.120TS	6L41844	04.211.032TS	6L44672	201.761TS	5L56141	202.210TS	6L16462
04.200.050TS	6L39076	04.210.122TS	6L29737	04.211.032TS	6L44659	201.761TS	6L33579	202.210TS	6L18456
04.200.050TS	6L39310	04.210.122TS	6L30106	04.211.034TS	5L60147	201.761TS	6L34624	202.210TS	6L18457
04.200.055TS	5L60141	04.210.122TS	6L30497	04.211.034TS	6L37351	201.762TS	6L28390	202.210TS	6L20271
04.200.055TS	6L18240	04.210.122TS	6L30710	04.211.034TS	6L44255	201.762TS	6L29354	202.210TS	6L20278
04.200.055TS	6L20123	04.210.122TS	6L40761	04.211.036TS	5L49948	201.762TS	6L18452	202.210TS	6L21119
04.200.055TS	6L28024	04.210.122TS	6L41845	04.211.036TS	6L21023	201.762TS	6L30712	202.210TS	6L21120
04.200.055TS	6L34086	04.210.122TS	6L42155	04.211.036TS	6L30963	201.762TS	6L31060	202.210TS	6L21121
04.200.055TS	6L35332	04.210.124TS	6L29732	04.211.036TS	6L30977	201.763TS	6L31491	202.210TS	6L21122
04.200.055TS	6L37475	04.210.124TS	6L30110	04.211.036TS	6L37134	201.764TS	6L21157	202.210TS	6L21123
04.200.055TS	6L37476	04.210.124TS	6L30505	04.211.036TS	6L44673	201.764TS	6L20268	202.210TS	6L21124
04.200.055TS	6L37477	04.210.126TS	6L35537	04.211.036TS	6L44660	201.764TS	6L27693	202.210TS	6L21125
04.200.055TS	6L37709	04.210.126TS	6L37131	04.211.038TS	6L39439	201.764TS	6L28391	202.210TS	6L21126
04.200.055TS	6L37710	04.210.128TS	6L42860	04.211.038TS	6L42502	201.764TS	6L29355	202.210TS	6L22738
04.200.055TS	6L38771	04.210.128TS	6L43999	04.211.040TS	5L52727	201.764TS	6L30713	202.210TS	6L27699
04.200.055TS	6L39311	04.211.012TS	6L21155	04.211.040TS	5L52597	201.766TS	5L57913	202.210TS	6L27703
04.200.055TS	6L39605	04.211.012TS	6L30960	04.211.040TS	6L30964	201.766TS	6L34625	202.210TS	6L27706
04.200.055TS	6L39606	04.211.012TS	6L30974	04.211.040TS	6L40753	201.768TS	5L52596	202.210TS	6L27617
04.210.112TS	6L27691	04.211.012TS	6L36828	04.211.040TS	6L44597	201.772TS	5L61673	202.210TS	6L26622
04.210.112TS	6L30105	04.211.012TS	6L44670	04.211.040TS	6L44661	201.772TS	6L30507	202.210TS	6L26623
04.210.112TS	6L35431	04.211.012TS	6L44653	04.211.042TS	5L47723	201.774TS	6L30508	202.210TS	6L26628
04.210.114TS	6L22712	04.211.014TS	6L30961	04.211.042TS	6L37350	201.776TS	5L53673	202.210TS	6L27694
04.210.114TS	6L27692	04.211.014TS	6L30975	04.211.042TS	6L37620	201.776TS	6L29356	202.210TS	6L26620
04.210.114TS	6L30107	04.211.014TS	6L36829	04.211.044TS	6L22713	201.776TS	6L31500	202.210TS	6L26621
04.210.114TS	6L30495	04.211.014TS	6L37132	04.211.044TS	6L28381	201.778TS	6L29357	202.210TS	6L26624
04.210.114TS	6L30707	04.211.014TS	6L44654	04.211.044TS	6L40747	201.780TS	5L61674	202.210TS	6L26625
04.210.114TS	6L31051	04.211.016TS	6L22737	04.211.044TS	6L44598	201.782TS	6L29358	202.210TS	6L26626
04.210.114TS	6L32161	04.211.016TS	6L25154	04.211.046TS	6L22714	201.784TS	6L41814	202.210TS	6L26627
04.210.114TS	6L33566	04.211.016TS	6L29343	04.211.046TS	6L31058	201.786TS	6L34484	202.210TS	6L26629
04.210.116TS	6L28379	04.211.018TS	6L21156	04.211.048TS	6L20124	201.788TS	6L29359	202.210TS	6L26630
04.210.116TS	6L30108	04.211.018TS	6L25155	04.211.048TS	6L28025	201.790TS	6L29360	202.210TS	6L26631
04.210.116TS	6L30496	04.211.018TS	6L28380	04.211.048TS	6L30111	201.886TS	5L49982	202.210TS	6L28055
04.210.116TS	6L30708	04.211.018TS	6L30711	04.211.052TS	6L35432	201.886TS	5L60148	202.210TS	6L28056
04.210.116TS	6L31052	04.211.018TS	6L31056	04.211.054TS	6L22715	201.886TS	5L61675	202.210TS	6L28059
04.210.116TS	6L32163	04.211.018TS	6L44655	04.211.056TS	6L25344	201.886TS	6L20269	202.210TS	6L28392
04.210.116TS	6L40755	04.211.020TS	5L53672	04.211.056TS	6L27616	201.886TS	6L20270	202.210TS	6L18455
04.210.116TS	6L42152	04.211.020TS	6L37619	04.211.056TS	6L40746	201.886TS	6L29361	202.210TS	6L18453
04.210.116TS	6L41842	04.211.022TS	6L44656	201.358TS	6L29345	201.886TS	6L31501	202.210TS	6L18454
04.210.118TS	5L65754	04.211.024TS	5L60582	201.359TS	6L29346	201.886TS	6L31502	202.210TS	6L29362
04.210.118TS	6L37831	04.211.024TS	5L65756	201.360TS	6L20266	201.886TS	6L31503	202.210TS	6L29363
04.210.118TS	6L38920	04.211.024TS	6L16460	201.361TS	5L49949	201.886TS	6L31504	202.210TS	6L29364
04.210.118TS	6L40819	04.211.026TS	6L35433	201.361TS	6L29347	201.886TS	6L31506	202.210TS	6L29365
04.210.118TS	6L42153	04.211.026TS	6L44595	201.362TS	4L56211	201.886TS	6L31508	202.210TS	6L29366
04.210.118TS	6L41843	04.211.026TS	6L44657	201.362TS	4L56101	202.210TS	5L52599	202.210TS	6L29387
04.210.118TS	6L42859	04.211.028TS	6L30962	201.363TS	5L49980	202.210TS	5L53674	202.210TS	6L29390
04.210.118TS	6L42501	04.211.028TS	6L30976	201.363TS	6L29348	202.210TS	5L54930	202.210TS	6L29367
04.210.118TS	6L44208	04.211.028TS	6L41846	201.366TS	6L29349	202.210TS	5L54932	202.210TS	6L29368
04.210.120TS	6L29342	04.211.030TS	5L60146	201.366TS	6L31059	202.210TS	5L54942	202.210TS	6L29369
04.210.120TS	6L29731	04.211.030TS	5L61671	201.368TS	6L29350	202.210TS	5L54944	202.210TS	6L29370

Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.
202.210TS	6L29371	202.874TS	6L28394	206.050TS	6L42059	212.101TS	4L56215	212.101TS	6L34856
202.210TS	6L29373	202.874TS	6L29735	206.055TS	6L33812	212.101TS	4L56103	212.101TS	6L34857
202.210TS	6L29335	202.874TS	6L31064	207.010TS	6L34098	212.101TS	4L56104	212.101TS	6L35104
202.210TS	6L29336	202.876TS	5L53675	207.012TS	6L34099	212.101TS	4L56122	212.101TS	6L35221
202.210TS	6L29337	202.876TS	6L18461	207.012TS	6L38866	212.101TS	4L56123	212.101TS	6L35222
202.210TS	6L30132	202.876TS	6L28395	207.014TS	6L33813	212.101TS	6L29344	212.101TS	6L34858
202.210TS	6L30130	202.876TS	6L29736	207.014TS	6L37797	212.101TS	4L56247	212.101TS	6L34860
202.210TS	6L30133	202.876TS	6L30511	212.101TS	5L47725	212.101TS	6L33830	212.101TS	6L34859
202.210TS	6L30135	202.876TS	6L31065	212.101TS	5L47745	212.101TS	6L33831	212.101TS	6L34862
202.210TS	6L30136	202.878TS	6L18462	212.101TS	5L49985	212.101TS	6L33832	212.101TS	6L35521
202.210TS	6L30137	202.878TS	6L28396	212.101TS	5L49986	212.101TS	6L33833	212.101TS	6L35520
202.210TS	6L30138	202.878TS	6L29374	212.101TS	5L49988	212.101TS	6L33834	212.101TS	6L35522
202.210TS	6L30112	202.878TS	6L31066	212.101TS	5L49940	212.101TS	6L34105	212.101TS	6L35101
202.210TS	6L30097	202.880TS	6L29375	212.101TS	5L52603	212.101TS	6L34106	212.101TS	6L35223
202.210TS	6L30098	202.880TS	6L31067	212.101TS	5L52146	212.101TS	6L34107	212.101TS	6L35224
202.210TS	6L30099	202.882TS	6L18463	212.101TS	5L52147	212.101TS	6L34109	212.101TS	6L35517
202.210TS	6L30100	202.882TS	6L29376	212.101TS	5L52728	212.101TS	6L33816	212.101TS	6L35518
202.210TS	6L30101	202.884TS	6L29377	212.101TS	5L49983	212.101TS	6L33817	212.101TS	6L35519
202.210TS	6L29733	202.884TS	6L33568	212.101TS	5L49987	212.101TS	6L33818	212.101TS	6L35524
202.210TS	6L30509	202.886TS	6L29378	212.101TS	5L52148	212.101TS	6L33820	212.101TS	6L35526
202.210TS	6L30131	202.888TS	6L29379	212.101TS	5L52136	212.101TS	6L33821	212.101TS	6L35527
202.210TS	6L30514	202.892TS	6L29380	212.101TS	5L52142	212.101TS	6L33822	212.101TS	6L35528
202.210TS	6L30515	202.894TS	6L29381	212.101TS	5L53684	212.101TS	6L33823	212.101TS	6L35639
202.210TS	6L30517	202.896TS	6L29382	212.101TS	5L53685	212.101TS	6L33824	212.101TS	6L35640
202.210TS	6L30518	202.900TS	6L29383	212.101TS	5L53677	212.101TS	6L33825	212.101TS	6L35641
202.210TS	6L30096	202.962TS	6L29384	212.101TS	5L53678	212.101TS	6L34380	212.101TS	6L35642
202.210TS	6L31061	202.963TS	6L30113	212.101TS	5L53679	212.101TS	6L34381	212.101TS	6L35643
202.210TS	6L31062	202.964TS	6L29385	212.101TS	5L53680	212.101TS	6L34374	212.101TS	6L35634
202.210TS	6L31063	202.965TS	6L30114	212.101TS	5L53681	212.101TS	6L34377	212.101TS	6L35635
202.210TS	6L30716	206.014TS	6L42057	212.101TS	5L53682	212.101TS	6L34378	212.101TS	6L35636
202.210TS	6L30714	206.016TS	5L53676	212.101TS	5L54927	212.101TS	6L34103	212.101TS	6L35637
202.210TS	6L31437	206.016TS	6L35103	212.101TS	5L55203	212.101TS	6L34104	212.101TS	6L35638
202.210TS	6L31510	206.016TS	6L39342	212.101TS	5L54928	212.101TS	6L34108	212.101TS	6L36945
202.210TS	6L31572	206.016TS	6L42058	212.101TS	5L54949	212.101TS	6L34641	212.101TS	6L36838
202.210TS	6L31453	206.018TS	5L54926	212.101TS	5L54957	212.101TS	6L34637	212.101TS	6L36831
202.210TS	6L33599	206.018TS	6L35102	212.101TS	5L56144	212.101TS	6L34379	212.101TS	6L37250
202.210TS	6L33600	206.018TS	6L39343	212.101TS	5L56145	212.101TS	6L34368	212.101TS	6L36940
202.210TS	6L33601	206.018TS	6L41722	212.101TS	5L57939	212.101TS	6L34369	212.101TS	6L36942
202.210TS	6L32165	206.020TS	6L34363	212.101TS	5L57941	212.101TS	6L34371	212.101TS	6L36947
202.210TS	6L29372	206.020TS	6L41723	212.101TS	5L60152	212.101TS	6L34372	212.101TS	6L37143
202.210TS	6L37251	206.022TS	6L34094	212.101TS	5L60157	212.101TS	6L34373	212.101TS	6L36988
202.210TS	6L39620	206.022TS	6L36612	212.101TS	5L60158	212.101TS	6L34376	212.101TS	6L36968
202.210TS	6L40674	206.026TS	6L34095	212.101TS	5L60159	212.101TS	6L34639	212.101TS	6L37484
202.210TS	6L42072	206.026TS	6L35629	212.101TS	5L61687	212.101TS	6L34640	212.101TS	6L37485
202.210TS	6L42735	206.028TS	6L36830	212.101TS	5L61680	212.101TS	6L34642	212.101TS	6L37486
202.210TS	6L42736	206.028TS	6L44674	212.101TS	5L61690	212.101TS	6L34643	212.101TS	6L37478
202.866TS	6L32166	206.028TS	6L44662	212.101TS	5L61691	212.101TS	6L34861	212.101TS	6L37479
202.868TS	6L20279	206.030TS	5L52600	212.101TS	5L61692	212.101TS	6L34863	212.101TS	6L37480
202.870TS	6L18458	206.032TS	6L37798	212.101TS	5L61693	212.101TS	6L34636	212.101TS	6L37481
202.870TS	6L30510	206.035TS	6L44269	212.101TS	5L60584	212.101TS	6L34638	212.101TS	6L37482
202.872TS	6L18459	206.035TS	6L44256	212.101TS	5L60151	212.101TS	6L34370	212.101TS	6L37794
202.872TS	6L28393	206.035TS	6L44663	212.101TS	5L65788	212.101TS	6L34365	212.101TS	6L37793
202.872TS	6L29734	206.036TS	6L36613	212.101TS	5L64098	212.101TS	6L34852	212.101TS	6L37803
202.874TS	6L18460	206.040TS	6L34096	212.101TS	6L28007	212.101TS	6L34855	212.101TS	6L37800
202.874TS	6L20280	206.040TS	6L35631	212.101TS	6L29394	212.101TS	6L34854	212.101TS	6L37711
202.874TS	6L27700	206.050TS	6L34097	212.101TS	6L29395	212.101TS	6L34853	212.101TS	6L37795

Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.
212.101TS	6L37792	212.101TS	6L41724	212.810TS	6L36625	401.356TS	6L31519	401.876TS	5L47726
212.101TS	6L37804	212.101TS	6L41738	212.810TS	6L36626	401.356TS	6L31520	401.876TS	6L31570
212.101TS	6L37802	212.101TS	6L41771	212.810TS	6L36627	401.356TS	6L31516	401.876TS	6L33594
212.101TS	6L39084	212.101TS	6L41772	212.810TS	6L37145	401.356TS	6L33569	401.877TS	5L47727
212.101TS	6L38774	212.101TS	6L41773	212.810TS	6L37256	401.356TS	6L33570	401.877TS	6L33595
212.101TS	6L38772	212.101TS	6L41774	212.810TS	6L37799	401.356TS	6L33571	401.878TS	6L30978
212.101TS	6L38873	212.101TS	6L41775	212.810TS	6L39348	401.356TS	6L33573	401.878TS	6L30965
212.101TS	6L38874	212.101TS	6L41776	212.810TS	6L39366	401.356TS	6L33574	401.879TS	6L31074
212.101TS	6L38878	212.101TS	6L41786	212.810TS	6L39367	401.356TS	6L33575	401.879TS	6L33596
212.101TS	6L38875	212.101TS	6L41787	212.810TS	6L39368	401.356TS	6L33576	401.881TS	6L44271
212.101TS	6L38876	212.101TS	6L41788	212.810TS	6L41793	401.356TS	6L32167	401.881TS	6L44258
212.101TS	6L39079	212.101TS	6L41761	212.810TS	6L41794	401.356TS	6L32168	401.882TS	6L44073
212.101TS	6L39345	212.101TS	6L41736	212.810TS	6L41795	401.356TS	6L36628	401.882TS	6L42732
212.101TS	6L39082	212.101TS	6L41737	212.810TS	6L42747	401.356TS	6L44270	401.882TS	6L43786
212.101TS	6L39083	212.101TS	6L41739	212.810TS	6L42728	401.356TS	6L44257	401.884TS	6L29747
212.101TS	6L39360	212.101TS	6L42384	212.810TS	6L42729	401.356TS	6L42730	401.884TS	6L30129
212.101TS	6L39362	212.101TS	6L42385	212.810TS	6L42742	401.356TS	6L43784	401.884TS	6L31571
212.101TS	6L39364	212.101TS	6L42386	212.810TS	6L42743	401.356TS	6L43785	401.884TS	6L33597
212.101TS	6L39352	212.101TS	6L42387	212.810TS	6L42745	401.368TS	6L30122	401.886TS	6L44272
212.101TS	6L39353	212.101TS	6L42388	212.810TS	6L42741	401.368TS	6L33572	401.886TS	6L44259
212.101TS	6L39355	212.101TS	6L42739	212.810TS	6L43805	401.757TS	6L31569	401.890TS	6L31075
212.101TS	6L39358	212.101TS	6L42377	212.810TS	6L43807	401.757TS	6L33578	401.890TS	6L33598
212.101TS	6L39312	212.101TS	6L43794	212.810TS	6L42744	401.757TS	6L33593	401.892TS	6L39666
212.101TS	6L39077	212.101TS	6L43795	212.810TS	6L42746	401.758TS	6L33577	401.894TS	6L29723
212.101TS	6L39078	212.101TS	6L43799	212.810TS	6L43801	401.758TS	6L42064	401.896TS	5L47733
212.101TS	6L39080	212.101TS	6L42727	212.810TS	6L43803	401.762TS	5L61688	401.896TS	6L21158
212.101TS	6L39346	212.101TS	6L42769	212.810TS	6L43804	401.762TS	6L16463	401.896TS	6L31076
212.101TS	6L39347	212.101TS	6L42738	212.810TS	6L43806	401.762TS	6L29742	401.896TS	6L42733
212.101TS	6L39361	212.101TS	6L42740	212.810TS	6L43808	401.762TS	6L30512	401.896TS	6L43787
212.101TS	6L39365	212.101TS	6L42720	212.810TS	6L43802	401.762TS	6L41745	401.898TS	5L46317
212.101TS	6L39354	212.101TS	6L43796	212.810TS	6L44301	401.764TS	5L53683	401.898TS	6L31085
212.101TS	6L39356	212.101TS	6L43797	212.810TS	6L44077	401.764TS	6L27701	401.898TS	6L36969
212.101TS	6L39357	212.101TS	6L43798	212.810TS	6L44262	401.764TS	6L29743	401.900TS	6L29722
212.101TS	6L39621	212.101TS	6L44076	212.810TS	6L44263	401.764TS	6L30124	401.900TS	6L42734
212.101TS	6L39622	212.101TS	6L44081	401.356TS	6L29386	401.764TS	6L30513	401.900TS	6L43788
212.101TS	6L39363	212.101TS	6L44665	401.356TS	6L29738	401.764TS	6L30718	402.216TS	5L54933
212.101TS	6L40811	212.810TS	5L52601	401.356TS	6L29739	401.764TS	6L41746	402.216TS	6L18464
212.101TS	6L40812	212.810TS	4L56126	401.356TS	6L29740	401.766TS	5L54929	402.216TS	6L21159
212.101TS	6L40679	212.810TS	6L29725	401.356TS	6L29741	401.766TS	6L27702	402.216TS	6L27704
212.101TS	6L40814	212.810TS	6L31527	401.356TS	6L30119	401.766TS	6L29744	402.216TS	6L28057
212.101TS	6L40815	212.810TS	6L31436	401.356TS	6L30120	401.766TS	6L30125	402.216TS	6L28397
212.101TS	6L40816	212.810TS	6L33835	401.356TS	6L30121	401.766TS	6L41747	402.216TS	6L29388
212.101TS	6L40817	212.810TS	6L33836	401.356TS	6L30115	401.768TS	5L61689	402.216TS	6L30134
212.101TS	6L40818	212.810TS	6L34110	401.356TS	6L30116	401.768TS	6L29745	402.216TS	6L30516
212.101TS	6L40675	212.810TS	6L34382	401.356TS	6L30117	401.768TS	6L30126	402.216TS	6L30715
212.101TS	6L40671	212.810TS	6L34627	401.356TS	6L30118	401.768TS	6L33580	402.216TS	6L31086
212.101TS	6L40673	212.810TS	6L35225	401.356TS	6L30123	401.768TS	6L41754	402.218TS	5L54934
212.101TS	6L41770	212.810TS	6L35226	401.356TS	6L31068	401.770TS	6L29746	402.218TS	6L27705
212.101TS	6L41777	212.810TS	6L35644	401.356TS	6L31069	401.770TS	6L30127	402.218TS	6L28058
212.101TS	6L40678	212.810TS	6L35645	401.356TS	6L31073	401.770TS	6L33588	402.218TS	6L28406
212.101TS	6L41789	212.810TS	6L35646	401.356TS	6L31070	401.772TS	6L30128	402.218TS	6L29389
212.101TS	6L41790	212.810TS	6L35647	401.356TS	6L31071	401.772TS	6L33589	402.866TS	6L27687
212.101TS	6L41791	212.810TS	6L35648	401.356TS	6L31072	401.776TS	6L36826	402.868TS	6L12426
212.101TS	6L41792	212.810TS	6L35649	401.356TS	6L31517	401.778TS	6L33591	402.868TS	6L42065
212.101TS	6L41760	212.810TS	6L35650	401.356TS	6L31521	401.780TS	6L33592	402.870TS	6L27708
212.101TS	6L41744	212.810TS	6L35651	401.356TS	6L31518	401.780TS	6L42731	402.870TS	6L33602

Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.	Part No.	Lot No.
402.872TS	6L34628	402.962TS	5L47736	406.045TS	6L44274	407.050TS	6L38872
402.874TS	5L65776	402.962TS	6L24947	406.045TS	6L44261	407.050TS	6L42063
402.874TS	6L16464	402.962TS	6L39081	406.050TS	6L36837	407.055TS	6L34367
402.874TS	6L18465	402.963TS	5L56142	406.050TS	6L37141	407.055TS	6L37801
402.874TS	6L20281	402.963TS	6L27618	406.055TS	6L36973	407.055TS	6L37796
402.874TS	6L27709	402.963TS	6L27716	407.010TS	6L33826		
402.874TS	6L28407	402.963TS	6L31574	407.010TS	6L36975		
402.874TS	6L29391	402.963TS	6L33615	407.012TS	6L34713		
402.874TS	6L30717	402.963TS	6L37253	407.012TS	6L42382		
402.874TS	6L31573	402.964TS	6L33616	407.012TS	6L44075		
402.874TS	6L31087	402.965TS	6L27619	407.014TS	6L34631		
402.874TS	6L28060	402.965TS	6L31088	407.016TS	6L34632		
402.874TS	6L32169	402.965TS	6L40672	407.016TS	6L36622		
402.874TS	6L33603	402.966TS	5L47737	407.018TS	5L47739		
402.876TS	5L47734	402.966TS	6L27620	407.018TS	6L33827		
402.876TS	6L28061	402.966TS	6L27717	407.018TS	6L36976		
402.876TS	6L30139	402.966TS	6L36971	407.018TS	6L42726		
402.876TS	6L30519	402.967TS	6L28319	407.020TS	6L34633		
402.878TS	6L27710	402.967TS	6L37254	407.020TS	6L36623		
402.878TS	6L30140	402.968TS	5L56143	407.022TS	6L33828		
402.878TS	6L33604	402.968TS	6L24948	407.022TS	6L34626		
402.880TS	6L16465	402.968TS	6L29396	407.022TS	6L36986		
402.880TS	6L36836	402.968TS	6L36629	407.024TS	6L34100		
402.882TS	5L47735	402.968TS	6L36972	407.024TS	6L34385		
402.882TS	6L27711	406.010TS	6L37255	407.024TS	6L38867		
402.882TS	6L30141	406.012TS	5L47738	407.026TS	5L47740		
402.884TS	6L27712	406.012TS	6L35540	407.026TS	6L33829		
402.884TS	6L30521	406.012TS	6L37137	407.026TS	6L33814		
402.886TS	6L27713	406.014TS	6L44666	407.026TS	6L36987		
402.886TS	6L33613	406.016TS	5L54945	407.028TS	5L52143		
402.888TS	6L27714	406.016TS	6L39350	407.028TS	6L33815		
402.888TS	6L42378	406.016TS	6L40676	407.028TS	6L34634		
402.890TS	6L27715	406.016TS	6L41755	407.030TS	5L52144		
402.890TS	6L33614	406.018TS	5L54947	407.030TS	5L61679		
402.892TS	6L29392	406.018TS	6L34629	407.030TS	6L34101		
402.892TS	6L36743	406.018TS	6L39351	407.030TS	6L34386		
402.892TS	6L38773	406.018TS	6L40677	407.030TS	6L38868		
402.892TS	6L44514	406.018TS	6L41756	407.030TS	6L42383		
402.894TS	6L29393	406.020TS	5L54948	407.035TS	6L34635		
402.894TS	6L37252	406.020TS	6L34630	407.035TS	6L35632		
402.894TS	6L43789	406.020TS	6L39359	407.035TS	6L42060		
402.896TS	6L22740	406.020TS	6L41757	407.040TS	5L47743		
402.896TS	4L55529	406.022TS	6L44260	407.040TS	6L35633		
402.896TS	6L30142	406.022TS	6L44273	407.040TS	6L42061		
402.896TS	6L42379	406.024TS	6L39505	407.040TS	6L44667		
402.896TS	6L42737	406.026TS	5L60585	407.045TS	6L34366		
402.896TS	6L43790	406.026TS	5L60578	407.045TS	6L35106		
402.898TS	6L22742	406.026TS	6L35541	407.045TS	6L36624		
402.898TS	6L36970	406.026TS	6L37139	407.045TS	6L38869		
402.898TS	6L42380	406.028TS	6L37483	407.045TS	6L39344		
402.898TS	6L43791	406.030TS	5L46327	407.045TS	6L41759		
402.900TS	6L22743	406.032TS	6L35525	407.045TS	6L42062		
402.900TS	4L55530	406.032TS	6L41758	407.050TS	5L47744		
402.900TS	6L36744	406.035TS	6L43793	407.050TS	6L34102		
402.900TS	6L42381	406.045TS	5L49984	407.050TS	6L34364		
402.900TS	6L43792	406.045TS	6L35542	407.050TS	6L35105		