

November 29, 2017

To: Risk Managers
cc: Chairman Medical Board and relevant Head of Departments

Subject: **URGENT MEDICAL DEVICE RECALL**

Affected Product: See Attachment 2 – Affected Product List




Figure 1: Knee Femoral Component, Hip Multipolar Shell and Hip CPT Stem with LDPE bag residue on implant

Zimmer Biomet is conducting a medical device recall for various highly polished hip and knee implants that were packaged in a low density polyethylene (LDPE) bag. The LDPE bag may become adhered to the highly polished surface, leaving residue or material from the LDPE bag on the implant after it is removed from the bag.

Beginning in July 2012 and extending into 2013, Zimmer Biomet implemented a phased change to package the highly polished implants with a new LDPE bag. Testing has shown that this new bag resolves the issue and prevents the bag from adhering to the implants. A similar recall was initiated in January 2016 to remove remaining inventory packaged in the old bag. Further evaluation identified additional lots packaged in the old bag that were not included in the January 2016 recall; therefore, resulting in this new recall for the additional lots packaged in the old bag.

Zimmer Biomet has performed an extensive evaluation of the potential risks associated with this type of event and has concluded that it is unlikely that adhesion of the LDPE bag would cause an adverse effect to either the patient or function of the implant. This conclusion was based on the following:

- The LDPE material is biocompatible, similar to ultra-high molecular weight polyethylene (UHMWPE).
- LDPE is softer than the two mating materials of UHMWPE and cobalt chrome molybdenum (CoCrMo), and therefore, is not expected to scratch either wear surface, which could increase wear rates and possibly lead to osteolysis
- If there were wear particulates generated from the adhered film/residue of a LDPE bag they would be expected to elicit similar biologic reactions as those from UHMWPE and are unlikely to increase the likelihood of peri-prosthetic osteolysis

Risks		
Describe immediate health consequences (injuries or illness) that may result from use of or exposure to the product issue.	Most Probable	Highest Severity
	<i>Extended surgery time less than 30 minutes</i>	<i>Extended surgery time greater than 30 minutes</i>
Describe long range health consequences (injuries or illness) that may result from use of or exposure to the product issue.	Most Probable	Highest Severity
	<i>No long range health consequences are expected.</i>	<i>No long range health consequences are expected.</i>

Our records indicate that you may have received one or more of the affected products. The affected units were distributed between May 2011 and September 2017.

Risk Manager Responsibilities:

1. Review this notification and ensure that affected personnel are aware of the contents.
2. If you have affected product at your facility, assist your Zimmer Biomet sales representative and quarantine all affected product. Your Zimmer Biomet sales representative will remove the affected product from your facility.
3. Complete **Attachment 1 – Certificate of Acknowledgement** and send to CorporateQuality.PostMarket@zimmerbiomet.com. This form must be returned even if you do not have affected products at your facility.
4. Retain a copy of the acknowledgement form with your recall records in the event of a compliance audit of your facility's documentation.
5. If you have further questions or concerns after reviewing this notice, please call customer service at 574-371-3071 between 8:00 am and 5:00pm EST, Monday through Friday. Calls received outside of call center operating hours will receive a voicemail prompt or be transferred to an on-call representative in the event of an emergency. Alternatively, your questions may be emailed to CorporateQuality.PostMarket@zimmerbiomet.com.

Other Information

This medical device recall was reported to the U.S. Food and Drug Administration and will be reported to other Competent Authorities, Notified Bodies, and Regulatory Authorities as required.

- Med Watch Reporting: Adverse reactions or quality problems experienced with the use of this product may be reported to the FDA's Med Watch Adverse Event Reporting program either online, by mail, or by fax.
- Online: www.fda.gov/medwatch/report.htm
- Mail: Use postage paid, pre-addressed form FDA 3500, available at: www.fda.gov/MedWatch/getforms.htm
- Fax: 1-800-FDA-0178


Under 21 CFR 803, manufacturers are also required to report any serious injuries where a product has contributed or may have contributed to the event. Please keep Zimmer Biomet informed of any adverse events associated with this product or any other Zimmer Biomet product by emailing Zimmer.PER@zimmerbiomet.com.

Please be aware that the names of user facilities notified are routinely provided to the Competent Authorities for audit purposes. This action is being taken with the knowledge of Health Canada and is in compliance with regulations set forth by Health Canada. Your urgent cooperation is needed.

The undersigned confirms that this notice has been delivered to the appropriate Regulatory Agencies.

Thank you for your assistance. We regret any inconvenience caused by this recall.

Sincerely,



Kevin W. Escapule
Post Market Surveillance & Regulatory Compliance Director



ATTACHMENT 1 Certificate of Acknowledgement

IMMEDIATE RESPONSE REQUIRED – TIME SENSITIVE ACTION NEEDED

Affected Product: Hip & Knee Implants packaged in LDPE bag

Field Action Reference: ZFA 2017-510

<p>Do you have affected product in your facility?</p> <p><input type="checkbox"/> Yes, we currently have one or more affected items in our facility.</p> <p><input type="checkbox"/> No, we currently have no affected items in our facility.</p>
--

By signing below, I acknowledge that the required actions have been taken in accordance with this recall notice.

Printed Name: _____ **Signature:** _____

Title: _____ **Telephone:** () _____ - _____ **Date:** ____/____/____

Facility Name: _____

Facility Address: _____

City: _____ **State:** _____ **ZIP:** _____

Note: This form must be returned to Zimmer Biomet before this action is closed for your account. It is important that you complete this form and email a copy to CorporateQuality.PostMarket@zimmerbiomet.com or fax to 574-372-4265.

ATTACHMENT 2 - Affected Product List

Affected Hip Products

BIPOLAR SHELL 48MM OD	VERSYS 10 INCH BEADED FC REV 13.5X250MM BWD LT	VERSYS 6 INCH BEADED FC 15X160MM STD BODY STD NECK	VERSYS 8 INCH BEADED FC 13X200MM STR STD BODY/NECK
Item Number: 00-5001-048-00	Item Number: 00-7843-013-11	Item Number: 00-7843-015-06	Item Number: 00-7843-013-08
Lot Number: 62202263	Lot Numbers: 61879448 61768680 77000811	Lot Numbers: 11004305 61729918	Lot Number: 62184774
CPT 12/14 COCR REVISION SIZE 2 180MM	VERSYS 10 INCH BEADED FC REV 13.5X250MM BWD RT	VERSYS 6 INCH BEADED FC STEM 12X160MM LM	VERSYS 8 INCH BEADED FC 14X200MM STR STD BODY/NECK
Item Number: 00-8114-002-18	Item Number: 00-7843-013-12	Item Number: 00-7843-012-36	Item Number: 00-7843-014-08
Lot Number: 61750651	Lot Numbers: 62029061 61852818 62116481 61937672 62143765	Lot Number: 62132189	Lot Numbers: 62184777 62158847 62184778
CPT 12/14 COCR SIZE 0 STD	VERSYS 10 INCH BEADED FC REV 15.0X250MM BWD LT	VERSYS 6 INCH BEADED FC STEM 13X160MM LM	VERSYS 8 INCH BEADED FC REV 13.5X200MM BWD LT
Item Number: 00-8114-000-00	Item Number: 00-7843-015-11	Item Number: 00-7843-013-36	Item Number: 00-7843-013-81
Lot Numbers: 61779346 61783708	Lot Numbers: 62097294 11006832 62116495 61981821	Lot Numbers: 62085292 62120673	Lot Numbers: 61987954 61815082 62116484 61956726
CPT 12/14 COCR SIZE 1 EXT	VERSYS 10 INCH BEADED FC REV 15.0X250MM BWD RT	VERSYS 6 INCH BEADED FC STEM 15X160MM LM	VERSYS 8 INCH BEADED FC 15X200MM STR STD BODY/NECK
Item Number: 00-8114-001-10	Item Number: 00-7843-015-12	Item Number: 00-7843-015-36	Item Number: 00-7843-015-08
Lot Number: 61767833	Lot Numbers: 62085294 11004177 62097244 61759259 62133149 61807097 77000821 61963730	Lot Numbers: 61695988 62089503	Lot Number: 62158851
CPT 12/14 COCR SIZE 2 EXT	VERSYS 6 INCH BEADED FC 12X160MM LM BODY EXT NECK	VERSYS 6 INCH BEADED FC STEM 17X160MM LM	VERSYS 8 INCH BEADED FC REV 13.5X200MM BWD RT
Item Number: 00-8114-002-10	Item Number: 00-7843-012-56	Item Number: 00-7843-017-36	Item Number: 00-7843-013-82
Lot Numbers: 62250377 61744467 62250378 61762055 62257221 61783714 62257223 62249024 62261916 62249025 62261917	Lot Numbers: 62089495 11007679 62097295	Lot Number: 62097290	Lot Numbers: 61937673 61759258 61951276 61815083 62120674 61829373 77000819
CPT 12/14 COCR SIZE 2 STD	VERSYS 6 INCH BEADED FC 13X160MM LM BODY EXT NECK	VERSYS 7.5 IN BEADED FC +10 CALCAR 13.5X180MM STR	VERSYS 8 INCH BEADED FC REV 15.0X200MM BWD LT
Item Number: 00-8114-002-00	Item Number: 00-7843-013-56	Item Number: 00-7843-013-07	Item Number: 00-7843-015-81
Lot Numbers: 62326105 61779347 62327662 62240744 62327663 62240745 62330760 62250379 62333960 62256011 62340054 62264280 62344682 62268572 62370963 62273048 62370965 62273049 62370966 62275085 62298278 62276441 62370975 62277490 62298279 62370968 62303495 62370969 62317907 62278684 62320168 62284463 62320169 62370971 62370976 62370972 62370978 62291922 62370979 62370973 62370990 62326104	Lot Numbers: 61993635 61987936 62011961	Lot Numbers: 61946025 61735805 61993658 61768679 62029060 61815080 62043103 61879441 77000599	Lot Numbers: 61968101 61781981 61968103
	VERSYS 6 INCH BEADED FC 13X160MM STD BODY STD NECK	VERSYS 7.5 IN BEADED FC +10 CALCAR 15.0X180MM STR	VERSYS 8 INCH BEADED FC REV 15.0X200MM BWD RT
	Item Number: 00-7843-013-06	Item Number: 00-7843-015-07	Item Number: 00-7843-015-82
	Lot Numbers: 61721714 62097225	Lot Number: 62135320	Lot Numbers: 61834826 61932605
	VERSYS 6 INCH BEADED FC 14X160MM STD BODY EXT NECK	VERSYS 7.5 IN BEADED FC +20 CALCAR 16.5X170MM STR	VERSYS 8 INCH BEADED FC REV 16.5X200MM BWD LT
	Item Number: 00-7843-014-26	Item Number: 00-7843-016-17	Item Number: 00-7843-016-81
	Lot Number: 61741806	Lot Numbers: 61879439 61993648	Lot Numbers: 61735809 61759260
	VERSYS 6 INCH BEADED FC 14X160MM STD BODY STD NECK	VERSYS 7.5 INCH BEADED FC REV 13.5X190MM STRAIGHT	VERSYS 8 INCH BEADED FC REV 16.5X200MM BWD RT
	Item Number: 00-7843-014-06	Item Number: 00-7843-013-18	Item Number: 00-7843-016-82
	Lot Numbers: 61886028 61752736 62120676	Lot Number: 62158840	Lot Numbers: 62029071 61775246 62029071N 61981820 77000820
		VERSYS 8 INCH BEADED FC 12X200MM STR STD BODY/NECK	
		Item Number: 00-7843-012-08	
		Lot Numbers: 62053685 62116471	

Affected Knee Products

CR POR FEM HATCP SIZE E/L	CR-FLEX GSF PRECOAT FEMORAL, SIZE D-RT	CR-FLEX GSF PRECOAT FEMORAL, SIZE E-LT	CR-FLEX GSF PRECOAT FEMORAL, SIZE E-RT	
Item Number: 65-5972-015-01 Lot Number: 62086845	Item Number: 00-5750-014-02 Lot Numbers: 61970681 61843131 62072514 61862167 62088280 61870724 61978269 61890017 62093271 61905264 62132082 61914484 62134707 62059374 62159269 61951689 62170344	Item Number: 00-5750-015-01 Lot Numbers: 62155453 62002081 62172750 62002082	Item Number: 00-5750-015-02 Lot Numbers: 62020823 62166658 62173170	
CR POR FEM HATCP SIZE F/L		CR-FLEX GSF PRECOAT FEMORAL, SIZE E-LT MINUS	CR-FLEX GSF PRECOAT FEMORAL, SIZE E-RT MINUS	
Item Number: 65-5972-016-01 Lot Numbers: 11007043 62105229		Item Number: 00-5750-015-05 Lot Numbers: 62122660 61829059 62133157 61979151 62170401 62003579 62075016 61829059 62122660 61979151 62133157 62003579 62170401 62075016	Item Number: 00-5750-015-06 Lot Numbers: 62060528 61797739 62091562 61830366 62122538 62003582 62141936	
CR-FLEX GSF PRECOAT FEMORAL, SIZE C-LT	CR-FLEX GSF PRECOAT FEMORAL, SIZE D-RT MINUS		CR-FLEX GSF PRECOAT FEMORAL, SIZE F-LT	
Item Number: 00-5750-013-01 Lot Numbers: 62003770 61951107 62088332 61990633 62168463	Item Number: 00-5750-014-06 Lot Numbers: 62135959 62121229 62170385		Item Number: 00-5750-016-01 Lot Numbers: 61996959 61794718 62005925 61796369 62007016 61804735 62013511 61823887 62022324 61834861 62022325 61843134 62034567 61879330 62044626 61886333 62063739 61890709 62072219 61905190 62080345 61905191 62092539 61914487 62093019 61929422 62112670 61936004 62119497 61946957 62125812 61955579 62137994 61962742 62144301 61962744 62152969 61976418 62155068	
CR-FLEX GSF PRECOAT FEMORAL, SIZE C-RT	CR-FLEX GSF PRECOAT FEMORAL, SIZE E-LT	CR-FLEX GSF PRECOAT FEMORAL, SIZE E-RT		
Item Number: 00-5750-013-02 Lot Numbers: 62093265 61970680 62166656 62004806 62093265 61970680 62166656 62004806	Item Number: 00-5750-015-01 Lot Numbers: 62018053 11400226 62018142 61797750 62028178 61804265 62029873 61809715 62029875 61810409 62035840 61819936 62042331 61820642 62049388 61829061 62055039 61834876 62060477 61842352 62064338 61847666 62086643 62011329 62172749 61863943 62086644 61870821 62091506 61878293 62091507 61889443 62106697 61895335 62116013 61902970 62121556 61905209 62121558 61908315 62132523 61909358 62132525 61916173 61916174 61921446 62140119 61931666 62148243 61941604 62152391 61952956 62170954 61958754 62170955 61963888 62170956 61966239 62170957 61966240 62170958 61979149 62172733 61986108 62172737 61996141 62172742	Item Number: 00-5750-015-02 Lot Numbers: 61876301 61794716 61879242 61810412 61885147 61811776 61886280 61820645 61890043 61824008 61890706 61834877 61904206 61838758 61910227 61847670 61910229 61872331 61862132 61910889 61996356 61918992 61996357 61923590 62011576 61936794 62011577 61950709 62011578 61950710 62011579 61960710 62060528 61964117 62071255 61970682 62075660 61973267 62079865 61986430 62087941 61986739 62092277 62094385 62022109 62106698 62028437 62107126 62033069 62117325 62036622 62121774 62043322 62132865 62050766 62132878 62166748 62141732 62171101 62141927 62171102 62147818 62171757 62152394 62171762 62156825 62171764		
CR-FLEX GSF PRECOAT FEMORAL, SIZE C-RT MINUS			CR-FLEX GSF PRECOAT FEMORAL, SIZE F-LT	
Item Number: 00-5750-013-06 Lot Number: 62120693 11400140 62004808 11400196 62029480 61827938 62072208			Item Number: 00-5750-016-01 Lot Numbers: 61982466 62171793	
CR-FLEX GSF PRECOAT FEMORAL, SIZE D-LT			CR-FLEX GSF PRECOAT FEMORAL, SIZE F-LT MINUS	
Item Number: 00-5750-014-01 Lot Numbers: 62087759 61844939 62093268 61879304 62121194 61905173 62134274 61929402 62159864 61951674 62170341 61964632 62170342 61972167 62170343			Item Number: 00-5750-016-05 Lot Numbers: 62093021 11006472 62100379 61836353 62126817 61970691	
CR-FLEX GSF PRECOAT FEMORAL, SIZE D-LT MINUS			CR-FLEX GSF PRECOAT FEMORAL, SIZE F-RT	
Item Number: 00-5750-014-05 Lot Numbers: 62121555 62166655			Item Number: 00-5750-016-02 Lot Numbers: 61968824 11400202 61968857 11400318 61976430 61796371 61988466	
CR-FLEX GSF PRECOAT FEMORAL, SIZE D-RT				
Item Number: 00-5750-014-02 Lot Numbers: 62022989 11400111 62011560 61791935 62170384 61791936 62171084 61813666 61990915				

CR-FLEX GSF PRECOAT FEMORAL, SIZE F-RT <u>Item Number: 00-5750-016-02</u> <u>Lot Numbers:</u> 61798627 61989668 61809722 61997351 61819934 62005963 62014079 61823891 62014080 61823892 62022326 61829840 62022327 61844041 62038034 61846258 62046244 61867611 62064361 61876284 62073119 61888046 62080827 61893643 62088862 61900478 62094391 61905269 62106555 61905271 62112673 61914490 62130983 61922259 62138093 61936008 62146939 61936677 62154050 61946960 62163077 61959246 62171794 61959248	CR-FLEX GSF PRECOAT FEMORAL, SIZE G-RT MINUS <u>Item Number: 00-5750-017-06</u> <u>Lot Numbers:</u> 62134277 62173183	CR-FLEX POR HATCP FEM C-L <u>Item Number: 65-5952-013-01</u> <u>Lot Number:</u> 62020175	CR-FLEX POR HATCP FEM C-R <u>Item Number: 65-5952-013-02</u> <u>Lot Numbers:</u> 62038825 62110571	CR-FLEX POR HATCP FEM D-L <u>Item Number: 65-5952-014-01</u> <u>Lot Numbers:</u> 61866624 11004583 61987911 61782590 62020176 61815515 62046099 61846120 62068159 61855468 62117121	CR-FLEX POR HATCP FEM D-R <u>Item Number: 65-5952-014-02</u> <u>Lot Numbers:</u> 62013405 61822350 62032575 61848475 62068160 61899497 62110572 61974938 62139310	CR-FLEX POR HATCP FEM E-L <u>Item Number: 65-5952-015-01</u> <u>Lot Numbers:</u> 61936668 11004349 61949553 11004533 61994913 11005729 62005727 11006459 62025668 11006641 62025669 61801431 62054539 61822351 62074885 61831562 62092422 61855469 62117122 61861694 62139311 61910510	CR-FLEX POR HATCP FEM E-R <u>Item Number: 65-5952-015-02</u> <u>Lot Numbers:</u> 11007056 11004148 61789012 11004348 61968653 61801432 62013406 61822358 62025670 61831563 62025671 61836562 62062054 61861695 62080588 61910511 62105374 61936669 62110573	CR-FLEX POR HATCP FEM F-R <u>Item Number: 65-5952-015-02</u> <u>Lot Numbers:</u> 11004504 61949554 62139312	CR-FLEX POR HATCP FEM F-L <u>Item Number: 65-5952-016-01</u> <u>Lot Numbers:</u> 61915979 11003498 61935649 11004357 61949555 11004435 61951537 11006442 62001383 61777546 62013407 61795935 62020217 61805663 62025673 61822352 62046102 61824508 62054540 61848477 62080589 61857406 62092423 61866625 62117123	CR-FLEX POR HATCP FEM F-R <u>Item Number: 65-5952-016-02</u> <u>Lot Numbers:</u> 61928637 11004351 61939608 11004432 61951538 11004505 62001384 61774128 62013408 61789013 62020218 61815516 62032576 61822353 62062055 61831564 62068161 61848478 62086793 61855470 62105375 61861696 62117124 61899498 62128162	CR-FLEX POR HATCP FEM G-L <u>Item Number: 65-5952-017-01</u> <u>Lot Numbers:</u> 61921289 11004091 61939610 11004205 61949556 11004582 61983388 11005728 61994915 61777548 62013409 61795936 62020181 61815517 62032577 61824509 62054542 61836564 62074887 61857409 62086794 61861697 62110574 61885032 62128163	CR-FLEX POR HATCP FEM G-R <u>Item Number: 65-5952-017-02</u> <u>Lot Numbers:</u> 61939611 11004506 61951540 11004666 61987913	CR-FLEX POR HATCP FEM G-R <u>Item Number: 65-5952-017-02</u> <u>Lot Numbers:</u> 62005728 11006037 62020179 61789014 62025674 61805664 62056871 61822360 62068162 61836565 62086795 61857410 62104152 61866626 62122759 61904113 62132130 61935650	CR-FLX POR HATCP FEM C-L MINUS <u>Item Number: 65-5952-013-05</u> <u>Lot Numbers:</u> 61819644 11004433 61855471	CR-FLX POR HATCP FEM C-R MINUS <u>Item Number: 65-5952-013-06</u> <u>Lot Number:</u> 62056868	CR-FLX POR HATCP FEM D-L MINUS <u>Item Number: 65-5952-014-05</u> <u>Lot Numbers:</u> 61987912 11003623 62104150 61782591	CR-FLX POR HATCP FEM D-R MINUS <u>Item Number: 65-5952-014-06</u> <u>Lot Numbers:</u> 62038826 61835232 62132128	CR-FLX POR HATCP FEM E-L MINUS <u>Item Number: 65-5952-015-05</u> <u>Lot Numbers:</u> 62020177 11004436 62046101 61782592 61857405 61822359 62122758 61939607	CR-FLX POR HATCP FEM E-R MINUS <u>Item Number: 65-5952-015-06</u> <u>Lot Numbers:</u> 61968654 11003505 62001382 11004261 62025672 61805662 62074886 61848476 62128161 61921288	CR-FLX POR HATCP FEM F-L MINUS <u>Item Number: 65-5952-016-05</u> <u>Lot Numbers:</u> 61928636 11004258 61951539 11004438 62020178 61801433 62056869 61836563 62092424 61857407 62139313
--	--	---	---	--	---	--	--	--	---	---	---	---	--	---	--	--	---	--	---	---

CR-FLX POR HATCP FEM F-R MINUS	LPS POROUS SRF HDN FEM SZ E-LT	LPS-FLEX GSF OPTION FEMORAL, SIZE D-LT	LPS-FLEX GSF OPTION FEMORAL, SIZE D-RT
Item Number: 65-5952-016-06 Lot Numbers: 61939609 11004434 62002454 61777547 62020219 61822354 62062056 61857408 62104151 61904112	Item Number: 00-5992-015-31 Lot Number: 61933109 LPS PRECOAT SRF HDN FEM SZ D-L Item Number: 00-5990-014-31 Lot Number: 61938548 LPS PRECOAT SRF HDN FEM SZ D-R Item Number: 00-5990-014-32 Lot Number: 61967973 LPS PRECOAT SRF HDN FEM SZ E-L Item Number: 00-5990-015-31 Lot Number: 61949343	Item Number: 00-5764-014-51 Lot Numbers: 62023806 61868887 62096024 61881032 62100482 61881034 62100483 61885761 62101706 61888051 62106707 61892383 62107956 61896058 62110811 61896059 62119557 61897301 62119558 61902968 62120493 61910188 62131261 61914468 62132092 61918704 62143628 61918705 62151117 61924472 62151121 61926153 62151252 61936817 62152410 61943203 62152412 61946974 62162780 61958271 62169040 61958748 62169058 61965573 62169059 61965574 62169060 61967192 62169061 61967194 62170076 61974066 62170079 61974068 62170080 61979793 62170082 61984251 62173155 61985016 62173156 61986228 62173167 61986231 62173168 62017272 62173193 62023333 62173194	Item Number: 00-5764-014-52 Lot Numbers: 62019511 61869869 62098784 61874487 62098834 61881657 62099264 61881659 62099786 61881660 62102049 61892387 62107957 61893247 62107995 61894545 62107999 61899648 62114385 61901228 62115960 61902966 62116005 61907467 62121784 61915329 62132886 61921477 62132887 61926156 62133207 61932311 62133209 61946459 62146949 61956807 62154048 61956808 62160342 61958269 62160892 61964453 62165941 61965505 62170397 61972316 62170948 61972317 62170951 61972754 62170952 61976340 62170953 61976342 62171096 61980592 62171097 61984199 62171098 61992067 62171099 61992069 62171584 62018156 62174422
CR-FLX POR HATCP FEM G-L MINUS Item Number: 65-5952-017-05 Lot Numbers: 62005729 61824510 62025675 61893908 62086012 61939612 62132131	LPS-FLEX GSF OPTION FEMORAL, SIZE C-LT Item Number: 00-5764-013-51 Lot Numbers: 62022335 61874481 62029870 61892716 62107937 61946510 62121782 61954552 62155687 61977566 62159335 62003780		
CR-FLX POR HATCP FEM G-R MINUS Item Number: 65-5952-017-06 Lot Numbers: 61861698 11003018 61915980 11003671 61987914 11004260 62020182 61795937 62038827 61822355 62105368 61846121	LPS-FLEX GSF OPTION FEMORAL, SIZE C-RT Item Number: 00-5764-013-52 Lot Numbers: 62023325 61800630 62030453 61808793 62096017 61819071 62102077 61885786 62127093 61892717 62137523 61914463 62168360 61949724 62169039 61954981 62172720 61977567 62172721 62003781 62173154		
Gender PFJ FEMORAL COMP, CEMENTED, SIZE 1-LT Item Number: 00-5926-011-01 Lot Number: 61777451 Gender PFJ FEMORAL COMP, CEMENTED, SIZE 2-LT Item Number: 00-5926-012-01 Lot Numbers: 11003321 61777452			
Gender PFJ FEMORAL COMP, CEMENTED, SIZE 3-LT Item Number: 00-5926-013-01 Lot Numbers: 61768086 11002963 61771794			
Gender PFJ FEMORAL COMP, CEMENTED, SIZE 3-RT Item Number: 00-5926-013-02 Lot Numbers: 61758943 61771796	LPS-FLEX GSF OPTION FEMORAL, SIZE D-LT Item Number: 00-5764-014-51 Lot Numbers: 62029858 11400112 62030266 61787880 62030268 61793522 62030269 61796450 62030578 61802094 62038007 61816554 62038008 61817772 62038009 61817773 62039076 61829076 62049410 61829078 62049412 61834863 62071907 61835849 62076334 61842366 62084449 61858497 62085549 61868886 62096020		
Gender PFJ FEMORAL COMP, CEMENTED, SIZE 4-RT Item Number: 00-5926-014-02 Lot Numbers: 62264027 62258127 62265282 62262925 62269962			
Gender PFJ FEMORAL COMP, CEMENTED, SIZE 5-RT Item Number: 00-5926-015-02 Lot Number: 61771799 LPS POROUS SRF HDN FEM SZ D-LT Item Number: 00-5992-014-31 Lot Number: 61933108		LPS-FLEX GSF OPTION FEMORAL, SIZE D-RT Item Number: 00-5764-014-52 Lot Numbers: 62023807 61796361 62024586 61803334 62026104 61804253 62030580 61809580 62032079 61816405 62039171 61819068 62039173 61821435 62040748 61822979 62046929 61832458 62047537 61833619 62048449 61837343 62063766 61837344 62064493 61842369 62069351 61854802 62077281 61859250 62077309 61868889 62085550 61868890 62085552	LPS-FLEX GSF OPTION FEMORAL, SIZE E-LT Item Number: 00-5764-015-51 Lot Numbers: 61819644 11006800 61986731 11400127 61986732 11400159 61990932 11400326 62013572 61778270 62014088 61779535 62015860 61779538 62015863 61780342 62023803 61780343 62023804 61780344 62024965 61793548 62026103 61793550 62026945 61817771 62030613 61819065 62030641 61819067 62039186 61819232 62039946 61821147 62040768 61827947 62044625 61827948 62050091

LPS-FLEX GSF OPTION FEMORAL, SIZE E-LT		LPS-FLEX GSF OPTION FEMORAL, SIZE E-LT		LPS-FLEX GSF OPTION FEMORAL, SIZE E-RT		LPS-FLEX GSF OPTION FEMORAL, SIZE F-LT	
<u>Item Number: 00-5764-015-51</u>		<u>Item Number: 00-5764-015-51</u>		<u>Item Number: 00-5764-015-52</u>		<u>Item Number: 00-5764-016-51</u>	
<u>Lot Numbers:</u>		<u>Lot Numbers:</u>		<u>Lot Numbers:</u>		<u>Lot Numbers:</u>	
61827949	62051061	61984216	62174389	61878281	61977573	61821411	61982477
61829859	62051206	61985033	62174418	61878285	62110820	61829141	62021945
61829861	62051207	61986425	62174419	61878286	62110822	61832279	62030615
61837319	62052132	61986426	62174420	61879238	62111573	61835848	62030617
61839696	62064405			61880111	62112040	61835848	62039992
61839697	62068557	LPS-FLEX GSF OPTION FEMORAL, SIZE E-RT		61881691	62119555	61855948	62050760
61839701	62072217	<u>Item Number: 00-5764-015-52</u>		61885151	62119556	61861376	62063780
61846271	62073121	<u>Lot Numbers:</u>	61980327	61885151	62123093	61866734	62077021
61846271	62080805	11400123	61980331	61889072	62123102	61880092	62097837
61846295	62084279	11400124	61980588	61889073	62123109	61884397	62097838
61855939	62084447	11400323	61980589	61896056	62129369	61889176	62111956
61856755	62092535	11400331	61980589	61897298	62129497	61889178	62112984
61866732	62096018	11400331	61987068	61900472	62129499	61906078	62118647
61866733	62096629	11400407	61987069	61900474	62129501	61906079	62129502
61872345	62099061	61780357	61988203	61900475	62129501	61912350	62129823
61876931	62099278	61780358	62003783	61904939	62129501	61912350	62129823
61880112	62099947	61780359	62003783	61905204	62142625	61912350	62133843
61880113	62099950	61780360	62005010	61910884	62143321	61926158	62141734
61881692	62103221	61780360	62013660	61916102	62147830	61933850	62154053
61883280	62103221	61781294	62014093	61916102	62148253	61948147	62154079
61883281	62107993	61782058	62017638	61916165	62148806	61956821	62157969
61883281	62107994	61782059	62017639	61916166	62149021	61958019	62169062
61883283	62109175	61782060	62018580	61916167	62152332	61971189	62169063
61884198	62114466	61782061	62027210	61917767	62152339	61971190	62170085
61885185	62117288	61785877	62028434	61920457	62160344	61977411	62170087
61885186	62117289	61785878	62028459	61920458	62160908	LPS-FLEX GSF OPTION FEMORAL, SIZE F-RT	
61902963	62121804	61794729	62028460	61927979	62163090	<u>Item Number: 00-5764-016-52</u>	
61902964	62122578	61794731	62033741	61927979	62163092	<u>Lot Numbers:</u>	
61904203	62127094	61806293	62033742	61929419	62163115	11400420	61964451
61906073	62127493	61806294	62036666	61932307	62164681	61793531	61977415
61908309	62128369	61814773	62038027	61935999	62171082	61794710	61977420
61915325	62128382	61817767	62039089	61936001	62171103	61794710	61977569
61915327	62132116	61817768	62040769	61943100	62171104	61798624	61984192
61915328	62132118	61819248	62040770	61947025	62171105	61807823	61990934
61916101	62132869	61819250	62042055	61949720	62171106	61816385	61990935
61918735	62146907	61819270	62042058	61952892	62171795	61816386	61997092
61924478	62146914	61820658	62051073	61954534	62171796	61817749	61997655
61926148	62152632	61820669	62052008	61954551	62171797	61819049	62005011
61933818	62152649	61820732	62052684	61958765	62171798	61832270	62022987
61933846	62160897	61821146	62055603	61960727	62171798	61834999	62030644
61949719	62162211	61822977	62063754	61960730	62171851	61835000	62036758
61950809	62169025	61824793	62064397	61965494	62171852	61857909	62040775
61952891	62169026	61832298	62074510	61965498	62172424	61857910	62052777
61953109	62169027	61832299	62077305	61971185	62172425	61863969	62072221
61958030	62169028	61833582	62077869	61971734	62173197	61878280	62077311
61958031	62169064	61837336	62085543	61971735	62173198	61879235	62085964
61958327	62169065	61845036	62085555	61971736	62174394	61880390	62098850
61958329	62170094	61853393	62085958	61976427	62174395	61890013	62098851
61963907	62170097	61853394	62094393	LPS-FLEX GSF OPTION FEMORAL, SIZE F-LT		61900476	62110243
61963908	62170102	61854830	62096027	<u>Item Number: 00-5764-016-51</u>		61900477	62111640
61964178	62170104	61866727	62097833	<u>Lot Numbers:</u>		61906080	62122673
61971133	62170402	61871541	62097834	61792793	61982478	61907464	62127562
61971181	62170403	61871542	62098073	61792794	61990937	61915334	62134718
61972190	62170404	61876293	62098802	61804255	61995694	61916171	62140669
61972194	62171080	61876294	62098847	61806271	61995696	61927986	62149048
61978305	62173192	61876295	62098848	61814754	61997354	61933852	62154294
61978545	62174388	61876929	62110241		62005967	61939198	62156828

LPS-FLEX GSF OPTION FEMORAL, SIZE F-RT <u>Item Number: 00-5764-016-52</u> <u>Lot Numbers:</u> 61964450 61949725 62163088 61954907 62168669 61955570 62168677 61960732 62168688 61964191 62168693	M/G UNI KNEE PRECOAT FEM, MED-LT/LAT-RT, 47MM A/P <u>Item Number: 00-5790-024-00</u> <u>Lot Numbers:</u> 62025572 61822278 62105235 61928510 62110658 61952565 77001780	M/G UNI KNEE PRECOAT FEM, MED-RT/LAT-LT, 50MM A/P <u>Item Number: 00-5790-027-00</u> <u>Lot Numbers:</u> 62045723 61782626 62068115 61824553 62092386 61921258 77001783 61952568	NATURAL-KNEE GSF NP FLEX FEM, SIZE 4-RT <u>Item Number: 00-5414-017-02</u> <u>Lot Numbers:</u> 62001895 61824684 62005771 61871418 62025844 61876224 62046213 61889363 62054634 61899535 62062145 61904179 62074937 61928675 62104630 61933394 62110727 61950164
LPS-FLEX GSF OPTION FEMORAL, SIZE G-LT <u>Item Number: 00-5764-017-51</u> <u>Lot Numbers:</u> 62034482 61807828 62087938 61835350 62117139 61907463 62126104 61944906 62136387 61967907 62171585 61996353	M/G UNI KNEE PRECOAT FEM, MED-LT/LAT-RT, 50MM A/P <u>Item Number: 00-5790-026-00</u> <u>Lot Numbers:</u> 61952567 11006199 62005615 61777601 62062058 61815430 62086067 61824552 77001782	M/G UNI KNEE PRECOAT FEM, MED-RT/LAT-LT, 52MM A/P <u>Item Number: 00-5790-029-00</u> <u>Lot Numbers:</u> 62005616 61805632 77001785	NATURAL-KNEE GSM POR FLEX FEM, SIZE 2-RT <u>Item Number: 00-5412-015-02</u> <u>Lot Numbers:</u> 61975655 11001440 61983536 11003149 61988106 11004422 61988107 11004477 61995014 11005840 62001897 11006300 62005775 61777756 62013575 61782977 62013586 61796047 62013587 61806562 62020383 61806563 62025800 61867212 62025801 61867213 62032711 61871424 62041160 61872035 62054603 61876226 62056958 61885240 62062122 61889370 62068182 61894436 62074905 61899511 62080694 61904154 62086883 61916029 62104205 61921361 62104646 61928677 62110690 61940207 62117168 61952635 62122805 61963807 62128264 61968948
LPS-FLEX GSF OPTION FEMORAL, SIZE G-RT <u>Item Number: 00-5764-017-52</u> <u>Lot Numbers:</u> 61982476 61806276 61996355 61835845 62067956 61907462 62090171 61917771 62136388 61951114 62156850 61967909 62171591	M/G UNI KNEE PRECOAT FEM, MED-LT/LAT-RT, 56MM A/P <u>Item Number: 00-5790-030-00</u> <u>Lot Numbers:</u> 61968455 61881087 77001786	M/G UNI KNEE PRECOAT FEM, MED-RT/LAT-LT, 56MM A/P <u>Item Number: 00-5790-031-00</u> <u>Lot Numbers:</u> 62002487 77001787	NATURAL-KNEE GSF NP FLEX FEM, SIZE 0-RT <u>Item Number: 00-5414-013-02</u> <u>Lot Numbers:</u> 62001882 11006436 62020365 61768287 62032695 61768287R 62046201 61871391 62054637 61876203 62054637R 61889334 62074928 61910523 62094873 61940029 62103499 61971424 62117187 61983528 62128275 61994980 77001023
LPS-FLEX GSF POROUS FEMORAL, SIZE F-RT <u>Item Number: 00-5762-016-52</u> <u>Lot Number:</u> 11005195	M/G UNI KNEE PRECOAT FEM, MED-LT/LAT-RT, 42MM A/P <u>Item Number: 00-5790-020-00</u> <u>Lot Numbers:</u> 61933229 11004267 61949771 11005144 62062057 61855504 77001776 61915949	M/G UNI KNEE PRECOAT FEM, MED-RT/LAT-LT, 44MM A/P <u>Item Number: 00-5790-021-00</u> <u>Lot Numbers:</u> 61933230 11005163 62013359 11005243 62086803 11006254 77001777	NATURAL-KNEE GSF NP FLEX FEM, SIZE 4-LT <u>Item Number: 00-5414-017-01</u> <u>Lot Numbers:</u> 61963823 11004925 61971475 11006490 61988083 61777741 61994991 61801560 62005770 61837057 62025843 61867209 62041147 61872013 62054638 61885276 62062144 61889362 62074934 61904178 62094875 61904178R 62117197 61910561 77001027 61928674 77001464 61952632
M/G UNI KNEE PRECOAT FEM, MED-LT/LAT-RT, 44MM A/P <u>Item Number: 00-5790-022-00</u> <u>Lot Numbers:</u> 61932359 11004689 61939861 11004850 61958404 11006797 61987879 61788933 61993535 61810168 62041075 61831670 62074833 61894364 77001778 61910487	M/G UNI KNEE PRECOAT FEM, MED-RT/LAT-LT, 47MM A/P <u>Item Number: 00-5790-023-00</u> <u>Lot Numbers:</u> 61928509 11003416 61933231 11003680 61949247 11003843 61993546 61795864 62020300 61810169 62056497 61889207 62074834 61904081 77001779	M/G UNI KNEE PRECOAT FEM, MED-RT/LAT-LT, 47MM A/P <u>Item Number: 00-5790-025-00</u> <u>Lot Numbers:</u> 61952566 11005169 62000808 11006839 62020301 61777599 62056498 61795865 62074835 61822279 62080646 61899451 62104129 61928511 77001781 61936678	NATURAL-KNEE GSM POR FLEX FEM, SIZE 3-LT <u>Item Number: 00-5412-016-01</u> <u>Lot Numbers:</u> 61950171 11003194 61952636 11003207 61952637 11004099 61958443 11004486 61958444 11004488 61963808 11004569 61963809 11004571 61968949 11004855 61968950 11004880 61968951 11004921 61975656
M/G UNI KNEE PRECOAT FEM, MED-LT/LAT-RT, 47MM A/P <u>Item Number: 00-5790-024-00</u> <u>Lot Numbers:</u> 62054459 61788934 62068114 61805631 62080645		NATURAL-KNEE GSF NP FLEX FEM, SIZE 4-RT <u>Item Number: 00-5414-017-02</u> <u>Lot Numbers:</u> 61950165 11005236 61950166 11005311 61952633 11006896 61958472 61774297 61994989	

NATURAL-KNEE GSM POR FLEX FEM, SIZE 3-LT		NATURAL-KNEE GSM POR FLEX FEM, SIZE 3-LT		NATURAL-KNEE GSM POR FLEX FEM, SIZE 4-LT		NATURAL-KNEE GSM POR FLEX FEM, SIZE 4-RT	
<u>Item Number: 00-5412-016-01</u>		<u>Item Number: 00-5412-016-01</u>		<u>Item Number: 00-5412-017-01</u>		<u>Item Number: 00-5412-017-02</u>	
<u>Lot Numbers:</u>	61975657	<u>Lot Numbers:</u>	62092462	<u>Lot Numbers:</u>	62025803	<u>Lot Numbers:</u>	61952643
11005070	61975658	61899549	62092463	61810282	62025804	61782981	61995024
11005071	61975659	61904155	62094870	61815751	62032716	61782982	61995025
11005190	61975660	61904156	62103504	61815752	62032717	61789829	62001907
11005766	61983537	61910532	62104180	61815753	62041168	61789841	62001908
11005970	61983538	61910533	62104181	61822644	62041169	61796056	62001909
11005971	61983550	61910534	62104207	61824652	62041170	61796057	62005782
11006008	61986587	61916030	62104207R	61824653	62046178	61796057R	62005783
11006100	61988108	61921362	62110691	61831774	62054608	61796058	62005784
11006234	61988109	61921363	62110692	61837034	62054609	61801509	62013582
11006309	61991990	61928678	62110699	61837035	62056964	61801510	62013583
11006457	61991994	61928679	62122806	61867303	62056965	61801511	62013584
11006982	61995015	61933415	62128265	61867303R	62056965R	61806572	62013593
11007026	61995016	61933416	62128266	61871431	62062130	61810247	62020394
61777751	61995017	61936720	62128267	61871433	62062131	61810283	62020395
61777757	62001898	61940208	62132476	61872040	62062132	61810284	62025805
61777758	62001899	NATURAL-KNEE GSM POR FLEX FEM, SIZE 4-LT		61876231	62068188	61822645	62025806
61782969	62001900	<u>Item Number: 00-5412-017-01</u>		61885245	62068189	61822646	62032718
61782978	62001901	<u>Lot Numbers:</u>	61936722	61885246	62074911	61824654	62032719
61789818	62005776	11003172	61936723	61885247	62074912	61824655	62032720
61789819	62005777	11003197	61940212	61889376	62080692	61837036	62041171
61789820	62013576	11003218	61940213	61894442	62080695	61837037	62041172
61789821	62013577	11003218	61940213	61894443	62086890	61862204	62046179
61796048	62013588	11004423	61950175	61899513	62092467	61867215	62046180
61796049	62013589	11004611	61952640	61899551	62092468	61871436	62054610
61796050	62020384	11004762	61952641	61904159	62092469	61871437	62054611
61796051	62020385	11004927	61958447	61910535	62104206	61872041	62056966
61806564	62020386	11004963	61958448	61910567	62104647	61872042	62056969
61806565	62020387	11004964	61963812	61916033	62110695	61876232	62062133
61810244	62032712	11005044	61963813	61916034	62110696	61885248	62062134
61810277	62032713	11005301	61968954	61916035	62117175	61885249	62068190
61810278	62041161	11005768	61968955	61921366	62122808	61885250	62074913
61822640	62041162	11006518	61975664	61921367	62128270	61889377	62074914
61824650	62041163	11006585	61975665	61928682		61894444	62080693
61831766	62041164	11006765	61975666	NATURAL-KNEE GSM POR FLEX FEM, SIZE 4-RT		61894445	62086891
61831767	62046173	11006815	61975667	<u>Item Number: 00-5412-017-02</u>		61899514	62086892
61831768	62046174	11007002	61983542	<u>Lot Numbers:</u>	61958449	61899552	62092470
61831769	62054604	61777761	61983543	11003153	61958450	61899552	62092471
61831770	62054605	61777762	61983544	11003217	61963814	61904160	62092472
61855584	62056959	61782970	61983552	11003217	61963814	61904161	62092472
61858029	62056960	61782971	61986590	11003238	61963815	61910568	62104183
61862203	62056961	61782980	61988112	11004424	61968956	61916036	62104183R
61862254	62062123	61789826	61991992	11004734	61968957	61916037	62104632
61867302	62062124	61789827	61991996	11004734	61968957	61921368	62104657
61871426	62062125	61789827	61991996	11004761	61968958	61921369	62110697
61871427	62068183	61789828	61995022	11004761	61968958	61921369	62110697
61872036	62074906	61789828	61995022	11005011	61975668	61928683	62110698
61872037	62074907	61789840	61995023	11005011	61975668	61933420	62117176
61876227	62074908	61796054	62001905	11005302	61975669	61936724	62122809
61876228	62080686	61796054	62001905	11005310	61975670	61940214	62122810
61885241	62080687	61796055	62001906	11005889	61975671	61950173	62128271
61885242	62080688	61801505	62005780	11006045	61983545	61950174	62128272
61889371	62080689	61801506	62005781	11006492	61983546	61950174	62128272
61889372	62086884	61801507	62005805	11006519	61983547	61952642	62132481
61889373	62086885	61801508	62013580	11006758	61983553	NATURAL-KNEE GSM POR FLEX FEM, SIZE 5-LT	
61894437	61894438	61801526	62013581	11007001	61986591	<u>Item Number: 00-5412-018-01</u>	
		61806570	62013592	11007081	61988113	<u>Lot Numbers:</u>	61986592
		61810246	62020392	61777752	61991993	11005767	61988114
		61810281	62020393	61777763	61991997		

NATURAL-KNEE GSM POR FLEX FEM, SIZE 5-LT		NEXGEN CR OPTION FEMORAL, SIZE D-LT		NEXGEN CR OPTION FEMORAL, SIZE E-LT		NEXGEN CR OPTION FEMORAL, SIZE E-RT	
<u>Item Number: 00-5412-018-01</u>		<u>Item Number: 00-5966-014-01</u>		<u>Item Number: 00-5966-015-01</u>		<u>Item Number: 00-5966-015-02</u>	
<u>Lot Numbers:</u>	61995026	<u>Lot Numbers:</u>	62033088	<u>Lot Numbers:</u>	61996942	<u>Lot Numbers:</u>	61991760
61777753	62001910	61883263	62038014	11400253	61996943	61782660	61998807
61782983	62001910R	61883264	62047985	11400295	61998627	61851727	61998808
61867216	62005785	61886303	62047986	11400336	61998629	61851728	62003609
61872043	77001458	61890020	62067952	61788743	61998630	61861361	62003615
61889378	77001459	61890021	62073104	61851725	62003109	61865915	62008744
61910536	77001459R	61890694	62080815	61851726	62003113	61865916	62008746
61958451	77001460	61896063	62085948	61865917	62008742	61870728	62008747
61958488	77001461	61916149	62092287	61874466	62008743	61874467	62020809
61971299	77001462	61918723	62110232	61881044	62020808	61879228	62028080
61975672	77001483	61918967	62110233	61881045	62033090	61881671	62028084
61975673		61924446	62119544	61885779	62034603	61881672	62028447
		61933809	62123030	61886270	62034604	61886323	62038814
		61943608	62133158	61899631	62038572	61886335	62039938
		61946984	62139972	61899632	62038575	61901223	62040758
<u>Item Number: 00-5412-018-02</u>		61971724	62139995	61901222	62038808	61901224	62040759
<u>Lot Numbers:</u>	61928684	61971725	62142356	61907445	62040757	61901225	62049389
11003166	61971308	61974122	62152617	61907446	62048441	61902951	62049390
11003184	61983548	61974123	62152618	61907447	62048442	61902952	62067349
11003184R	61988115	61978521	62164318	61908303	62066772	61902953	62067350
11004266	61991998	61978522	62173182	61908304	62074210	61908284	62067903
11004519	61995027	61988452		61909335	62079856	61908285	62074224
11005930	62001911	NEXGEN CR OPTION FEMORAL, SIZE D-RT		61917749	62080306	61908286	62080307
61777754	62001911R	<u>Item Number: 00-5966-014-02</u>		61917750	62086629	61909339	62080308
61782984	62005786	<u>Lot Numbers:</u>	61996348	61917751	62086630	61909480	62080813
61801512	62117177	61851720	61996349	61918721	62092284	61917752	62080814
61806574	77001452	61874459	62003768	61918722	62092524	61917753	62086631
61810285	77001453	61879226	62003769	61918970	62098791	61917754	62086632
61850060	77001454	61883265	62022323	61918971	62098792	61918719	62092525
61867304	77001455	61883266	62028039	61918972	62105730	61918720	62092526
61872044	77001456	61890696	62028045	61920436	62110795	61920437	62098968
61889379	77001457	61893232	62033089	61924448	62110809	61920438	62099266
61910537		61893234	62038015	61924449	62110810	61920439	62099267
NEXGEN CR OPTION FEMORAL, SIZE C-LT		61902949	62047525	61937799	62117283	61920440	62105731
<u>Item Number: 00-5966-013-01</u>		61902949	62047525	61942841	62123034	61924450	62105732
<u>Lot Numbers:</u>	61997338	61916150	62073105	61943188	62128363	61926126	62106560
61893236	61997339	61918968	62081425	61946975	62128375	61939182	62106566
61939720	62164866	61918969	62081426	61946980	62131027	61939184	62117328
61976414	62173181	61924447	62085949	61964947	62137079	61939185	62117329
NEXGEN CR OPTION FEMORAL, SIZE C-RT		61933810	62098790	61964949	62146937	61946981	62128376
<u>Item Number: 00-5966-013-02</u>		61937797	62105729	61965576	62155683	61946982	62129349
<u>Lot Numbers:</u>	61997340	61947190	62110794	61977306	62156809	61961808	62131207
61893622	61998618	61947191	62120689	61977316	62159281	61961809	62131208
61939721	62142586	61971726	62123033	61979126	62164823	61962865	62137083
61976415	62172419	61971727	62133828	61979130	62171734	61965577	62142595
NEXGEN CR OPTION FEMORAL, SIZE D-LT		61971728	62140115	61979774	62171789	61965934	62146938
<u>Item Number: 00-5966-014-01</u>		61976416	62141175	61979776	62174412	61965935	62156810
<u>Lot Numbers:</u>	61989037	61977294	62147815	61987062	62174413	61977550	62157840
11400307	61996346	61978523	62152642	61987063		61977551	62157958
61782655	61996347	61979109	62164670	NEXGEN CR OPTION FEMORAL, SIZE E-RT		61977552	62157959
61850922	62002067	61979116	62164672	<u>Item Number: 00-5966-015-02</u>		61982743	62164865
61855960	62002068	61989039	62164822	<u>Lot Numbers:</u>	61988191	61982744	62171845
61873439	62022322			11400252	61988194	61983003	62171846
				11400352	61991759		

NEXGEN CR OPTION FEMORAL, SIZE F-LT		NEXGEN CR OPTION FEMORAL, SIZE F-RT		NEXGEN CR OPTION FEMORAL, SIZE G-RT		NEXGEN CR PRECOAT FEMORAL, SIZE G-LT	
<u>Item Number: 00-5966-016-01</u>		<u>Item Number:00-5966-016-02</u>		<u>Item Number: 00-5966-017-02</u>		<u>Item Number: 00-5970-017-01</u>	
<u>Lot Numbers:</u>	61999559	<u>Lot Numbers:</u>	61999561	<u>Lot Numbers:</u>	62004999	<u>Lot Number:</u>	11003049
11400415	61999560	61930159	62080786	61884189	62017262	NEXGEN CR PRECOAT FEMORAL, SIZE H-LT	
61782661	62014391	61944889	62080787	61885172	62017263	<u>Item Number: 00-5970-018-01</u>	
61851729	62014392	61944890	62086610	61890763	62027200	<u>Lot Numbers:</u>	11003323
61881674	62025349	61944891	62086611	61890765	62028075	61758971	
61881675	62025350	61946488	62092527	61909343	62029523	NEXGEN CR PRECOAT FEMORAL, SIZE H-RT	
61886337	62025351	61957853	62093279	61909344	62029524	<u>Item Number: 00-5970-018-02</u>	
61890515	62028448	61959003	62098796	61913198	62040756	<u>Lot Numbers:</u>	11003322
61892006	62028449	61970683	62098837	61918727	62046898	11003330	
61892007	62029475	61970684	62106708	61923573	62046900	NEXGEN CR-FLEX OPTION FEMORAL, SIZE C-LT	
61895868	62039939	61970685	62106709	61923574	62067912	<u>Item Number: 00-5956-013-01</u>	
61895869	62042302	61972170	62118583	61923575	62086636	<u>Lot Numbers:</u>	62075014
61899630	62046231	61972171	62121797	61931673	62107212	62117143	
61912359	62046232	61972172	62121798	61936774	62122565	62056264	
61912360	62046233	61972740	62137517	61957855	62131282	NEXGEN CR-FLEX OPTION FEMORAL, SIZE C-RT	
61913193	62046890	61980577	62137537	61962870	62131283	<u>Item Number: 00-5956-013-02</u>	
61917760	62046891	61980578	62140666	61962871	62137758	<u>Lot Numbers:</u>	62086765
61918724	62063280	61980579	62141726	61962872	62140645	62117144	
61923577	62067905	61982463	62142358	61964945	62152311	NEXGEN CR-FLEX OPTION FEMORAL, SIZE C-RT MINUS	
61924451	62086608	61988454	62151140	61782666	62002039	<u>Item Number: 00-5956-013-06</u>	
61930158	62086609	61988455	62152272	61879329	62004998	<u>Lot Numbers:</u>	62124799
61935980	62098794	61992056	62159329	61964946		NEXGEN CR-FLEX OPTION FEMORAL, SIZE D-LT	
61944886	62098795	61992057	62168049	NEXGEN CR POROUS FEMORAL, SIZE E-LT		<u>Item Number: 00-5972-015-01</u>	
61944887	62106701	61926130	62074271	<u>Item Number: 00-5972-015-01</u>		<u>Lot Number:</u>	
61962873	62106702	NEXGEN CR OPTION FEMORAL, SIZE G-LT		<u>Lot Number:</u>		61777593	
61962874	62117331	<u>Item Number: 00-5966-017-01</u>		NEXGEN CR POROUS FEMORAL, SIZE H-RT		<u>Item Number: 00-5972-018-02</u>	
61966241	62121775	<u>Lot Numbers:</u>	61989042	<u>Item Number: 00-5972-018-02</u>		<u>Lot Number:</u>	
61966242	62127470	11400356	61992058	<u>Lot Number:</u>		4002192	
61972741	62131252	61782664	61992059	NEXGEN CR PRECOAT FEMORAL, SIZE B-RT		<u>Item Number: 00-5970-012-02</u>	
61977556	62137515	61866711	62000974	<u>Item Number: 00-5970-012-02</u>		<u>Lot Number:</u>	
61977557	62140660	61875442	62000975	<u>Lot Number:</u>		61878272	
61988453	62151136	61884187	62017250	NEXGEN CR PRECOAT FEMORAL, SIZE C-RT		<u>Item Number: 00-5970-013-02</u>	
61991776	62152287	61890027	62017261	<u>Item Number: 00-5970-013-02</u>		<u>Lot Numbers:</u>	
61991778	62159285	61890755	62027198	<u>Lot Numbers:</u>		61777576	
NEXGEN CR OPTION FEMORAL, SIZE F-RT		61908292	62027199	NEXGEN CR PRECOAT FEMORAL, SIZE D-RT		<u>Item Number: 00-5970-014-02</u>	
<u>Item Number:00-5966-016-02</u>		61909340	62029521	<u>Item Number: 00-5970-014-02</u>		<u>Lot Numbers:</u>	
<u>Lot Numbers:</u>	62003771	61913196	62047527	<u>Lot Numbers:</u>		61777576	
11400303	62003772	61918726	62080816	NEXGEN CR PRECOAT FEMORAL, SIZE E-RT		<u>Item Number: 00-5970-015-02</u>	
61782663	62014393	61921453	62086633	<u>Item Number: 00-5970-015-02</u>		<u>Lot Numbers:</u>	
61788749	62014394	61921454	62086634	<u>Lot Numbers:</u>		62132399	
61883268	62025176	61922262	62107194	NEXGEN CR PRECOAT FEMORAL, SIZE F-RT		<u>Item Number: 00-5970-016-02</u>	
61884184	62025177	61931672	62107198	<u>Item Number: 00-5970-016-02</u>		<u>Lot Numbers:</u>	
61888037	62025348	61957854	62128377	<u>Lot Numbers:</u>		11400192	
61892369	62028450	61962314	62131268	<u>Lot Numbers:</u>		62122787	
61892371	62039940	61962316	62140649	<u>Lot Numbers:</u>		62132399	
61895870	62042313	61962317	62140653	NEXGEN CR PRECOAT FEMORAL, SIZE G-RT		<u>Item Number: 00-5966-017-02</u>	
61895871	62042323	61972742	62159330	<u>Item Number: 00-5966-017-02</u>		<u>Lot Numbers:</u>	
61895872	62043145	61978274	62172420	<u>Lot Numbers:</u>		61973242	
61895875	62046215	61988456	62174409	NEXGEN CR OPTION FEMORAL, SIZE G-RT		<u>Item Number: 00-5966-017-02</u>	
61899629	62046216	NEXGEN CR OPTION FEMORAL, SIZE G-RT		<u>Item Number: 00-5966-017-02</u>		<u>Lot Numbers:</u>	
61913194	62046230	<u>Item Number: 00-5966-017-02</u>		<u>Lot Numbers:</u>		11006417	
61918725	62059377	<u>Lot Numbers:</u>		<u>Lot Numbers:</u>		61989045	
61924453	62062230	11400382		NEXGEN CR PRECOAT FEMORAL, SIZE F-RT		<u>Item Number: 00-5970-016-02</u>	
61924454	62067908	62000976		<u>Item Number: 00-5970-016-02</u>		<u>Lot Numbers:</u>	
61924455	62074267			<u>Lot Numbers:</u>		61777581	
				<u>Lot Numbers:</u>		62135936	
				<u>Lot Numbers:</u>		61777581	
				<u>Lot Numbers:</u>		62135936	

<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE E-LT</p> <p><u>Item Number: 00-5956-015-01</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>11005209</td><td>61972737</td></tr> <tr><td>11400132</td><td>61982306</td></tr> <tr><td>61803565</td><td>61986420</td></tr> <tr><td>61821415</td><td>61986702</td></tr> <tr><td>61821416</td><td>62024955</td></tr> <tr><td>61826928</td><td>62046227</td></tr> <tr><td>61835330</td><td>62076994</td></tr> <tr><td>61854818</td><td>62076999</td></tr> <tr><td>61854818</td><td>62084437</td></tr> <tr><td>61858473</td><td>62109654</td></tr> <tr><td>61867606</td><td>62110605</td></tr> <tr><td>61869874</td><td>62119541</td></tr> <tr><td>61876279</td><td>62127087</td></tr> <tr><td>61895313</td><td>62138241</td></tr> <tr><td>61915309</td><td>62144132</td></tr> <tr><td>61926128</td><td>62152926</td></tr> <tr><td>61934960</td><td>62162199</td></tr> <tr><td>61958022</td><td>62172418</td></tr> <tr><td>61962868</td><td></td></tr> </table>	11005209	61972737	11400132	61982306	61803565	61986420	61821415	61986702	61821416	62024955	61826928	62046227	61835330	62076994	61854818	62076999	61854818	62084437	61858473	62109654	61867606	62110605	61869874	62119541	61876279	62127087	61895313	62138241	61915309	62144132	61926128	62152926	61934960	62162199	61958022	62172418	61962868		<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE F-LT</p> <p><u>Item Number: 00-5956-016-01</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>61989682</td><td>61810703</td></tr> <tr><td>62030345</td><td>61826933</td></tr> <tr><td>62030593</td><td>61833599</td></tr> <tr><td>62054911</td><td>61864881</td></tr> <tr><td>62083126</td><td>61895864</td></tr> <tr><td>62083127</td><td>61907448</td></tr> <tr><td>62109161</td><td>61915313</td></tr> <tr><td>62114098</td><td>61932289</td></tr> <tr><td>62117146</td><td>61932291</td></tr> <tr><td>62126728</td><td>61971163</td></tr> <tr><td>62154274</td><td>61978525</td></tr> <tr><td>62163103</td><td>61986178</td></tr> </table>	61989682	61810703	62030345	61826933	62030593	61833599	62054911	61864881	62083126	61895864	62083127	61907448	62109161	61915313	62114098	61932289	62117146	61932291	62126728	61971163	62154274	61978525	62163103	61986178	<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE G-LT MINUS</p> <p><u>Item Number: 00-5956-017-05</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62126102</td></tr> <tr><td>62134747</td></tr> </table>	62126102	62134747	<p>NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE C-RT</p> <p><u>Item Number: 00-5950-013-02</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62030313</td></tr> </table>	62030313																																									
11005209	61972737																																																																																																												
11400132	61982306																																																																																																												
61803565	61986420																																																																																																												
61821415	61986702																																																																																																												
61821416	62024955																																																																																																												
61826928	62046227																																																																																																												
61835330	62076994																																																																																																												
61854818	62076999																																																																																																												
61854818	62084437																																																																																																												
61858473	62109654																																																																																																												
61867606	62110605																																																																																																												
61869874	62119541																																																																																																												
61876279	62127087																																																																																																												
61895313	62138241																																																																																																												
61915309	62144132																																																																																																												
61926128	62152926																																																																																																												
61934960	62162199																																																																																																												
61958022	62172418																																																																																																												
61962868																																																																																																													
61989682	61810703																																																																																																												
62030345	61826933																																																																																																												
62030593	61833599																																																																																																												
62054911	61864881																																																																																																												
62083126	61895864																																																																																																												
62083127	61907448																																																																																																												
62109161	61915313																																																																																																												
62114098	61932289																																																																																																												
62117146	61932291																																																																																																												
62126728	61971163																																																																																																												
62154274	61978525																																																																																																												
62163103	61986178																																																																																																												
62126102																																																																																																													
62134747																																																																																																													
62030313																																																																																																													
<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE E-LT MINUS</p> <p><u>Item Number: 00-5956-015-05</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62122982</td><td>62091540</td></tr> <tr><td>62127088</td><td>62119542</td></tr> <tr><td>62138258</td><td></td></tr> </table>	62122982	62091540	62127088	62119542	62138258		<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE F-LT MINUS</p> <p><u>Item Number: 00-5956-016-05</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62138681</td><td>62121616</td></tr> <tr><td>62162202</td><td>62127085</td></tr> <tr><td>62170903</td><td></td></tr> </table>	62138681	62121616	62162202	62127085	62170903		<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE G-RT</p> <p><u>Item Number: 00-5956-017-02</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62030636</td></tr> <tr><td>61803569</td><td>62051197</td></tr> <tr><td>61815712</td><td>62073095</td></tr> <tr><td>61855965</td><td>62096010</td></tr> <tr><td>61885778</td><td>62114025</td></tr> <tr><td>61901226</td><td>62117149</td></tr> <tr><td>61929404</td><td>62126736</td></tr> <tr><td>61971165</td><td>62143640</td></tr> <tr><td>61991756</td><td>62161216</td></tr> <tr><td>62007004</td><td>62170057</td></tr> </table>	62030636	61803569	62051197	61815712	62073095	61855965	62096010	61885778	62114025	61901226	62117149	61929404	62126736	61971165	62143640	61991756	62161216	62007004	62170057	<p>NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE C-RT MINUS</p> <p><u>Item Number: 00-5950-013-06</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>11003339</td><td>62171088</td></tr> </table>	11003339	62171088																																																																									
62122982	62091540																																																																																																												
62127088	62119542																																																																																																												
62138258																																																																																																													
62138681	62121616																																																																																																												
62162202	62127085																																																																																																												
62170903																																																																																																													
62030636																																																																																																													
61803569	62051197																																																																																																												
61815712	62073095																																																																																																												
61855965	62096010																																																																																																												
61885778	62114025																																																																																																												
61901226	62117149																																																																																																												
61929404	62126736																																																																																																												
61971165	62143640																																																																																																												
61991756	62161216																																																																																																												
62007004	62170057																																																																																																												
11003339	62171088																																																																																																												
<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE E-RT</p> <p><u>Item Number: 00-5956-015-02</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>61982462</td><td>11400116</td></tr> <tr><td>61982742</td><td>61803567</td></tr> <tr><td>61988148</td><td>61822947</td></tr> <tr><td>62002066</td><td>61822948</td></tr> <tr><td>62012658</td><td>61835393</td></tr> <tr><td>62022038</td><td>61842348</td></tr> <tr><td>62028028</td><td>61849045</td></tr> <tr><td>62047517</td><td>61859222</td></tr> <tr><td>62072210</td><td>61864878</td></tr> <tr><td>62078811</td><td>61864879</td></tr> <tr><td>62085937</td><td>61869876</td></tr> <tr><td>62098284</td><td>61876938</td></tr> <tr><td>62110230</td><td>61897287</td></tr> <tr><td>62127119</td><td>61901218</td></tr> <tr><td>62127468</td><td>61916147</td></tr> <tr><td>62144290</td><td>61927965</td></tr> <tr><td>62154043</td><td>61929400</td></tr> <tr><td>62162764</td><td>61963686</td></tr> <tr><td>62169024</td><td>61973239</td></tr> <tr><td>62171732</td><td>61973241</td></tr> </table>	61982462	11400116	61982742	61803567	61988148	61822947	62002066	61822948	62012658	61835393	62022038	61842348	62028028	61849045	62047517	61859222	62072210	61864878	62078811	61864879	62085937	61869876	62098284	61876938	62110230	61897287	62127119	61901218	62127468	61916147	62144290	61927965	62154043	61929400	62162764	61963686	62169024	61973239	62171732	61973241	<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE F-RT</p> <p><u>Item Number: 00-5956-016-02</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>61986713</td><td>61811794</td></tr> <tr><td>62026095</td><td>61817755</td></tr> <tr><td>62043310</td><td>61827940</td></tr> <tr><td>62077294</td><td>61865912</td></tr> <tr><td>62083130</td><td>61873444</td></tr> <tr><td>62100468</td><td>61880101</td></tr> <tr><td>62107984</td><td>61894398</td></tr> <tr><td>62115961</td><td>61910871</td></tr> <tr><td>62117147</td><td>61921452</td></tr> <tr><td>62126729</td><td>61964438</td></tr> <tr><td>62145821</td><td>61978526</td></tr> <tr><td>62156807</td><td>61986709</td></tr> <tr><td>62162770</td><td></td></tr> </table>	61986713	61811794	62026095	61817755	62043310	61827940	62077294	61865912	62083130	61873444	62100468	61880101	62107984	61894398	62115961	61910871	62117147	61921452	62126729	61964438	62145821	61978526	62156807	61986709	62162770		<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE G-RT MINUS</p> <p><u>Item Number: 00-5956-017-06</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62160793</td><td>62079857</td></tr> <tr><td>62169051</td><td></td></tr> </table>	62160793	62079857	62169051		<p>NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE D-LT</p> <p><u>Item Number: 00-5950-014-01</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>61999556</td><td>11400334</td></tr> <tr><td>62005884</td><td>61817750</td></tr> <tr><td>62011561</td><td>61827935</td></tr> <tr><td>62022990</td><td>61832285</td></tr> <tr><td>62030319</td><td>61835326</td></tr> <tr><td>62039177</td><td>61857891</td></tr> <tr><td>62051191</td><td>61857892</td></tr> <tr><td>62062520</td><td>61900460</td></tr> <tr><td>62092262</td><td>61900462</td></tr> <tr><td>62097716</td><td>61906384</td></tr> <tr><td>62100467</td><td>61921448</td></tr> <tr><td>62107134</td><td>61933834</td></tr> <tr><td>62124802</td><td>61940956</td></tr> <tr><td>62131186</td><td>61964633</td></tr> <tr><td>62142614</td><td>61974381</td></tr> <tr><td>62160335</td><td>61980308</td></tr> <tr><td>62166742</td><td>61992035</td></tr> <tr><td>62171089</td><td>61992036</td></tr> </table>	61999556	11400334	62005884	61817750	62011561	61827935	62022990	61832285	62030319	61835326	62039177	61857891	62051191	61857892	62062520	61900460	62092262	61900462	62097716	61906384	62100467	61921448	62107134	61933834	62124802	61940956	62131186	61964633	62142614	61974381	62160335	61980308	62166742	61992035	62171089	61992036
61982462	11400116																																																																																																												
61982742	61803567																																																																																																												
61988148	61822947																																																																																																												
62002066	61822948																																																																																																												
62012658	61835393																																																																																																												
62022038	61842348																																																																																																												
62028028	61849045																																																																																																												
62047517	61859222																																																																																																												
62072210	61864878																																																																																																												
62078811	61864879																																																																																																												
62085937	61869876																																																																																																												
62098284	61876938																																																																																																												
62110230	61897287																																																																																																												
62127119	61901218																																																																																																												
62127468	61916147																																																																																																												
62144290	61927965																																																																																																												
62154043	61929400																																																																																																												
62162764	61963686																																																																																																												
62169024	61973239																																																																																																												
62171732	61973241																																																																																																												
61986713	61811794																																																																																																												
62026095	61817755																																																																																																												
62043310	61827940																																																																																																												
62077294	61865912																																																																																																												
62083130	61873444																																																																																																												
62100468	61880101																																																																																																												
62107984	61894398																																																																																																												
62115961	61910871																																																																																																												
62117147	61921452																																																																																																												
62126729	61964438																																																																																																												
62145821	61978526																																																																																																												
62156807	61986709																																																																																																												
62162770																																																																																																													
62160793	62079857																																																																																																												
62169051																																																																																																													
61999556	11400334																																																																																																												
62005884	61817750																																																																																																												
62011561	61827935																																																																																																												
62022990	61832285																																																																																																												
62030319	61835326																																																																																																												
62039177	61857891																																																																																																												
62051191	61857892																																																																																																												
62062520	61900460																																																																																																												
62092262	61900462																																																																																																												
62097716	61906384																																																																																																												
62100467	61921448																																																																																																												
62107134	61933834																																																																																																												
62124802	61940956																																																																																																												
62131186	61964633																																																																																																												
62142614	61974381																																																																																																												
62160335	61980308																																																																																																												
62166742	61992035																																																																																																												
62171089	61992036																																																																																																												
<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE E-RT MINUS</p> <p><u>Item Number: 00-5956-015-06</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62085936</td><td>62127118</td></tr> </table>	62085936	62127118	<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE F-RT MINUS</p> <p><u>Item Number: 00-5956-016-06</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62127086</td><td>62075017</td></tr> <tr><td>62169047</td><td>62107985</td></tr> </table>	62127086	62075017	62169047	62107985	<p>NEXGEN CR-FLEX POROUS FEMORAL, SIZE F-LT</p> <p><u>Item Number: 00-5952-016-01</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>61774124</td><td>61782595</td></tr> </table>	61774124	61782595	<p>NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE D-LT MINUS</p> <p><u>Item Number: 00-5950-014-05</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62058212</td><td>61835327</td></tr> <tr><td>62122662</td><td></td></tr> </table>	62058212	61835327	62122662																																																																																															
62085936	62127118																																																																																																												
62127086	62075017																																																																																																												
62169047	62107985																																																																																																												
61774124	61782595																																																																																																												
62058212	61835327																																																																																																												
62122662																																																																																																													
	<p>NEXGEN CR-FLEX OPTION FEMORAL, SIZE G-LT</p> <p><u>Item Number: 00-5956-017-01</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>61996944</td><td>61797740</td></tr> <tr><td>62018673</td><td>61813670</td></tr> <tr><td>62052670</td><td>61814765</td></tr> <tr><td>62078825</td><td>61837325</td></tr> <tr><td>62084440</td><td>61876918</td></tr> <tr><td>62086771</td><td>61903674</td></tr> <tr><td>62112977</td><td>61918975</td></tr> <tr><td>62125785</td><td>61946490</td></tr> <tr><td>62143177</td><td>61971164</td></tr> <tr><td>62160337</td><td>61973245</td></tr> </table>	61996944	61797740	62018673	61813670	62052670	61814765	62078825	61837325	62084440	61876918	62086771	61903674	62112977	61918975	62125785	61946490	62143177	61971164	62160337	61973245	<p>NEXGEN CR-FLEX POROUS FEMORAL, SIZE F-RT MINUS</p> <p><u>Item Number: 00-5952-016-06</u></p> <p><u>Lot Number:</u></p> <table border="1"> <tr><td>11004288</td></tr> </table>	11004288	<p>NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE D-RT</p> <p><u>Item Number: 00-5950-014-02</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62007000</td><td>11400337</td></tr> <tr><td>62011562</td><td>61827936</td></tr> <tr><td>62019491</td><td>61830369</td></tr> <tr><td>62030320</td><td>61844039</td></tr> <tr><td>62040699</td><td>61847663</td></tr> <tr><td>62058769</td><td>61857895</td></tr> <tr><td>62068569</td><td>61857896</td></tr> <tr><td>62092285</td><td>61900463</td></tr> <tr><td>62097717</td><td>61901179</td></tr> <tr><td>62097718</td><td>61922252</td></tr> <tr><td>62119167</td><td>61939801</td></tr> <tr><td>62126098</td><td>61939802</td></tr> <tr><td>62126099</td><td>61964638</td></tr> <tr><td>62137073</td><td>61980573</td></tr> <tr><td>62148759</td><td>61980574</td></tr> <tr><td>62149010</td><td>61992037</td></tr> <tr><td>62167498</td><td>61999557</td></tr> <tr><td>62171787</td><td></td></tr> </table>	62007000	11400337	62011562	61827936	62019491	61830369	62030320	61844039	62040699	61847663	62058769	61857895	62068569	61857896	62092285	61900463	62097717	61901179	62097718	61922252	62119167	61939801	62126098	61939802	62126099	61964638	62137073	61980573	62148759	61980574	62149010	61992037	62167498	61999557	62171787																																																		
61996944	61797740																																																																																																												
62018673	61813670																																																																																																												
62052670	61814765																																																																																																												
62078825	61837325																																																																																																												
62084440	61876918																																																																																																												
62086771	61903674																																																																																																												
62112977	61918975																																																																																																												
62125785	61946490																																																																																																												
62143177	61971164																																																																																																												
62160337	61973245																																																																																																												
11004288																																																																																																													
62007000	11400337																																																																																																												
62011562	61827936																																																																																																												
62019491	61830369																																																																																																												
62030320	61844039																																																																																																												
62040699	61847663																																																																																																												
62058769	61857895																																																																																																												
62068569	61857896																																																																																																												
62092285	61900463																																																																																																												
62097717	61901179																																																																																																												
62097718	61922252																																																																																																												
62119167	61939801																																																																																																												
62126098	61939802																																																																																																												
62126099	61964638																																																																																																												
62137073	61980573																																																																																																												
62148759	61980574																																																																																																												
62149010	61992037																																																																																																												
62167498	61999557																																																																																																												
62171787																																																																																																													
		<p>NEXGEN CR-FLEX POROUS FEMORAL, SIZE G-LT</p> <p><u>Item Number: 00-5952-017-01</u></p> <p><u>Lot Number:</u></p> <table border="1"> <tr><td>11005095</td></tr> </table>	11005095																																																																																																										
11005095																																																																																																													
		<p>NEXGEN CR-FLEX PRECOAT FEMORAL SIZE E-LT MINUS</p> <p><u>Item Number: 00-5950-015-05</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62030633</td><td>61823885</td></tr> <tr><td>62054070</td><td>61858477</td></tr> <tr><td>62091541</td><td>61876281</td></tr> <tr><td>62138242</td><td>61953097</td></tr> <tr><td>62172417</td><td>61999987</td></tr> </table>	62030633	61823885	62054070	61858477	62091541	61876281	62138242	61953097	62172417	61999987																																																																																																	
62030633	61823885																																																																																																												
62054070	61858477																																																																																																												
62091541	61876281																																																																																																												
62138242	61953097																																																																																																												
62172417	61999987																																																																																																												
		<p>NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE B-RT</p> <p><u>Item Number: 00-5950-012-02</u></p> <p><u>Lot Number:</u></p> <table border="1"> <tr><td>61777545</td></tr> </table>	61777545																																																																																																										
61777545																																																																																																													
		<p>NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE C-LT</p> <p><u>Item Number: 00-5950-013-01</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62092032</td><td>61854815</td></tr> <tr><td>62171087</td><td></td></tr> </table>	62092032	61854815	62171087																																																																																																								
62092032	61854815																																																																																																												
62171087																																																																																																													
		<p>NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE C-RT</p> <p><u>Item Number: 00-5950-013-02</u></p> <p><u>Lot Numbers:</u></p> <table border="1"> <tr><td>62005882</td></tr> </table>	62005882																																																																																																										
62005882																																																																																																													

NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE G-LT		NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE G-RT		NEXGEN KNEE LPS OPTION FEMORAL, SIZE C-LT		NEXGEN KNEE LPS OPTION FEMORAL, SIZE D-LT	
<u>Item Number: 00-5950-017-01</u>		<u>Item Number: 00-5950-017-02</u>		<u>Item Number: 00-5996-013-01</u>		<u>Item Number: 00-5996-014-01</u>	
<u>Lot Numbers:</u>	61999571	<u>Lot Numbers:</u>	61990919	<u>Lot Numbers:</u>	62170068	<u>Lot Numbers:</u>	61979137
11400135	62005000	11400232	62005953	62008974	62172687	61855971	62058223
11400139	62005001	11400312	62005954	62051204		61868871	62066791
11400156	62012670	61800618	62005955	NEXGEN KNEE LPS OPTION FEMORAL, SIZE C-RT		61868872	62066792
11400282	62018146	61800619	62006741	<u>Item Number: 00-5996-013-02</u>		61868873	62066793
11400313	62018672	61806282	62006743	<u>Lot Numbers:</u>	62003779	61871529	62066794
11400317	62034607	61815714	62019996	11400172	62004852	61871530	62067351
11400354	62034610	61827942	62026806	61782578	62009310	61886272	62074350
11400357	62038017	61829062	62030596	61850932	62009311	61886273	62074353
61800616	62038018	61832308	62036588	61864885	62080341	61886274	62076157
61800617	62052668	61835933	62036609	61922268	62091041	61886326	62076322
61814766	62052669	61841279	62051196	61924458	62100473	61886327	62076323
61814767	62060542	61844942	62054065	61974046	62109165	61890029	62078827
61815711	62065390	61870726	62060479	61974131	62112981	61890030	62085957
61829065	62066548	61880104	62060480	61974132	62124859	61890031	62092269
61829066	62078822	61898497	62067915	61974133	62134305	61895324	62092270
61834859	62078824	61898837	62073089	61977563	62152623	61895325	62092271
61840473	62084439	61899628	62086359	61978530	62161218	61895326	62097831
61844902	62092283	61902954	62086360	61989687	62170070	61895327	62097832
61854822	62095239	61912358	62091543	61990926		61896041	62098067
61860221	62102785	61917758	62096011	NEXGEN KNEE LPS OPTION FEMORAL, SIZE D-LT		61896042	62098068
61899627	62109668	61920447	62096621	<u>Item Number: 00-5996-014-01</u>		61910215	62098069
61903675	62112976	61920448	62105718	<u>Lot Numbers:</u>	61979138	61910216	62102068
61912355	62120702	61927969	62110814	61782547	61979139	61914494	62102069
61916155	62124820	61934965	62113002	61782562	61990929	61914495	62102070
61918974	62125144	61934966	62113003	61804267	61992039	61914496	62109168
61934963	62134731	61952883	62126115	61804268	61992040	61914497	62109669
61936776	62143173	61952884	62126732	61808647	61992041	61915316	62109670
61936778	62154276	61956827	62132540	61808651	61997344	61918706	62110799
61946988	62160283	61961812	62135963	61809730	61997345	61918707	62110815
61950760	62170059	61962312	62143620	61809731	61998635	61918709	62110816
61951696	62170392	61974125	62149040	61809731	61998635	61918710	62110817
61960712	62170393	61974126	62155001	61820657	62004855	61918711	62111851
61960713	11400238	61982246	62160794	61820668	62004858	61922269	62113006
61960714	11400440	61982464	62170053	61820730	62004861	61922270	62113007
61973243	61833598	61990918	62170055	61829136	62004861	61922271	62113008
61973244	61937290	NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE G-RT MINUS		61829136	62005005	61924459	62114442
61980313	61960715	<u>Item Number: 00-5950-017-06</u>		61835341	62010337	61924460	62114446
61984231	61999572	<u>Lot Numbers:</u>	62030646	61835342	62010338	61924461	62114449
61990637	62036734	61816389	62067918	61835835	62010339	61926136	62118556
61996360	62061873	61917756	62134748	61835836	62011349	61930162	62118557
61996361	62121230	61951701	62160792	61843606	62018147	61931676	62124313
61996945	62126101	62000986	62169050	61844576	62029528	61935988	62124315
NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE G-LT MINUS		NEXGEN KNEE LPS OPTION FEMORAL, SIZE C-LT		61844577	62029529	61935989	62124316
<u>Item Number: 00-5950-017-05</u>		<u>Item Number: 00-5996-013-01</u>		61844578	62040744	61935990	62124840
<u>Lot Numbers:</u>	61999572	<u>Lot Numbers:</u>	62091551	61844906	62040764	61935991	62131254
11400238	62036734	61782577	62100471	61844907	62040765	61936798	62131255
11400440	62061873	61863954	62107986	61848266	62050086	61941590	62134284
61833598	62121230	61924457	62107987	61848267	62050087	61941591	62134285
61937290	62126101	61977561	62118558	61848268	62050088	61964442	62134287
61960715		61977562	62124601	61848269	62050285	61964443	62141180
NEXGEN CR-FLEX PRECOAT FEMORAL, SIZE G-RT		61977562	62124601	61848270	62052112	61971730	62141727
<u>Item Number: 00-5950-017-02</u>		62003777	62134289	61853388	62052113	61971731	62142361
<u>Lot Numbers:</u>	62000987	62003778	62151257	61853389	62052114	61971732	62150038
		62008972	62161217	61853390	62058204	61972175	62150039
				61855970	62058222	61972176	62152624

NEXGEN KNEE LPS OPTION FEMORAL, SIZE E-LT		NEXGEN KNEE LPS OPTION FEMORAL, SIZE E-RT		NEXGEN KNEE LPS OPTION FEMORAL, SIZE E-RT		NEXGEN KNEE LPS OPTION FEMORAL, SIZE E-RT	
<u>Item Number: 00-5996-015-01</u>		<u>Item Number: 00-5996-015-02</u>		<u>Item Number: 00-5996-015-02</u>		<u>Item Number: 00-5996-015-02</u>	
<u>Lot Numbers:</u>	61982319	<u>Lot Numbers:</u>	61985032	<u>Lot Numbers:</u>	62055037	<u>Lot Numbers:</u>	62111952
61916156	62107117	11005547	61986205	61847709	62055038	61929409	62116024
61916157	62107946	11400108	61986207	61847710	62056776	61929410	62116026
61916158	62107947	11400143	61986211	61850229	62056777	61929411	62116027
61916175	62111953	11400220	61993954	61850230	62061653	61929412	62117279
61917761	62111994	11400433	61993957	61850231	62061654	61930160	62117280
61918957	62111995	61782568	61993960	61850232	62061655	61934968	62117281
61920450	62114026	61797751	61995688	61850933	62061880	61934969	62117282
61920451	62114027	61797752	61995689	61850934	62062554	61934970	62121232
61920452	62114099	61804739	62000000	61851731	62062556	61934971	62121233
61920453	62114101	61804740	62000001	61851732	62064894	61937807	62121777
61920454	62114451	61804741	62000004	61851733	62064895	61937808	62121778
61926141	62114452	61804742	62000977	61852861	62065327	61938903	62128379
61927970	62114455	61806283	62000978	61852862	62065340	61938904	62128380
61927971	62119549	61806284	62000979	61890015	62065351	61938905	62129225
61927972	62119550	61813678	62002044	61890016	62065354	61939186	62129227
61927973	62124865	61813679	62002045	61890691	62066543	61946494	62133179
61932298	62124874	61814768	62007007	61890692	62066544	61946495	62134714
61933648	62125505	61814769	62007008	61890702	62073111	61958544	62134749
61933658	62125506	61814770	62012110	61890703	62073112	61958545	62134750
61933664	62125507	61815716	62012115	61890704	62074299	61958546	62134751
61937297	62132867	61815717	62012122	61899622	62081869	61958763	62144186
61937762	62137990	61815718	62012127	61899623	62083135	61958764	62144187
61937763	62137991	61815719	62012663	61899624	62083137	61960722	62149529
61937764	62138243	61819956	62012664	61899636	62087936	61964972	62149546
61937804	62142619	61820646	62013563	61899637	62088885	61964973	62149547
61937805	62142620	61820647	62022048	61899638	62088888	61965580	62150036
61944892	62149475	61820648	62022049	61900454	62090402	61965581	62152934
61944893	62149476	61820714	62022328	61900455	62091020	61965582	62152935
61958028	62149526	61823898	62023800	61900456	62091476	61966246	62152936
61958029	62152622	61824782	62023801	61900457	62091485	61966247	62154044
61963902	62152644	61824783	62023802	61900464	62091548	61971171	62160287
61963903	62152645	61826936	62024962	61900465	62091550	61971178	62160340
61963904	62157961	61834872	62024963	61905293	62096714	61971180	62160341
61963905	62159292	61834873	62028455	61905294	62096715	61976422	62160801
61964138	62159293	61834874	62028456	61906071	62097828	61977366	62160802
61964145	62164669	61835824	62032745	61906389	62097830	61977374	62166651
61964444	62164819	61835825	62033091	61906390	62098064	61977380	62166652
61964445	62164821	61835826	62033092	61907454	62098065	61982471	62166740
61968859	62164858	61838756	62033714	61909349	62103217	61982472	62166741
61968868	62164863	61846251	62033715	61909350	62105738	61982473	62167495
61968871	62164864	61846252	62034612	61910208	62105740	61982474	62167496
61970686	62165853	61846253	62042145	61917762	62105741	61982781	62167497
61970687	62165856	61846263	62042147	61917763	62106540	61982793	62171095
61974134	62165931	61846264	62042151	61917764	62106543	61983009	62171593
61976419	62168303	61846265	62042153	61917765	62106545	61983010	62171594
61976420	62168304	61846266	62042285	61917766	62106691	61983011	62171595
61976421	62170395	61846267	62043136	61918729	62106692	61985030	62171596
61979784	62170396	61847657	62048259	61918730	62106693	61985031	62173187
61979785	62170910	61847658	62048265	61918731	62106694	NEXGEN KNEE LPS OPTION FEMORAL, SIZE F-LT	
61980584	62172676	61847659	62048436	61920455	62111494		
61980585	62172678	61847660	62048437	61921455	62111495	<u>Lot Numbers:</u>	61980320
61980586	62172681	61847677	62048438	61921456	62111496	11400216	61980580
61980587	62172682	61847678	62050078	61921457	62111498	11400223	61980581
61982315	62172683	61847679	62055036		61921458		

NEXGEN KNEE LPS OPTION FEMORAL, SIZE F-LT		NEXGEN KNEE LPS OPTION FEMORAL, SIZE F-LT		NEXGEN KNEE LPS OPTION FEMORAL, SIZE F-RT		NEXGEN KNEE LPS OPTION FEMORAL, SIZE F-RT	
<u>Item Number: 00-5996-016-01</u>		<u>Item Number: 00-5996-016-01</u>		<u>Item Number: 00-5996-016-02</u>		<u>Item Number: 00-5996-016-02</u>	
<u>Lot Numbers:</u>	61914456	<u>Lot Numbers:</u>	62107991	<u>Lot Numbers:</u>	61937810	<u>Lot Numbers:</u>	62148244
11400246	61986215	61914457	62117332	61844591	62046237	61937811	62155685
11400248	61986217	61914458	62117333	61844592	62046238	61937812	62155686
11400279	61990279	61914459	62117334	61846268	62052672	61938951	62160285
11400292	61990280	61918732	62117335	61846269	62052673	61946501	62166661
11400503	61990281	61918733	62121206	61846270	62052674	61946502	62166662
61782551	61990282	61918734	62121234	61849052	62052675	61980582	62166746
61792785	62002070	61918734	62121234	61849053	62052675	61980583	62166747
61794725	62002071	61918952	62126117	61849054	62052678	61982307	62167511
61794726	62008279	61922266	62126118	61849055	62052722	61982308	62167512
61794727	62008280	61922267	62126745	61849056	62059387	61982309	62167513
61804272	62013564	61923578	62131023	61849056	62059388	61982467	62168311
61804273	62014085	61923579	62137109	61849057	62059389	61986727	62168313
61804274	62014086	61926142	62137520	61850223	62060481	NEXGEN KNEE LPS OPTION FEMORAL, SIZE G-LT	
61811786	62014086	61926143	62144316	61851706	62060482		
61811787	62019501	61926144	62144317	61851707	62067953	<u>Lot Numbers:</u>	61982468
61811787	62028457	61926145	62148245	61851708	62067954	11400269	61982469
61817762	62029546	61926146	62156811	61851709	62068099	11400338	61997346
61817763	62029547	61932299	62156812	61856761	62068312	11400386	62015858
61817764	62030605	61932300	62159294	61856762	62068313	61796376	62015859
61820742	62030606	61932301	62159331	61864890	62068314	61796377	62016808
61821428	62036145	61936805	62164359	61864891	62070836	61804745	62026099
61821429	62036146	61936806	62164675	61889065	62070837	61804746	62026100
61835343	62036147	61936807	62164676	61889066	62078832	61808890	62027205
61835344	62036148	61936808	62164831	61889067	62080339	61808891	62034616
61835345	62047530	61936809	62164832	61896050	62084443	61808892	62034617
61835346	62047531	61936810	62164871	61896051	62084444	61835839	62043318
61839733	62050778	61944902	62164872	61896052	62092529	61842361	62043319
61844593	62051199	61979783	62168354	61897294	62092530	61850225	62050083
61844912	62051202	61980319	62168355	61901177	62092531	61850226	62050084
61844913	62052108	NEXGEN KNEE LPS OPTION FEMORAL, SIZE F-RT		61901178	62098840	61852840	62054089
61844914	62052109	<u>Item Number: 00-5996-016-02</u>		61907457	62098841	61852841	62065621
61844952	62058220	<u>Lot Numbers:</u>	61986728	61907458	62099270	61870738	62074427
61848296	62058221	11005506	61990643	61910876	62099271	61883275	62079860
61848297	62066551	11400217	61990646	61910877	62099272	61883276	62086376
61848298	62066552	11400222	61990648	61914460	62099910	61889068	62093706
61848299	62067353	11400222	62002072	61914461	62107992	61889077	62093707
61848301	62067354	11400278	62002073	61915322	62109172	61889078	62099048
61850937	62072524	61794724	62002073	61915323	62109173	61912362	62099049
61850938	62072525	61804275	62008281	61918953	62109174	61915324	62099273
61851704	62076158	61804276	62008282	61918954	62118555	61921459	62109747
61851705	62076326	61804277	62014087	61918955	62121802	61930164	62113005
61855974	62076329	61804743	62014427	61918956	62121803	61931682	62117336
61863958	62077301	61804744	62014429	61923580	62127174	61935993	62122866
61863959	62078810	61810418	62014431	61923581	62127175	61935994	62131270
61894499	62083131	61810419	62014432	61923582	62127176	61944862	62145407
61894501	62083133	61835347	62015857	61924456	62131204	61962880	62148771
61894502	62091885	61835348	62019503	61927975	62131205	61962881	62154075
61895328	62092096	61835841	62019504	61927976	62131206	61962882	62161225
61895329	62092097	61835842	62029869	61927977	62131217	61966911	62165233
61900466	62100477	61841282	62032090	61929413	62131218	61966920	62165937
61900467	62101634	61841283	62032869	61929414	62138245	61966928	62169037
61907455	62102071	61841284	62032879	61933842	62138246	61977558	
61907456	62107935	61844052	62039083	61933843	62142599		
61910213	62107936	61844053	62039084	61937765	62142601		
61910874	62107952	61844590	62039085	61937809	62147822		

NEXGEN KNEE LPS OPTION FEMORAL, SIZE G-RT		NEXGEN LCCK FEMORAL, SIZE E-LT		NEXGEN LPS OPTION FEMORAL, SIZE E-LT		NEXGEN LPS PRECOAT FEMORAL, SIZE E-RT	
<u>Item Number: 00-5996-017-02</u>		<u>Item Number: 00-5994-015-91</u>		<u>Item Number: 00-5996-015-51</u>		<u>Item Number: 00-5990-015-02</u>	
<u>Lot Numbers:</u>	61978301	<u>Lot Number:</u>	11003022	<u>Lot Numbers:</u>	62166597	<u>Lot Numbers:</u>	11003324
11400221	61983005	NEXGEN LCCK FEMORAL, SIZE E-RT		62121199	62170916	61776965	
11400300	62018152	<u>Item Number: 00-5994-015-92</u>		62133177	62171094	NEXGEN LPS PRECOAT FEMORAL, SIZE G-LT	
61797756	62025181	<u>Lot Numbers:</u>	61771826	62138260	62173164	<u>Item Number: 00-5990-017-91</u>	
61798634	62025182	61771822	61771827	NEXGEN LPS OPTION FEMORAL, SIZE E-RT		<u>Lot Number:</u> 61771841	
61806288	62025372	61771825	61771828	<u>Item Number: 00-5996-015-52</u>		NEXGEN LPS-FLEX OPTION FEMORAL, C-LT	
61806289	62025376	NEXGEN LCCK FEMORAL, SIZE F-LT		<u>Lot Numbers:</u>	62171831	<u>Item Number: 00-5964-013-01</u>	
61806290	62043153	<u>Item Number: 00-5994-016-91</u>		62121801	62171832	<u>Lot Number:</u> 77001358	
61835840	62043154	<u>Lot Numbers:</u>	61768099	62135973	62172688	NEXGEN LPS-FLEX OPTION FEMORAL, SIZE C-LT	
61842362	62048444	11005670	61771829	62136386	62172690	<u>Item Number: 00-5964-013-51</u>	
61850928	62049391	NEXGEN LCCK FEMORAL, SIZE F-RT		62160803	62172691	<u>Lot Numbers:</u>	
61850929	62055051	<u>Item Number: 00-5994-016-92</u>		62171830	62009298		
61850930	62060537	<u>Lot Numbers:</u>	61768104	NEXGEN LPS OPTION FEMORAL, SIZE F-LT		61808871	62018048
61852845	62067352	61771832	61771832	<u>Item Number: 00-5996-016-51</u>	61876900		62080785
61884197	62074433	NEXGEN LPS OPTION FEMORAL, SIZE B-LT		<u>Lot Numbers:</u>	62137521	61890747	62120585
61889079	62074436	<u>Item Number: 00-5996-012-01</u>		61792784	62137541	61890748	62160331
61889080	62080817	<u>Lot Number:</u>	11004400	62121779	61990909		62170030
61894563	62086993	NEXGEN LPS OPTION FEMORAL, SIZE B-RT		NEXGEN LPS OPTION FEMORAL, SIZE F-RT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE C-RT	
61894564	62099934	<u>Item Number: 00-5996-012-02</u>		<u>Item Number: 00-5996-016-52</u>	<u>Item Number: 00-5964-013-52</u>		
61898502	62099938	<u>Lot Number:</u>	62068093	<u>Lot Numbers:</u>	62138262	<u>Lot Numbers:</u>	61999924
61913202	62100474	62131269	62068093	61782556	62139469	11400373	62010044
61916162	62110801	NEXGEN LPS OPTION FEMORAL, SIZE C-RT		62122572	62170933	61858467	62083118
61922264	62110818	<u>Item Number: 00-5996-013-51</u>		NEXGEN LPS OPTION FEMORAL, SIZE G-LT		61876901	62107109
61929417	62110819	<u>Lot Numbers:</u>	62118559	<u>Item Number: 00-5996-017-51</u>	<u>Item Number: 00-5996-017-51</u>		61876902
61930163	62117338	61782571	62134293	<u>Lot Numbers:</u>	62169055	61992032	62131121
61934977	62117339	NEXGEN LPS OPTION FEMORAL, SIZE C-RT		62122868	62173166	NEXGEN LPS-FLEX OPTION FEMORAL, SIZE D-LT	
61934978	62123889	<u>Item Number: 00-5996-013-52</u>		62131271	<u>Item Number: 00-5964-014-51</u>		
61939188	62132112	<u>Lot Numbers:</u>	62118560	NEXGEN LPS OPTION FEMORAL, SIZE G-RT		<u>Lot Numbers:</u>	61992698
61939189	62145408	62131269	62172686	<u>Item Number: 00-5996-017-52</u>	<u>Item Number: 00-5996-017-52</u>		61794708
61941599	62148247	NEXGEN LPS OPTION FEMORAL, SIZE D-LT		<u>Lot Numbers:</u>	62169057	61796364	61992701
61962883	62154285	<u>Item Number: 00-5996-014-51</u>		62123891	62170066	61804727	61999966
61962884	62161596	<u>Lot Numbers:</u>	62134711	62140118	62173191	61809705	61999969
61967904	62166737	62124846	62169053	NEXGEN LPS OPTION FEMORAL, SIZE H-LT		61809707	62000965
61967905	62169056	62131280	62169054	<u>Item Number: 00-5996-018-51</u>	<u>Item Number: 00-5996-018-51</u>		61809745
61968858	62173190	NEXGEN LPS OPTION FEMORAL, SIZE D-RT		<u>Lot Number:</u>	61777248	61810398	62010054
NEXGEN KNEE LPS OPTION FEMORAL, SIZE H-RT		NEXGEN LPS OPTION FEMORAL, SIZE D-RT		NEXGEN LPS OPTION FEMORAL, SIZE H-LT		61810398	62010322
<u>Item Number: 00-5996-018-92</u>		<u>Item Number: 00-5996-014-51</u>		<u>Item Number: 00-5996-018-51</u>		61813688	62016669
<u>Lot Number:</u>	61836677	<u>Lot Numbers:</u>	62134711	<u>Lot Number:</u>	61777248	61813689	62016674
NEXGEN LCCK FEMORAL, SIZE C-LT		<u>Lot Numbers:</u>	62124846	NEXGEN LPS POROUS FEMORAL, SIZE F-RT		61826922	62023783
<u>Item Number: 00-5994-013-91</u>		62131280	62169054	<u>Item Number: 00-5992-016-02</u>		61826923	62023784
<u>Lot Number:</u>	61758275	NEXGEN LPS OPTION FEMORAL, SIZE D-RT		<u>Lot Number:</u>	61782589	61827924	62024952
NEXGEN LCCK FEMORAL, SIZE C-RT		<u>Item Number: 00-5996-014-52</u>		NEXGEN LPS POROUS FEMORAL, SIZE G-LT		61827925	62026044
<u>Item Number: 00-5994-013-92</u>		<u>Lot Numbers:</u>	62119552	<u>Item Number: 00-5992-017-91</u>		61827926	62026045
<u>Lot Number:</u>	61758283	11002937	62137519	<u>Lot Number:</u>	61928631	61827927	62026089
NEXGEN LCCK FEMORAL, SIZE D-RT		11003004	62171791	NEXGEN LPS PRECOAT FEMORAL, SIZE D-RT		61829824	62034568
<u>Item Number: 00-5994-014-92</u>		11003326	62171792	<u>Item Number: 00-5990-014-02</u>		61839693	62047521
<u>Lot Numbers:</u>	11006826	61782574	NEXGEN LPS OPTION FEMORAL, SIZE E-LT		<u>Lot Number:</u>	61839694	62051185
11004374	61768114	<u>Item Number: 00-5996-015-51</u>		<u>Item Number: 00-5996-015-51</u>		61840456	62051186
11004966		<u>Lot Numbers:</u>	62170913	<u>Lot Number:</u>		61846248	62051187
		62120709	62170915	61776963		61846249	62074155
						61847652	62077286
						61859212	62077837

NEXGEN LPS-FLEX OPTION FEMORAL, SIZE D-LT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE D-RT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE D-RT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE E-LT	
<u>Item Number: 00-5964-014-51</u>		<u>Item Number: 00-5964-014-52</u>		<u>Item Number: 00-5964-014-52</u>		<u>Item Number: 00-5964-015-51</u>	
<u>Lot Numbers:</u>	61992051	<u>Lot Numbers:</u>	61954540	<u>Lot Numbers:</u>	62154067	<u>Lot Numbers:</u>	61996950
61859214	62079851	61812820	62018047	61965929	62160333	61870713	62059370
61859215	62084425	61814746	62022925	61965930	62161196	61879221	62063717
61866692	62088867	61814747	62022927	61973058	62163070	61881047	62064893
61873428	62089663	61814749	62027189	61973228	62163071	61881048	62065289
61876904	62096004	61814750	62028442	61973229	62163097	61881676	62069327
61876905	62096005	61814751	62028443	61979766	62164813	61881678	62070226
61884182	62096006	61815695	62035834	61979767	62173159	61886313	62070227
61885169	62100460	61829046	62036141	61979768	62173160	61886314	62070820
61886258	62100461	61832272	62041496	61980293	62173162	61886315	62078048
61886259	62102053	61832273	62048233	61980296	62173175	61888026	62079812
61890010	62102695	61832274	62048236	61988446	62173176	61889083	62081419
61890011	62109151	61839689	62052657	61989660	62173178	61889084	62081865
61890690	62109649	61839691	62052658	61989661	62173179	61889085	62086334
61905149	62109650	61847653	62052659	61989662	62173180	61889086	62086604
61905256	62112541	61847654	62076318	61993903		61896039	62086605
61914474	62112685	61847655	62077284	NEXGEN LPS-FLEX OPTION			
61914475	62112972	61860203	62077839	FEMORAL, SIZE E-LT			
61914478	62118561	61860204	62079815	<u>Item Number: 00-5964-015-51</u>			
61914479	62118562	61860205	62084433	<u>Lot Numbers:</u>	61996951	61902941	62096035
61915292	62124294	61860206	62085524	11400142	61996952	61902942	62096720
61915293	62130970	61866693	62085525	11400242	61999573	61904182	62096721
61918961	62130973	61869843	62085941	61820639	62000988	61904183	62096722
61930141	62130974	61870709	62096703	61820640	62000989	61904184	62098857
61931654	62134259	61874450	62098052	61820641	62002051	61904915	62098858
61939789	62134260	61878262	62098053	61821139	62003650	61904916	62098859
61951105	62141164	61878263	62099804	61829067	62003776	61904917	62099077
61951646	62141723	61883260	62099848	61829068	62009308	61906382	62099080
61951647	62142299	61884178	62099860	61829069	62015855	61907433	62099260
61951648	62146620	61886312	62100456	61829070	62015856	61907434	62102872
61964434	62151772	61890007	62103206	61829071	62016807	61918714	62105756
61965568	62152639	61890008	62105721	61829072	62016807	61918715	62105758
61965569	62152897	61890008	62105721	61829846	62017264	61918716	62105759
61965927	62152904	61894531	62105722	61835354	62017265	61918717	62105760
61965928	62161209	61895306	62105723	61835356	62017266	61918958	62113022
61968882	62163072	61895308	62112971	61835357	62017597	61918959	62113023
61972732	62163099	61896038	62112988	61840459	62017599	61918959	62113024
61978248	62170387	61901209	62112989	61840460	62017600	61922244	62113024
61986694	62171766	61901210	62112994	61844225	62027453	61922245	62114104
61987054	62171767	61901211	62114093	61844227	62027465	61923559	62114105
61987055	62171782	61902939	62114094	61844228	62028235	61923560	62114499
61988120	62171783	61906064	62114367	61844885	62028236	61926121	62114503
61988444	62171833	61906381	62119153	61850913	62028237	61927955	62114505
61992034	62171834	61915294	62124787	61850914	62029263	61933492	62119528
		61915295	62125121	61850914	62029263	61933492	62119528
		61920422	62126096	61856735	62029266	61933493	62119529
		61931652	62130977	61856736	62029856	61936376	62119530
		61931653	62131185	61857881	62029877	61943463	62119532
		61936489	62131201	61857882	62030480	61943489	62124608
		61939172	62134704	61858464	62046036	61943490	62124886
		61939175	62134745	61861342	62046037	61947006	62124888
		61949694	62143136	61861343	62046235	61947020	62125119
		61951104	62143139	61861344	62046236	61947181	62131029
		61954537	62148749	61864866	62052711	61947182	62131192
		61954538	62154036	61867593	62054075	61947183	62131211
		61954539	62154066	61870710	62054076	61952862	62131212
				61870711	62059369	61952863	62134264
						61952864	62134265

NEXGEN LPS-FLEX OPTION FEMORAL, SIZE E-LT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE E-RT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE E-RT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE F-LT	
<u>Item Number: 00-5964-015-51</u>		<u>Item Number: 00-5964-015-52</u>		<u>Item Number: 00-5964-015-52</u>		<u>Item Number: 00-5964-016-51</u>	
<u>Lot Numbers:</u>	62142580	<u>Lot Numbers:</u>	61924437	<u>Lot Numbers:</u>	62116037	<u>Lot Numbers:</u>	61990661
61960816	62143148	61844929	62029558	61926116	62116038	11005907	61994834
61962852	62143154	61851713	62029559	61926117	62117292	11400319	61994838
61972318	62143616	61853363	62030482	61926119	62117294	61783049	61994843
61974115	62143617	61853364	62030620	61926120	62118407	61783050	61996971
61974117	62150842	61853365	62030622	61933827	62121213	61783051	62000983
61976406	62150849	61854839	62030623	61933829	62121214	61784877	62000984
61976407	62151225	61855955	62043314	61934948	62121217	61784878	62002023
61976408	62151230	61856738	62046906	61934950	62121240	61815795	62002063
61976409	62152196	61856739	62047528	61939791	62121241	61817742	62002077
61977235	62152215	61862123	62048276	61940348	62125513	61817744	62002079
61982271	62152615	61864868	62048277	61943491	62125515	61817746	62008258
61982282	62160278	61872323	62054077	61946455	62126107	61821407	62010325
61982451	62160279	61873425	62054650	61946456	62132072	61821408	62010327
61982452	62160332	61873426	62054651	61948134	62132076	61821409	62011224
61982453	62161211	61880089	62061868	61949691	62132094	61822939	62018137
61983002	62166751	61880090	62061869	61949692	62132740	61822963	62020797
61987053	62167520	61881030	62061870	61953034	62132872	61822964	62020798
61987073	62168057	61881679	62062236	61953036	62135958	61822965	62024953
61988467	62168350	61881680	62063719	61953037	62136379	61822966	62024954
61988468	62168352	61881681	62064892	61970671	62136380	61822966	62026092
61995684	62169036	61883269	62072145	61970672	62137052	61832304	62026092
61995685	62169044	61883270	62072149	61970673	62144123	61833584	62027191
61995686	62169045	61883270	62072149	61970673	62144123	61833585	62030583
61996350	62169046	61888027	62072202	61971150	62144288	61833588	62030585
61996948		61888028	62072203	61971151	62144312	61835820	62033833
		61888029	62072507	61972319	62144929	61841269	62036524
		61888030	62087012	61972744	62149361	61842336	62041511
		61897280	62087943	61972746	62149362	61844230	62041520
		61897281	62087944	61973249	62149363	61844231	62042334
		61898483	62087946	61973250	62149518	61848276	62042335
		61904185	62087947	61973251	62149533	61859217	62048451
		61904186	62088863	61973252	62149551	61859218	62049418
		61904187	62091538	61973253	62152640	61862126	62049420
		61905145	62091539	61974043	62152911	61862128	62050095
		61905146	62098077	61980344	62152916	61872326	62050096
		61905147	62098078	61982821	62161204	61876276	62061688
		61907436	62098571	61982822	62161222	61876907	62062240
		61907437	62098572	61984203	62162207	61881037	62062242
		61907438	62098804	61984205	62162229	61881038	62070236
		61908278	62098805	61985043	62164877	61885165	62070240
		61917741	62100195	61986758	62164878	61885166	62072204
		61917742	62100213	61987052	62165873	61885167	62077290
		61917743	62100484	61990662	62165945	61885168	62077843
		61917778	62100485	61990663	62166665	61890036	62077844
		61918712	62107942	61990664	62170033	61890759	62082375
		61918713	62107944	61991999	62170035	61890760	62083081
		61920424	62107960	61992001	62170036	61892380	62086987
		61920425	62109184	61992070	62170037	61892382	62087934
		61920426	62109186	61992780	62170388	61899641	62092003
		61920427	62109751	61992781	62170389	61899645	62092257
		61920428	62109754	61992992	62170390	61899646	62092259
		61921432	62110246	61993977	62170391	61900470	62092260
		61921433	62110796			61901204	62099275
		61921434	62116035			61901206	62099276
		61921435	62116036			61901207	62099941

NEXGEN LPS-FLEX OPTION FEMORAL, SIZE F-LT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE F-RT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE F-RT		NEXGEN LPS-FLEX OPTION FEMORAL, SIZE G-LT	
<u>Item Number: 00-5964-016-51</u>		<u>Item Number: 00-5964-016-52</u>		<u>Item Number: 00-5964-016-52</u>		<u>Item Number: 00-5964-017-51</u>	
<u>Lot Numbers:</u>	61990657	<u>Lot Numbers:</u>	61943183	<u>Lot Numbers:</u>	62121772	<u>Lot Numbers:</u>	62006726
61901412	62100478	61821125	62014396	61952877	62121794	61909456	62096003
61902961	62100480	61823875	62014397	61953086	62126683	61910862	62096697
61909333	62100481	61823877	62018587	61953089	62127080	61912346	62107981
61909467	62101679	61830375	62023296	61953093	62132450	61912347	62115243
61909471	62101704	61834855	62023797	61953095	62132454	61918718	62115255
61916144	62102075	61834856	62023798	61953096	62132862	61923565	62115256
61917744	62112032	61837315	62023799	61960708	62132876	61938744	62117299
61917746	62112982	61843117	62024561	61962876	62138238	61939177	62117301
61917747	62112983	61843119	62024570	61962879	62143738	61939178	62122655
61917748	62113013	61844583	62024960	61965932	62144287	61943517	62126686
61918739	62115992	61844584	62024961	61965933	62145398	61949697	62131222
61920464	62117311	61855951	62033086	61971720	62146935	61966236	62131246
61923587	62117313	61855952	62034601	61971721	62146956	61966237	62132055
61926122	62117314	61857886	62039091	61971723	62148238	61970679	62141170
61926125	62117315	61857887	62044629	61980301	62148258	61971153	62142287
61931662	62120595	61857888	62044630	61980302	62154272	61976411	62148755
61931663	62124289	61869844	62046964	61980569	62154957	61977546	62154269
61936787	62124759	61869846	62048290	61980570	62154958	61982454	62154994
61943495	62125260	61876908	62048291	61980571	62156819	61988126	62159252
61943496	62131291	61876909	62056792	61980572	62157853	61990911	62168117
61946954	62132080	61878268	62057299	61982291	62166739	61990912	62169033
61951655	62132098	61883271	62057306	61989035	62167493	61990913	62169034
61952869	62132101	61883272	62058232	61989665	62171785	61992052	62170490
61959249	62137510	61884190	62059372	61996344	62171837	62000967	62170491
61959250	62143614	61884192	62059398	61996345	62171838	62000968	62170492
61959252	62143654	61885768	62067920	61997333	62171839	NEXGEN LPS-FLEX OPTION FEMORAL, SIZE G-RT	
61960706	62143734	61888047	62072206	61997336	62171840		
61960707	62150066	61888048	62072511	61997337	62172413	<u>Item Number: 00-5964-017-52</u>	
61970688	62150068	61895303	62077288	61998614		<u>Lot Numbers:</u>	61998800
61970690	62150087	61895304	62078816	NEXGEN LPS-FLEX OPTION FEMORAL, SIZE G-LT		11400311	62003765
61971186	62154039	61895305	62078817	<u>Item Number: 00-5964-017-51</u>		11400341	62004792
61971187	62154068	61901214	62083121	<u>Lot Numbers:</u>		61807817	62004794
61971739	62154069	61901215	62083123	61823882	62012394	61820636	62004993
61971740	62159325	61902943	62083124	61823883	62012654	61822968	62014073
61971741	62165929	61902944	62083125	61829831	62017256	61823880	62014074
61971742	62165930	61902945	62098280	61842343	62018139	61830363	62014387
61976432	62166647	61902946	62098667	61844892	62018140	61843123	62019486
61978516	62166649	61902947	62098670	61861346	62028013	61844893	62019487
61985037	62170884	61904188	62098786	61865905	62029429	61844931	62019489
61986240	62170890	61904189	62098787	61874455	62035837	61860211	62029441
61986243	62170893	61904190	62098788	61874456	62042284	61866703	62029446
61986245	62171085	61904191	62098789	61881661	62043132	61875434	62038010
61986733	62171086	61915298	62102055	61881663	62044610	61875435	62043281
61987070	62171726	61915301	62102056	61881663	62054831	61883261	62043283
61988463	62171730	61915302	62102057	61886260	62058757	61883262	62046217
NEXGEN LPS-FLEX OPTION FEMORAL, SIZE F-RT		61922249	62102712	61886261	62059371	61886317	62055031
<u>Item Number: 00-5964-016-52</u>		61922250	62102714	61886262	62068311	61886318	62055523
<u>Lot Numbers:</u>	61998615	61923561	62114023	61892362	62073097	61892359	62058211
11006368	62006438	61927959	62114334	61893239	62074152	61892360	62068306
61816395	62006449	61932281	62114346	61898485	62085944	61899616	62074167
61816396	62006465	61937792	62115228	61898487	62085945	61899617	62076154
61820650	62014083	61940950	62119151	61905152	62086606	61899618	62086981
61820651	62014084	61943181	62119152	61905155	62092518	61901208	62086983
		61943182	62119534	61905169	62093708	61905263	62098049

NEXGEN LPS-FLEX OPTION FEMORAL, SIZE G-RT <u>Item Number: 00-5964-017-52</u> <u>Lot Numbers:</u> 61998799 62098051 61909460 62098278 61910185 62106732 61910186 62106733 61910189 62115851 61910193 62115852 61924441 62117296 61935976 62117297 61935978 62122860 61939181 62127394 61939795 62132078 61946486 62132108 61949695 62138255 61954543 62140623 61955558 62140628 61966806 62149037 61966807 62149365 61966808 62155191 61978256 62156806 61978264 62159249 61982458 62168118 61986416 62169042 61986695 62169043 61992981 62170028 61993919 62170029 61993921 62170029	NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE C-RT <u>Item Number: 00-5960-013-52</u> <u>Lot Numbers:</u> 62120588 62130963 NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE D-LT <u>Item Number: 00-5960-014-51</u> <u>Lot Numbers:</u> 62120591 11006176 62131197 11006333 62134257 11400330 62170032 NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE D-RT <u>Item Number: 00-5960-014-52</u> <u>Lot Numbers:</u> 62131248 11006544 62134703 62012423 62171835 62061891 62174384 62121191 NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE E-LT <u>Item Number: 00-5960-015-51</u> <u>Lot Numbers:</u> 62134706 62121224 62168365 62131267 NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE E-RT <u>Item Number: 00-5960-015-52</u> <u>Lot Numbers:</u> 62121770 61641772 62131277 62121768 62135957 NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE F-LT <u>Item Number: 00-5960-016-51</u> <u>Lot Numbers:</u> 62121792 11400431 62131249 61784879 62131284 61793525 62171784 62121790 NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE F-RT <u>Item Number: 00-5960-016-52</u> <u>Lot Numbers:</u> 62122521 11005960 62133156 61794705 62137533 NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE G-LT <u>Item Number: 00-5960-017-51</u> <u>Lot Numbers:</u> 62131221 114691 62137511 62122653 62168367 NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE G-RT <u>Item Number: 00-5960-017-52</u> <u>Lot Numbers:</u> 62128361 62122857 62132861	NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE H-LT <u>Item Number: 00-5960-018-51</u> <u>Lot Numbers:</u> 11005166 11004155 11005910 NEXGEN LPS-FLEX PRECOAT FEMORAL SIZE H-RT <u>Item Number: 00-5960-018-52</u> <u>Lot Numbers:</u> 11005164 11004157 11005165 NKII NP FEMORAL LT, SIZE 0 <u>Item Number: 6307-00-008</u> <u>Lot Number:</u> 62125578 NKII NP FEMORAL LT, SIZE 2 <u>Item Number: 6307-00-020</u> <u>Lot Numbers:</u> 61975784 62098206 NKII NP FEMORAL LT, SIZE 3 <u>Item Number: 6307-00-030</u> <u>Lot Numbers:</u> 11004845 62142198 NKII NP FEMORAL RT, SIZE 0 <u>Item Number: 6307-00-009</u> <u>Lot Number:</u> 62107229 NKII NP FEMORAL RT, SIZE 2 <u>Item Number: 6307-00-021</u> <u>Lot Numbers:</u> 62178731 11004878 62180483 11007656 62181538 11007666 62183228 11007677 62185696 11007719 62188142 11007831 62189950 11007922 62190389 11007925 62191204 11007944 62192798 11008018 62193587 62155319 62194526 62160819 62196210 62162357 62200791 62164084 62203902 62172373 62206609 62173770 62207438 62174770 62208953 62176798 62210107 62177132 62210574 62178184	PRCT MIS CR MOB PLT SZ 6 <u>Item Number: 00-5930-060-01</u> <u>Lot Number:</u> 61774234 PRCT MIS CR MOB PLT SZ 7 <u>Item Number: 00-5930-070-01</u> <u>Lot Number:</u> 62139076 PRCT MIS CR MOB PLT SZ 7 <u>Item Number: 00-5930-070-02</u> <u>Lot Number:</u> 61780002 PRCT MIS CR MOB PLT SZ 8 <u>Item Number: 00-5930-080-01</u> <u>Lot Numbers:</u> 61774229 61780004 PRCT MIS CR MOB PLT SZ 8 <u>Item Number: 00-5930-080-02</u> <u>Lot Numbers:</u> 61774216 61774215 61785662 ZUK HIGH FLEX PRECOAT FEM SZ B,(RT/M,LT/L) <u>Item Number: 00-5842-012-02</u> <u>Lot Numbers:</u> 61774093 61777606 ZUK HIGH FLEX PRECOAT FEM SZ C,(LT/M,RT/L) <u>Item Number: 00-5842-013-01</u> <u>Lot Number:</u> 61774095 ZUK HIGH FLEX PRECOAT FEM SZ C,(RT/M,LT/L) <u>Item Number: 00-5842-013-02</u> <u>Lot Number:</u> 61777611 ZUK HIGH FLEX PRECOAT FEM SZ D,(LT/M,RT/L) <u>Item Number: 00-5842-014-01</u> <u>Lot Numbers:</u> 11002921 61777613 ZUK HIGH FLEX PRECOAT FEM SZ D,(RT/M,LT/L) <u>Item Number: 00-5842-014-02</u> <u>Lot Numbers:</u> 61777615 61777616 ZUK HIGH FLEX PRECOAT FEM SZ E,(LT/M,RT/L) <u>Item Number: 00-5842-015-01</u> <u>Lot Number:</u> 61777617 ZUK HIGH FLEX PRECOAT FEM SZ F,(RT/M,LT/L) <u>Item Number: 00-5842-016-02</u> <u>Lot Numbers:</u> 11002915 61756843
		PRCT MIS CR MOB PLT SZ 4 <u>Item Number: 00-5930-040-01</u> <u>Lot Numbers:</u> 61774225 11005139 61774226 61774223 61780003 PRCT MIS CR MOB PLT SZ 5 <u>Item Number: 00-5930-050-01</u> <u>Lot Numbers:</u> 11004031 61774236	