

URGENT MEDICAL DEVICE RECALL NOTIFICATION
LFIT™ Anatomic CoCr V40™ Femoral Heads

Oct 13th, 2016

FSCA identifier: Product Field Action RA2016-028

Type of Action: Field Safety Corrective Action: **Recall**

Description: LFIT™ Anatomic CoCr V40™ Femoral Heads

Catalog #: 6260-9-236, 6260-9-240, 6260-9-244, 6260-9-340, 6260-9-344, 6260-9-440, 6260-9-444

License #: DE0003747 and DE0008048

Lot Code: 30647803, 34354901, 82YMEA, 9K0MHA, E32MKA, MHNX13, MHNX16, MJE70V
See attached list for globally affected lot numbers

Product owner: Howmedica Osteonics Corp., 325 Corporate Drive, Mahwah, NJ 07430, USA

License holder: Stryker Singapore Pte Ltd, 108 Pasir Panjang Road #03-04 Golden Agri Plaza Singapore 118535

Dear Hospital Representative,

Stryker has initiated in Singapore a voluntary medical device recall for the following Femoral Heads.

Catalog Number	Head Diameter	Offset
6260-9-236	36mm	+5
6260-9-240	40mm	+4
6260-9-244	44mm	+4
6260-9-340	40mm	+8
6260-9-440	40mm	+12
6260-9-344	44mm	+8
6260-9-444	44mm	+12

The intent of this letter is to describe all potential hazards associated with the below noted issue, and any risk mitigation factors associated with the use of the product.

Our records indicate that you have received the above referenced product. It is Stryker's responsibility as the manufacturer to ensure that customers who may have received these affected products also receive this important communication.

Reason for the Voluntary Recall:

Stryker has received higher than expected complaints of taper lock failure for specific lots of the affected sizes of LFIT™ Anatomic CoCr V40™ Femoral Heads manufactured prior to 2011.

Potential Hazards:

The potential hazards may include:

- Disassociation of femoral head from hip stem
- Fractured hip stem trunnion
- Excessive metallic debris
- Insufficient ROM
- Insufficient soft tissue tension
- Noise
- Loss of implant: bone fixation strength
- Excessive wear debris (polymeric)

- Implant construct with a shortened neck length

The aforementioned potential hazards may result in one or more of the following potential patient harms:

- User annoyance
- Loss of mobility
- Pain requiring revision
- Inflammatory response
- Adverse local tissue reaction
- Dislocation
- Joint instability
- Revision to alleviate hazardous situation
- Pain associated with implant loosening
- Periprosthetic fracture
- Leg length discrepancy

Follow up:

Implanted patients with LFIT™ Anatomic CoCr V40™ Femoral Heads as described above should continue to be followed per the normal protocol established by his/her surgeon.

Required actions:

- Please inform users of this Urgent Medical Device Recall Notification and forward this notice to all those individuals who need to be aware within your organization.
- Complete and sign the enclosed Business Reply Form and ASEAN.PMS@stryker.com or hand to your Stryker sales representative.

Our records indicate that you have received the above referenced product. It is our responsibility to ensure that customers who may have received this affected product also receive this important communication. **Please assist us in meeting our regulatory obligation by signing and returning the attached Business Reply Form within 5 days of receipt of this letter.**


Reporting of Adverse Events

Please report any adverse events or product quality problems associated with this device to Stryker. Healthcare Professionals may also report any suspected adverse events associated with these devices to the Vigilance and Compliance Branch, Health Products Regulation Group, HSA at Tel: 6866 3538, Fax: 6478 9069, or report online at www.hsa.gov.sg/ae online.

For patient questions, Stryker has established a dedicated call center at +1 978-338-3143.

We regret any inconvenience this action may cause you and if you have any questions, feel free to contact your Stryker representative or the signatory below.

Sincerely,



Sara Jato
QA Specialist
Stryker ASEAN
sara.jato@stryker.com
ASEAN.PMS@stryker.com
Tel. +84 (0)8 3827 5399 Ext.12

URGENT MEDICAL DEVICE RECALL NOTIFICATION
LFIT™ Anatomic CoCr V40™ Femoral Heads

Oct 13th, 2016

FSCA identifier: Product Field Action RA2016-028

Type of Action: Field Safety Corrective Action: **Recall**

Description: LFIT™ Anatomic CoCr V40™ Femoral Heads

Catalog #: 6260-9-236, 6260-9-240, 6260-9-244, 6260-9-340, 6260-9-344, 6260-9-440, 6260-9-444

License #: DE0003747 and DE0008048

Lot Code: 30647803, 34354901, 82YMEA, 9K0MHA, E32MKA, MHNX13, MHNX16, MJE70V
See attached list for globally affected lot numbers

Product owner: Howmedica Osteonics Corp., 325 Corporate Drive, Mahwah, NJ 07430, USA

License holder: Stryker Singapore Pte Ltd, 108 Pasir Panjang Road #03-04 Golden Agri Plaza
Singapore 118535

I have received the product recall letter from Stryker Orthopaedics dated Oct 13th, 2016 stating that the company has initiated a voluntary lot-specific recall of the above referenced products.

Hospital Name

Date

Hospital Representative
(Name and signature)

Hospital stamp

PLEASE COMPLETE THIS FORM WITHIN 5 BUSINESS DAYS AND RETURN IT VIA EMAIL
(ASEAN.PMS@stryker.com) OR HAND TO YOUR STRYKER SALES REPRESENTATIVE

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-244	000MHD	6260-9-244	0KXMNAX1	6260-9-240	1RTMLD	6260-9-240	2WJMMD	6260-9-240	LP7MLE
6260-9-244	001MHD	6260-9-244	0KYMNA	6260-9-240	1THMLD	6260-9-244	2WKMAE	6260-9-240	LP8MLE
6260-9-244	002MHD	6260-9-244	0KYMNAX1	6260-9-244	1X6MRA	6260-9-240	2WYMMD	6260-9-240	LP9MLE
6260-9-244	004MHD	6260-9-244	0KYMNAX2	6260-9-244	1X6MRAX1	6260-9-244	2X2MAE	6260-9-240	LPAMED
6260-9-340	004MLD	6260-9-244	0M1MNA	6260-9-244	1X7MRA	6260-9-240	30MMMD	6260-9-240	LRWMTA
6260-9-244	005MHD	6260-9-444	0NLMAD	6260-9-244	1X8MRA	6260-9-240	30MMMDX1	6260-9-240	LRXMTA
6260-9-340	005MLD	6260-9-444	0NNMAD	6260-9-244	1X8MRAX1	6260-9-240	30NMMD	6260-9-240	LRYMTA
6260-9-340	006MLD	6260-9-444	0NNMAD	6260-9-240	1X9MMD	6260-9-244	31RMJD	6260-9-240	LRYMTAX1
6260-9-340	00AMPE	6260-9-340	0PHMLD	6260-9-244	1X9MRA	6260-9-244	31TMJD	6260-9-240	LT0MTA
6260-9-340	00DMJE	6260-9-340	0THMLD	6260-9-244	1X9MRAX1	6260-9-240	31TMTD	6260-9-240	LT1MTA
6260-9-240	00JMKE	6260-9-244	0VOMLA	6260-9-244	1X9MRAX1X1	6260-9-244	31VMJD	6260-9-340	LT3MEE
6260-9-240	00RMLA	6260-9-244	0V8MAE	6260-9-244	1XAMRA	6260-9-240	31VMTD	6260-9-340	LT5MEE
6260-9-340	00TMLA	6260-9-340	0Y5MHE	6260-9-244	1XAMRAX1	6260-9-240	31WMTD	6260-9-340	LT6MEE
6260-9-240	017MKE	6260-9-340	0Y9MHE	6260-9-240	1XDMMD	6260-9-240	321MJD	6260-9-340	LT7MEE
6260-9-340	018MPE	6260-9-340	0Y9MHX1	6260-9-244	20JMJD	6260-9-240	323MJD	6260-9-240	LWVMND
6260-9-240	019MKE	6260-9-244	0Y9MLA	6260-9-244	20JMJD	6260-9-240	353MNE	6260-9-240	LWWMND
6260-9-340	01AMPE	6260-9-240	0YAMHE	6260-9-244	20KMJD	6260-9-240	354MNE	6260-9-240	LWXMND
6260-9-340	01DMJE	6260-9-240	0YDMHE	6260-9-240	20KMMD	6260-9-240	355MNE	6260-9-240	LWYMND
6260-9-340	01DMPE	6260-9-244	0YDMLA	6260-9-244	20LMJD	6260-9-240	357MNE	6260-9-240	LX3MND
6260-9-240	02NMKE	6260-9-240	0YEMHE	6260-9-244	20MMJD	6260-9-240	376MRD	6260-9-240	LX4MND
6260-9-240	02PMKE	6260-9-240	0YHMHE	6260-9-240	211MLE	6260-9-240	37RMMD	6260-9-340	LX9MKD
6260-9-340	02VMRA	6260-9-340	111MPE	6260-9-240	22AMLE	6260-9-240	396MRD	6260-9-340	LY4MLE
6260-9-340	02WMRA	6260-9-244	115MHD	6260-9-240	22WMLE	6260-9-240	3A9MMD	6260-9-340	LY5MLE
6260-9-240	038MRA	6260-9-244	115MPE	6260-9-240	23WMLE	6260-9-240	3J9MMD	6260-9-340	LY6MLE
6260-9-240	03AMRA	6260-9-244	116MHD	6260-9-240	23XMMD	6260-9-240	3JNMMD	6260-9-340	LY7MLE
6260-9-240	03AMRAX1	6260-9-340	116MPE	6260-9-240	241MLE	6260-9-240	3LDMTD	6260-9-240	LY9MED
6260-9-240	03AMRAX2	6260-9-244	117MHD	6260-9-240	242MLE	6260-9-240	3LEMTD	6260-9-240	LYHMED
6260-9-240	03EMKE	6260-9-244	118MHD	6260-9-244	250MHD	6260-9-240	3LHMTD	6260-9-240	M1AMED
6260-9-240	03LMKE	6260-9-340	11DMPE	6260-9-244	251MHD	6260-9-240	3P1MRA	6260-9-240	M1LMED
6260-9-240	03XMEE	6260-9-340	11JMPE	6260-9-244	252MHD	6260-9-240	3P1MRAX1	6260-9-240	M1NMED
6260-9-240	04PMRA	6260-9-244	11KMPE	6260-9-244	253MHD	6260-9-240	3P1MRAX2	6260-9-240	M2HMED
6260-9-244	05NMAD	6260-9-240	12DMLE	6260-9-244	254MHD	6260-9-240	3P5MRA	6260-9-240	M30MNA
6260-9-344	06PMAD	6260-9-240	12EMLE	6260-9-240	26RMNE	6260-9-340	3P7MRA	6260-9-240	M30MNAX1
6260-9-344	06RMAD	6260-9-240	12HMLE	6260-9-240	27MMPA	6260-9-340	3P9MRA	6260-9-240	M37MPA
6260-9-344	06TMAD	6260-9-340	12JMKE	6260-9-340	27NMMPA	6260-9-340	3P9MRAX1	6260-9-240	M37MPAX1
6260-9-340	06TMLD	6260-9-240	12JMLE	6260-9-340	27PMPA	6260-9-340	3PAMRA	6260-9-340	M3RMPA
6260-9-344	06VMAD	6260-9-240	12NMMLD	6260-9-340	27RMPA	6260-9-240	3RHMMD	6260-9-340	M3RMPAX1
6260-9-340	06WMMLD	6260-9-340	14WMKE	6260-9-240	285MJE	6260-9-240	3TRMRD	6260-9-340	M3WMPA
6260-9-244	071MPE	6260-9-244	14YMDE	6260-9-240	286MJE	6260-9-244	3W7MHE	6260-9-240	M4DMPA
6260-9-244	073MPE	6260-9-244	150MDE	6260-9-340	289MJE	6260-9-240	3WKMMMD	6260-9-240	M4DMPAX1
6260-9-244	074MPE	6260-9-244	151MDE	6260-9-340	293MJE	6260-9-240	3YDMJD	6260-9-340	M4LMPA
6260-9-240	08NMKE	6260-9-244	152MDE	6260-9-340	296MJE	6260-9-240	3YEMJD	6260-9-444	M4RMAD
6260-9-340	08PMJE	6260-9-244	161MME	6260-9-244	297MJE	6260-9-240	3YHMJD	6260-9-340	M5NMTA
6260-9-340	08PMJEX1	6260-9-240	18NMDE	6260-9-244	298MJE	6260-9-240	3YJMJD	6260-9-340	M5PMTA
6260-9-340	0A9MHE	6260-9-344	1AKMAD	6260-9-244	299MJE	6260-9-240	3YLMJD	6260-9-340	M5TMTA
6260-9-340	0AHMJE	6260-9-240	1EPMLE	6260-9-244	2A1MJE	6260-9-240	3YPMJD	6260-9-444	M6EMAD
6260-9-340	0DLMKE	6260-9-240	1ETMLE	6260-9-240	2EMMMD	6260-9-240	3YVMJD	6260-9-444	M6EMADX1
6260-9-240	0DMMKE	6260-9-244	1H1MAE	6260-9-240	2LTMMMD	6260-9-240	400MJD	6260-9-240	M6PMED
6260-9-340	0DNMKE	6260-9-244	1J0MAE	6260-9-240	2LTMMDX1	6260-9-240	41WMAE	6260-9-240	M6PMRA
6260-9-340	0HJMHE	6260-9-244	1J2MAE	6260-9-240	2N1MMD	6260-9-240	41WMMMD	6260-9-240	M6RMED
6260-9-244	0J3MNA	6260-9-240	1KTMDE	6260-9-244	2TEMMAE	6260-9-340	42DMMD	6260-9-240	M6TMED
6260-9-244	0J3MNAX1	6260-9-240	1KVMDE	6260-9-244	2THMAE	6260-9-340	42JMJD	6260-9-240	M6XMRA
6260-9-244	0JWMNA	6260-9-240	1KWMDE	6260-9-244	2TJMAE	6260-9-340	42KMJD	6260-9-240	M72MED
6260-9-244	0KXMNA	6260-9-240	1P2MDE	6260-9-244	2TKMAE	6260-9-340	42LMJD	6260-9-240	M73MRA
6260-9-340	43NMJD	6260-9-340	55RMNDX1	6260-9-244	66YMKE	6260-9-340	72RMPA	6260-9-240	M75MRA
6260-9-340	43PMJD	6260-9-244	55RMRA	6260-9-244	675MMA	6260-9-340	73XMDE	6260-9-240	M7RMRA
6260-9-340	43RMJD	6260-9-340	55WMPD	6260-9-244	676MMA	6260-9-244	743MMA	6260-9-240	M7VMRA
6260-9-240	43YMTD	6260-9-340	55YMPD	6260-9-244	677MMA	6260-9-244	744MMA	6260-9-240	M7XMRA
6260-9-240	441MTD	6260-9-340	55YMPDX1	6260-9-244	677MMAX1	6260-9-244	746MMA	6260-9-240	M7XMRAX1
6260-9-240	44HMRD	6260-9-340	560MPD	6260-9-244	678MMA	6260-9-244	747MMA	6260-9-240	M81MRA

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-244	44NMHE	6260-9-340	560MPDX1	6260-9-244	678MMAX2	6260-9-244	747MMAX1	6260-9-240	M81MRAX1
6260-9-240	453MRD	6260-9-340	56HMRD	6260-9-244	67XMKE	6260-9-240	74YMNE	6260-9-240	M84MRA
6260-9-240	464MRD	6260-9-340	57YMND	6260-9-244	68NMKE	6260-9-240	752MNE	6260-9-240	M85MRA
6260-9-244	47LMHE	6260-9-340	59EMND	6260-9-240	694MED	6260-9-240	779MNE	6260-9-444	MD5MAD
6260-9-244	49RMLE	6260-9-340	59EMNDX1	6260-9-240	695MED	6260-9-240	77EMPA	6260-9-444	MD5MADX1
6260-9-244	49TMLE	6260-9-340	5A4MRD	6260-9-240	696MED	6260-9-240	77HMPA	6260-9-444	MDDMAD
6260-9-244	49VMLE	6260-9-244	5JAMKE	6260-9-240	69PMPA	6260-9-240	77HMPAX1	6260-9-444	MDEMAD
6260-9-244	49XMLE	6260-9-340	5JXMRD	6260-9-240	69RMPA	6260-9-240	77JMPA	6260-9-444	MDHMDAD
6260-9-244	49XMLEX1	6260-9-240	5K8MPA	6260-9-240	69TMPA	6260-9-240	791MTA	6260-9-444	MDJMDAD
6260-9-240	4ALMED	6260-9-240	5K9MPA	6260-9-240	69TMPAX1	6260-9-240	791MTAX1	6260-9-240	MEPDR8
6260-9-240	4D8MED	6260-9-240	5KAMPA	6260-9-340	6DYMPA	6260-9-240	792MTA	6260-9-240	MEPDT4
6260-9-240	4DAMED	6260-9-240	5KAMPAX1	6260-9-340	6E0MPA	6260-9-240	793MTA	6260-9-240	MEPDT4X1
6260-9-340	4DDMED	6260-9-244	5L4MKE	6260-9-340	6E1MPA	6260-9-240	794MPD	6260-9-240	MEPDXE
6260-9-340	4DEMED	6260-9-340	5LTMEE	6260-9-340	6E2MPA	6260-9-240	794MTA	6260-9-240	MEPDY8
6260-9-340	4DHMED	6260-9-340	5LVMEE	6260-9-240	6J8MAD	6260-9-240	795MTA	6260-9-240	MEPDY9
6260-9-340	4DJMED	6260-9-340	5LWMEE	6260-9-240	6J8MADX1	6260-9-240	7E4MPD	6260-9-240	MEPE04
6260-9-240	4DVMPD	6260-9-340	5LXMEE	6260-9-240	6JAMAD	6260-9-240	7H3MKD	6260-9-240	MEPEOK
6260-9-240	4DVMPDX1	6260-9-240	5M3MPA	6260-9-340	6JTMDE	6260-9-244	7HLMMA	6260-9-244	MEPHD4
6260-9-244	4DWMHE	6260-9-240	5M3MPAX1	6260-9-340	6K1MKE	6260-9-244	7HMMMA	6260-9-244	MEPHLJ
6260-9-240	4DWMPD	6260-9-240	5M4MPA	6260-9-340	6K3MPA	6260-9-240	7HTMKD	6260-9-244	MEPHM5
6260-9-240	4DXMPD	6260-9-240	5M4MPAX1	6260-9-240	6RLMPD	6260-9-244	7JYMMA	6260-9-244	MEPHM5X1
6260-9-240	4HHMAE	6260-9-240	5M5MPA	6260-9-340	6T3MAE	6260-9-244	7K0MMA	6260-9-244	MEPHT5
6260-9-340	4NJMND	6260-9-240	5M3MND	6260-9-340	6T8MAE	6260-9-244	7K1MMA	6260-9-240	MEPJA0
6260-9-340	4NJMNDX1	6260-9-240	5MKMND	6260-9-340	6TMMAE	6260-9-244	7K1MMAX1	6260-9-240	MEPJD6
6260-9-240	4PKMAE	6260-9-340	5MVMRD	6260-9-244	6VPMRA	6260-9-244	7K2MMA	6260-9-240	MEPJEW
6260-9-240	4V6MAE	6260-9-240	5N3MND	6260-9-244	6VPMRAX1	6260-9-244	7K2MMAX1	6260-9-240	MEPJEWX1
6260-9-340	4XKMND	6260-9-240	5N4MND	6260-9-244	6VRMRA	6260-9-244	7K3MMA	6260-9-240	MEPJ8
6260-9-340	4XKMNDX1	6260-9-240	5N6MND	6260-9-244	6VRMRAX1	6260-9-244	7K4MMA	6260-9-240	MEPJM1
6260-9-344	513MTD	6260-9-240	5N7MND	6260-9-244	6VWMRA	6260-9-244	7K5MMA	6260-9-240	MEPJX0
6260-9-344	514MTD	6260-9-244	5RMMKE	6260-9-244	6VWMRAX1	6260-9-244	7K6MMA	6260-9-240	MEPJX1
6260-9-344	515MTD	6260-9-244	5V1MKE	6260-9-240	6W9MNE	6260-9-240	7L7MPD	6260-9-240	MEPL42
6260-9-344	516MTD	6260-9-244	5V2MKE	6260-9-240	6WNNME	6260-9-244	7N8MMA	6260-9-240	MEPL57
6260-9-340	529MND	6260-9-244	5V3MKE	6260-9-240	6WPMME	6260-9-244	7N9MMA	6260-9-240	MEPLA8
6260-9-340	529MNDX1	6260-9-240	617MKD	6260-9-240	6WRMME	6260-9-244	7N9MMAX1	6260-9-240	MEPLAJ
6260-9-444	52YMTD	6260-9-244	61WMNE	6260-9-240	6WTMME	6260-9-244	7NHMMA	6260-9-240	MEPLKE
6260-9-444	530MTD	6260-9-244	638MNE	6260-9-340	6YDMTA	6260-9-244	7NLMMMA	6260-9-240	MEPLKEX1
6260-9-340	53XMPD	6260-9-244	639MNE	6260-9-340	6YMMTA	6260-9-244	7NMMMA	6260-9-340	MEPLH
6260-9-340	54VMPD	6260-9-244	63AMNE	6260-9-340	6YNMTA	6260-9-244	7NPMMA	6260-9-340	MEPLH
6260-9-340	54WMPD	6260-9-244	63DMNE	6260-9-340	6YPMTA	6260-9-244	7NPMMA	6260-9-340	MEPLH
6260-9-340	54WMPDX1	6260-9-240	63WMTA	6260-9-340	6YRMTA	6260-9-240	7NTMKD	6260-9-340	MEPLM0
6260-9-340	54XMND	6260-9-240	63WMTAX1	6260-9-240	707MPD	6260-9-240	7RNMLD	6260-9-340	MEPLMP
6260-9-340	54XMNDX1	6260-9-240	63XMTA	6260-9-240	709MPD	6260-9-240	7RNMLDX1	6260-9-240	MEPR2P
6260-9-244	557MKE	6260-9-240	640MTA	6260-9-240	709MPDX1	6260-9-240	7RVMLD	6260-9-240	MEPR2V
6260-9-244	55NMRA	6260-9-240	641MTA	6260-9-340	72JMDE	6260-9-240	7RVMLDX1	6260-9-240	MEPR3L
6260-9-244	55PMRA	6260-9-240	642MTA	6260-9-340	72KMDE	6260-9-240	7RVMLDX1X1	6260-9-240	MEPR3LX1
6260-9-244	55PMRAX1	6260-9-340	664MDE	6260-9-340	72LMDE	6260-9-240	7V0MPD	6260-9-240	MEPR3Y
6260-9-244	55PMRAX1X1	6260-9-244	66PMKE	6260-9-340	72MMDE	6260-9-240	7V1MPD	6260-9-340	MEPT7N
6260-9-340	55RMND	6260-9-244	66XMKE	6260-9-340	72PMPA	6260-9-240	7V2MPD	6260-9-340	MEPT84
6260-9-240	7V3MPD	6260-9-340	8K4MKDX1	6260-9-244	946MND	6260-9-244	9NNMKD	6260-9-340	MEPT8T
6260-9-240	7X8MMD	6260-9-340	8K5MKD	6260-9-244	950MPA	6260-9-244	9NPMKD	6260-9-340	MEPT8V
6260-9-240	7X9MMD	6260-9-340	8K6MKD	6260-9-244	959MND	6260-9-240	9P0MKD	6260-9-340	MEPT94
6260-9-240	7X9MPA	6260-9-340	8K7MKD	6260-9-240	95HMDD	6260-9-240	9P1MKD	6260-9-240	MEPVNY
6260-9-240	7XAMMD	6260-9-240	8LJMEE	6260-9-240	95HMDDX1	6260-9-244	9PAMKD	6260-9-240	MEPV7
6260-9-240	7XAMPA	6260-9-240	8P4MEE	6260-9-240	95JMDD	6260-9-244	9PAMKDX1	6260-9-240	MEPVR6
6260-9-340	7XDMPA	6260-9-240	8P6MEE	6260-9-240	95KMDD	6260-9-340	9R2MDD	6260-9-240	MEPVW0
6260-9-340	7XEMPA	6260-9-240	8P7MEE	6260-9-240	95LMDD	6260-9-340	9R2MDDX1	6260-9-240	MEPVW0X1
6260-9-340	7XEMPAX1	6260-9-240	8P8MEE	6260-9-240	95MMDD	6260-9-340	9R8MDD	6260-9-240	MER4KA
6260-9-240	7XVMMD	6260-9-240	8P9MEE	6260-9-240	95NMDD	6260-9-340	9RHMD	6260-9-244	MER4LR
6260-9-340	7XXMMD	6260-9-244	8P9MKD	6260-9-340	98DMME	6260-9-340	9RWMD	6260-9-244	MER4RL
6260-9-240	7Y6MME	6260-9-240	8PAMEE	6260-9-244	99RLA	6260-9-244	9V2MRA	6260-9-244	MER4RLX1

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-240	7YJMKD	6260-9-244	8PAMKD	6260-9-240	9A8MJE	6260-9-244	9V3MRA	6260-9-240	MER4RT
6260-9-340	80LMMD	6260-9-244	8PAMKDX1	6260-9-240	9ADMJE	6260-9-244	9V4MRA	6260-9-244	MER4V2
6260-9-340	816MMD	6260-9-340	8PDMKD	6260-9-240	9AKMJE	6260-9-240	9V8MKD	6260-9-244	MER4WP
6260-9-340	818MTA	6260-9-340	8PEMKD	6260-9-240	9ANMJE	6260-9-340	9X4MTA	6260-9-240	MER4WR
6260-9-340	81AMTA	6260-9-240	8R1MAE	6260-9-244	9D8MPD	6260-9-340	9X4MTAX1	6260-9-240	MER4WT
6260-9-240	81HMME	6260-9-240	8R2MAE	6260-9-244	9D8MPDX1	6260-9-340	9X6MTA	6260-9-240	MER4WW
6260-9-240	81VMME	6260-9-240	8R3MAE	6260-9-340	9EKMMMA	6260-9-340	9X6MTAX1	6260-9-240	MER4WWX1
6260-9-240	81WMME	6260-9-240	8R4MAE	6260-9-240	9ELMMA	6260-9-340	9X7MTA	6260-9-240	MER4WX
6260-9-240	81XMME	6260-9-340	8RTMKD	6260-9-240	9ELMMAX1	6260-9-340	9X7MTAX1	6260-9-244	MER4YY
6260-9-340	82JMMD	6260-9-340	8RVMKD	6260-9-340	9EMMMA	6260-9-244	A0RMPD	6260-9-240	MER585
6260-9-340	82JMMDX1	6260-9-244	8RWMKD	6260-9-340	9ENMMA	6260-9-244	A0WMPD	6260-9-240	MER5EP
6260-9-240	82NMKD	6260-9-244	8RWMKDX1	6260-9-244	9H1MLA	6260-9-244	A11MKD	6260-9-240	MER5J2
6260-9-340	83LMTA	6260-9-340	8T1MKD	6260-9-244	9H1MLAX1	6260-9-244	A11MKDX1	6260-9-240	MER738
6260-9-340	83MMTA	6260-9-340	8T3MAE	6260-9-244	9H2MLA	6260-9-240	A15MPA	6260-9-240	MER739
6260-9-240	856MJE	6260-9-240	8WVMRA	6260-9-240	9H2MMA	6260-9-240	A16MPA	6260-9-340	MER7HJ
6260-9-240	8AEMKD	6260-9-240	8WVMRAX1	6260-9-240	9HDMMA	6260-9-240	A17MPA	6260-9-340	MER7HP
6260-9-340	8AHMMA	6260-9-240	8WVMRAX2	6260-9-244	9HHMPD	6260-9-244	A1HMPD	6260-9-340	MER7HX
6260-9-340	8AKMMA	6260-9-240	8WVMRA	6260-9-244	9HHMPDX1	6260-9-240	A2MMPA	6260-9-340	MER7J3
6260-9-340	8D6MMA	6260-9-340	8WXMMD	6260-9-240	9J2MMA	6260-9-240	A2MMPAX1	6260-9-244	MERA9M
6260-9-340	8D6MMAX1	6260-9-340	8WYMMD	6260-9-240	9J3MMA	6260-9-240	A2NMPA	6260-9-244	MERAAP
6260-9-240	8DDMJE	6260-9-340	8X0MMD	6260-9-240	9J5MMA	6260-9-240	A2PMPA	6260-9-244	MERAAT
6260-9-340	8DMMMA	6260-9-340	8X0MMDX1	6260-9-240	9J7MMA	6260-9-240	A2RMPA	6260-9-244	MERAAW
6260-9-340	8DMMMAX1	6260-9-340	8X4MMD	6260-9-240	9J8MMA	6260-9-240	A2RMPAX1	6260-9-244	MERAAWX1
6260-9-340	8DNMMA	6260-9-340	8X5MMD	6260-9-240	9J8MMAX1	6260-9-244	A2RMPD	6260-9-244	MERAAAX
6260-9-340	8DPMMA	6260-9-340	8X5MMDX1	6260-9-340	9JRMMA	6260-9-240	A2TMPA	6260-9-244	MERAAAX1
6260-9-340	8DRMMA	6260-9-240	8XWMEE	6260-9-340	9JRMMAAX1	6260-9-340	A2VMPA	6260-9-240	MERAAV
6260-9-340	8DTMMA	6260-9-340	8Y5MMD	6260-9-340	9JTMMMA	6260-9-340	A2VMPAX1	6260-9-240	MERAD1
6260-9-340	8DVMMMA	6260-9-240	927MLD	6260-9-340	9JTMMAAX1	6260-9-340	A2VMPAX2	6260-9-240	MERAD2
6260-9-340	8DWMMMA	6260-9-240	927MLDX1	6260-9-340	9JTMMAAX2	6260-9-240	A2WMMA	6260-9-240	MERAD4
6260-9-340	8DXMMA	6260-9-240	928MLD	6260-9-340	9JVMMA	6260-9-340	A2WMPA	6260-9-244	MERAW4
6260-9-240	8E0MLD	6260-9-340	929MME	6260-9-340	9M6MMA	6260-9-340	A2XMPA	6260-9-244	MERAW5
6260-9-244	8HPMPA	6260-9-240	931MAD	6260-9-340	9M6MMAX1	6260-9-340	A2YMPA	6260-9-240	MERD8A
6260-9-240	8J1MJE	6260-9-240	932MAD	6260-9-340	9M7MMA	6260-9-240	A4YMKD	6260-9-240	MERD8M
6260-9-240	8J2MJE	6260-9-240	933MAD	6260-9-340	9M8MMA	6260-9-240	A56MMA	6260-9-240	MERD8N
6260-9-240	8J2MJEX1	6260-9-240	933MADX1	6260-9-340	9M9MMA	6260-9-240	A56MMAX1	6260-9-240	MERDE9
6260-9-340	8JEMRA	6260-9-240	935MAD	6260-9-340	9MAMMA	6260-9-240	A5EMLD	6260-9-340	MERJ71
6260-9-340	8JEMRAX1	6260-9-240	935MADX1	6260-9-340	9MDMMA	6260-9-240	A5JMLD	6260-9-340	MERJTK
6260-9-244	8K0MKD	6260-9-240	936MAD	6260-9-340	9MEMMA	6260-9-240	A5JMLDX1	6260-9-340	MERJTL
6260-9-244	8K1MKD	6260-9-240	937MAD	6260-9-340	9MHMMA	6260-9-240	A62MME	6260-9-340	MERJTM
6260-9-244	8K2MKD	6260-9-340	943MME	6260-9-340	9MJMMA	6260-9-244	A6KMPD	6260-9-340	MERJTV
6260-9-244	8K3MKD	6260-9-340	944MME	6260-9-340	9MKMMA	6260-9-244	A70MPD	6260-9-240	MERL3D
6260-9-340	8K4MKD	6260-9-340	945MME	6260-9-244	9NMMKD	6260-9-240	A82MJE	6260-9-240	MERL3E
6260-9-240	A83MJE	6260-9-244	D33MLE	6260-9-240	DY2MME	6260-9-240	EH8MMA	6260-9-240	MET6PN
6260-9-240	A83MJEX1	6260-9-244	D35MLE	6260-9-240	DY3MME	6260-9-240	EH9MMA	6260-9-240	MET6PP
6260-9-240	A85MJE	6260-9-240	D3NMMA	6260-9-240	DYEMKD	6260-9-240	EHAMMA	6260-9-240	MET6YK
6260-9-240	A8AMKD	6260-9-240	D3RMLD	6260-9-240	DYHMKD	6260-9-240	EJ0MMA	6260-9-240	MET6YL
6260-9-240	A8AMKDX1	6260-9-240	D3VMLD	6260-9-240	DYHMKDX1	6260-9-240	EJ2MMA	6260-9-240	MET729
6260-9-240	A8AMMA	6260-9-240	D4KMMA	6260-9-240	DYJMKD	6260-9-240	EK4MMA	6260-9-240	MET729X1
6260-9-244	A8DMKD	6260-9-344	D97MTA	6260-9-240	DYVME	6260-9-240	EKAMMA	6260-9-340	MET9HD
6260-9-240	A8KMMA	6260-9-444	D98MTA	6260-9-240	E0KMKD	6260-9-340	EKJMMMA	6260-9-340	MET9KP
6260-9-240	A8KMMAX1	6260-9-444	D98MTAX1	6260-9-240	E0KMKDX1	6260-9-340	EKKMMA	6260-9-340	MET9MV
6260-9-240	A8KMMAX2	6260-9-244	DALMME	6260-9-240	E0LMKD	6260-9-240	ELJMKD	6260-9-340	MET9N2
6260-9-244	AA7MTA	6260-9-244	DAPMME	6260-9-244	E0PMKD	6260-9-240	EN9MKD	6260-9-340	MET9V6
6260-9-244	AA9MTA	6260-9-244	DARMME	6260-9-244	E0PMKDX1	6260-9-340	EP5MMA	6260-9-340	MET9VE
6260-9-244	AA9MTAX1	6260-9-340	DEHMTD	6260-9-240	E0RMKD	6260-9-344	ERAMAD	6260-9-340	MET9VEX1
6260-9-244	AAAMTA	6260-9-340	DEJMTD	6260-9-244	E10MKD	6260-9-344	ERDMAD	6260-9-340	MET9VW
6260-9-244	AADMTA	6260-9-340	DEKMTD	6260-9-244	E10MKDX1	6260-9-344	EREMAD	6260-9-340	MET9VY
6260-9-244	AAHMTA	6260-9-240	DHJMRD	6260-9-240	E12MKD	6260-9-340	ERLMMA	6260-9-240	MET9W2
6260-9-244	AAHMTAX1	6260-9-240	DJ8MDE	6260-9-240	E16MRD	6260-9-240	ET7MME	6260-9-240	MET9WA
6260-9-240	AAJMJE	6260-9-240	DJ9MDE	6260-9-240	E17MRD	6260-9-240	ETWMLE	6260-9-240	MET9WAX1

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-240	ADKMHE	6260-9-240	DJAMDE	6260-9-240	E18MRD	6260-9-340	EV2MPA	6260-9-240	MET9XX
6260-9-240	ADMMHE	6260-9-240	DJLMMA	6260-9-240	E18MRDX1	6260-9-340	EV2MPAX1	6260-9-240	MET9Y0
6260-9-240	AH4MKD	6260-9-240	DJRMMA	6260-9-240	E1XMEE	6260-9-240	EV4MLE	6260-9-240	METE02
6260-9-240	AH4MKDX1	6260-9-240	DJRMMAX1	6260-9-240	E1XMKD	6260-9-240	EV5MPA	6260-9-240	METE05
6260-9-244	AHAMKD	6260-9-240	DL5MJJE	6260-9-240	E1XMKDX1	6260-9-240	EV6MLE	6260-9-240	METE45
6260-9-240	AHDMJE	6260-9-240	DL5MRD	6260-9-240	E1YMHJE	6260-9-240	EWYMLE	6260-9-240	METE45X1
6260-9-244	AJVMKD	6260-9-240	DL6MRD	6260-9-240	E23MEE	6260-9-240	EX8MHD	6260-9-240	METE4J
6260-9-240	AJWMKD	6260-9-240	DLTMND	6260-9-240	E35MHD	6260-9-240	EX8MHDX1	6260-9-240	METE4JX1
6260-9-244	AK0MAD	6260-9-240	DLVMND	6260-9-240	E3JMND	6260-9-240	EX9MHD	6260-9-240	METE5J
6260-9-244	AK0MADX1	6260-9-240	DMLMJJE	6260-9-240	E3NMMA	6260-9-240	EXDMLA	6260-9-240	METE5T
6260-9-244	AK1MAD	6260-9-240	DMTMND	6260-9-240	E3NMMAAX1	6260-9-240	EXKMHD	6260-9-240	METE64
6260-9-244	AK2MAD	6260-9-240	DMVMND	6260-9-240	E3RMMA	6260-9-240	EXLMHD	6260-9-340	METHE9
6260-9-244	AK3MAD	6260-9-240	DMVMNDX1	6260-9-240	E3RMMAAX1	6260-9-240	EXMMHD	6260-9-340	METHEA
6260-9-240	AM1MKD	6260-9-244	DN5MKD	6260-9-240	E43MHE	6260-9-240	EYMLA	6260-9-340	METHED
6260-9-240	AMPMMA	6260-9-240	DNVMJE	6260-9-240	E4EMKD	6260-9-240	EYRMLA	6260-9-340	METHEE
6260-9-240	AMPMMAX1	6260-9-240	DNVMJEX1	6260-9-240	E5LMHE	6260-9-240	H0DMLA	6260-9-340	METHEEX1
6260-9-240	ANYMDE	6260-9-240	DNXMDE	6260-9-340	E5RMTD	6260-9-340	H0DMMA	6260-9-340	METHJL
6260-9-240	AV5MMA	6260-9-240	DP4MJJE	6260-9-340	E5VMTD	6260-9-240	H0EMLA	6260-9-340	METHJM
6260-9-240	AV5MMAX1	6260-9-240	DP7MMA	6260-9-340	E5WMTD	6260-9-340	H0EMMA	6260-9-340	METHJMX1
6260-9-240	AV9MLD	6260-9-240	DPJMJJE	6260-9-340	E6AMLE	6260-9-340	H0JMMA	6260-9-340	METHJN
6260-9-240	AVDMLD	6260-9-340	DVVMHE	6260-9-340	E6DMLE	6260-9-240	H0KMMA	6260-9-340	METHJP
6260-9-240	AVHMLD	6260-9-340	DVWMHE	6260-9-244	E6LMJE	6260-9-240	H2WMME	6260-9-340	METHJPX1
6260-9-240	AVJMLD	6260-9-340	DVXMHE	6260-9-244	E6MMJE	6260-9-240	H3EMME	6260-9-240	METKMR
6260-9-244	AWTMRA	6260-9-240	DVXMKD	6260-9-244	E6NMJE	6260-9-240	H3HMME	6260-9-240	METKMV
6260-9-244	AWTMRAX1	6260-9-240	DVXMKDX1	6260-9-244	E6PMJE	6260-9-240	H68MHD	6260-9-240	METKMW
6260-9-244	AWVMRA	6260-9-340	DVYMHE	6260-9-240	E96MHE	6260-9-240	H6HMHD	6260-9-240	METKMX
6260-9-244	AWWMRA	6260-9-240	DVYMKD	6260-9-240	ED4MLE	6260-9-240	H7HMME	6260-9-240	METKYN
6260-9-244	AWWMRAX1	6260-9-240	DW0MKD	6260-9-240	ED5MLE	6260-9-340	H7RMRD	6260-9-240	METKYP
6260-9-244	AWXMRA	6260-9-344	DW1MAD	6260-9-240	ED6MLE	6260-9-340	H7RMRDX1	6260-9-240	METKYR
6260-9-244	AWXMRAX1	6260-9-344	DW2MAD	6260-9-240	ED7MLE	6260-9-240	H92MNE	6260-9-240	METKYT
6260-9-240	AY2MMA	6260-9-344	DW2MADX1	6260-9-240	EDNMEE	6260-9-240	H9HMHD	6260-9-240	METLLE
6260-9-340	D0AMLE	6260-9-344	DW2MADX1X1	6260-9-240	EEKMKD	6260-9-240	H9HMHDX1	6260-9-240	METLM2
6260-9-340	D2JMLE	6260-9-344	DW3MAD	6260-9-244	EH4MME	6260-9-340	H9PMRD	6260-9-240	METLM3
6260-9-244	D2WMLE	6260-9-240	DXDMEE	6260-9-244	EH5MME	6260-9-240	HAHMNE	6260-9-240	METLMA
6260-9-244	D2XMME	6260-9-240	DY0MME	6260-9-244	EH6MME	6260-9-240	HD7MNE	6260-9-240	METLME
6260-9-244	D2YMLE	6260-9-240	DY1MME	6260-9-244	EH7MME	6260-9-240	HD8MNE	6260-9-240	METLML
6260-9-240	HDNMHD	6260-9-240	J12MLA	6260-9-344	K8YMA	6260-9-340	KY4MJJE	6260-9-240	METLYP
6260-9-340	HH5MPD	6260-9-240	J13MHE	6260-9-344	K90MAD	6260-9-340	L11MHD	6260-9-240	METM1M
6260-9-240	HJ2MLD	6260-9-240	J14MHE	6260-9-344	K91MAD	6260-9-244	L1AMEE	6260-9-240	METM27
6260-9-240	HJKMRA	6260-9-240	J1TMLA	6260-9-344	K92MAD	6260-9-244	L1EMEE	6260-9-244	MHA2AW
6260-9-240	HJKMRAX1	6260-9-240	J1WMHE	6260-9-240	K92MKE	6260-9-244	L1JMEE	6260-9-244	MHA2AX
6260-9-240	HJMMRA	6260-9-240	J28MHE	6260-9-344	K93MAD	6260-9-244	L1KMEE	6260-9-244	MHA2AY
6260-9-240	HJMMRAX1	6260-9-240	J2HMLA	6260-9-240	KA0MKE	6260-9-340	L34MHD	6260-9-244	MHA2D0
6260-9-340	HJPMRA	6260-9-340	J3PMME	6260-9-240	KA0MKEX1	6260-9-240	L4XMPA	6260-9-244	MHA2D3
6260-9-340	HJPMRAX1	6260-9-240	J3PMRD	6260-9-240	KA8MKE	6260-9-240	L4YMPA	6260-9-244	MHA2D5
6260-9-340	HJWMRA	6260-9-340	J3TMME	6260-9-240	KD7MDD	6260-9-240	L50MPA	6260-9-244	MHA2D6
6260-9-444	HK1MPA	6260-9-340	J3VMME	6260-9-240	KD8MDD	6260-9-240	L50MPAX1	6260-9-244	MHA2D7
6260-9-444	HK1MPAX1	6260-9-340	J3WMME	6260-9-240	KDAMKE	6260-9-240	L51MPA	6260-9-244	MHA3D0
6260-9-240	HKNMAE	6260-9-340	J4VMAE	6260-9-340	KDKMHD	6260-9-340	L53MHD	6260-9-244	MHA3D9
6260-9-240	HL8MHD	6260-9-340	J4WMAE	6260-9-240	KEEMHE	6260-9-344	L6XMAD	6260-9-244	MHA3DY
6260-9-340	HL8MLD	6260-9-340	J4XMAE	6260-9-240	KEHMHE	6260-9-340	L8MMHD	6260-9-244	MHA3H0
6260-9-340	HL9MLD	6260-9-340	J5JMKE	6260-9-240	KEJMHE	6260-9-240	L9YMPA	6260-9-240	MHAD2H
6260-9-240	HMDMHD	6260-9-340	J64MKE	6260-9-340	KEKMHE	6260-9-340	LAEMPA	6260-9-240	MHAD2K
6260-9-240	HMDMHDX1	6260-9-340	J6YMKE	6260-9-340	KELMHE	6260-9-340	LAEMPAX1	6260-9-240	MHAD2L
6260-9-340	HMLMRD	6260-9-340	J7KMKE	6260-9-340	KEMMHE	6260-9-340	LAHMPA	6260-9-240	MHAD2M
6260-9-240	HNHMAE	6260-9-340	J7LMKE	6260-9-340	KENMHE	6260-9-240	LAWMMA	6260-9-240	MHADWJ
6260-9-240	HP8MHD	6260-9-340	J7MMKE	6260-9-240	KH2MKE	6260-9-240	LAXMMA	6260-9-240	MHADWJX1
6260-9-240	HP8MHDX1	6260-9-340	J7NMKE	6260-9-340	KHEMHD	6260-9-444	LD1MAD	6260-9-240	MHADWL
6260-9-340	HPVMRD	6260-9-240	J8WMRD	6260-9-340	KHEMHDX1	6260-9-444	LD8MAD	6260-9-340	MHADYA
6260-9-340	HR2MPD	6260-9-240	J9YMRD	6260-9-240	KKNMKE	6260-9-344	LDHMAD	6260-9-340	MHADYD

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-340	HR4MPD	6260-9-240	JA9MRD	6260-9-240	KL4MDE	6260-9-344	LDHMADX1	6260-9-340	MHAE0Y
6260-9-340	HR5MPD	6260-9-340	JAEMKE	6260-9-240	KL5MDE	6260-9-344	LDKMAD	6260-9-340	MHAE11
6260-9-340	HR5MPDX1	6260-9-340	JAHMKE	6260-9-240	KL6MDE	6260-9-344	LDLMAD	6260-9-340	MHAE77
6260-9-340	HR7MPD	6260-9-240	JD4MKD	6260-9-240	KL7MDE	6260-9-344	LDMMAD	6260-9-340	MHAEEH
6260-9-240	HRXMAE	6260-9-240	JD7MKD	6260-9-240	KM3MKE	6260-9-344	LDNMAD	6260-9-340	MHAE EJ
6260-9-240	HRYMHD	6260-9-240	JDNMTA	6260-9-340	KMEMHD	6260-9-344	LDPMAD	6260-9-340	MHAE EK
6260-9-244	HV7MRA	6260-9-240	JDRMTA	6260-9-240	KMNMD D	6260-9-444	LE2MAD	6260-9-340	MHAE EKX1
6260-9-240	HVJMRA	6260-9-240	JE0MTA	6260-9-240	KMPMD D	6260-9-444	LE2MADX1	6260-9-340	MHAE EN
6260-9-240	HVKMRA	6260-9-240	JE1MTA	6260-9-240	KMRMD D	6260-9-240	LE6MHD	6260-9-340	MHAE EY
6260-9-340	HVLMRA	6260-9-340	JH4MPD	6260-9-240	KMTMD D	6260-9-240	LE7MHD	6260-9-340	MHAE H1
6260-9-340	HVMMRA	6260-9-340	JK1MKE	6260-9-240	KMVMDD	6260-9-240	LE8MHD	6260-9-340	MHAE H2
6260-9-240	HW8MAE	6260-9-240	JK3MRD	6260-9-240	KMWMD D	6260-9-240	LE8MHD X1	6260-9-340	MHAE H3
6260-9-340	HW9MRD	6260-9-340	JK6MKE	6260-9-240	KMWMDX1	6260-9-240	LJLMMA	6260-9-340	MHAHDN
6260-9-244	HWDMD	6260-9-240	JLMMRD	6260-9-240	KMXMD D	6260-9-240	LJMMMA	6260-9-340	MHAHDNX1
6260-9-244	HWDMDX1	6260-9-340	JPHMME	6260-9-240	KMYMD D	6260-9-240	LJMMMAX1	6260-9-240	MHALR6
6260-9-244	HWHMAD	6260-9-340	JPRMME	6260-9-340	KPTMKD	6260-9-240	LJPMMA	6260-9-240	MHALYM
6260-9-244	HWJMAD	6260-9-340	JPWMME	6260-9-340	KPTMKDX1	6260-9-240	LJRMMA	6260-9-240	MHALYN
6260-9-244	HWJMADX1	6260-9-340	JPXMME	6260-9-240	KRRMMA	6260-9-240	LJTMMA	6260-9-240	MHALYP
6260-9-244	HWLMAD	6260-9-240	JW2MHE	6260-9-240	KV9MME	6260-9-240	LJTMMAX1	6260-9-240	MHALYPX1
6260-9-244	HWMMAD	6260-9-240	K1HMKE	6260-9-240	KVAMME	6260-9-240	LJWMMA	6260-9-240	MHAM22
6260-9-244	HWMMADX1	6260-9-240	K43MKE	6260-9-240	KVDMME	6260-9-240	LK0MMA	6260-9-240	MHAM23
6260-9-244	HWPMD	6260-9-240	K4KMKE	6260-9-240	KVEMME	6260-9-340	LMDMMA	6260-9-240	MHAM2T
6260-9-340	J01MAE	6260-9-340	K6JMTA	6260-9-340	KVVMHD	6260-9-340	LMEMMA	6260-9-240	MHAM94
6260-9-240	J05MLA	6260-9-340	K6KMTA	6260-9-240	KVXMMA	6260-9-340	LMHMMA	6260-9-240	MHAMEE
6260-9-240	J06MLA	6260-9-340	K7WMJE	6260-9-240	KVXMMAX1	6260-9-340	LMRMMA	6260-9-240	MHAMEEX1
6260-9-240	J07MLA	6260-9-340	K7XMJE	6260-9-240	KVYMMA	6260-9-340	LMWMMA	6260-9-240	MHAMH0
6260-9-240	J07MLAX1	6260-9-340	K7YMJE	6260-9-244	KWNMHE	6260-9-340	LMWMMAX1	6260-9-240	MHAM P3
6260-9-240	J07MLAX1X1	6260-9-340	K80MJE	6260-9-244	KWRMHE	6260-9-340	LYYMMA	6260-9-240	MHAMVV
6260-9-240	J08MLA	6260-9-240	K84MKE	6260-9-244	KWTMHE	6260-9-340	LN8MMA	6260-9-240	MHAN37
6260-9-244	J0DMPA	6260-9-344	K8XMAD	6260-9-244	KWVMHE	6260-9-240	LP6MLE	6260-9-240	MHD2HM
6260-9-244	MHD671	6260-9-340	MHDV0W	6260-9-244	MHE254X1	6260-9-340	MHEMH9X1	6260-9-240	MHH2RV
6260-9-244	MHD671X1	6260-9-340	MHDV0X	6260-9-244	MHE259	6260-9-240	MHEMHJ	6260-9-240	MHH2RVX1
6260-9-244	MHD67M	6260-9-340	MHDV0XX1	6260-9-244	MHE26P	6260-9-240	MHEMHJX1	6260-9-240	MHH2X2
6260-9-244	MHD68N	6260-9-340	MHDV8M	6260-9-244	MHE2AH	6260-9-240	MHEMHK	6260-9-240	MHH2X5
6260-9-244	MHD69E	6260-9-340	MHDV8MX1	6260-9-244	MHE2E0	6260-9-240	MHEMHN	6260-9-240	MHH2XR
6260-9-244	MHD69T	6260-9-340	MHDVEX	6260-9-244	MHE2KH	6260-9-240	MHEMHNX1	6260-9-240	MHH2XRX1
6260-9-244	MHD69X	6260-9-340	MHDVHY	6260-9-240	MHE97W	6260-9-240	MHEMJW	6260-9-240	MHH2XT
6260-9-244	MHD6AR	6260-9-340	MHDVJ2	6260-9-240	MHE999	6260-9-240	MHEMJWX1	6260-9-244	MHH5HD
6260-9-244	MHD6MD	6260-9-340	MHDVL4	6260-9-240	MHE9HP	6260-9-244	MHEMKR	6260-9-244	MHH5JA
6260-9-244	MHD6RJ	6260-9-340	MHDVL4X1	6260-9-240	MHE9JD	6260-9-244	MHEMKRX1	6260-9-244	MHH5M2
6260-9-244	MHD6TM	6260-9-340	MHDVLA	6260-9-240	MHE9K9	6260-9-244	MHEMKRX1X1	6260-9-244	MHH5PP
6260-9-244	MHD6TR	6260-9-240	MHDVM7	6260-9-240	MHEL52	6260-9-240	MHEMKY	6260-9-244	MHH5W1
6260-9-244	MHD6TRX1	6260-9-240	MHDVRY	6260-9-240	MHEL52X1	6260-9-240	MHEMM6	6260-9-244	MHH62E
6260-9-244	MHD6VA	6260-9-240	MHDVRYX1	6260-9-240	MHEL8L	6260-9-244	MHEMMD	6260-9-244	MHH659
6260-9-244	MHD6VD	6260-9-240	MHDVT1	6260-9-240	MHELLD	6260-9-240	MHEMNN	6260-9-244	MHH65H
6260-9-244	MHD6VE	6260-9-240	MHDVT2	6260-9-240	MHELLV	6260-9-340	MHERWN	6260-9-244	MHH65R
6260-9-244	MHD6VH	6260-9-240	MHDVT5	6260-9-240	MHELLW	6260-9-340	MHERWP	6260-9-240	MHH662
6260-9-244	MHD7JJ	6260-9-240	MHDVT6	6260-9-240	MHELLWX1	6260-9-340	MHERWR	6260-9-240	MHH680
6260-9-240	MHDMM9	6260-9-240	MHDWWD	6260-9-240	MHEL PW	6260-9-340	MHERWT	6260-9-240	MHH682
6260-9-240	MHDMMA	6260-9-240	MHDWX2	6260-9-340	MHEM89	6260-9-240	MHET5J	6260-9-240	MHH689
6260-9-240	MHDMMD	6260-9-240	MHDWXK	6260-9-340	MHEM89X1	6260-9-240	MHET99	6260-9-240	MHH68D
6260-9-240	MHDMME	6260-9-240	MHDY0A	6260-9-340	MHEM9P	6260-9-240	MHETH9	6260-9-240	MHH707
6260-9-240	MHDNPH	6260-9-240	MHDY0H	6260-9-340	MHEMAE	6260-9-240	MHETH9X1	6260-9-240	MHH74N
6260-9-240	MHDNPJ	6260-9-240	MHE0N9	6260-9-244	MHEMAH	6260-9-240	MHETHJ	6260-9-240	MHH74R
6260-9-240	MHDNPK	6260-9-240	MHE0W5	6260-9-340	MHEMD1	6260-9-240	MHETHJX1	6260-9-240	MHHLW0
6260-9-240	MHDNPL	6260-9-240	MHE0W8	6260-9-340	MHEMD1X1	6260-9-240	MHETLM	6260-9-240	MHHLYN
6260-9-240	MHDP9N	6260-9-240	MHE107	6260-9-244	MHEMD3	6260-9-240	MHETLP	6260-9-240	MHHM2E
6260-9-240	MHDP9P	6260-9-240	MHE108	6260-9-244	MHEMD3X1	6260-9-240	MHETMP	6260-9-240	MHHM5V
6260-9-240	MHDPLL	6260-9-240	MHE116	6260-9-340	MHEMD6	6260-9-240	MHETMPX1	6260-9-240	MHHMA7
6260-9-240	MHDPLN	6260-9-240	MHE117	6260-9-340	MHEMD6X1	6260-9-240	MHETMR	6260-9-240	MHHMAK

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-240	MHDPT3	6260-9-244	MHE22P	6260-9-244	MHEMDD	6260-9-240	MHH1X0	6260-9-240	MHHN0V
6260-9-240	MHDR96	6260-9-244	MHE22Y	6260-9-244	MHEMDDX1	6260-9-240	MHH28V	6260-9-240	MHHN11
6260-9-340	MHDR9Y	6260-9-244	MHE23D	6260-9-340	MHEMEH	6260-9-240	MHH2DK	6260-9-240	MHHN15
6260-9-340	MHDRA8	6260-9-244	MHE23E	6260-9-340	MHEMEHX1	6260-9-240	MHH2E0	6260-9-240	MHHN4N
6260-9-240	MHDRAH	6260-9-244	MHE254	6260-9-340	MHEMH9	6260-9-240	MHH2RE	6260-9-240	MHHN4P
6260-9-240	MHHN4R	6260-9-240	MHK12A	6260-9-240	MHKPRN	6260-9-244	MHLW4K	6260-9-240	MHMAAK
6260-9-240	MHHN4T	6260-9-240	MHK12D	6260-9-244	MHKPT2	6260-9-244	MHLW76	6260-9-240	MHMAAL
6260-9-240	MHHNJL	6260-9-240	MHK1K7	6260-9-244	MHKPV3	6260-9-344	MHLX8Y	6260-9-240	MHMAAM
6260-9-240	MHHNJN	6260-9-240	MHK1RX	6260-9-244	MHKPVE	6260-9-344	MHLXAP	6260-9-240	MHMAAN
6260-9-240	MHHNLM	6260-9-240	MHK1V0	6260-9-244	MHKPVT	6260-9-344	MHLXD9	6260-9-240	MHMAAP
6260-9-240	MHHNLN	6260-9-240	MHK1VD	6260-9-240	MHKPY7	6260-9-344	MHLXDA	6260-9-240	MHMAAR
6260-9-240	MHHP9N	6260-9-240	MHK1VL	6260-9-240	MHKPY8	6260-9-344	MHLXDD	6260-9-240	MHMAAT
6260-9-240	MHHP9T	6260-9-240	MHK1VV	6260-9-240	MHKPY9	6260-9-344	MHLXDE	6260-9-240	MHMAATX1
6260-9-240	MHHP9V	6260-9-240	MHK1W6	6260-9-240	MHKPYA	6260-9-340	MHLXL2	6260-9-240	MHMAAV
6260-9-240	MHHP9W	6260-9-244	MHK5E7	6260-9-244	MHKR01	6260-9-340	MHLXL3	6260-9-240	MHMAAW
6260-9-240	MHHPMV	6260-9-244	MHK5J1	6260-9-244	MHKR29	6260-9-340	MHLXL4	6260-9-240	MHMAAX
6260-9-340	MHJOWY	6260-9-244	MHK5M6	6260-9-244	MHKR3X	6260-9-340	MHLXL5	6260-9-240	MHMAAY
6260-9-340	MHJ16J	6260-9-244	MHK5N8	6260-9-244	MHKR6E	6260-9-340	MHLY2P	6260-9-240	MHMAD0
6260-9-340	MHJ16W	6260-9-244	MHK5RR	6260-9-240	MHKR87	6260-9-340	MHLY2R	6260-9-240	MHMAD1
6260-9-340	MHJ172	6260-9-244	MHK61D	6260-9-340	MHKRJR	6260-9-340	MHLY2T	6260-9-340	MHMD86
6260-9-340	MHJ19T	6260-9-244	MHK61J	6260-9-340	MHKRYK	6260-9-340	MHLY30	6260-9-340	MHMD87
6260-9-340	MHJ19W	6260-9-244	MHK61K	6260-9-340	MHKRYX	6260-9-240	MHLYD3	6260-9-340	MHMD88
6260-9-340	MHJ19X	6260-9-240	MHKLYR	6260-9-340	MHKT3E	6260-9-240	MHM001	6260-9-340	MHMD89
6260-9-340	MHJ19Y	6260-9-240	MHKM96	6260-9-340	MHKT3H	6260-9-240	MHM003	6260-9-340	MHMDYK
6260-9-340	MHJ26L	6260-9-240	MHKMDR	6260-9-340	MHKT3J	6260-9-240	MHM00H	6260-9-340	MHMDYKX1
6260-9-340	MHJ27J	6260-9-240	MHKMMJ	6260-9-340	MHKT3K	6260-9-240	MHM00P	6260-9-340	MHMDYL
6260-9-340	MHJ27W	6260-9-240	MHKMMM	6260-9-340	MHKT3M	6260-9-240	MHM00X	6260-9-340	MHME08
6260-9-340	MHJ282	6260-9-240	MHKMMN	6260-9-340	MHKT3P	6260-9-240	MHM00Y	6260-9-340	MHME09
6260-9-340	MHJ283	6260-9-240	MHKMN7	6260-9-340	MHKTRJ	6260-9-240	MHM010	6260-9-340	MHMMAA
6260-9-244	MHJ2NL	6260-9-240	MHKMN8	6260-9-340	MHKTRN	6260-9-240	MHM04N	6260-9-340	MHMMMJ
6260-9-244	MHJ2RE	6260-9-240	MHKMNK	6260-9-340	MHKV40	6260-9-240	MHM04P	6260-9-340	MHMMN2
6260-9-244	MHJ2VJ	6260-9-240	MHKMNL	6260-9-340	MHKV41	6260-9-240	MHM04R	6260-9-340	MHMMN3
6260-9-244	MHJ31N	6260-9-240	MHKMNX	6260-9-340	MHKV42	6260-9-240	MHM8YX	6260-9-340	MHMMNH
6260-9-244	MHJ31R	6260-9-240	MHKMP3	6260-9-340	MHKV7L	6260-9-240	MHM8YY	6260-9-340	MHMMP6
6260-9-244	MHJ359	6260-9-240	MHKMXN	6260-9-340	MHKV7N	6260-9-240	MHM900	6260-9-340	MHMMP7
6260-9-244	MHJ360	6260-9-240	MHKMXR	6260-9-340	MHKV7P	6260-9-240	MHM99X	6260-9-340	MHMPYV
6260-9-244	MHJ38Y	6260-9-240	MHKMXX	6260-9-340	MHKVMA	6260-9-240	MHM99Y	6260-9-340	MHMR0N
6260-9-244	MHJ390	6260-9-240	MHKMXY	6260-9-340	MHKVMD	6260-9-240	MHM9A0	6260-9-340	MHMR0R
6260-9-240	MHJJKJ	6260-9-240	MHKN5E	6260-9-340	MHKVT4	6260-9-240	MHM9DM	6260-9-340	MHMR0T
6260-9-240	MHJJK7X	6260-9-240	MHKN5K	6260-9-340	MHL00E	6260-9-240	MHM9DN	6260-9-340	MHMR1K
6260-9-240	MHJK88	6260-9-240	MHKN5M	6260-9-340	MHL06H	6260-9-240	MHM9E5	6260-9-244	MHN115
6260-9-240	MHJK8J	6260-9-240	MHKNED	6260-9-340	MHL06J	6260-9-240	MHM9E6	6260-9-244	MHN169
6260-9-240	MHJK8T	6260-9-240	MHKNPR	6260-9-340	MHL0KK	6260-9-240	MHM9HY	6260-9-244	MHN19M
6260-9-240	MHJKT0	6260-9-240	MHKNPV	6260-9-340	MHL0KN	6260-9-240	MHM9HYX1	6260-9-244	MHN19N
6260-9-240	MHJKV3	6260-9-240	MHKNR3	6260-9-340	MHLJA4	6260-9-240	MHM9J0	6260-9-244	MHN19NX1
6260-9-240	MHJP2J	6260-9-240	MHKNRD	6260-9-340	MHLJA5	6260-9-240	MHM9J2	6260-9-244	MHN1K7
6260-9-240	MHJP92	6260-9-240	MHKNRE	6260-9-340	MHLJA6	6260-9-240	MHM9J3	6260-9-244	MHN1NX
6260-9-240	MHJP99	6260-9-240	MHKNRT	6260-9-340	MHLJA7	6260-9-240	MHM9J4	6260-9-244	MHN1TA
6260-9-240	MHJP9J	6260-9-240	MHKNTD	6260-9-340	MHLJRD	6260-9-240	MHM9J5	6260-9-244	MHN1TY
6260-9-240	MHJPM0	6260-9-240	MHKNYR	6260-9-340	MHLJRE	6260-9-240	MHM9J6	6260-9-244	MHN1V2
6260-9-240	MHJPRO	6260-9-240	MHKPOP	6260-9-340	MHLJRH	6260-9-240	MHM9J7	6260-9-244	MHN1V3
6260-9-240	MHJPTX	6260-9-240	MHKPOR	6260-9-340	MHLK16	6260-9-240	MHM9J8	6260-9-244	MHN1WJ
6260-9-240	MHJPTY	6260-9-240	MHKPOT	6260-9-340	MHLK1T	6260-9-240	MHM9J9	6260-9-244	MHN1WL
6260-9-240	MHJPYY	6260-9-240	MHKPOV	6260-9-340	MHLK26	6260-9-240	MHM9JA	6260-9-244	MHN1WN
6260-9-240	MHJR1R	6260-9-240	MHKP9L	6260-9-340	MHLK3K	6260-9-240	MHMA9Y	6260-9-244	MHN1WNX1
6260-9-240	MHJR36	6260-9-240	MHKP9M	6260-9-244	MHLV9X	6260-9-240	MHMAA1	6260-9-244	MHN2T2
6260-9-240	MHJR3K	6260-9-240	MHKPDE	6260-9-244	MHLVD4	6260-9-240	MHMAA4	6260-9-244	MHN2T2X1
6260-9-240	MHJR3R	6260-9-240	MHKPJH	6260-9-244	MHLVD4X1	6260-9-240	MHMAA6	6260-9-244	MHN2WT
6260-9-240	MHJR3Y	6260-9-240	MHKPM9	6260-9-244	MHLVD5	6260-9-240	MHMAA7	6260-9-244	MHN2WV
6260-9-240	MHJR40	6260-9-240	MHKPMA	6260-9-244	MHLVD6	6260-9-240	MHMAA7X1	6260-9-244	MHN2WVW

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-240	MHJR4D	6260-9-240	MHKPPE	6260-9-244	MHLVD7	6260-9-240	MHMAA8	6260-9-244	MHN2X6
6260-9-340	MHK0RN	6260-9-240	MHKPR0	6260-9-244	MHLVME	6260-9-240	MHMAAA	6260-9-244	MHN2X6X1
6260-9-340	MHK0RW	6260-9-240	MHKPR7	6260-9-244	MHLVMJ	6260-9-240	MHMAAE	6260-9-244	MHN2X7
6260-9-240	MHK127	6260-9-240	MHKPRE	6260-9-244	MHLVMK	6260-9-240	MHMAAH	6260-9-244	MHN2X8
6260-9-240	MHK128	6260-9-240	MHKPRH	6260-9-244	MHLVML	6260-9-240	MHMAAJ	6260-9-340	MHN2Y7
6260-9-240	MHK129	6260-9-240	MHKPRM	6260-9-244	MHLW4H	6260-9-240	MHMAAJX1	6260-9-340	MHN2YT
6260-9-340	MHN2YV	6260-9-340	MHNX1H	6260-9-240	MHR13D	6260-9-340	MHRMYD	6260-9-240	MJA7KD
6260-9-340	MHN2YW	6260-9-240	MHNXMD	6260-9-240	MHR13J	6260-9-240	MHRN1V	6260-9-240	MJA7KJ
6260-9-340	MHN33Y	6260-9-240	MHNXMK	6260-9-240	MHR275	6260-9-240	MHRN1W	6260-9-240	MJA7KK
6260-9-240	MHN94N	6260-9-240	MHNXXM	6260-9-240	MHR27P	6260-9-240	MHRN1X	6260-9-240	MJA7KL
6260-9-240	MHN95K	6260-9-240	MHNXN3	6260-9-240	MHR27R	6260-9-240	MHRN1Y	6260-9-340	MJA7XR
6260-9-240	MHN95L	6260-9-240	MHNXPH	6260-9-240	MHR27T	6260-9-240	MHT7HV	6260-9-340	MJA7XT
6260-9-240	MHN95M	6260-9-240	MHNXPJ	6260-9-240	MHR27V	6260-9-240	MHT7X0	6260-9-340	MJA7XV
6260-9-240	MHN98X	6260-9-240	MHNXPK	6260-9-340	MHR2N8	6260-9-240	MHT897	6260-9-340	MJA7XW
6260-9-240	MHN98Y	6260-9-240	MHNXPL	6260-9-340	MHR2ND	6260-9-240	MHT89Y	6260-9-240	MJA7YN
6260-9-240	MHN990	6260-9-240	MHNY44	6260-9-340	MHR2NL	6260-9-240	MHT8AK	6260-9-240	MJA7YP
6260-9-240	MHN9AT	6260-9-240	MHNY4T	6260-9-340	MHR2NR	6260-9-240	MHT8AR	6260-9-340	MJA82H
6260-9-240	MHN9EA	6260-9-240	MHNY4W	6260-9-240	MHR33H	6260-9-240	MHT8JM	6260-9-340	MJA82K
6260-9-240	MHN9ED	6260-9-240	MHNY4Y	6260-9-240	MHR33J	6260-9-240	MHT8K9	6260-9-340	MJA82L
6260-9-240	MHN9EE	6260-9-244	MHP427	6260-9-240	MHR33K	6260-9-240	MHT8L3	6260-9-340	MJA82M
6260-9-240	MHN9EH	6260-9-244	MHP428	6260-9-240	MHR33L	6260-9-240	MHT8LJ	6260-9-340	MJA82N
6260-9-240	MHN9MD	6260-9-244	MHP42A	6260-9-240	MHR3LH	6260-9-240	MHT957	6260-9-240	MJA83P
6260-9-240	MHN9ME	6260-9-244	MHP42D	6260-9-240	MHR4VT	6260-9-240	MHT96W	6260-9-244	MJA865
6260-9-240	MHN9MH	6260-9-244	MHP5Y8	6260-9-240	MHR4VV	6260-9-240	MHT975	6260-9-244	MJA87M
6260-9-240	MHN9MK	6260-9-244	MHP679	6260-9-240	MHR4VW	6260-9-240	MHT99K	6260-9-244	MJA8JE
6260-9-240	MHNA1V	6260-9-244	MHP807	6260-9-240	MHR4VY	6260-9-240	MHT99M	6260-9-244	MJA8MR
6260-9-240	MHNALE	6260-9-244	MHP808	6260-9-244	MHR7LL	6260-9-240	MHT9A9	6260-9-244	MJA8VJ
6260-9-240	MHNALH	6260-9-244	MHP80A	6260-9-244	MHR7LV	6260-9-240	MHT9AV	6260-9-244	MJA9H4
6260-9-240	MHNALJ	6260-9-244	MHP80D	6260-9-244	MHR7LW	6260-9-240	MHT9AW	6260-9-244	MJA9JH
6260-9-240	MHNALK	6260-9-244	MHP80E	6260-9-244	MHR7LX	6260-9-340	MHTH6E	6260-9-244	MJA9K8
6260-9-240	MHNAV0	6260-9-244	MHP80H	6260-9-244	MHR7M0	6260-9-340	MHTH9W	6260-9-244	MJA9NL
6260-9-240	MHNAV1	6260-9-244	MHP80K	6260-9-244	MHR7M1	6260-9-340	MHTHN2	6260-9-244	MJA9RR
6260-9-240	MHNAV2	6260-9-244	MHP80L	6260-9-244	MHR7M2	6260-9-340	MHTHNL	6260-9-244	MJA9XK
6260-9-240	MHNAV3	6260-9-244	MHP80M	6260-9-244	MHR7M4	6260-9-340	MHTHVJ	6260-9-244	MJAA08
6260-9-240	MHNDY	6260-9-244	MHP80P	6260-9-340	MHREL3	6260-9-340	MHTHYT	6260-9-244	MJAAAD1
6260-9-240	MHNDE1	6260-9-240	MHPJJ2	6260-9-340	MHREL5	6260-9-340	MHTJ2K	6260-9-244	MJAAAD3
6260-9-240	MHNDE4	6260-9-240	MHPJJ3	6260-9-340	MHREL7	6260-9-340	MHTJ2L	6260-9-244	MJAAE7
6260-9-240	MHNDE9	6260-9-240	MHPJJ4	6260-9-340	MHREL8	6260-9-340	MHTJ42	6260-9-244	MJAAE8
6260-9-240	MHNDVL	6260-9-240	MHPJJ6	6260-9-240	MHRH8H	6260-9-340	MHTJ43	6260-9-244	MJAAE9
6260-9-240	MHNDVM	6260-9-240	MHPJJN	6260-9-240	MHRH8J	6260-9-340	MHTJ44	6260-9-244	MJAAH8
6260-9-240	MHNDX8	6260-9-240	MHPJW8	6260-9-240	MHRH8K	6260-9-340	MHTJ48	6260-9-240	MJAN0V
6260-9-240	MHNE0J	6260-9-240	MHPJW9	6260-9-240	MHRH8L	6260-9-340	MHTJ49	6260-9-240	MJAN7N
6260-9-240	MHNE0K	6260-9-240	MHPJW9X1	6260-9-240	MHRHME	6260-9-240	MHTJ6T	6260-9-240	MJAN94
6260-9-240	MHNE0L	6260-9-240	MHPJWD	6260-9-240	MHRHMH	6260-9-240	MHTJ86	6260-9-240	MJAN95
6260-9-240	MHNE0M	6260-9-240	MHPJWE	6260-9-240	MHRHMJ	6260-9-240	MHTJ95	6260-9-240	MJAN9X
6260-9-240	MHNE7N	6260-9-240	MHPKML	6260-9-240	MHRHPA	6260-9-240	MHTJER	6260-9-240	MJAN9Y
6260-9-240	MHNE7P	6260-9-240	MHPKMM	6260-9-344	MHRL7D	6260-9-240	MHTJHT	6260-9-240	MJANAA
6260-9-240	MHNE7V	6260-9-240	MHPKMN	6260-9-344	MHRL9R	6260-9-244	MHTKAH	6260-9-240	MJANAD
6260-9-240	MHNEMW	6260-9-240	MHPKMP	6260-9-344	MHRLN0	6260-9-244	MHTKKW	6260-9-240	MJANAX
6260-9-240	MHNEMWX1	6260-9-240	MHPKMR	6260-9-344	MHRLNX	6260-9-244	MHTKLW	6260-9-240	MJAND7
6260-9-240	MHNEMY	6260-9-240	MHPKMT	6260-9-344	MHRLNY	6260-9-244	MHTKNO	6260-9-240	MJANE1
6260-9-340	MHNH2P	6260-9-240	MHPKMV	6260-9-344	MHRLP0	6260-9-244	MHTKPW	6260-9-240	MJANH0
6260-9-340	MHNH3A	6260-9-240	MHPKMW	6260-9-344	MHRLP1	6260-9-244	MHTKTW	6260-9-240	MJANHA
6260-9-240	MHNH58	6260-9-340	MHPRVE	6260-9-344	MHRLP2	6260-9-244	MHTKV9	6260-9-240	MJANV4
6260-9-340	MHNH78	6260-9-340	MHPRVH	6260-9-344	MHRLP3	6260-9-244	MHTKVR	6260-9-240	MJAP17
6260-9-340	MHNH79	6260-9-340	MHPRVJ	6260-9-344	MHRLP4	6260-9-244	MHTL8M	6260-9-240	MJARON
6260-9-340	MHNH8A	6260-9-340	MHPRVK	6260-9-344	MHRLP5	6260-9-244	MHTLPW	6260-9-240	MJAR4N
6260-9-340	MHNW50	6260-9-340	MHPT78	6260-9-344	MHRMT9	6260-9-244	MHTLPX	6260-9-240	MJARDT
6260-9-340	MHNW52	6260-9-340	MHPT79	6260-9-344	MHRMTL	6260-9-244	MHTLPY	6260-9-240	MJD7T5
6260-9-340	MHNW54	6260-9-340	MHPT7A	6260-9-344	MHRMVK	6260-9-244	MHTLYE	6260-9-240	MJD7Y5

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-244	MJEKT0	6260-9-240	MJJEVW	6260-9-240	MJKAYA	6260-9-244	MJKRPT	6260-9-340	MJLE6T
6260-9-244	MJEKXK	6260-9-244	MJJEW8	6260-9-340	MJKAYD	6260-9-244	MJKRPV	6260-9-340	MJLE6V
6260-9-240	MJEMWW	6260-9-244	MJJEW9	6260-9-340	MJKAYE	6260-9-244	MJKTPY	6260-9-240	MJLEPM
6260-9-240	MJEMWY	6260-9-244	MJJEWA	6260-9-340	MJKAYH	6260-9-244	MJKTR0	6260-9-240	MJLH3L
6260-9-240	MJEMX0	6260-9-244	MJJEDW	6260-9-240	MJKD22	6260-9-244	MJKTR1	6260-9-240	MJLHE4
6260-9-240	MJEMX2	6260-9-244	MJJH4L	6260-9-240	MJKD53	6260-9-244	MJKTR2	6260-9-240	MJLHLY
6260-9-240	MJEMYL	6260-9-244	MJJH4M	6260-9-240	MJKD95	6260-9-240	MJL08R	6260-9-240	MJLHX8
6260-9-240	MJEMYP	6260-9-244	MJJHA1	6260-9-240	MJKDAX	6260-9-240	MJL246	6260-9-240	MJLK8N
6260-9-240	MJLKJW	6260-9-244	MJMN98	6260-9-240	MJNKR1	6260-9-240	MJPD9X	6260-9-240	MJRXAA
6260-9-340	MJLLR2	6260-9-340	MJMN9X	6260-9-340	MJNKY3	6260-9-240	MJPE02	6260-9-240	MJRYTH
6260-9-340	MJLLX2	6260-9-340	MJMN9V	6260-9-340	MJNL8Y	6260-9-244	MJPE1L	6260-9-340	MJT01D
6260-9-240	MJLMA5	6260-9-344	MJMNYP	6260-9-240	MJNLVE	6260-9-244	MJPE3K	6260-9-244	MJT071
6260-9-240	MJLME0	6260-9-240	MJMP6J	6260-9-340	MJNM6V	6260-9-340	MJPE5E	6260-9-240	MJT0ND
6260-9-240	MJLMX9	6260-9-240	MJMP77	6260-9-240	MJNMY4	6260-9-244	MJPH9W	6260-9-240	MJT118
6260-9-340	MJLNPP	6260-9-240	MJMR15	6260-9-240	MJNNDN	6260-9-240	MJPJ6P	6260-9-240	MJT13L
6260-9-240	MJLP2T	6260-9-240	MJMR1W	6260-9-240	MJNNMT	6260-9-244	MJPJNN	6260-9-244	MJT143
6260-9-240	MJLP7W	6260-9-340	MJMR9X	6260-9-240	MJNNRN	6260-9-340	MJPK30	6260-9-340	MJT22E
6260-9-240	MJLPET	6260-9-240	MJMRD0	6260-9-340	MJNP6N	6260-9-240	MJPK97	6260-9-344	MJT2AH
6260-9-340	MJLPLE	6260-9-240	MJMRT2	6260-9-240	MJNP71	6260-9-244	MJPKEA	6260-9-244	MJT2NH
6260-9-240	MJLR9L	6260-9-240	MJMT7M	6260-9-240	MJNPH9	6260-9-240	MJPKHK	6260-9-340	MJT2W8
6260-9-340	MJLRWM	6260-9-340	MJMV1Y	6260-9-340	MJNPT4	6260-9-244	MJPKM3	6260-9-244	MJT4DM
6260-9-240	MJLRWW	6260-9-244	MJMV23	6260-9-240	MJNR4J	6260-9-240	MJPKPR	6260-9-244	MJT4H5
6260-9-240	MJLTON	6260-9-244	MJMVJ4	6260-9-240	MJNR7E	6260-9-340	MJPL95	6260-9-244	MJT4NN
6260-9-340	MJLT85	6260-9-244	MJMW3R	6260-9-240	MJNRE9	6260-9-244	MJPLD7	6260-9-244	MJT5D1
6260-9-340	MJLTT1	6260-9-344	MJMWE8	6260-9-240	MJNRJ3	6260-9-240	MJPM2M	6260-9-340	MJT61R
6260-9-240	MJLV01	6260-9-340	MJMWP4	6260-9-340	MJNRN2	6260-9-240	MJPMY1	6260-9-244	MJT62X
6260-9-240	MJLV4X	6260-9-340	MJMX4D	6260-9-240	MJNT90	6260-9-240	MJPNXW	6260-9-444	MJT678
6260-9-340	MJLWMR	6260-9-240	MJMY1X	6260-9-240	MJNTAK	6260-9-240	MJPP65	6260-9-340	MJT6LP
6260-9-240	MJLWV2	6260-9-244	MJMYDJ	6260-9-240	MJNTEY	6260-9-240	MJPP8A	6260-9-240	MJT6RK
6260-9-340	MJLWX3	6260-9-340	MJMYN7	6260-9-240	MJNTMP	6260-9-244	MJPPVT	6260-9-240	MJT6YY
6260-9-240	MJLWXW	6260-9-344	MJMYPE	6260-9-340	MJNV7X	6260-9-244	MJPR33	6260-9-244	MJT71L
6260-9-340	MJLY01	6260-9-340	MJMYWT	6260-9-240	MJNVJ5	6260-9-240	MJPR50	6260-9-244	MJT78W
6260-9-240	MJLY4R	6260-9-340	MJN008	6260-9-240	MJNW4W	6260-9-240	MJPRMD	6260-9-340	MJT7ER
6260-9-240	MJLYAA	6260-9-340	MJN0AA	6260-9-244	MJNWDW	6260-9-244	MJPTEW	6260-9-340	MJT7P2
6260-9-340	MJLYHM	6260-9-240	MJN0EE	6260-9-240	MJNWE5	6260-9-240	MJPV5J	6260-9-340	MJT88L
6260-9-340	MJLYJR	6260-9-240	MJN259	6260-9-240	MJNWX2	6260-9-340	MJPVH8	6260-9-240	MJT88V
6260-9-240	MJLYT6	6260-9-240	MJN2WR	6260-9-240	MJNWKD	6260-9-244	MJPVPK	6260-9-340	MJT8RJ
6260-9-340	MJM07M	6260-9-240	MJN3JM	6260-9-240	MJNWW5	6260-9-240	MJPWPY	6260-9-240	MJT95E
6260-9-340	MJM0EW	6260-9-240	MJN3R5	6260-9-240	MJNX1T	6260-9-240	MJPX9D	6260-9-444	MJT9KD
6260-9-244	MJM1NY	6260-9-240	MJN3VY	6260-9-240	MJNX7A	6260-9-240	MJPXH3	6260-9-244	MJT9PW
6260-9-240	MJM24V	6260-9-240	MJN4HV	6260-9-240	MJNX9Y	6260-9-240	MJPXJP	6260-9-340	MJT9R1
6260-9-340	MJM35W	6260-9-244	MJN5E7	6260-9-240	MJNXNM	6260-9-244	MJPY01	6260-9-340	MJT9YM
6260-9-244	MJM37J	6260-9-340	MJN5J4	6260-9-340	MJNYWY	6260-9-240	MJPY7R	6260-9-240	MJTANL
6260-9-340	MJM5LE	6260-9-340	MJN5W6	6260-9-240	MJP13Y	6260-9-240	MJPYDL	6260-9-444	MJTAPX
6260-9-340	MJM5P2	6260-9-240	MJN607	6260-9-340	MJP1EX	6260-9-244	MJPYLM	6260-9-340	MJTAV9
6260-9-340	MJM5TP	6260-9-240	MJN72A	6260-9-244	MJP1V2	6260-9-240	MJR3M8	6260-9-244	MJTAVJ
6260-9-344	MJM61X	6260-9-340	MJN7ME	6260-9-244	MJP23E	6260-9-240	MJR43J	6260-9-240	MJTAVP
6260-9-240	MJM69W	6260-9-340	MJN7NJ	6260-9-340	MJP2RW	6260-9-240	MJR49T	6260-9-240	MJTDHM
6260-9-340	MJM6JT	6260-9-240	MJN7VD	6260-9-244	MJP39P	6260-9-240	MJR610	6260-9-240	MJTDKL
6260-9-244	MJM6WN	6260-9-244	MJN8K4	6260-9-244	MJP3K1	6260-9-240	MJR62A	6260-9-244	MJTE4W
6260-9-240	MJM84N	6260-9-340	MJN8ME	6260-9-244	MJP3KN	6260-9-240	MJR76M	6260-9-340	MJTEXR
6260-9-244	MJM9DK	6260-9-240	MJN9L2	6260-9-240	MJP3LV	6260-9-340	MJR7PA	6260-9-240	MJTHT4
6260-9-244	MJM9J1	6260-9-340	MJN9M4	6260-9-240	MJP3N8	6260-9-240	MJRAY4	6260-9-240	MJTHW8
6260-9-240	MJM9XW	6260-9-340	MJN9VY	6260-9-340	MJP3V0	6260-9-240	MJRDRR	6260-9-340	MJTHWK
6260-9-244	MJMA24	6260-9-240	MJNA64	6260-9-240	MJP4AX	6260-9-240	MJRE9M	6260-9-240	MJTJ1E
6260-9-340	MJMAA1	6260-9-240	MJNEV4	6260-9-240	MJP4TA	6260-9-240	MJRH62	6260-9-340	MJTJXM
6260-9-340	MJMADP	6260-9-240	MJNH07	6260-9-240	MJP5J6	6260-9-240	MJRJ33	6260-9-244	MJTJYW
6260-9-240	MJMD14	6260-9-240	MJNH2X	6260-9-240	MJP63V	6260-9-240	MJRKR5	6260-9-244	MJTJOL
6260-9-240	MJMDLY	6260-9-240	MJNH7M	6260-9-244	MJP6L6	6260-9-240	MJRNL6	6260-9-240	MJTJTP
6260-9-244	MJMEM5	6260-9-240	MJNHAK	6260-9-240	MJP730	6260-9-240	MJRP5N	6260-9-340	MJTJTM

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-340	MJMEX5	6260-9-240	MJNHEW	6260-9-240	MJP76D	6260-9-240	MJRR16	6260-9-340	MJTKW3
6260-9-340	MJMH7K	6260-9-340	MJNHL7	6260-9-244	MJP8LR	6260-9-240	MJRRJM	6260-9-244	MJTL4K
6260-9-340	MJMH6X	6260-9-240	MJNJ38	6260-9-240	MJP8M4	6260-9-240	MJRT8E	6260-9-244	MJTLAP
6260-9-240	MJMJJ0	6260-9-240	MJNJRK	6260-9-240	MJP8TM	6260-9-340	MJRTT1	6260-9-244	MJTLT8
6260-9-340	MJMKJ0	6260-9-244	MJNJVE	6260-9-240	MJP8WN	6260-9-240	MJRV5V	6260-9-244	MJTLWK
6260-9-240	MJML17	6260-9-340	MJNK6K	6260-9-240	MJP8XJ	6260-9-240	MJRV9A	6260-9-340	MJTLYA
6260-9-344	MJML4T	6260-9-240	MJNKEN	6260-9-244	MJPA1D	6260-9-240	MJRVAJ	6260-9-340	MJTMET
6260-9-240	MJMLJ5	6260-9-240	MJNKJ3	6260-9-340	MJPAHW	6260-9-340	MJRVTO	6260-9-244	MJTMME
6260-9-340	MJMM1W	6260-9-340	MJNKP8	6260-9-240	MJPAVJ	6260-9-240	MJRW3M	6260-9-244	MJTM3T
6260-9-340	MJTND4	6260-9-240	N1JMDDX1	6260-9-244	P8VMPA	6260-9-240	RYRMPD	6260-9-240	WDMME
6260-9-340	MJTNE1	6260-9-240	N1KMDD	6260-9-244	P8VMPAX1	6260-9-240	RYTMPD	6260-9-240	WDMME
6260-9-444	MJTNMW	6260-9-340	N4JMPA	6260-9-244	P8VMPAX2	6260-9-240	TOHMPD	6260-9-340	WE4MLE
6260-9-340	MJTNP6	6260-9-340	N4KMPA	6260-9-444	P8XMAD	6260-9-240	TOJME	6260-9-340	WE5MLE
6260-9-244	MJTNNW	6260-9-340	N4LMPA	6260-9-444	P95MAD	6260-9-240	TOPMPD	6260-9-340	WE6MLE
6260-9-444	MJTP17	6260-9-244	N61MMA	6260-9-244	PARMLA	6260-9-240	TORMPD	6260-9-340	WE8MLE
6260-9-340	MJTP5R	6260-9-244	N63MMA	6260-9-244	PARMLAX1	6260-9-240	TOTMPD	6260-9-240	WJOMME
6260-9-240	MJTR5P	6260-9-244	N65MMA	6260-9-244	PATMLA	6260-9-240	T19MAE	6260-9-240	WK5MME
6260-9-444	MJTT9M	6260-9-244	N65MMAX1	6260-9-244	PAVMLA	6260-9-240	T1AMAE	6260-9-240	WL6MME
6260-9-340	MJTTMJ	6260-9-244	N65MMAX2	6260-9-244	PAVMLAX1	6260-9-240	T1DMAE	6260-9-240	WLYMME
6260-9-240	MJTNN1	6260-9-244	N66MMA	6260-9-244	PAWMLA	6260-9-240	T1EMAE	6260-9-340	WN2MME
6260-9-344	MJTVRM	6260-9-244	N67MMA	6260-9-244	PAXMLA	6260-9-240	T4DMAE	6260-9-340	WN3MME
6260-9-240	MJTW3H	6260-9-244	N68MMA	6260-9-244	PAYMLA	6260-9-244	T4VMDE	6260-9-244	WNEMMA
6260-9-344	MJTW42	6260-9-244	N68MMAX1	6260-9-244	PD0MLA	6260-9-244	T4WMDE	6260-9-244	WNHMMMA
6260-9-240	MJTWEX	6260-9-444	N71MAD	6260-9-240	PKAMDE	6260-9-244	T4XMDE	6260-9-244	WNJMMMA
6260-9-240	MJTX9J	6260-9-444	N72MAD	6260-9-244	PN7MDE	6260-9-244	T4YMDE	6260-9-240	WNXMME
6260-9-344	MJTXE1	6260-9-244	NH3MTA	6260-9-244	PN8MDE	6260-9-340	T69MNE	6260-9-244	WP7MMA
6260-9-244	MJTXHK	6260-9-244	NH3MTAX1	6260-9-244	PNAMDE	6260-9-340	T6HMNE	6260-9-244	WP7MMAX1
6260-9-340	MJTXHV	6260-9-244	NH9MTA	6260-9-244	PNDMDE	6260-9-340	T9MMLE	6260-9-244	WPJMMMA
6260-9-344	MJTY72	6260-9-244	NHHMTA	6260-9-240	PNEMDE	6260-9-340	T9PMLE	6260-9-244	WPNMMMA
6260-9-240	MJTYEM	6260-9-244	NHHMTAX1	6260-9-240	PNHMDE	6260-9-340	T9RMLE	6260-9-244	WPNMMAX1
6260-9-240	MJTYP2	6260-9-244	NJ2MMD	6260-9-240	PNJMDE	6260-9-340	T9TMLE	6260-9-340	X37MKA
6260-9-240	MJVMKE	6260-9-244	NJYMPD	6260-9-240	PNKMDE	6260-9-340	TADMNE	6260-9-340	X37MKAX1
6260-9-240	MK7MKE	6260-9-444	NK5MAD	6260-9-444	PNYMAD	6260-9-340	TAJMNE	6260-9-340	X3HMKA
6260-9-240	MK8MKE	6260-9-444	NK6MAD	6260-9-340	PT7MMA	6260-9-244	TMHMMD	6260-9-340	X3HMKAX1
6260-9-240	MK9MKE	6260-9-244	NK6MPD	6260-9-340	PT8MMA	6260-9-244	TRNMMD	6260-9-340	X3KMKA
6260-9-240	MKA1A1	6260-9-444	NK8MAD	6260-9-340	PT8MMAX1	6260-9-240	TT6MPA	6260-9-340	X3KMKAX1
6260-9-240	MKA1T2	6260-9-344	NKHMPA	6260-9-244	PVXMMMA	6260-9-240	TT7MPA	6260-9-340	X3MMKA
6260-9-240	MKA367	6260-9-344	NKHMPAX1	6260-9-244	PVYMMMA	6260-9-240	TT8MPA	6260-9-340	X3WMKA
6260-9-240	MKA4Y9	6260-9-240	NR6MDD	6260-9-244	PVYMMAX1	6260-9-340	TTAMPA	6260-9-244	X84MLE
6260-9-240	MKA6D4	6260-9-240	NRAMDD	6260-9-244	PW0MMA	6260-9-340	TTEMPA	6260-9-244	X85MLE
6260-9-240	MKA7V5	6260-9-240	NRAMDDX1	6260-9-244	PW0MMAX1	6260-9-340	TTHMPA	6260-9-340	XAVMDE
6260-9-240	MKA82V	6260-9-240	NRNMDD	6260-9-244	PW1MMA	6260-9-244	TX8MMD	6260-9-340	XAXMDE
6260-9-240	MKA90V	6260-9-244	NVHMLA	6260-9-244	PW2MMA	6260-9-244	VOPMLE	6260-9-244	XDVMLE
6260-9-240	MKAAJW	6260-9-244	NVMMLA	6260-9-240	R1TMMA	6260-9-244	VOTMLE	6260-9-244	XDXMLE
6260-9-240	MKAH1A	6260-9-244	NVNMLA	6260-9-240	R1TMMAX1	6260-9-244	VOVMLE	6260-9-240	XL2MNE
6260-9-240	MKAHMA	6260-9-244	P0PMPD	6260-9-240	R1TMMAX2	6260-9-244	V0VMLE	6260-9-240	XNXMHE
6260-9-240	MKAHYR	6260-9-244	P0RMPD	6260-9-240	R1VMMA	6260-9-244	V2DMMD	6260-9-340	XPJMKA
6260-9-240	MKAKP5	6260-9-240	P19MLE	6260-9-240	R1WMMA	6260-9-244	V3YMED	6260-9-340	XR0MKA
6260-9-240	MKAY76	6260-9-240	P24MLE	6260-9-240	RANMDD	6260-9-244	V42MED	6260-9-340	XR0MKAX1
6260-9-240	MKAY8R	6260-9-240	P25MLE	6260-9-240	RAPMDD	6260-9-244	V46MED	6260-9-240	XR5MKA
6260-9-240	MPNMEE	6260-9-240	P26MLE	6260-9-340	RL5MMA	6260-9-244	V4DMED	6260-9-340	XW2MKA
6260-9-240	MPRMEE	6260-9-444	P2MMAD	6260-9-340	RL5MMAX1	6260-9-240	VJ8MEE	6260-9-340	XW3MKA
6260-9-240	MPTMEE	6260-9-340	P2XMDD	6260-9-340	RPAMMA	6260-9-240	VJ9MEE	6260-9-240	XWKMHE
6260-9-240	MPVMEE	6260-9-340	P30MDD	6260-9-340	RPDMMA	6260-9-240	VJAMEE	6260-9-240	XWKMHX1
6260-9-244	MW4MRD	6260-9-244	P3RMLA	6260-9-340	RPDMMAX1	6260-9-240	VJDMEE	6260-9-240	XX3MHE
6260-9-244	MWHMRD	6260-9-244	P3TMLA	6260-9-340	RPEMMA	6260-9-240	VVYMEE	6260-9-240	XY0MNE
6260-9-340	MY3MED	6260-9-244	P3VMLA	6260-9-340	RPEMMAX1	6260-9-240	VYWMME	6260-9-240	XY7MNE
6260-9-244	MY3MRD	6260-9-244	P3WMLA	6260-9-340	RPHMMA	6260-9-240	VYXMME	6260-9-240	Y23MHE
6260-9-340	MY4MED	6260-9-244	P3WMLAX1	6260-9-340	RPLMMA	6260-9-240	VYXMME	6260-9-240	Y64MNE
6260-9-340	MY4MEDX1	6260-9-244	P3XMLA	6260-9-340	RPNMMA	6260-9-240	W00MME	6260-9-240	Y66MNE

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-236	11991501	6260-9-236	20983001	6260-9-236	29209902	6260-9-236	108MPD	6260-9-236	20142803X1X1
6260-9-236	13692402	6260-9-236	20983002	6260-9-236	29273301	6260-9-236	109MPD	6260-9-236	20142803X2
6260-9-236	14973301	6260-9-236	20983003	6260-9-236	30032402	6260-9-236	12LMTD	6260-9-236	20142804X1
6260-9-236	15186903	6260-9-236	20983004	6260-9-236	31337901	6260-9-236	12WMEE	6260-9-236	20142804X3
6260-9-236	19806601	6260-9-236	21035301	6260-9-236	32693902	6260-9-236	12XMEE	6260-9-236	20142805X1
6260-9-236	19865901	6260-9-236	21035302	6260-9-236	33187601	6260-9-236	12YMEE	6260-9-236	20142805X2
6260-9-236	19865902	6260-9-236	21035303	6260-9-236	34354901	6260-9-236	153MDE	6260-9-236	20368701X1
6260-9-236	19865903	6260-9-236	21035304	6260-9-236	34711701	6260-9-236	154MDE	6260-9-236	20368701X2
6260-9-236	19865904	6260-9-236	21035305	6260-9-236	35921901	6260-9-236	155MDE	6260-9-236	20368702X1
6260-9-236	19865905	6260-9-236	21164801	6260-9-236	36141505	6260-9-236	17KMDE	6260-9-236	20368703X1
6260-9-236	19880501	6260-9-236	21164802	6260-9-236	36923401	6260-9-236	17KMDEX1	6260-9-236	20368704X1
6260-9-236	19880502	6260-9-236	21164803	6260-9-236	37063304	6260-9-236	18PMDE	6260-9-236	20368705X1
6260-9-236	19880503	6260-9-236	21164804	6260-9-236	37194801	6260-9-236	19806602X1	6260-9-236	20368706X1
6260-9-236	19880504	6260-9-236	21220201	6260-9-236	37262101	6260-9-236	19865904X1	6260-9-236	20368707X1
6260-9-236	19880505	6260-9-236	21220202	6260-9-236	37381201	6260-9-236	19865904X2	6260-9-236	20368708X1
6260-9-236	19880506	6260-9-236	21220203	6260-9-236	37381203	6260-9-236	19865905X1	6260-9-236	20783304X1
6260-9-236	20113201	6260-9-236	21220204	6260-9-236	37381204	6260-9-236	19880501X1	6260-9-236	20783305X1
6260-9-236	20113202	6260-9-236	21220205	6260-9-236	37633001	6260-9-236	19880501X1X1	6260-9-236	20783306X1
6260-9-236	20113203	6260-9-236	21220206	6260-9-236	37633003	6260-9-236	19880502X1X1	6260-9-236	20904201X1
6260-9-236	20113204	6260-9-236	21263501	6260-9-236	37633004	6260-9-236	19880503X1	6260-9-236	20934304X1
6260-9-236	20142801	6260-9-236	21263502	6260-9-236	37633006	6260-9-236	19880504X1	6260-9-236	20AMMA
6260-9-236	20142802	6260-9-236	21263503	6260-9-236	38005403	6260-9-236	19880504X2	6260-9-236	20DMMA
6260-9-236	20142803	6260-9-236	21263504	6260-9-236	38005404	6260-9-236	19880505X1	6260-9-236	20HMMA
6260-9-236	20142804	6260-9-236	21263505	6260-9-236	38005405	6260-9-236	19880505X1X1	6260-9-236	21164803X1
6260-9-236	20142805	6260-9-236	21297001	6260-9-236	38271201	6260-9-236	19880505X2	6260-9-236	21220201X1
6260-9-236	20368701	6260-9-236	21297002	6260-9-236	38271202	6260-9-236	19880506X1	6260-9-236	21220206X1
6260-9-236	20368702	6260-9-236	21297003	6260-9-236	38271203	6260-9-236	1AMMEE	6260-9-236	21263502X1
6260-9-236	20368703	6260-9-236	21297004	6260-9-236	38271204	6260-9-236	1EP02H	6260-9-236	21297002X1
6260-9-236	20368704	6260-9-236	21297005	6260-9-236	38271205	6260-9-236	1EP07W	6260-9-236	21297003X1
6260-9-236	20368705	6260-9-236	21318701	6260-9-236	00MMKE	6260-9-236	1EP07XX1	6260-9-236	21297005X1
6260-9-236	20368706	6260-9-236	21318703	6260-9-236	00PMKE	6260-9-236	1EP07XX2	6260-9-236	21318702X1
6260-9-236	20368707	6260-9-236	21318704	6260-9-236	00TMKE	6260-9-236	1EP082	6260-9-236	21318706X1
6260-9-236	20368708	6260-9-236	21318705	6260-9-236	01RMTD	6260-9-236	1EP083	6260-9-236	21667802X1
6260-9-236	20637501	6260-9-236	21318706	6260-9-236	01TMTD	6260-9-236	1EP085	6260-9-236	245MLE
6260-9-236	20637502	6260-9-236	21451801	6260-9-236	01WMTD	6260-9-236	1H4MEE	6260-9-236	25RMNA
6260-9-236	20637503	6260-9-236	21451802	6260-9-236	03KMLE	6260-9-236	1H4MEEX1	6260-9-236	25RMNAX1
6260-9-236	20637504	6260-9-236	21451803	6260-9-236	0K8MTD	6260-9-236	1HN01H	6260-9-236	25VMNA
6260-9-236	20749503	6260-9-236	21451804	6260-9-236	0KMMTD	6260-9-236	1HN01J	6260-9-236	25VMNAX1
6260-9-236	20769701	6260-9-236	21451805	6260-9-236	0KPMLE	6260-9-236	1HN01R	6260-9-236	25VMNA
6260-9-236	20769702	6260-9-236	21451806	6260-9-236	0L2MTD	6260-9-236	1HN025	6260-9-236	261MHD
6260-9-236	20769703	6260-9-236	21667801	6260-9-236	0L5MMD	6260-9-236	1HN026	6260-9-236	26KMNE
6260-9-236	20769704	6260-9-236	21667802	6260-9-236	0N5MEE	6260-9-236	1K0MED	6260-9-236	27DMNE
6260-9-236	20769705	6260-9-236	21667803	6260-9-236	0N8MEE	6260-9-236	1KVMED	6260-9-236	27WMLE
6260-9-236	20783301	6260-9-236	21667804	6260-9-236	0NAMEE	6260-9-236	1KVMEDX1	6260-9-236	296MHD
6260-9-236	20783302	6260-9-236	21667805	6260-9-236	0R5MPA	6260-9-236	1KXMDE	6260-9-236	296MHDX1
6260-9-236	20783303	6260-9-236	21726401	6260-9-236	0R5MPAX1	6260-9-236	1LRMED	6260-9-236	299MKD
6260-9-236	20783304	6260-9-236	21781902	6260-9-236	0RDMPA	6260-9-236	1N2MHE	6260-9-236	29RMKD
6260-9-236	20783305	6260-9-236	22268801	6260-9-236	0RDMPAX2	6260-9-236	1N3MHE	6260-9-236	29RMKDX1
6260-9-236	20783306	6260-9-236	22268802	6260-9-236	0RJMPA	6260-9-236	1N4MHE	6260-9-236	29TMKD
6260-9-236	20783307	6260-9-236	22268804	6260-9-236	0RJMPAX1	6260-9-236	1VMMNA	6260-9-236	2A6MEE
6260-9-236	20783308	6260-9-236	22268805	6260-9-236	0RLMPA	6260-9-236	1VNMNA	6260-9-236	2ATMNE
6260-9-236	20904201	6260-9-236	22268806	6260-9-236	0V1MMD	6260-9-236	20113201X1	6260-9-236	2D2MNE
6260-9-236	20904202	6260-9-236	22358202	6260-9-236	0V2MMD	6260-9-236	20113202X1	6260-9-236	2D4MNE
6260-9-236	20904203	6260-9-236	22358204	6260-9-236	0W8MMD	6260-9-236	20113202X1X1	6260-9-236	2D4MNEX1
6260-9-236	20904204	6260-9-236	22358206	6260-9-236	0W9MMD	6260-9-236	20113203X1	6260-9-236	2EHMLE
6260-9-236	20934301	6260-9-236	22379903	6260-9-236	0WJMMD	6260-9-236	20113204X1	6260-9-236	2ETMKD
6260-9-236	20934302	6260-9-236	22379906	6260-9-236	0WKMMMD	6260-9-236	20142801X1	6260-9-236	2LPMJE
6260-9-236	20934303	6260-9-236	24269001	6260-9-236	0Y0MMD	6260-9-236	20142802X1	6260-9-236	2LRMJE
6260-9-236	20934304	6260-9-236	25406902	6260-9-236	106MPD	6260-9-236	20142802X2	6260-9-236	2LTMJE
6260-9-236	20950001	6260-9-236	28846502	6260-9-236	107MPD	6260-9-236	20142803X1	6260-9-236	2M1MLE

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-236	2NEMTD	6260-9-236	4M0MHE	6260-9-236	73PMDE	6260-9-236	90MMEE	6260-9-236	A1EMPD
6260-9-236	2R5MHD	6260-9-236	4M1MHE	6260-9-236	742MDE	6260-9-236	90NMEE	6260-9-236	A3XMHE
6260-9-236	2WEMND	6260-9-236	4VLMKD	6260-9-236	743MDE	6260-9-236	90NMRD	6260-9-236	A5DMPD
6260-9-236	2WNMTD	6260-9-236	4VMMKD	6260-9-236	75DMPA	6260-9-236	90PMEE	6260-9-236	A5KMME
6260-9-236	30VMHD	6260-9-236	4VNMKD	6260-9-236	75DMPAX2	6260-9-236	91NMRD	6260-9-236	A5LMME
6260-9-236	315MTD	6260-9-236	4VPMKD	6260-9-236	75EMPA	6260-9-236	91RMJE	6260-9-236	A6PMPD
6260-9-236	33LMTD	6260-9-236	4WRMKD	6260-9-236	75EMPAX1	6260-9-236	91TMJE	6260-9-236	A90MJE
6260-9-236	33MMLE	6260-9-236	4WTMKD	6260-9-236	75JMPA	6260-9-236	91VMJE	6260-9-236	A9EMJE
6260-9-236	33MMTD	6260-9-236	51VMJE	6260-9-236	7L4MED	6260-9-236	932MJA	6260-9-236	AAAMJE
6260-9-236	344MEA	6260-9-236	51WMJE	6260-9-236	7LLMLA	6260-9-236	93NMDA	6260-9-236	AADMJE
6260-9-236	345MEA	6260-9-236	51XMJE	6260-9-236	7R4MED	6260-9-236	93NMDAS1	6260-9-236	ADWMRD
6260-9-236	346MEA	6260-9-236	58XMLD	6260-9-236	7V9MHHE	6260-9-236	93PMDA	6260-9-236	AE3MHA
6260-9-236	34711702D	6260-9-236	597MLD	6260-9-236	7VDMHE	6260-9-236	93PMDAS1	6260-9-236	AE4MHA
6260-9-236	347MEA	6260-9-236	599MLD	6260-9-236	7VEMHE	6260-9-236	93RMDA	6260-9-236	AH2MEE
6260-9-236	348MEA	6260-9-236	59WMLD	6260-9-236	7W3MLA	6260-9-236	93RMDAS1	6260-9-236	AH8MRD
6260-9-236	353MHD	6260-9-236	5DMMLD	6260-9-236	7WHMHE	6260-9-236	93TMDA	6260-9-236	AH8MRDX1
6260-9-236	365MRD	6260-9-236	5ELMLD	6260-9-236	7WJMHE	6260-9-236	93TMDAS1	6260-9-236	AK4MRD
6260-9-236	377MND	6260-9-236	5KJMHD	6260-9-236	7WKMHHE	6260-9-236	957MND	6260-9-236	ALAMEE
6260-9-236	38449901C	6260-9-236	5KJMHDX1	6260-9-236	7X3MHHE	6260-9-236	958MND	6260-9-236	ALDMEE
6260-9-236	38449902X	6260-9-236	5L1MHD	6260-9-236	7X4MHHE	6260-9-236	95AMND	6260-9-236	AXKMPD
6260-9-236	38449903X	6260-9-236	5MPMHA	6260-9-236	7X5MHHE	6260-9-236	96JMJA	6260-9-236	AXKMPDX1
6260-9-236	38449904X	6260-9-236	5MRMHA	6260-9-236	7YPMCA	6260-9-236	96KMJA	6260-9-236	AYLMHA
6260-9-236	38AMHD	6260-9-236	5MWVMA	6260-9-236	7YVMCA	6260-9-236	974MJA	6260-9-236	AYMMHA
6260-9-236	38EMHD	6260-9-236	5MYMHA	6260-9-236	7YVMCAS1	6260-9-236	97XMME	6260-9-236	D07MED
6260-9-236	38HMHD	6260-9-236	5R2MHE	6260-9-236	7YWMCA	6260-9-236	97YMME	6260-9-236	D63MED
6260-9-236	394MLE	6260-9-236	5R3MHE	6260-9-236	7YYMCA	6260-9-236	980MME	6260-9-236	D70MEE
6260-9-236	3AYMLE	6260-9-236	5R4MHE	6260-9-236	82YMEA	6260-9-236	9A8MJA	6260-9-236	D71MEE
6260-9-236	3RXMMA	6260-9-236	5R7MDD	6260-9-236	832MEA	6260-9-236	9AVMRD	6260-9-236	D72MEE
6260-9-236	3RXMMAX1	6260-9-236	5R8MDD	6260-9-236	833MEA	6260-9-236	9H4MRD	6260-9-236	D8PMTA
6260-9-236	3RYMMA	6260-9-236	5R9MDD	6260-9-236	846MME	6260-9-236	9JYMHA	6260-9-236	D8RMTA
6260-9-236	3T0MMA	6260-9-236	5T0MJE	6260-9-236	847MME	6260-9-236	9K0MHA	6260-9-236	D8VMTA
6260-9-236	3T1MMA	6260-9-236	5T5MJE	6260-9-236	84HMED	6260-9-236	9K1MPE	6260-9-236	D8WMTA
6260-9-236	3T3MMA	6260-9-236	668MME	6260-9-236	877MED	6260-9-236	9K2MPE	6260-9-236	D95MHA
6260-9-236	3T4MMA	6260-9-236	691MND	6260-9-236	877MME	6260-9-236	9KJMJA	6260-9-236	D97MHA
6260-9-236	3T4MND	6260-9-236	694MND	6260-9-236	87VMED	6260-9-236	9L8MPE	6260-9-236	DA3MRD
6260-9-236	3TPMRD	6260-9-236	696MND	6260-9-236	89EMDE	6260-9-236	9M5MND	6260-9-236	DDMMLA
6260-9-236	3V1MKD	6260-9-236	6AKMPA	6260-9-236	89HMDE	6260-9-236	9PJMDE	6260-9-236	DDTMLA
6260-9-236	3WYMKD	6260-9-236	6ANMPA	6260-9-236	89JMDE	6260-9-236	9PLMDE	6260-9-236	DH3MED
6260-9-236	3YTMND	6260-9-236	6ANMPAX1	6260-9-236	8AAMED	6260-9-236	9PLMDX1	6260-9-236	DKJMHA
6260-9-236	40NMRD	6260-9-236	6APMPA	6260-9-236	8DDMRD	6260-9-236	9PMMDE	6260-9-236	DL4MJE
6260-9-236	40RMRD	6260-9-236	6EJMLE	6260-9-236	8ENMEA	6260-9-236	9R4MJD	6260-9-236	DL6MED
6260-9-236	412MKD	6260-9-236	6EXMND	6260-9-236	8L2MJD	6260-9-236	9R5MJD	6260-9-236	DMRMJE
6260-9-236	431MRD	6260-9-236	6EYMND	6260-9-236	8L3MJD	6260-9-236	9RDMJD	6260-9-236	DNXMED
6260-9-236	432MRD	6260-9-236	6H0MND	6260-9-236	8L4MJD	6260-9-236	9RHMJD	6260-9-236	DP2MKA
6260-9-236	434MRD	6260-9-236	6HHMND	6260-9-236	8L5MJD	6260-9-236	9RJMJJD	6260-9-236	DPHMJE
6260-9-236	44MMND	6260-9-236	6HLMKA	6260-9-236	8L6MJD	6260-9-236	9RKMJD	6260-9-236	DR2MKA
6260-9-236	44PMLD	6260-9-236	6L6MME	6260-9-236	8LJMED	6260-9-236	9RLMJJD	6260-9-236	DR3MKA
6260-9-236	455MLD	6260-9-236	6MVMLD	6260-9-236	8LLMED	6260-9-236	9RMMJD	6260-9-236	E04MRD
6260-9-236	477MKD	6260-9-236	6N5MPD	6260-9-236	8LMMED	6260-9-236	9TLMJA	6260-9-236	E18MKE
6260-9-236	47TMLD	6260-9-236	6N5MPDX1	6260-9-236	8M1MRD	6260-9-236	9V8MJA	6260-9-236	E1WMKA
6260-9-236	48JMJE	6260-9-236	6N6MPD	6260-9-236	8M1MRDX1	6260-9-236	9VAMHE	6260-9-236	E22MHA
6260-9-236	48JMJEX1	6260-9-236	6N7MPD	6260-9-236	8THMRD	6260-9-236	9WAMEE	6260-9-236	E23MHA
6260-9-236	48LMJE	6260-9-236	6NHMKA	6260-9-236	8WNMLA	6260-9-236	9WAMHA	6260-9-236	E2NMKA
6260-9-236	4AEMMD	6260-9-236	6NJMKA	6260-9-236	8WPMLA	6260-9-236	9WDMHA	6260-9-236	E30MKA
6260-9-236	4ARMJE	6260-9-236	6NKMKA	6260-9-236	8WRMLA	6260-9-236	9Y4MEE	6260-9-236	E32MKA
6260-9-236	4D0MKD	6260-9-236	6NLMKA	6260-9-236	8WTMLA	6260-9-236	9Y7MEE	6260-9-236	E36MKA
6260-9-236	4EVMKD	6260-9-236	6NMMKA	6260-9-236	8X4MND	6260-9-236	A0PMME	6260-9-236	E3AMKA
6260-9-236	4H0MMD	6260-9-236	6NMMLD	6260-9-236	8X9MND	6260-9-236	A11MME	6260-9-236	E3AMKAX1
6260-9-236	4KWMJE	6260-9-236	6NWMLD	6260-9-236	8X9MNDX1	6260-9-236	A17MAE	6260-9-236	E3PMAE
6260-9-236	4LYMHE	6260-9-236	6P2MLD	6260-9-236	8XHMND	6260-9-236	A18MME	6260-9-236	E3RMAE

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-236	E79MMD	6260-9-236	JT7MAE	6260-9-236	LNMMLE	6260-9-236	MEPH7K	6260-9-236	MET55N
6260-9-236	E7MMMD	6260-9-236	JT8MAE	6260-9-236	LNRMLE	6260-9-236	MEPJXJ	6260-9-236	MET55R
6260-9-236	E82MMD	6260-9-236	JTEMPA	6260-9-236	LNXML	6260-9-236	MEPK7Y	6260-9-236	MET5J9
6260-9-236	E90MMD	6260-9-236	JTEMPAX1	6260-9-236	LP6MAE	6260-9-236	MEPKLY	6260-9-236	MET7KN
6260-9-236	E90MMDX1	6260-9-236	JTHMPA	6260-9-236	LPDMKD	6260-9-236	MEPNXY	6260-9-236	MET7N0
6260-9-236	EA4MAE	6260-9-236	JTJMPA	6260-9-236	LPHMKD	6260-9-236	MEPNY6	6260-9-236	MET7N3
6260-9-236	EA5MAE	6260-9-236	K41MEE	6260-9-236	LPJMKD	6260-9-236	MEPNYE	6260-9-236	MET7N4
6260-9-236	EA6MAE	6260-9-236	K46MEE	6260-9-236	LPLMLE	6260-9-236	MEPVNH	6260-9-236	MET7X4
6260-9-236	EA6MKE	6260-9-236	K47MEE	6260-9-236	LPNMND	6260-9-236	MEPVP2	6260-9-236	MET803
6260-9-236	EA7MKE	6260-9-236	K5DMJA	6260-9-236	LRTMND	6260-9-236	MEPVPH	6260-9-236	MET804
6260-9-236	ED8MMD	6260-9-236	K66MEE	6260-9-236	LT2MHE	6260-9-236	MEPVXW	6260-9-236	MET806
6260-9-236	ED8MMDX1	6260-9-236	K90MKE	6260-9-236	LV2MHE	6260-9-236	MEPW8T	6260-9-236	MET8LR
6260-9-236	EEKMME	6260-9-236	K91MKE	6260-9-236	LV3MHE	6260-9-236	MEPY08	6260-9-236	MET8MH
6260-9-236	EJ1MME	6260-9-236	K93MJA	6260-9-236	LV4MHE	6260-9-236	MER1DJ	6260-9-236	MET8NP
6260-9-236	EL1MME	6260-9-236	K96MPD	6260-9-236	LV5MHE	6260-9-236	MER1DJX1	6260-9-236	MET8NX
6260-9-236	EV5MKA	6260-9-236	KD6MEE	6260-9-236	LV6MHE	6260-9-236	MER1EM	6260-9-236	MET8PE
6260-9-236	EWPMKA	6260-9-236	KD9MDD	6260-9-236	LW0MDD	6260-9-236	MER1JK	6260-9-236	MET93X
6260-9-236	EX4MKA	6260-9-236	KD9MEE	6260-9-236	LW0MDDX1	6260-9-236	MER1LK	6260-9-236	MET96E
6260-9-236	EX6MKA	6260-9-236	KJ0MJA	6260-9-236	LW1MDD	6260-9-236	MER1VD	6260-9-236	MET9E8
6260-9-236	EX7MKA	6260-9-236	KM5MDD	6260-9-236	LW1MDDX1	6260-9-236	MER237X1	6260-9-236	MET9HK
6260-9-236	EX7MKAX1	6260-9-236	KM6MDD	6260-9-236	LW2MDD	6260-9-236	MER243	6260-9-236	MET9L3
6260-9-236	EX9MKA	6260-9-236	KM7MDD	6260-9-236	LW2MDDX1	6260-9-236	MER243X1	6260-9-236	MET9N0
6260-9-236	EXAMKA	6260-9-236	KM7MDDX1	6260-9-236	LW3MDD	6260-9-236	MER243X2	6260-9-236	MET9VK
6260-9-236	EXHMDA	6260-9-236	KT4MKD	6260-9-236	LW5MKD	6260-9-236	MER268	6260-9-236	METHE6
6260-9-236	EXJMMDA	6260-9-236	KVXML	6260-9-236	LW6MKD	6260-9-236	MER268X1	6260-9-236	METHE7
6260-9-236	EXKMMDA	6260-9-236	KW2MTD	6260-9-236	LW9MNE	6260-9-236	MER26R	6260-9-236	METLK7
6260-9-236	EXLMMDA	6260-9-236	KWJMJA	6260-9-236	LWDMNE	6260-9-236	MER27A	6260-9-236	METLK8
6260-9-236	H0RMLE	6260-9-236	KWMMTD	6260-9-236	LWEMNE	6260-9-236	MER27J	6260-9-236	METLMK
6260-9-236	H22MLE	6260-9-236	KYEMPD	6260-9-236	LWYMKD	6260-9-236	MER29D	6260-9-236	METLYR
6260-9-236	H24MLE	6260-9-236	KYWMPD	6260-9-236	LX7MKD	6260-9-236	MER3DH	6260-9-236	METM4P
6260-9-236	H5JMPD	6260-9-236	L1EMTD	6260-9-236	M00MHE	6260-9-236	MER3DT	6260-9-236	METM5N
6260-9-236	H9LMPD	6260-9-236	L1HMTD	6260-9-236	M01MHE	6260-9-236	MER3KN	6260-9-236	METM7T
6260-9-236	H9MMPD	6260-9-236	L1JMTD	6260-9-236	M02MHE	6260-9-236	MER3P1	6260-9-236	METMET
6260-9-236	H9NMPD	6260-9-236	L1NMJA	6260-9-236	M0NMNE	6260-9-236	MER3V0	6260-9-236	METMK7
6260-9-236	HA2MLD	6260-9-236	L22MND	6260-9-236	M23MHD	6260-9-236	MER43W	6260-9-236	MHAA9D
6260-9-236	HA7MLD	6260-9-236	L2VMLE	6260-9-236	M35MNE	6260-9-236	MER822	6260-9-236	MHAAMP
6260-9-236	HAEMNE	6260-9-236	L4MMKD	6260-9-236	M3JMN	6260-9-236	MER87K	6260-9-236	MHAAN9
6260-9-236	HD1MPD	6260-9-236	L55MND	6260-9-236	M3KMNE	6260-9-236	MER9DT	6260-9-236	MHAAPE
6260-9-236	HD2MPD	6260-9-236	L67MMD	6260-9-236	M5HMDE	6260-9-236	MER9JN	6260-9-236	MHAA6
6260-9-236	HD9MNE	6260-9-236	L67MMDX1	6260-9-236	M61MKE	6260-9-236	MER9R5	6260-9-236	MHAD31
6260-9-236	HDEMPD	6260-9-236	L73MME	6260-9-236	M7EMHD	6260-9-236	MER9V6	6260-9-236	MHADA8
6260-9-236	HEPMLD	6260-9-236	L7KMKE	6260-9-236	MA1MKE	6260-9-236	MERA15	6260-9-236	MHADKJ
6260-9-236	HERMLD	6260-9-236	L9JMAE	6260-9-236	MANMDE	6260-9-236	MERA9A	6260-9-236	MHADKL
6260-9-236	HHEMLD	6260-9-236	L9RMJA	6260-9-236	MEP9N9	6260-9-236	MERAA7	6260-9-236	MHADY5
6260-9-236	HVHMND	6260-9-236	LA3MJA	6260-9-236	MEP9TT	6260-9-236	MERAA7X1	6260-9-236	MHADY6
6260-9-236	J0TMEE	6260-9-236	LAPMAD	6260-9-236	MEPA6A	6260-9-236	MERAA7X3	6260-9-236	MHAE59
6260-9-236	J3EMND	6260-9-236	LAPMADX1	6260-9-236	MEPA6AX1	6260-9-236	MERAA9	6260-9-236	MHAE9E
6260-9-236	J3WMND	6260-9-236	LARMAD	6260-9-236	MEPAJL	6260-9-236	MERD73	6260-9-236	MHALR9
6260-9-236	J52MND	6260-9-236	LATMAD	6260-9-236	MEPAP5	6260-9-236	MERD73X1	6260-9-236	MHALYD
6260-9-236	J53MND	6260-9-236	LAVMAD	6260-9-236	MEPATR	6260-9-236	MERD8E	6260-9-236	MHAM0X
6260-9-236	J57MTD	6260-9-236	LAXMAD	6260-9-236	MEPAWA	6260-9-236	MET12M	6260-9-236	MHAM1J
6260-9-236	JA2MTD	6260-9-236	LH5MHD	6260-9-236	MEPAWE	6260-9-236	MET12P	6260-9-236	MHAM27
6260-9-236	JA2MTDX1	6260-9-236	LJ2MHD	6260-9-236	MEPD8T	6260-9-236	MET12T	6260-9-236	MHAM2D
6260-9-236	JE2MTD	6260-9-236	LJEMLE	6260-9-236	MEPDDT	6260-9-236	MET24Y	6260-9-236	MHAM6M
6260-9-236	JENMTD	6260-9-236	LLJMME	6260-9-236	MEPE35	6260-9-236	MET251	6260-9-236	MHAMH9
6260-9-236	JETMTD	6260-9-236	LLTMAE	6260-9-236	MEPEJD	6260-9-236	MET252	6260-9-236	MHAMH9X1
6260-9-236	JEXMTD	6260-9-236	LLTMLE	6260-9-236	MEPH55	6260-9-236	MET47M	6260-9-236	MHAMVX
6260-9-236	JHKMLE	6260-9-236	LM2MME	6260-9-236	MEPH62	6260-9-236	MET540	6260-9-236	MHAN09
6260-9-236	JHLMLE	6260-9-236	LM4MME	6260-9-236	MEPH6M	6260-9-236	MET541	6260-9-236	MHAN38
6260-9-236	JT6MAE	6260-9-236	LMJMHD	6260-9-236	MEPH75	6260-9-236	MET542	6260-9-236	MHAN3J

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-236	MHAN3X	6260-9-236	MHE724	6260-9-236	MHHNDW	6260-9-236	MHKT0R	6260-9-236	MHLM6X
6260-9-236	MHAN4V	6260-9-236	MHE725	6260-9-236	MHHNJV	6260-9-236	MHKT0T	6260-9-236	MHLM73
6260-9-236	MHAN5N	6260-9-236	MHE7AL	6260-9-236	MHHNM7	6260-9-236	MHKT4H	6260-9-236	MHLM7P
6260-9-236	MHAN72	6260-9-236	MHE7ED	6260-9-236	MHHNVX	6260-9-236	MHKT7R	6260-9-236	MHLM7T
6260-9-236	MHANH9	6260-9-236	MHE7EE	6260-9-236	MHHP07	6260-9-236	MHKT08	6260-9-236	MHLM8E
6260-9-236	MHANLN	6260-9-236	MHEKKM	6260-9-236	MHHP2M	6260-9-236	MHKT0H	6260-9-236	MHLM8H
6260-9-236	MHANYD	6260-9-236	MHEKKMX1	6260-9-236	MHHP2N	6260-9-236	MHKT0N	6260-9-236	MHLM8V
6260-9-236	MHAP0M	6260-9-236	MHEKVV	6260-9-236	MHHP35	6260-9-236	MHKT05	6260-9-236	MHLM8X
6260-9-236	MHAP80	6260-9-236	MHEL1R	6260-9-236	MHHP8W	6260-9-236	MHKT0L	6260-9-236	MHLM8Y
6260-9-236	MHAP81	6260-9-236	MHEL40	6260-9-236	MHHPDT	6260-9-236	MHKT0M	6260-9-236	MHLM9M
6260-9-236	MHAP82	6260-9-236	MHEL40X1	6260-9-236	MHHPHE	6260-9-236	MHKT0V	6260-9-236	MHLM9N
6260-9-236	MHAPDY	6260-9-236	MHEL53	6260-9-236	MHHPNR	6260-9-236	MHKT0W	6260-9-236	MHLM9R
6260-9-236	MHD351	6260-9-236	MHEL8A	6260-9-236	MHHPRH	6260-9-236	MHKT0Y	6260-9-236	MHLM9T
6260-9-236	MHD3ET	6260-9-236	MHELDE	6260-9-236	MHHPV8	6260-9-236	MHKT0Z	6260-9-236	MHLM9V
6260-9-236	MHD3J1	6260-9-236	MHELM	6260-9-236	MHHPW5	6260-9-236	MHKT0A	6260-9-236	MHLM9W
6260-9-236	MHD3J9	6260-9-236	MHELNJ	6260-9-236	MHHPW6	6260-9-236	MHKT0B	6260-9-236	MHLM9X
6260-9-236	MHD3JV	6260-9-236	MHEP3T	6260-9-236	MHHPXP	6260-9-236	MHKT0C	6260-9-236	MHLM9Y
6260-9-236	MHD3K2	6260-9-236	MHEPD3	6260-9-236	MHHPY6	6260-9-236	MHKT0D	6260-9-236	MHLM9Z
6260-9-236	MHD3PL	6260-9-236	MHEPJY	6260-9-236	MHHPY7	6260-9-236	MHKT0E	6260-9-236	MHLM9A
6260-9-236	MHD40M	6260-9-236	MHEPK0	6260-9-236	MHHPYD	6260-9-236	MHKT0F	6260-9-236	MHLM9B
6260-9-236	MHD47R	6260-9-236	MHEPK0X1	6260-9-236	MHHPYM	6260-9-236	MHKT0G	6260-9-236	MHLM9C
6260-9-236	MHD4MD	6260-9-236	MHEPKH	6260-9-236	MHHR5J	6260-9-236	MHKT0H	6260-9-236	MHLM9D
6260-9-236	MHD4MH	6260-9-236	MHEPP4	6260-9-236	MHHRDA	6260-9-236	MHKT0I	6260-9-236	MHLM9E
6260-9-236	MHD4XW	6260-9-236	MHEPTV	6260-9-236	MHJK2X	6260-9-236	MHKT0J	6260-9-236	MHLM9F
6260-9-236	MHDHND	6260-9-236	MHEPVD	6260-9-236	MHJK7H	6260-9-236	MHKT0K	6260-9-236	MHLM9G
6260-9-236	MHDHRV	6260-9-236	MHEPVDX1	6260-9-236	MHJK7W	6260-9-236	MHKT0L	6260-9-236	MHLM9H
6260-9-236	MHDHTJ	6260-9-236	MHER07	6260-9-236	MHJK84	6260-9-236	MHKT0M	6260-9-236	MHLM9I
6260-9-236	MHDHTJX2	6260-9-236	MHERJW	6260-9-236	MHJK8W	6260-9-236	MHKT0N	6260-9-236	MHLM9J
6260-9-236	MHDHV8	6260-9-236	MHERK5	6260-9-236	MHJKJ9	6260-9-236	MHKT0O	6260-9-236	MHLM9K
6260-9-236	MHDHXY	6260-9-236	MHERVV	6260-9-236	MHJKJA	6260-9-236	MHKT0P	6260-9-236	MHLM9L
6260-9-236	MHDHY6	6260-9-236	MHERX8	6260-9-236	MHJL94	6260-9-236	MHKT0Q	6260-9-236	MHLM9M
6260-9-236	MHDJ28	6260-9-236	MHET6H	6260-9-236	MHJL9J	6260-9-236	MHKT0R	6260-9-236	MHLM9N
6260-9-236	MHDJ42	6260-9-236	MHH1X2	6260-9-236	MHJN0L	6260-9-236	MHKT0S	6260-9-236	MHLM9O
6260-9-236	MHDJEH	6260-9-236	MHH1X2X1	6260-9-236	MHJN58	6260-9-236	MHKT0T	6260-9-236	MHLM9P
6260-9-236	MHDJJN	6260-9-236	MHH256	6260-9-236	MHJN6W	6260-9-236	MHKT0U	6260-9-236	MHLM9Q
6260-9-236	MHDJMY	6260-9-236	MHH28R	6260-9-236	MHJN7J	6260-9-236	MHKT0V	6260-9-236	MHLM9R
6260-9-236	MHDK16	6260-9-236	MHH28RX1	6260-9-236	MHJN7V	6260-9-236	MHKT0W	6260-9-236	MHLM9S
6260-9-236	MHDK16X1	6260-9-236	MHH2DJ	6260-9-236	MHJNAR	6260-9-236	MHKT0X	6260-9-236	MHLM9T
6260-9-236	MHDK3A	6260-9-236	MHH2E1	6260-9-236	MHJNEV	6260-9-236	MHKT0Y	6260-9-236	MHLM9U
6260-9-236	MHDK3N	6260-9-236	MHH2WY	6260-9-236	MHJNLP	6260-9-236	MHKT0Z	6260-9-236	MHLM9V
6260-9-236	MHDK4H	6260-9-236	MHH2X0	6260-9-236	MHJNNR	6260-9-236	MHKT0A	6260-9-236	MHLM9W
6260-9-236	MHDK50	6260-9-236	MHH2XP	6260-9-236	MHJNWV	6260-9-236	MHKT0B	6260-9-236	MHLM9X
6260-9-236	MHDK84	6260-9-236	MHH782	6260-9-236	MHJP80	6260-9-236	MHKT0C	6260-9-236	MHLM9Y
6260-9-236	MHDKH3	6260-9-236	MHH7EA	6260-9-236	MHJP8L	6260-9-236	MHKT0D	6260-9-236	MHLM9Z
6260-9-236	MHDKH7	6260-9-236	MHH7MH	6260-9-236	MHJP93	6260-9-236	MHKT0E	6260-9-236	MHLM9A
6260-9-236	MHDKH7X1	6260-9-236	MHH7TA	6260-9-236	MHJP9L	6260-9-236	MHKT0F	6260-9-236	MHLM9B
6260-9-236	MHDKXT	6260-9-236	MHH7XW	6260-9-236	MHJPKT	6260-9-236	MHKT0G	6260-9-236	MHLM9C
6260-9-236	MHDL24	6260-9-236	MHH812	6260-9-236	MHJPM	6260-9-236	MHKT0H	6260-9-236	MHLM9D
6260-9-236	MHDL4X	6260-9-236	MHH813	6260-9-236	MHJPPX	6260-9-236	MHKT0I	6260-9-236	MHLM9E
6260-9-236	MHDL62	6260-9-236	MHH82T	6260-9-236	MHJPTV	6260-9-236	MHKT0J	6260-9-236	MHLM9F
6260-9-236	MHDL6L	6260-9-236	MHH83E	6260-9-236	MHJPYX	6260-9-236	MHKT0K	6260-9-236	MHLM9G
6260-9-236	MHDL6LX1	6260-9-236	MHH8K6	6260-9-236	MHJR1T	6260-9-236	MHKT0L	6260-9-236	MHLM9H
6260-9-236	MHDL6M	6260-9-236	MHH8PW	6260-9-236	MHJR34	6260-9-236	MHKT0M	6260-9-236	MHLM9I
6260-9-236	MHE6LD	6260-9-236	MHH903	6260-9-236	MHJR3L	6260-9-236	MHKT0N	6260-9-236	MHLM9J
6260-9-236	MHE6MD	6260-9-236	MHH904	6260-9-236	MHKRRK	6260-9-236	MHKT0O	6260-9-236	MHLM9K
6260-9-236	MHE6ME	6260-9-236	MHH97T	6260-9-236	MHKRRM	6260-9-236	MHKT0P	6260-9-236	MHLM9L
6260-9-236	MHE6MEX1	6260-9-236	MHH97W	6260-9-236	MHKRW0	6260-9-236	MHKT0Q	6260-9-236	MHLM9M
6260-9-236	MHE6P9	6260-9-236	MHH994	6260-9-236	MHKRW1	6260-9-236	MHKT0R	6260-9-236	MHLM9N
6260-9-236	MHE6PA	6260-9-236	MHHN4M	6260-9-236	MHKRY7	6260-9-236	MHKT0S	6260-9-236	MHLM9O
6260-9-236	MHE722	6260-9-236	MHHN7P	6260-9-236	MHKT0D	6260-9-236	MHKT0T	6260-9-236	MHLM9P

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-236	MHNJ7K	6260-9-236	MHPV93	6260-9-236	MHT576	6260-9-236	MJAV08	6260-9-236	MJEL2V
6260-9-236	MHNJ7T	6260-9-236	MHPV94	6260-9-236	MHT5JW	6260-9-236	MJAV09	6260-9-236	MJEL5Y
6260-9-236	MHNJ7V	6260-9-236	MHPVAA	6260-9-236	MHT5N6	6260-9-236	MJAV0A	6260-9-236	MJEN4D
6260-9-236	MHNJAE	6260-9-236	MHPVAD	6260-9-236	MHT5T5	6260-9-236	MJAV1J	6260-9-236	MJEP9N
6260-9-236	MHNJD1	6260-9-236	MHPVD4	6260-9-236	MHT622	6260-9-236	MJAV1K	6260-9-236	MJEPE4
6260-9-236	MHNJD2	6260-9-236	MHPVEA	6260-9-236	MHT62W	6260-9-236	MJAV1L	6260-9-236	MJEPK2
6260-9-236	MHNJD3	6260-9-236	MHPVHV	6260-9-236	MHT641	6260-9-236	MJAV1M	6260-9-236	MJEPKX
6260-9-236	MHNJJ7	6260-9-236	MHPVM1	6260-9-236	MHT642	6260-9-236	MJAV1N	6260-9-236	MJERLN
6260-9-236	MHNJJ8	6260-9-236	MHPVT4	6260-9-236	MHTH9X	6260-9-236	MJAV1P	6260-9-236	MJERVE
6260-9-236	MHNJJ9	6260-9-236	MHPW99	6260-9-236	MHTHYR	6260-9-236	MJAW0P	6260-9-236	MJH5P3
6260-9-236	MHNKEE	6260-9-236	MHPWJ4	6260-9-236	MHTJ45	6260-9-236	MJD47P	6260-9-236	MJH5P5
6260-9-236	MHNKHE	6260-9-236	MHPWL1	6260-9-236	MHTJ4A	6260-9-236	MJD47R	6260-9-236	MJH8K8
6260-9-236	MHNVRM	6260-9-236	MHPWLL	6260-9-236	MHTJ4D	6260-9-236	MJD47T	6260-9-236	MJH8N4
6260-9-236	MHNVRN	6260-9-236	MHPWLM	6260-9-236	MHTJ6V	6260-9-236	MJD4RT	6260-9-236	MJH99P
6260-9-236	MHNVRP	6260-9-236	MHPWMD	6260-9-236	MHTJEK	6260-9-236	MJD4VL	6260-9-236	MJH99R
6260-9-236	MHNW4R	6260-9-236	MHPWMR	6260-9-236	MHTJHR	6260-9-236	MJD4VM	6260-9-236	MJH99T
6260-9-236	MHNW4W	6260-9-236	MHPWPL	6260-9-236	MHTJMN	6260-9-236	MJD4VN	6260-9-236	MJH9K4
6260-9-236	MHNW4X	6260-9-236	MHPWVD	6260-9-236	MHTPLM	6260-9-236	MJD4WA	6260-9-236	MJH9K5
6260-9-236	MHNW5P	6260-9-236	MHPWXJ	6260-9-236	MHTPPWA	6260-9-236	MJD4WD	6260-9-236	MJH9K6
6260-9-236	MHNWD1	6260-9-236	MHPX1H	6260-9-236	MHTR1H	6260-9-236	MJD4WE	6260-9-236	MJHELT
6260-9-236	MHNWRH	6260-9-236	MHPX40	6260-9-236	MHTR26	6260-9-236	MJD4WH	6260-9-236	MJHEP2
6260-9-236	MHNWYH	6260-9-236	MHRH6E	6260-9-236	MHTR7Y	6260-9-236	MJD4WK	6260-9-236	MJHERL
6260-9-236	MHNX10	6260-9-236	MHRH75	6260-9-236	MHTRJ4	6260-9-236	MJD4WM	6260-9-236	MJHERV
6260-9-236	MHNX13	6260-9-236	MHRH7E	6260-9-236	MHTRJ7	6260-9-236	MJD5XX	6260-9-236	MJHEVD
6260-9-236	MHNX16	6260-9-236	MHRH84	6260-9-236	MHTRKH	6260-9-236	MJD60H	6260-9-236	MJHEVE
6260-9-236	MHNX18	6260-9-236	MHRH92	6260-9-236	MHTRKK	6260-9-236	MJD63R	6260-9-236	MJHEVH
6260-9-236	MHNX4V	6260-9-236	MHRHH3	6260-9-236	MHTRY6	6260-9-236	MJD64R	6260-9-236	MJHEYH
6260-9-236	MHNXLV	6260-9-236	MHRHJW	6260-9-236	MJA017	6260-9-236	MJD668	6260-9-236	MJHH18
6260-9-236	MHNXM0	6260-9-236	MHRHNH	6260-9-236	MJA03T	6260-9-236	MJD66V	6260-9-236	MJHH5R
6260-9-236	MHNXM2	6260-9-236	MHRHPJ	6260-9-236	MJA0AH	6260-9-236	MJD66W	6260-9-236	MJHH91
6260-9-236	MHNXM6	6260-9-236	MHRHXV	6260-9-236	MJA12R	6260-9-236	MJD66X	6260-9-236	MJHH9P
6260-9-236	MHNY22	6260-9-236	MHRJ80	6260-9-236	MJA15W	6260-9-236	MJD670	6260-9-236	MJHJ03
6260-9-236	MHP01Y	6260-9-236	MHRJ8R	6260-9-236	MJA17J	6260-9-236	MJD671	6260-9-236	MJHJ15
6260-9-236	MHP021	6260-9-236	MHRJA4	6260-9-236	MJA1AA	6260-9-236	MJD672	6260-9-236	MJHJ32
6260-9-236	MHP022	6260-9-236	MHRJAK	6260-9-236	MJA1AE	6260-9-236	MJD6DN	6260-9-236	MJHJHN
6260-9-236	MHP023	6260-9-236	MHRJD3	6260-9-236	MJA1DK	6260-9-236	MJD6DP	6260-9-236	MJHJJJ
6260-9-236	MHP024	6260-9-236	MHRJD4	6260-9-236	MJA1JT	6260-9-236	MJD6DR	6260-9-236	MJHRH7
6260-9-236	MHP026	6260-9-236	MHRJKM	6260-9-236	MJA1TL	6260-9-236	MJE6HT	6260-9-236	MJHRH8
6260-9-236	MHP02M	6260-9-236	MHRJPJ	6260-9-236	MJA1W5	6260-9-236	MJE6R0	6260-9-236	MJHRH9
6260-9-236	MHP0YP	6260-9-236	MHRJXA	6260-9-236	MJA237	6260-9-236	MJE700	6260-9-236	MJHRJR
6260-9-236	MHP0YR	6260-9-236	MHRJYJ	6260-9-236	MJA278	6260-9-236	MJE70V	6260-9-236	MJHRL9
6260-9-236	MHP0YT	6260-9-236	MHRK20	6260-9-236	MJA2HY	6260-9-236	MJEA4L	6260-9-236	MJHRLD
6260-9-236	MHP11R	6260-9-236	MHRKK5	6260-9-236	MJA2KY	6260-9-236	MJEAV7	6260-9-236	MJHRLE
6260-9-236	MHP11T	6260-9-236	MHRKKR	6260-9-236	MJA2MP	6260-9-236	MJEAW0	6260-9-236	MJHRLR
6260-9-236	MHP11TX1	6260-9-236	MHT19H	6260-9-236	MJADJ2	6260-9-236	MJED0N	6260-9-236	MJHRMT
6260-9-236	MHP11V	6260-9-236	MHT1HL	6260-9-236	MJADLN	6260-9-236	MJED4Y	6260-9-236	MJHRVM
6260-9-236	MHP11X	6260-9-236	MHT1L7	6260-9-236	MJADMA	6260-9-236	MJED78	6260-9-236	MJHT1X
6260-9-236	MHP11Y	6260-9-236	MHT1R1	6260-9-236	MJADT3	6260-9-236	MJED97	6260-9-236	MJHT2T
6260-9-236	MHP120	6260-9-236	MHT23V	6260-9-236	MJADT4	6260-9-236	MJEDTR	6260-9-236	MJHT6J
6260-9-236	MHPLD8	6260-9-236	MHT23W	6260-9-236	MJADVP	6260-9-236	MJEDVE	6260-9-236	MJHTK0
6260-9-236	MHPLDT	6260-9-236	MHT23Y	6260-9-236	MJAH4J	6260-9-236	MJEDVM	6260-9-236	MJHTK1
6260-9-236	MHPLJ3	6260-9-236	MHT254	6260-9-236	MJAHD0	6260-9-236	MJEDVV	6260-9-236	MJHTK2
6260-9-236	MHPLR7	6260-9-236	MHT25J	6260-9-236	MJAHLN	6260-9-236	MJEDW7	6260-9-236	MJHTTT
6260-9-236	MHPRAK	6260-9-236	MHT4L7	6260-9-236	MJAHTP	6260-9-236	MJEDX5	6260-9-236	MJJ00L
6260-9-236	MHPRTP	6260-9-236	MHT4P6	6260-9-236	MJATDY	6260-9-236	MJEDXN	6260-9-236	MJJ00M
6260-9-236	MHPRVT	6260-9-236	MHT4PK	6260-9-236	MJATE0	6260-9-236	MJEH91	6260-9-236	MJJ00R
6260-9-236	MHPTEP	6260-9-236	MHT4PW	6260-9-236	MJATE1	6260-9-236	MJEH92	6260-9-236	MJJ0T4
6260-9-236	MHPTP2	6260-9-236	MHT4RK	6260-9-236	MJATE2	6260-9-236	MJEKP7	6260-9-236	MJJ0XX
6260-9-236	MHPTYV	6260-9-236	MHT513	6260-9-236	MJATE3	6260-9-236	MJEKP8	6260-9-236	MJJ0XY
6260-9-236	MHPV8P	6260-9-236	MHT542	6260-9-236	MJATE4	6260-9-236	MJEKP9	6260-9-236	MJJ0Y0

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-236	MJJ12H	6260-9-236	MRJMEE	6260-9-236	PWEMJE	6260-9-236	W1EMJA	6260-9-440	MJL91E
6260-9-236	MJJ12L	6260-9-236	MRKMEE	6260-9-236	PWHMJE	6260-9-236	W39MME	6260-9-440	MJL912
6260-9-236	MJJ152	6260-9-236	MRLMEE	6260-9-236	PWJMJE	6260-9-236	W3DMME	6260-9-440	MJL8Y9
6260-9-236	MJJ19M	6260-9-236	MTLMDE	6260-9-236	PY8MKE	6260-9-236	W3EMME	6260-9-440	MJJTD1
6260-9-236	MJJ902	6260-9-236	MTMMDE	6260-9-236	R1KMKE	6260-9-236	W3TMMMD	6260-9-440	MJL8YD
6260-9-236	MJJ91T	6260-9-236	MV5MKE	6260-9-236	R3DMKE	6260-9-236	W6JMEE	6260-9-440	MJL91D
6260-9-236	MJJ95K	6260-9-236	MVAMRD	6260-9-236	RD3MEA	6260-9-236	W6KMEE	6260-9-440	MJKA70
6260-9-236	MJJ984	6260-9-236	MVTMEA	6260-9-236	RD5MEA	6260-9-236	W6LMEE	6260-9-440	MJL91H
6260-9-236	MJJ9DP	6260-9-236	MVVMEA	6260-9-236	RD6MEA	6260-9-236	W7PMJA	6260-9-440	MJL91J
6260-9-236	MJJ9W6	6260-9-236	MWXMED	6260-9-236	RH2MLE	6260-9-236	W8EMJA	6260-9-440	MJM381
6260-9-236	MJJ9XW	6260-9-236	NO1MRD	6260-9-236	RJMMNA	6260-9-236	W8NMJE	6260-9-440	MJM67X
6260-9-236	MJJA4V	6260-9-236	NOYMKD	6260-9-236	RJMMNAX1	6260-9-236	W9DMJE	6260-9-440	MJMEA5
6260-9-236	MJJA7D	6260-9-236	N3HMNE	6260-9-236	RJRMNA	6260-9-236	W9EMJE	6260-9-240	MJD835
6260-9-236	MJJA9T	6260-9-236	N3YMAA	6260-9-236	RJTMNA	6260-9-236	W9HMJE	6260-9-240	MJD852
6260-9-236	MJJAAN	6260-9-236	N4EMKD	6260-9-236	RJTMNAX1	6260-9-236	WH2MHA	6260-9-240	MJD8AV
6260-9-236	MJJAYN	6260-9-236	N4KMPD	6260-9-236	RP5MLE	6260-9-236	WHMMHA	6260-9-240	MJD90E
6260-9-236	MJJAYP	6260-9-236	N4LMPD	6260-9-236	RV0MLE	6260-9-236	WHXMHA	6260-9-240	MJD90H
6260-9-236	MJJAYR	6260-9-236	N4MMPD	6260-9-236	RXMMEA	6260-9-236	WV6MJE	6260-9-240	MJD90J
6260-9-236	MJJDA4	6260-9-236	N55MKE	6260-9-236	RXXMEA	6260-9-236	WVNMEE	6260-9-240	MJD99H
6260-9-236	MJJMWY	6260-9-236	NAJMNE	6260-9-236	RXYMHD	6260-9-236	WVPMEE	6260-9-340	MJL2JX
6260-9-236	MJJMX0	6260-9-236	NAPMNE	6260-9-236	T12MEA	6260-9-236	WVRMEE	6260-9-240	MJL2NJ
6260-9-236	MJJMX1	6260-9-236	NAXMNE	6260-9-236	T15MEA	6260-9-236	WXVMMMA	6260-9-240	MJL39H
6260-9-236	MJJN7Y	6260-9-236	NEAMMA	6260-9-236	T3AMHD	6260-9-236	WXVMMAX1	6260-9-240	MJL3AL
6260-9-236	MJJN82	6260-9-236	NEDMMA	6260-9-236	T68MDE	6260-9-236	X6YMJE	6260-9-240	MJL3V4
6260-9-236	MJJN83	6260-9-236	NH3MKE	6260-9-236	TEKMDE	6260-9-236	XAKMJE	6260-9-340	MJL468
6260-9-236	MJJNVA	6260-9-236	NHMMMA	6260-9-236	THRMHE	6260-9-236	XD0MJE	6260-9-240	MJL4LL
6260-9-236	MJJNVD	6260-9-236	NMKMEE	6260-9-236	THTMHE	6260-9-236	XDKMJE	6260-9-240	MJL6A4
6260-9-236	MJJNVH	6260-9-236	NN3MEA	6260-9-236	THVMHEX1	6260-9-236	XE9MJE	6260-9-340	MJL6HL
6260-9-236	MJKDD0	6260-9-236	NN7MEA	6260-9-236	TK3MHD	6260-9-236	XELMJE	6260-9-340	MJL6VJ
6260-9-236	MJKDV5	6260-9-236	NNJMJE	6260-9-236	TNMMHD	6260-9-236	XV2MHE	6260-9-240	MJL7M4
6260-9-236	MJKDV8	6260-9-236	NVJMKD	6260-9-236	TPXMHD	6260-9-236	Y2NMHE	6260-9-444	MJL7NV
6260-9-236	MJKE2L	6260-9-236	NVJMKDX1	6260-9-236	TT0MMA	6260-9-236	Y9EMHE	6260-9-444	MJL7NW
6260-9-236	MJKE2M	6260-9-236	NXNMEE	6260-9-236	TTDMMA	6260-9-236	YJDMHE	6260-9-444	MJL7NX
6260-9-236	MJKE2N	6260-9-236	NXPMEE	6260-9-236	TTEMMA	6260-9-236	YJEMHE	6260-9-444	MJL7NY
6260-9-236	MJKEEV	6260-9-236	NXRMEE	6260-9-236	TTRMHD	6260-9-236	YJWMJE	6260-9-240	MJL7V5
6260-9-236	MJKEEW	6260-9-236	P3RMEA	6260-9-236	TTVMHD	6260-9-236	YRDMJE	6260-9-240	MJL7XH
6260-9-236	MJKEEY	6260-9-236	P3TMEA	6260-9-236	TVXMMMA	6260-9-440	MHKLTH	6260-9-444	MJL8KW
6260-9-236	MJKH4W	6260-9-236	P8NMLE	6260-9-236	TVYMMMA	6260-9-440	MHN2A0	6260-9-444	MJL8L0
6260-9-236	MJKHM4	6260-9-236	PE1MED	6260-9-236	V0WMHE	6260-9-440	MHKJAM	6260-9-444	MJL8L2
6260-9-236	MJKHM6	6260-9-236	PE3MED	6260-9-236	V0XMHE	6260-9-440	MHNY53	6260-9-444	MJL8L4
6260-9-236	MJKHPX	6260-9-236	PE4MED	6260-9-236	V0YMHE	6260-9-440	MHN2X0	6260-9-240	MJL8MK
6260-9-236	MJKRNN	6260-9-236	PE5MED	6260-9-236	V0YMHEX1	6260-9-440	MHN2WX	6260-9-444	MJL8Y6
6260-9-236	MJKRT7	6260-9-236	PHMMED	6260-9-236	V1DMDE	6260-9-440	MHN2A1	6260-9-444	MJL90A
6260-9-236	MJKT2V	6260-9-236	PHNMED	6260-9-236	V2MMDE	6260-9-440	MHN2A2	6260-9-444	MJL90D
6260-9-236	MJKTPT	6260-9-236	PHPMED	6260-9-236	V2NMDE	6260-9-440	MHTK9K	6260-9-444	MJL90E
6260-9-236	MJKV11	6260-9-236	PHRMED	6260-9-236	V4YMHA	6260-9-440	MHTJP0	6260-9-444	MJL911
6260-9-236	MJKVMV	6260-9-236	PLLME	6260-9-236	VATMEA	6260-9-440	MHTKA1	6260-9-444	MJL916
6260-9-236	MJKVWJ	6260-9-236	PLRMME	6260-9-236	VAVMEA	6260-9-440	MHTK9L	6260-9-444	MJL917
6260-9-236	MJKW33	6260-9-236	PM2MME	6260-9-236	VAVMEA	6260-9-244	MHTLYK	6260-9-444	MJL918
6260-9-236	MKXMPD	6260-9-236	PM6MAE	6260-9-236	VAXMEA	6260-9-244	MHTM2X	6260-9-444	MJL919
6260-9-236	ML7MNE	6260-9-236	PM7MAE	6260-9-236	VAYMEA	6260-9-244	MHTM42	6260-9-240	MJL98N
6260-9-236	MLLMNE	6260-9-236	PM8MAE	6260-9-236	VE1MDE	6260-9-240	MJA73E	6260-9-344	MJL9H2
6260-9-236	MLNMNE	6260-9-236	PNXMKE	6260-9-236	VE1MDEX1	6260-9-240	MJA73J	6260-9-344	MJL9H3
6260-9-236	MM9MKE	6260-9-236	PP0MKE	6260-9-236	VE3MJA	6260-9-240	MJA73L	6260-9-344	MJL9H4
6260-9-236	MN0MKE	6260-9-236	PP1MKE	6260-9-236	VJ3MJA	6260-9-240	MJA73M	6260-9-344	MJL9H5
6260-9-236	MNJMDE	6260-9-236	PRHMMD	6260-9-236	VJTMJA	6260-9-340	MJKDD1	6260-9-340	MJL9LD
6260-9-236	MNKMDE	6260-9-236	PRHMMDX1	6260-9-236	VJTMJAX1	6260-9-240	MJKDLH	6260-9-240	MJL9MR
6260-9-236	MPEMPD	6260-9-236	PRJMMD	6260-9-236	VM1MJA	6260-9-240	MJKDV1	6260-9-344	MJL9T3
6260-9-236	MPVMPD	6260-9-236	PRKMMD	6260-9-236	VV9MMA	6260-9-240	MJKDV2	6260-9-344	MJL9T5
6260-9-236	MR1MKE	6260-9-236	PWDMJE	6260-9-236	VVJMMA	6260-9-240	MJKDV3	6260-9-344	MJL9T8

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-340	MHNW55	6260-9-340	MHPT7D	6260-9-344	MHRMVL	6260-9-240	MJKDV4	6260-9-344	MJL9T9
6260-9-340	MHNWAY	6260-9-340	MHPT7T	6260-9-344	MHRMVY	6260-9-240	MJKE0J	6260-9-340	MJLAWVW
6260-9-340	MHNWK0	6260-9-340	MHPT94	6260-9-344	MHRMXH	6260-9-240	MJKE0K	6260-9-340	MJLDAE
6260-9-340	MHNWYJ	6260-9-340	MHPT95	6260-9-340	MHRMY7	6260-9-240	MJKE0L	6260-9-340	MJLDAH
6260-9-340	MHNX1A	6260-9-340	MHPT96	6260-9-340	MHRMY8	6260-9-240	MJKE0M	6260-9-340	MJLDAJ
6260-9-340	MHNX1D	6260-9-240	MHR0VV	6260-9-340	MHRMY9	6260-9-240	MJKE2P	6260-9-340	MJLDAK
6260-9-340	MHNX1E	6260-9-240	MHR137	6260-9-340	MHRMYA	6260-9-240	MJKE2R	6260-9-340	MJLDRR
6260-9-240	MJD9LK	6260-9-240	MJEMYT	6260-9-244	MJJHA3	6260-9-240	MJKE2T	6260-9-340	MJLDT9
6260-9-240	MJDA9V	6260-9-240	MJEMYW	6260-9-240	MJJHDL	6260-9-240	MJKE2V	6260-9-240	MJLDVJ
6260-9-240	MJDAA1	6260-9-240	MJER2Y	6260-9-240	MJJHHH	6260-9-340	MJKLKL	6260-9-340	MJLDW2
6260-9-240	MJDAA2	6260-9-240	MJERA6	6260-9-240	MJJHJW	6260-9-340	MJKLKN	6260-9-340	MJLE6J
6260-9-240	MJDAA3	6260-9-240	MJERA8	6260-9-240	MJJJJA	6260-9-340	MJKLKY	6260-9-340	MJLE6K
6260-9-240	MJDAHx	6260-9-240	MJERA9	6260-9-240	MJJJJP	6260-9-340	MJKLL4	6260-9-240	Y67MNE
6260-9-240	MJDAHx	6260-9-240	MJERAA	6260-9-240	MJJJR2	6260-9-240	MJKN45	6260-9-240	Y7HMHE
6260-9-240	MJDAJ0	6260-9-240	MJERJL	6260-9-244	MJJJV1	6260-9-240	MJKN47	6260-9-240	Y80MNE
6260-9-240	MJDAJ1	6260-9-240	MJERJM	6260-9-244	MJJJV2	6260-9-240	MJKN4J	6260-9-240	Y8AMNE
6260-9-240	MJDAMN	6260-9-240	MJERJN	6260-9-244	MJJJV3	6260-9-240	MJKN4K	6260-9-240	Y8EMNE
6260-9-340	MJDJP3	6260-9-240	MJERJP	6260-9-240	MJJJVW	6260-9-240	MJKN98	6260-9-240	YJ6MHE
6260-9-340	MJDJPR	6260-9-240	MJHE8V	6260-9-240	MJJK15	6260-9-240	MJKN99	6260-9-340	YJJMKD
6260-9-340	MJDJR6	6260-9-240	MJHEDV	6260-9-240	MJJK1N	6260-9-240	MJKN9A	6260-9-240	YK5MHE
6260-9-340	MJDJWA	6260-9-240	MJHEDX	6260-9-340	MJJTVK	6260-9-240	MJKN9D	6260-9-240	YK6MHE
6260-9-340	MJDK0J	6260-9-240	MJHEE0	6260-9-340	MJJV2E	6260-9-240	MJKND5	6260-9-240	YK7MHE
6260-9-340	MJDK7A	6260-9-240	MJHEE3	6260-9-340	MJJV2H	6260-9-240	MJKND6	6260-9-440	X2AMLE
6260-9-340	MJDK7V	6260-9-240	MJHETW	6260-9-340	MJJV2J	6260-9-240	MJKND7	6260-9-440	X9XMLE
6260-9-340	MJDKDP	6260-9-240	MJHETX	6260-9-340	MJJV2K	6260-9-240	MJKND8	6260-9-440	X29MLE
6260-9-340	MJDKEA	6260-9-240	MJHETY	6260-9-340	MJJV7W	6260-9-240	MJKNK3	6260-9-440	JVLMME
6260-9-340	MJDL00	6260-9-240	MJHEV1	6260-9-340	MJJV7X	6260-9-240	MJKNK4	6260-9-440	KKRMLE
6260-9-340	MJDL0D	6260-9-340	MJHEVX	6260-9-340	MJJV80	6260-9-240	MJKNK5	6260-9-440	JW4MME
6260-9-340	MJDL3A	6260-9-340	MJHEVY	6260-9-340	MJJV82	6260-9-240	MJKNK6	6260-9-440	418MLE
6260-9-340	MJDL4D	6260-9-340	MJHEYV	6260-9-340	MJJVA4	6260-9-240	MJKNKE	6260-9-440	JVHMME
6260-9-340	MJDL80	6260-9-240	MJHKK8	6260-9-340	MJJVA5	6260-9-240	MJKNMK	6260-9-440	JW1MME
6260-9-340	MJDLEA	6260-9-240	MJHKNW	6260-9-340	MJJVA6	6260-9-240	MJKNVX	6260-9-440	450MLE
6260-9-340	MJDLL9	6260-9-240	MJHKPT	6260-9-340	MJJVE0	6260-9-240	MJKNYA	6260-9-440	38MMNE
6260-9-340	MJDLP3	6260-9-240	MJHKPY	6260-9-340	MJJVE1	6260-9-240	MJKP09	6260-9-440	39EMNE
6260-9-240	MJE6HP	6260-9-240	MJHKYK	6260-9-340	MJJVE2	6260-9-240	MJKP1Y	6260-9-440	391MNE
6260-9-240	MJE6X6	6260-9-240	MJHKYL	6260-9-340	MJJVE3	6260-9-240	MJKP7D	6260-9-440	453MLEX1
6260-9-240	MJE6YY	6260-9-240	MJHKYN	6260-9-240	MJK2E1	6260-9-340	MJKPM8	6260-9-440	4JYMNE
6260-9-240	MJE70D	6260-9-240	MJHKYP	6260-9-240	MJK2E2	6260-9-340	MJKPN4	6260-9-440	JWLMME
6260-9-240	MJE70W	6260-9-340	MJHL26	6260-9-240	MJK2E3	6260-9-340	MJKPP4	6260-9-440	MEPMW8
6260-9-240	MJE7AP	6260-9-340	MJHL27	6260-9-240	MJK2E4	6260-9-340	MJKPT1	6260-9-440	MER5MN
6260-9-240	MJE7AR	6260-9-340	MJHL4P	6260-9-240	MJK3W4	6260-9-340	MJKPYE	6260-9-440	MEPMWA
6260-9-240	MJE7AT	6260-9-340	MJHL5X	6260-9-240	MJK40L	6260-9-340	MJKR0W	6260-9-440	MER5R1
6260-9-240	MJE7AV	6260-9-340	MJHL86	6260-9-240	MJK40M	6260-9-244	MJKRKR	6260-9-440	MER5NH
6260-9-240	MJE7JE	6260-9-340	MJHLHT	6260-9-240	MJK40N	6260-9-244	MJKRKT	6260-9-440	MEPMWAX1
6260-9-240	MJE7JH	6260-9-340	MJHLJW	6260-9-240	MJK40P	6260-9-244	MJKRKV	6260-9-440	MER5R2
6260-9-240	MJE7JJ	6260-9-340	MJHLKT	6260-9-240	MJK5XN	6260-9-244	MJKRKW	6260-9-440	MEPMW9
6260-9-240	MJE7JK	6260-9-240	MJHLM9	6260-9-240	MJK61R	6260-9-244	MJKRPP	6260-9-440	MEPMW8X1
6260-9-244	MJEHT4	6260-9-240	MJHLP7	6260-9-240	MJK63T	6260-9-244	MJKRPR	6260-9-440	MET8NW
6260-9-244	MJEHWH	6260-9-240	MJHLPP	6260-9-240	MJK66N	6260-9-244	W3VMMA	6260-9-440	X29MLEX1
6260-9-244	MJEJ4E	6260-9-240	MJHM28	6260-9-240	MJK6YK	6260-9-244	W4MMA	6260-9-440	MET909
6260-9-244	MJEJKM	6260-9-240	MJHM2V	6260-9-240	MJK70A	6260-9-244	W4DMED	6260-9-440	MET908
6260-9-244	MJEJLJ	6260-9-240	MJHM32	6260-9-240	MJK79E	6260-9-244	W4JMED	6260-9-440	MET90D
6260-9-244	MJEJMK	6260-9-240	MJHM58	6260-9-240	MJK79J	6260-9-240	W8EMLE	6260-9-440	MET8N9
6260-9-244	MJEJPH	6260-9-240	MJHM98	6260-9-240	MJK79K	6260-9-240	W8RMLE	6260-9-440	MET8N9X1
6260-9-244	MJEJR6	6260-9-244	MJJE8D	6260-9-240	MJK7TV	6260-9-240	W92MLE	6260-9-440	MEPMW7
6260-9-244	MJEJX4	6260-9-244	MJJE8E	6260-9-240	MJK7YN	6260-9-240	W98MJE	6260-9-440	MHA87E
6260-9-244	MJEK0M	6260-9-244	MJJE8H	6260-9-240	MJK7YR	6260-9-240	W99MLE	6260-9-440	MHAHRA
6260-9-244	MJEK2P	6260-9-244	MJJE8J	6260-9-240	MJK7YV	6260-9-240	WAXMJE	6260-9-440	MHA87K
6260-9-244	MJEKPV	6260-9-240	MJJEL8	6260-9-340	MJKAXJ	6260-9-440	H2NMND	6260-9-440	MHAHE2
6260-9-244	MJEKR3	6260-9-240	MJJER3	6260-9-240	MJKAXX	6260-9-440	H2JMND	6260-9-440	MHA87L

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-340	MY6MED	6260-9-244	P5LMLA	6260-9-340	RPNMMAX1	6260-9-440	MHJMPD	6260-9-440	MHAHL2
6260-9-340	MY6MEDX1	6260-9-244	P6WMMMD	6260-9-240	RRNMNA	6260-9-440	MENMPD	6260-9-440	MHAHL1
6260-9-244	N02MRD	6260-9-340	P7TMDD	6260-9-240	RRNMNAX1	6260-9-440	MEPMPD	6260-9-440	MHAHK4
6260-9-240	N19MDD	6260-9-340	P7VMDD	6260-9-240	RRRMNA	6260-9-440	MHKMPD	6260-9-440	MHA87N
6260-9-240	N1AMDD	6260-9-340	P7VMDDX1	6260-9-240	RRRMNAX1	6260-9-440	JWNMRD	6260-9-440	MHA87KX1
6260-9-240	N1DMDD	6260-9-340	P7WMDD	6260-9-244	RT3MME	6260-9-440	N56MPD	6260-9-440	MHA87NX1
6260-9-240	N1EMDD	6260-9-340	P7XMDD	6260-9-244	RT9MME	6260-9-440	JRTMRD	6260-9-440	MHA87EX1
6260-9-240	N1EMDDX1	6260-9-444	P8RMAD	6260-9-244	RTAMME	6260-9-440	JWXMRD	6260-9-440	MHDMMN
6260-9-240	N1HMDD	6260-9-244	P8RMPA	6260-9-244	RTDMME	6260-9-440	JTAMRD	6260-9-440	MHAHK8
6260-9-240	N1JMDD	6260-9-244	P8RMPAX1	6260-9-240	RX0MDE	6260-9-440	JTMMRD	6260-9-440	MHE149
6260-9-240	YRAMJE	6260-9-440	L3XMLAX1	6260-9-440	PVWMPAX1	6260-9-440	JVWMRD	6260-9-440	MHDT8X
6260-9-440	L6VMLA	6260-9-440	N6KMLA	6260-9-440	M2MMPAX1	6260-9-440	JVWMRDX1	6260-9-440	MHE15L
6260-9-440	EPDMLA	6260-9-440	MLRMLAX1	6260-9-440	PVXMPA	6260-9-440	94TMMDX1	6260-9-440	MHDMMM
6260-9-440	EPAMLA	6260-9-440	M3HMLAX1	6260-9-440	N4DMPAX1	6260-9-440	58AMAE	6260-9-440	MHDMMML
6260-9-440	L6TMLA	6260-9-440	N6KMLAX1	6260-9-440	LT4MDD	6260-9-440	74MMAE	6260-9-440	MHE16D
6260-9-440	DVWMLA	6260-9-440	L6DMLAX1	6260-9-440	JVHMTA	6260-9-440	667MAE	6260-9-440	MHERVE
6260-9-440	L3XMLA	6260-9-440	N71MLAX1	6260-9-440	RALMDD	6260-9-440	66EMAE	6260-9-440	MHDMMMP
6260-9-440	KAKMLA	6260-9-440	NV4MMAX1	6260-9-440	JVAMTA	6260-9-440	VH0MAE	6260-9-440	MHDTDK
6260-9-440	L4WMLA	6260-9-440	NV4MMA	6260-9-440	A1RMTA	6260-9-440	VEYMAE	6260-9-440	MHDMMMPX1
6260-9-440	DVVMLA	6260-9-440	NV5MMA	6260-9-440	4DLMED	6260-9-440	VEXMAE	6260-9-440	MHE120
6260-9-440	K0TMLA	6260-9-440	NV6MMA	6260-9-440	4DKMED	6260-9-440	VEWMAE	6260-9-440	MHDTDL
6260-9-440	L6LMLA	6260-9-440	NTVMMA	6260-9-440	JVDMTA	6260-9-440	9R3MJDX1	6260-9-440	MHDTDM
6260-9-440	RVLMMMA	6260-9-440	RRKMNA	6260-9-440	JR9MTA	6260-9-440	LM3MDE	6260-9-440	MHERJVX1
6260-9-440	LWKMLA	6260-9-440	NV4MMAX2	6260-9-440	M82MTA	6260-9-440	LM6MDE	6260-9-440	MJMMT9
6260-9-440	LWLMLA	6260-9-440	6KKMPA	6260-9-440	M86MTA	6260-9-440	M7RMDE	6260-9-440	MJNW37
6260-9-440	L6WMLA	6260-9-440	RVJMMA	6260-9-440	M83MTA	6260-9-440	M7NMDE	6260-9-440	MJNH5T
6260-9-440	NV3MMA	6260-9-440	RRMMNA	6260-9-440	LVVMDD	6260-9-440	AL2MEE	6260-9-440	MJNX81
6260-9-440	L65MLA	6260-9-440	6KLMPA	6260-9-440	M84MTA	6260-9-440	ARNMEE	6260-9-440	MJN3T9
6260-9-440	L4VMLA	6260-9-440	3D2MPA	6260-9-440	M85MTA	6260-9-440	D11MEE	6260-9-440	MJMWVPV
6260-9-440	L4XMLA	6260-9-440	6KMMPA	6260-9-440	JVHMTAX1	6260-9-440	AEWMEE	6260-9-440	MJP111
6260-9-440	KAHMLA	6260-9-440	DRYMPA	6260-9-440	4DMMED	6260-9-440	M3DMEE	6260-9-440	MJPY8A
6260-9-440	L3TMLA	6260-9-440	A35MPA	6260-9-440	RAMMDD	6260-9-440	M3JMEE	6260-9-440	MJPVJ6
6260-9-440	N6YMLA	6260-9-440	EV3MPA	6260-9-440	LANMHD	6260-9-440	MH9MEE	6260-9-440	MJPD LX
6260-9-440	N6LMLA	6260-9-440	A30MPA	6260-9-440	LD1MHD	6260-9-440	M28MEE	6260-9-440	MJPHYOY
6260-9-440	KAJMLA	6260-9-440	A31MPA	6260-9-440	LD5MHD	6260-9-440	4EWMHE	6260-9-440	MJMN34
6260-9-440	L6DMLA	6260-9-440	A33MPA	6260-9-440	LATMHD	6260-9-440	4ETMHE	6260-9-440	MJM71Y
6260-9-440	L6RMLAX1	6260-9-440	E2XMPA	6260-9-440	9HJMJ	6260-9-440	4ENMHE	6260-9-440	MJM4L5
6260-9-440	L65MLAX1	6260-9-440	6KLMPAX1	6260-9-440	9H4MJDX1	6260-9-440	LNLMMHE	6260-9-440	MJP3RT
6260-9-440	LWJMLAX2	6260-9-440	DTOMPA	6260-9-440	9EXMJ	6260-9-440	LP1MHE	6260-9-440	MJNK85
6260-9-440	KAHMLAX1	6260-9-440	DRWMPA	6260-9-440	9R1MJDX1	6260-9-440	4EVMHE	6260-9-440	MJMOM8
6260-9-440	EP9MLA	6260-9-440	N4DMPA	6260-9-440	9HAMJD	6260-9-440	LP2MHE	6260-9-440	MJME2E
6260-9-440	M3NMLA	6260-9-440	460MRA	6260-9-440	9H8MJD	6260-9-440	J91MJE	6260-9-440	MJMA4N
6260-9-440	L3WMLA	6260-9-440	45YMRA	6260-9-440	9R1MJD	6260-9-440	J90MJE	6260-9-440	MJMTON
6260-9-440	NW7MMA	6260-9-440	M2DMPA	6260-9-440	98XMKD	6260-9-440	J8YMJE	6260-9-440	MJMKLD
6260-9-440	NTTMMA	6260-9-440	8JHMRA	6260-9-440	9H4MJD	6260-9-440	LP3MHE	6260-9-440	MJRE4V
6260-9-440	M4TMLA	6260-9-440	3D2MPAX1	6260-9-440	9R2MJD	6260-9-440	233MKE	6260-9-440	MJRV4M
6260-9-440	L3RMLA	6260-9-440	M2MMPA	6260-9-440	9AXMKD	6260-9-440	1VYMKE	6260-9-440	MJMNK2
6260-9-440	HW2MLA	6260-9-440	6KMMPAX1	6260-9-440	9NVMKD	6260-9-440	1HJMKE	6260-9-440	MJM56E
6260-9-440	N70MLA	6260-9-440	8JJMRA	6260-9-440	8T5MKD	6260-9-440	J81MJE	6260-9-440	MJTE81
6260-9-440	L3VMLA	6260-9-440	PVWMPA	6260-9-440	9NRMKD	6260-9-440	JHTMJ	6260-9-440	MJTM3E
6260-9-440	NW9MMA	6260-9-440	PVYMPA	6260-9-440	8V8MKD	6260-9-440	1VVMKE	6260-9-440	MJTHDD
6260-9-440	N6JMLA	6260-9-440	M40MPA	6260-9-440	9NTMKD	6260-9-440	454MLE	6260-9-440	MJMDDEV
6260-9-440	L6RMLA	6260-9-440	A33MPAX1	6260-9-440	9R3MJD	6260-9-440	453MLE	6260-9-440	MJMED1
6260-9-440	LWJMLA	6260-9-440	A30MPAX1	6260-9-440	9M6MKD	6260-9-440	1VTMKE	6260-9-440	MJT85Y
6260-9-440	N71MLA	6260-9-440	PVXMPAX1	6260-9-440	9JTMKD	6260-9-440	MNVMLE	6260-9-440	MJTY6R
6260-9-440	MLRMLA	6260-9-440	HMAMRA	6260-9-440	9NVMKDX1	6260-9-440	MT9MLE	6260-9-440	MJNJT8
6260-9-440	MN1MLA	6260-9-440	M40MPAX1X1	6260-9-440	9NRMKDX1	6260-9-440	MT8MLE	6260-9-440	MJNA6X
6260-9-440	M3VMLA	6260-9-440	M3XMPA	6260-9-440	93RMKD	6260-9-440	X1WMLE	6260-9-440	MJNVP6
6260-9-440	L2RMLA	6260-9-440	45YMRAX1	6260-9-440	8Y8MMD	6260-9-440	X23MLE	6260-9-440	MJNK3R
6260-9-440	MMAMLA	6260-9-440	EV3MPAX1	6260-9-440	94TMMD	6260-9-440	MJPK9D	6260-9-440	MJMJRM

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-440	M3HMLA	6260-9-440	MJRMPA	6260-9-440	9H2MMD	6260-9-440	MHKJER	6260-9-440	MHRKMP
6260-9-440	N6HMLA	6260-9-440	RRMMNAX1	6260-9-440	9JTMKDX1	6260-9-440	MHKJA4	6260-9-440	MHRKMX
6260-9-440	MLXMLA	6260-9-440	9X8MTA	6260-9-440	9KVMMMD	6260-9-440	MHKK9D	6260-9-440	MHRKMR
6260-9-440	M4RMLA	6260-9-440	83NMTA	6260-9-440	8Y7MMD	6260-9-440	MHKJA2	6260-9-440	MHTJJN
6260-9-440	MLWMLA	6260-9-440	89TMTA	6260-9-440	9P7MMD	6260-9-440	MHKLPK	6260-9-440	MHTJMR
6260-9-440	M3TMLA	6260-9-440	M40MPAX1	6260-9-440	97PMMD	6260-9-440	MHKKEY	6260-9-440	MJH81K
6260-9-440	M3PMLA	6260-9-440	HVNMRA	6260-9-440	93JMMD	6260-9-440	MHKKK3	6260-9-440	MHRKMPV
6260-9-440	M3RMLA	6260-9-440	032MRA	6260-9-440	H67MND	6260-9-440	MHKJ8H	6260-9-440	MJHR1J
6260-9-440	M3TMLAX1	6260-9-440	N44MPA	6260-9-440	H2EMND	6260-9-440	MHKJEP	6260-9-440	MHTK92
6260-9-440	M4VMLA	6260-9-440	02YMRA	6260-9-440	H6TMND	6260-9-440	MHKJET	6260-9-440	MJHPRH
6260-9-440	M4TMLAX1	6260-9-440	45XMRA	6260-9-440	H2KMND	6260-9-440	MHKKM2	6260-9-440	MJH81L
6260-9-440	MHENKJX1	6260-9-440	MHN2WY	6260-9-440	MJM7AN	6260-9-440	MHKJAL	6260-9-440	MJJKNX
6260-9-440	MHERJV	6260-9-440	MHNV8E	6260-9-440	MJMPAW	6260-9-440	MHKK4V	6260-9-440	MJKLXD
6260-9-440	MHJ566	6260-9-440	MHNV8H	6260-9-440	MJML4N	6260-9-440	MHKL9X1	6260-9-440	MJKA4K
6260-9-440	MHJ57L	6260-9-440	MHN2X1	6260-9-440	MJM296	6260-9-440	MHKK0T	6260-9-440	MJHR41
6260-9-440	MHERVEX1	6260-9-440	MHNV8D	6260-9-440	MJM3PA	6260-9-440	MHKKMD	6260-9-440	MJJTRV
6260-9-440	MHDT78	6260-9-440	MHKL7	6260-9-440	MJMLE1	6260-9-440	MHKKJP	6260-9-440	MJKLXL
6260-9-440	MHDT7H	6260-9-440	MHNY50	6260-9-440	MJM1AL	6260-9-440	MHKJ9Y	6260-9-440	MJHPWV
6260-9-440	MHERK4	6260-9-440	MHNV8J	6260-9-440	MJMHKE	6260-9-440	MHKJAJ	6260-9-440	MJH81N
6260-9-440	MHE15LX1	6260-9-440	MJE7MW	6260-9-440	MJMJTT	6260-9-440	MHKLTV	6260-9-440	MJH81H
6260-9-440	MHENKJ	6260-9-440	MHNY54	6260-9-440	MJM5MK	6260-9-440	MHKJA3	6260-9-440	MJKLXH
6260-9-440	MHJ4YT	6260-9-440	MHTKA3	6260-9-440	MJMXTD	6260-9-440	MHKKJ8	6260-9-440	MJK9T6
6260-9-440	MHJ4YR	6260-9-440	MHNY52	6260-9-440	MJM127	6260-9-440	MHKK8J	6260-9-440	MJK9XT
6260-9-440	MHKKHE	6260-9-440	MHTK93	6260-9-440	MJNW6H	6260-9-440	MHKKK7	6260-9-440	MJKA0D
6260-9-440	MHKJXK	6260-9-440	MJJKTW	6260-9-440	MJJTDK	6260-9-440	MHKL6P	6260-9-440	MJKLY6
6260-9-440	MHKJEV	6260-9-440	MJJTJW	6260-9-440	MJK94T	6260-9-440	MHKKJ9	6260-9-440	MJKLXJ
6260-9-440	MHN2A3	6260-9-440	MJJTTK	6260-9-440	MJKA62	6260-9-440	MHKJ51	6260-9-440	MHKKN8
6260-9-440	MHKKRV	6260-9-440	MJL8Y8	6260-9-440	MJKA3R	6260-9-440	MHKLTV	6260-9-440	MHKKMW
6260-9-440	MHKL9	6260-9-440	MJL8YA	6260-9-440	MJK9VD	6260-9-440	MHKKMV		
6260-9-236	MHMR0P	6260-9-236	92546803	6260-9-236	15186902	6260-9-244	MJP9KV	6260-9-240	MJT6M1
6260-9-236	MHMT46X1	6260-9-236	93418702	6260-9-236	14946701	6260-9-440	MJPD65	6260-9-440	MJTR4V
6260-9-236	MHNJA5	6260-9-236	93418701	6260-9-236	86376102A	6260-9-340	MJPE16	6260-9-240	MJTWXA
6260-9-236	MHRKJ3	6260-9-236	93418703	6260-9-236	16462701	6260-9-340	MJPHA7	6260-9-340	MJTH2W
6260-9-236	MHTH6D	6260-9-236	94500002	6260-9-236	16551401	6260-9-340	MJPHHX	6260-9-240	MJTTJL
6260-9-236	MHTJ88	6260-9-236	94500001	6260-9-236	33637802	6260-9-440	MJPHR3	6260-9-444	MJTDNP
6260-9-236	MJA0RR	6260-9-236	93881203	6260-9-236	15629501	6260-9-440	MJPJ9W	6260-9-244	MJT4HE
6260-9-236	MJAE2T	6260-9-236	95178301	6260-9-236	16650801	6260-9-244	MJPJNW	6260-9-244	MJTXEY
6260-9-440	MJM3WV	6260-9-244	MJT4W7	6260-9-444	MJTH2J	6260-9-244	MJPK4W	6260-9-244	MJT9RP
6260-9-440	MJM9YJ	6260-9-340	MJTHPM	6260-9-240	MJTR6K	6260-9-340	MJPKA2	6260-9-340	MJTHMN
6260-9-440	MJMHMH	6260-9-444	MJT5KT	6260-9-244	MJT8JK	6260-9-340	MJPPVR	6260-9-244	MJTWK7
6260-9-440	MJN3V4	6260-9-240	MJTY84	6260-9-244	MJTNNW	6260-9-340	MJPR03	6260-9-340	MJT2AX
6260-9-440	MJNJ6J	6260-9-244	MJTJ12	6260-9-244	MJT7V3	6260-9-440	MJPR59	6260-9-440	MJTMN6
6260-9-244	MJP2WM	6260-9-244	MJT2HA	6260-9-340	MJTJ8H	6260-9-440	MJPR6P	6260-9-244	MJT6NJ
6260-9-340	MJP3TT	6260-9-440	MJT99A	6260-9-340	MJTXMY	6260-9-340	MJPRYY	6260-9-340	MJTW1P
6260-9-244	MJP49V	6260-9-244	MJT7DL	6260-9-444	MJT1AT	6260-9-440	MJPT47	6260-9-340	MJTJD7
6260-9-340	MJP691	6260-9-444	MJTXP1	6260-9-444	MJT9XM	6260-9-244	MJPWJE	6260-9-244	MJT4MH
6260-9-340	MJP76E	6260-9-440	MJT9RH	6260-9-244	MJTK68	6260-9-440	MJPHYX1	6260-9-340	MJTV6X
6260-9-440	MJP79L	6260-9-440	MJTKYX	6260-9-444	MJTEA3	6260-9-240	MJR1A7	6260-9-440	MJTVW8
6260-9-340	MJP83H	6260-9-240	MJT9Y8	6260-9-444	MJT1JH	6260-9-240	MJR423	6260-9-340	MJT1LA
6260-9-340	MJP8H4	6260-9-244	MJTJLY	6260-9-240	MJTK89	6260-9-240	MJR5L4	6260-9-244	MJT3V0
6260-9-444	MJTPD2	6260-9-444	MJTKY4	6260-9-444	MJTV9A	6260-9-240	MJR668	6260-9-244	MJTLEA
6260-9-340	MJT54P	6260-9-340	MJTK9L	6260-9-244	MJTX74	6260-9-440	MJRHPY	6260-9-240	MJTWLA
6260-9-240	MKAH3V	6260-9-240	MKA7WH	6260-9-244	MJTRJJ	6260-9-240	MJRHX4	6260-9-244	MJTLD3
6260-9-444	MJT480	6260-9-244	MJTWET	6260-9-444	MJTH8V	6260-9-340	MJRL6W	6260-9-440	MJTYPK
6260-9-240	MKAY2H	6260-9-244	MJT077	6260-9-340	MJTTPPT	6260-9-440	MJR8N8V	6260-9-236	32797902
6260-9-340	MJTNR8	6260-9-444	MJTYKP	6260-9-244	MJTVNL	6260-9-240	MJRNNW	6260-9-236	33637801
6260-9-240	MKAELX	6260-9-240	MJTHE2	6260-9-340	MJT0T5	6260-9-240	MJRVPL	6260-9-236	33187602
6260-9-444	MJTD3R	6260-9-340	MJT09H	6260-9-340	MJT8ED	6260-9-444	MJT08N	6260-9-236	32729101
6260-9-444	MJTVJK	6260-9-244	MJTY2M	6260-9-240	MJT810	6260-9-340	MJT5R1	6260-9-236	31337903
6260-9-444	MJTLW5	6260-9-240	MJTWVL	6260-9-240	MJT4W0	6260-9-444	MJTD9R	6260-9-236	37633002

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-444	MJTENA	6260-9-440	MJT65Y	6260-9-240	MJTA7J	6260-9-444	MJTEWK	6260-9-236	33187603
6260-9-444	MJTTEX	6260-9-444	MJTR71	6260-9-244	MJTKJ1	6260-9-240	MJTH41	6260-9-236	36141502
6260-9-444	MJT9PM	6260-9-240	MJTA44	6260-9-444	MJTML2	6260-9-340	MJTHX3	6260-9-236	33994401
6260-9-340	MJTP06	6260-9-240	MJTPRA	6260-9-244	MJTND4	6260-9-240	MJTJ62	6260-9-236	37381202
6260-9-444	MJTL7M	6260-9-240	MJTOWE	6260-9-340	MJT8D9	6260-9-440	MJTRYW	6260-9-236	37633005
6260-9-444	MJTAMH	6260-9-240	MJTEJN	6260-9-244	MJT044	6260-9-444	MJTVJ5	6260-9-236	22721102A
6260-9-240	MJT78H	6260-9-240	MJTH8L	6260-9-340	MJTP88	6260-9-440	MJTY2L	6260-9-236	27399002A
6260-9-340	MJTWJM	6260-9-444	MJT9WY	6260-9-240	MJT4T8	6260-9-444	MJTY83	6260-9-236	MER237
6260-9-236	30264601	6260-9-236	37063302	6260-9-236	96419702	6260-9-236	34711703	6260-9-236	34015802
6260-9-236	30501401	6260-9-236	36141501	6260-9-236	96419703F	6260-9-236	34711704	6260-9-236	25290802
6260-9-236	30647802	6260-9-236	35921902	6260-9-236	33637803	6260-9-236	34015801	6260-9-236	25290801
6260-9-236	30887201	6260-9-236	36498402	6260-9-236	33994402	6260-9-236	35086004	6260-9-236	90405401C
6260-9-236	30887202	6260-9-236	36141503	6260-9-236	32715601	6260-9-236	35086003	6260-9-236	25447902
6260-9-236	31337902	6260-9-236	37194802	6260-9-236	33994404	6260-9-236	35086001	6260-9-236	24806301
6260-9-236	32147901	6260-9-236	36498401	6260-9-236	34305004	6260-9-236	35086006	6260-9-236	26205701
6260-9-236	32147902	6260-9-236	36498403	6260-9-236	34305001	6260-9-236	35725603	6260-9-236	24173502
6260-9-236	32147903	6260-9-236	36923402	6260-9-236	33994403	6260-9-236	35725601	6260-9-236	24282601
6260-9-236	32147904	6260-9-236	16667001A	6260-9-236	33637804	6260-9-236	35725602	6260-9-236	25344101
6260-9-236	32693901	6260-9-236	37063303	6260-9-236	34305003	6260-9-236	36141504	6260-9-236	24806302
6260-9-236	32715602	6260-9-236	37262102	6260-9-236	34356401	6260-9-236	35725604	6260-9-236	25344102
6260-9-236	32797901	6260-9-236	37063301	6260-9-236	34711702	6260-9-236	35827401	6260-9-236	26286102
6260-9-236	35086002	6260-9-240	MJT06D	6260-9-240	MJTLKX	6260-9-236	3T4MMAX1	6260-9-236	27399001
6260-9-236	38005401	6260-9-236	34305002	6260-9-236	34711705	6260-9-236	35086005	6260-9-236	97848701A
6260-9-236	MER29DX1	6260-9-236	MJERMX	6260-9-236	96419701	6260-9-236	18165203	6260-9-236	25406901
6260-9-236	MER3DHX1	6260-9-236	MJENNM	6260-9-236	96419704	6260-9-236	17159201	6260-9-236	24867501
6260-9-236	MER3KNX1	6260-9-236	MJERPY	6260-9-236	95501601	6260-9-236	96419703D	6260-9-236	26286103
6260-9-236	MET24YX1	6260-9-236	MJJ12J	6260-9-236	98883601	6260-9-236	18748602	6260-9-236	26147801
6260-9-236	MET47MX1	6260-9-236	MHRKM1	6260-9-236	94194801	6260-9-236	16575701	6260-9-236	25447901
6260-9-236	MET542X1	6260-9-236	MJHR1M	6260-9-236	96837501	6260-9-236	18299001	6260-9-236	27320603
6260-9-236	MET5J9X1	6260-9-236	MJERVM	6260-9-236	95501602	6260-9-236	16667001	6260-9-236	26147802
6260-9-236	METLMKX1	6260-9-236	MJD66Y	6260-9-236	95178302	6260-9-236	16575702	6260-9-236	26205702
6260-9-236	METLMKX2	6260-9-236	MHRJDE	6260-9-236	96419703	6260-9-236	18165201	6260-9-236	26286101
6260-9-236	METLYRX1	6260-9-236	MJEPP9	6260-9-236	96419702	6260-9-236	17159202	6260-9-236	28340303
6260-9-236	METM29	6260-9-236	MJH5P4	6260-9-236	96837502	6260-9-236	96419703	6260-9-236	26205701C
6260-9-236	MHAAN9X1	6260-9-236	MJJJDHY	6260-9-236	9999	6260-9-236	19648801	6260-9-236	28340302
6260-9-236	MHADY5X1	6260-9-236	86333301	6260-9-236	12614602	6260-9-236	20892102	6260-9-236	28238401
6260-9-236	MHAM27X1	6260-9-236	86376104	6260-9-236	99370201	6260-9-236	18165202	6260-9-236	27257801
6260-9-236	MHAM8W	6260-9-236	MJENEK	6260-9-236	96904401	6260-9-236	18748601	6260-9-236	28281801
6260-9-236	MHAN38X1	6260-9-236	86466501	6260-9-236	10927701	6260-9-236	19648802	6260-9-236	28340301
6260-9-236	MHAN5NX1	6260-9-236	86376101	6260-9-236	97848702	6260-9-236	19780501	6260-9-236	27320601
6260-9-236	MHANEJ	6260-9-236	MJHH3R	6260-9-236	97848701	6260-9-236	18800801	6260-9-236	27399002
6260-9-236	MHD4MJ	6260-9-236	MJJDA9	6260-9-236	99143401	6260-9-236	93881203	6260-9-236	27320602
6260-9-236	MHDHTJX1	6260-9-236	86149602	6260-9-236	10888601	6260-9-236	20749501	6260-9-236	27399004
6260-9-236	MHDJMYX1	6260-9-236	86149601	6260-9-236	10353501	6260-9-236	18800802	6260-9-236	28340305
6260-9-236	MHDL4XX1	6260-9-236	86376105	6260-9-236	11135102	6260-9-236	19842001	6260-9-236	27399003
6260-9-236	MHDL62X1	6260-9-236	86149605A	6260-9-236	10176901	6260-9-236	19716301	6260-9-236	28340304
6260-9-236	MHE6LDX1	6260-9-236	86149603	6260-9-236	12849703A	6260-9-236	21688204	6260-9-236	28846501
6260-9-236	MHE6P9X1	6260-9-236	86466502	6260-9-236	11135101	6260-9-236	19842002	6260-9-236	28287701
6260-9-236	MHE7EEX1	6260-9-236	90405401	6260-9-236	14814401	6260-9-236	22663501	6260-9-236	28856601
6260-9-236	MHEL53X1	6260-9-236	89839901	6260-9-236	12849703	6260-9-236	21688201	6260-9-236	29209903
6260-9-236	MHEL54X1	6260-9-236	86333303	6260-9-236	12614601	6260-9-236	20892105	6260-9-236	30076901
6260-9-236	MHEL8AX1	6260-9-236	86333302	6260-9-236	12849701A	6260-9-236	20892104	6260-9-236	30032401
6260-9-236	MHELMMX1	6260-9-236	86149604	6260-9-236	12849701	6260-9-236	20892103	6260-9-236	30032403
6260-9-236	MHEPD3X1	6260-9-236	90405402	6260-9-236	12849704	6260-9-236	92535001	6260-9-236	30264602
6260-9-236	MHERK5X1	6260-9-236	89445901	6260-9-236	12656001	6260-9-236	20749502	6260-9-236	067MLE
6260-9-236	MHETOX	6260-9-236	86376103	6260-9-236	99143402	6260-9-236	19780502	6260-9-236	JOXMEE
6260-9-236	MHETOXX1	6260-9-236	86149605	6260-9-236	12614603	6260-9-236	20892101	6260-9-236	XAYMJE
6260-9-236	MHH256X1	6260-9-236	91124301	6260-9-236	11991502	6260-9-236	21688203	6260-9-236	MER87J
6260-9-236	MHH2WYX1	6260-9-236	86376102	6260-9-236	12849702	6260-9-236	20749504	6260-9-236	MEPVPJ
6260-9-236	MHH904X1	6260-9-236	89839902	6260-9-236	13692401	6260-9-236	21762401	6260-9-236	MERAA8

Affected Lots for RA2016-028 LFIT Anatomic CoCr V40 Femoral Heads

Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers	Item Numbers	Lot Numbers
6260-9-236	MHJKKE	6260-9-236	89445902	6260-9-236	14934501	6260-9-236	22818101	6260-9-236	VOKMME
6260-9-236	MHJL5A	6260-9-236	91798001	6260-9-236	15029101	6260-9-236	22663503	6260-9-236	MEP9PH
6260-9-236	MHJL8N	6260-9-236	91798101	6260-9-236	97848701	6260-9-236	22721101	6260-9-236	1EP07X
6260-9-236	MHLM6W	6260-9-236	91798002	6260-9-236	95501601	6260-9-236	22663502	6260-9-236	MEPERN
6260-9-236	MHMLLO	6260-9-236	89839903	6260-9-236	15186901	6260-9-236	22721102	6260-9-236	MEPVP2X1
6260-9-236	MHMMM9	6260-9-236	91798102	6260-9-236	15341601	6260-9-236	96837501	6260-9-236	NAWMNE
6260-9-236	MHMP40	6260-9-236	92546801	6260-9-236	16551402	6260-9-236	23524902	6260-9-236	MERD36
6260-9-236	MHMP6E	6260-9-236	92546802	6260-9-236	15759701	6260-9-236	22721104	6260-9-236	METHE5
6260-9-236	23472101	6260-9-236	30076902	6260-9-344	MJTPXN	6260-9-236	HA7MLDX1	6260-9-344	MJT18L
6260-9-236	23472102	6260-9-236	30501403	6260-9-344	MJTLPN	6260-9-236	PRJMMDX1	6260-9-344	MJTX58
6260-9-236	22721103	6260-9-236	29209901	6260-9-344	MJT9XH	6260-9-236	4AEMMDX1	6260-9-344	MJTOXH
6260-9-236	23524901	6260-9-236	30501402	6260-9-344	MJTWW5	6260-9-236	8DDMRDX1	6260-9-344	MJTVWR
6260-9-236	21688202	6260-9-236	24867601	6260-9-344	MJTAR8	6260-9-236	21164802X1	6260-9-236	19880502X1
6260-9-236	24173503	6260-9-236	30264603	6260-9-344	MJT585	6260-9-236	KM6MDDX1	6260-9-344	MJTH5J
6260-9-236	23472103	6260-9-236	30647803	6260-9-344	MJTP04	6260-9-236	T6WMHD	6260-9-236	19806601X1
6260-9-236	24173501	6260-9-236	30647801	6260-9-344	MJT3J6	6260-9-236	22379906X1	6260-9-236	20142804X1X1
6260-9-240	MKAJJ3	6260-9-236	0R5MPAX3	6260-9-344	MJT2H2	6260-9-236	H9LMPDX1	6260-9-344	MJTT4J
6260-9-444	MJTT36	6260-9-236	75JMPAX1	6260-9-344	MJT8M3	6260-9-236	21451805X1	6260-9-344	MJT4TP
6260-9-240	MKA0WM	6260-9-236	20142803X2X2	6260-9-344	MJT52A	6260-9-236	691MNDX1	6260-9-344	MJT130
6260-9-240	MKAREX	6260-9-236	20983003X1	6260-9-344	MJTPJX	6260-9-236	E3MMAE	6260-9-236	3RYMMAX1
6260-9-240	MKAX3R	6260-9-236	D8VMTAX1	6260-9-344	MJT4Y3	6260-9-236	1K0MEDX1	6260-9-344	MJTKNM
6260-9-344	MJTH6V	6260-9-236	D95MAE	6260-9-344	MJT36M	6260-9-236	J0PMEE	6260-9-236	25WMNAX1
6260-9-344	MJT102	6260-9-236	9R5MJDX1	6260-9-344	MJTTPK	6260-9-236	PRKMMDX1	6260-9-236	J1MMND
6260-9-236	RJMMNAX2	6260-9-236	MER1EMX1	6260-9-344	MJT91K	6260-9-236	Y2NMHEX1	6260-9-236	THVMHE
6260-9-344	MJMD2A	6260-9-236	21318701X1	6260-9-344	MJTWTO	6260-9-236	OW8MMDX1	6260-9-236	2J3MHD
6260-9-344	MJT10T	6260-9-236	9REMJD	6260-9-236	20113203X2	6260-9-236	MET7N3X1	6260-9-236	21035305X1
6260-9-344	MJTL85	6260-9-236	8XHMNDX1	6260-9-344	MJTVA1	6260-9-236	JENMTDX1	6260-9-236	LAVMADX1
6260-9-344	MJTV6R	6260-9-236	21220205X1	6260-9-236	ORDMPAX1	6260-9-236	MER1JKX1	6260-9-236	21318702
6260-9-344	MJTVMT	6260-9-240	MKAN00	6260-9-344	MJT7YX	6260-9-344	MJM90X	6260-9-344	MJT5E7
6260-9-344	MJT1PW								