

21.06.2016

PQA 2016-04: IHD / ID-DiaPanel and ID-DiaPanel-P

PRODUCT DETAILS:

Product Name	Catalog No	Lot No	Expiry date
ID-DiaPanel	004114	45161.85.1	01.08.2016
ID-DiaPanel-P	004214	45171.85.1	01.08.2016

DESCRIPTION OF THE ISSUE:

It has come to our notice that cell 7 of the above-mentioned panel is incorrectly declared as being Kp(a+) KEL:3 on the antigen table of the references 45161.85.x, 45171.85.x. This cell is Kp(a-) KEL:-3. All the other information on the table has been rechecked and confirmed as being correct.

IMPLEMENTATION INSTRUCTIONS:

A corrected antigen table is attached with the following new references 45161.85.z, 45171.85.z. Please ensure that all customers who have received this panel are informed and that they replace the antigen tables provided with the kit with the corrected antigen table attached.

If you have any additional questions on this matter, we kindly ask you to get in touch with our Complaint Handling department (email address: slabor_cressier@bio-rad.com) quoting this PQN reference.

With best regards,

Bio-Rad

Diana [REDACTED]
Site Complaint Management Representative
Immunohematology Division

[REDACTED]
Jérémy Poropane
Acting Director RA, EMEA, CDG
Immunohematology Division

[REDACTED]
[REDACTED]
[REDACTED] Manager
Immunohematology Division