

## Product Notification

27 Dec 2018

**Product: Panocell®-16**

**LOT** 41531

 2018-12-14

**REF** 0002332

**UDI (DI) 10888234000242**

**Dear Valued Customer:**  
**(cc Chairman Medical Board and relevant Head of Departments)**

You received the product referenced above. Immucor is issuing a field notification regarding an antigen typing for vial #4 of the affected products.

**Issue:**

Immucor is issuing this product notification regarding the Js<sup>a</sup> antigen typing for donor D2070, vial #4 of lots 34457 and 41531. Additional DNA testing of donor D2070 has confirmed the donor as Js(a+).

**Impact to Results:**

Since donor D2070 is Js(a+), the red blood cells in vial 4 will be reactive with specimens that may contain anti-Js<sup>a</sup>; Js<sup>a</sup> (KEL6) is a low-prevalence antigen in the Kell blood group system. You can evaluate the impact on any previous testing results using these lots per your procedures.

**Required Action:**

No further action is required.

Please acknowledge receipt of this notification by **returning a completed response form by January 3, 2019**, to Cherry, [REDACTED]

If you have questions or require additional information, please contact us at [REDACTED] (Cherry).

We apologize for the inconvenience this issue may have caused your laboratory.

Sincerely,



Cherry Teh  
Business Executive



## Product Notification

### RESPONSE FORM

I acknowledge that our facility is aware of the Product Notification FA-NOR-18-12
Name:
Position:
Facility:
Account number:
Address:
Telephone:
Stamp and Sign:
Date: