



URGENT MEDICAL DEVICE RECALL

Commercial name of the affected product:

Pipeline™ embolization device (Pipeline™ Classic)
Alligator™ retrieval device
Marathon™ microcatheter
UltraFlow™ HPC microcatheter
Xcelerator™ hydrophilic guidewire

FSCA-identifier:

FCA-2016-01

Type of action:

Recall

Date: 06 October 2016

Dear Healthcare Professional and Risk Manager

Medtronic has identified the potential for an issue with a specific set of lot numbers of Medtronic Neurovascular products as noted below where the PTFE (polytetrafluoroethylene) coating could delaminate and separate from the delivery wire. As a result, we are recalling these products from the specific lots as listed in Appendix 1 of this letter. This issue affects lots of the following Medtronic products manufactured between July 2014 and September 2016:

1. Pipeline™ embolization device
2. Alligator™ retrieval device
3. Marathon™ microcatheter stylets
4. UltraFlow™ HPC microcatheter stylets
5. Xcelerator™ hydrophilic guidewire

Note: for the Pipeline™ embolization device, this issue only affects the delivery wire. The implant itself is not affected.

Description of the problem:

Delamination and detachment of the PTFE coating material may lead to PTFE coating in the blood stream. PTFE in the blood stream, based on its specific size and quantity, may be safe or could lead to thromboembolic complications or irreversible injuries including but not limited to: intracranial edema, peripheral edema, incomplete treatment, infection, local inflammatory response, systemic inflammatory response, ischemic stroke, lysis/necrosis, neurological deficit, organ impairment, shock, space occupying lesion, thrombosis, and hemorrhage.

Through September 27, 2016, Medtronic has received a total of five (5) reports for the products that are potentially affected by this issue. Medtronic has received no reports of serious injuries or patient death as a result of this issue.

Recommendations on action to be taken by the user:

For affected product that has been used, no action is necessary and patients should continue to be managed in accordance with your standard patient management protocol.

For affected product that has not been used, Medtronic requests that you take immediate action. Our records show that your facility has received one or more lots of the affected products listed on page 1. Consequently, Medtronic requests that you immediately take the following actions:

1. Remove and quarantine all unused affected products in your inventory.
2. Return the potentially affected products to Medtronic. Your Medtronic representative can assist in facilitating the return of product as necessary. If replacement product is needed, your Medtronic representative can assist you with identifying suitable replacement product.
3. Complete the attached Customer Confirmation Certificate **and contact your local Sales representative.**

Medtronic has taken the necessary steps to prevent future shipment of the potentially affected product.

Transmission of this Field Safety Notice:

This field safety notice needs to be passed on to all those who need to be aware within your organization or to any organization where the potentially affected devices have been transferred.


Please transfer this notice to other organizations on which this action has an impact.

Regulatory notification:

Medtronic is communicating this information to the appropriate regulatory agency in your country.

We apologize for the inconvenience this may cause you; please be assured that patient safety and product quality remains our primary concern. Should you have any questions, please contact your Medtronic representative.

Sincerely,



Senior Regional Business Director - NV

Appendix 1 – List of Affected Lots (cont.):

Pipeline Classic Catalog Numbers

FA-71300-25V09	FA-71350-25V16	FA-71400-25V10	FA-71425-35V13	FA-71475-30V16
FA-71300-25V12	FA-71350-25V17	FA-71400-25V11	FA-71425-35V15	FA-71475-30V18
FA-71300-25V13	FA-71350-25V18	FA-71400-25V12	FA-71425-35V18	FA-71475-35V09
FA-71300-25V15	FA-71350-25V19	FA-71400-25V13	FA-71450-25V09	FA-71475-35V10
FA-71300-25V16	FA-71350-30V09	FA-71400-25V15	FA-71450-25V10	FA-71475-35V11
FA-71300-25V17	FA-71350-30V10	FA-71400-25V16	FA-71450-25V11	FA-71475-35V12
FA-71300-25V18	FA-71350-30V12	FA-71400-25V18	FA-71450-25V12	FA-71475-35V13
FA-71300-25V19	FA-71350-30V13	FA-71400-25V19	FA-71450-25V13	FA-71475-35V15
FA-71300-30V09	FA-71350-30V15	FA-71400-30V09	FA-71450-25V15	FA-71475-35V18
FA-71300-30V12	FA-71350-30V16	FA-71400-30V10	FA-71450-25V16	FA-71475-35V19
FA-71300-30V13	FA-71350-30V18	FA-71400-30V11	FA-71450-25V18	FA-71500-25V09
FA-71300-30V15	FA-71350-35V09	FA-71400-30V12	FA-71450-25V19	FA-71500-25V10
FA-71300-30V16	FA-71350-35V10	FA-71400-30V13	FA-71450-30V09	FA-71500-25V11
FA-71300-30V19	FA-71350-35V12	FA-71400-30V15	FA-71450-30V10	FA-71500-25V12
FA-71300-35V09	FA-71350-35V13	FA-71400-30V16	FA-71450-30V11	FA-71500-25V13
FA-71300-35V12	FA-71350-35V15	FA-71400-30V18	FA-71450-30V12	FA-71500-25V15
FA-71300-35V13	FA-71350-35V18	FA-71400-35V09	FA-71450-30V13	FA-71500-25V16
FA-71300-35V15	FA-71375-25V09	FA-71400-35V13	FA-71450-30V15	FA-71500-25V18
FA-71325-25V09	FA-71375-25V10	FA-71400-35V15	FA-71450-30V16	FA-71500-25V19
FA-71325-25V10	FA-71375-25V11	FA-71400-35V18	FA-71450-30V18	FA-71500-30V09
FA-71325-25V12	FA-71375-25V12	FA-71425-25V09	FA-71450-30V19	FA-71500-30V10
FA-71325-25V13	FA-71375-25V13	FA-71425-25V10	FA-71450-35V09	FA-71500-30V11
FA-71325-25V15	FA-71375-25V15	FA-71425-25V11	FA-71450-35V10	FA-71500-30V12
FA-71325-25V18	FA-71375-25V16	FA-71425-25V12	FA-71450-35V12	FA-71500-30V13
FA-71325-25V19	FA-71375-25V17	FA-71425-25V13	FA-71450-35V13	FA-71500-30V15
FA-71325-30V09	FA-71375-25V18	FA-71425-25V15	FA-71450-35V15	FA-71500-30V16
FA-71325-30V10	FA-71375-25V19	FA-71425-25V16	FA-71450-35V16	FA-71500-30V17
FA-71325-30V12	FA-71375-30V09	FA-71425-25V18	FA-71475-25V09	FA-71500-30V18
FA-71325-30V13	FA-71375-30V10	FA-71425-25V19	FA-71475-25V10	FA-71500-35V09
FA-71325-30V15	FA-71375-30V12	FA-71425-30V09	FA-71475-25V11	FA-71500-35V10
FA-71325-30V16	FA-71375-30V13	FA-71425-30V10	FA-71475-25V12	FA-71500-35V11
FA-71325-30V18	FA-71375-30V15	FA-71425-30V11	FA-71475-25V13	FA-71500-35V12
FA-71325-30V19	FA-71375-30V16	FA-71425-30V12	FA-71475-25V15	FA-71500-35V13
FA-71325-35V09	FA-71375-30V18	FA-71425-30V13	FA-71475-25V16	FA-71500-35V15

*Appendix 1 – List of Affected Lots (cont.):
Pipeline Classic Catalog Numbers*

FA-71325-35V12	FA-71375-35V09	FA-71425-30V15	FA-71475-25V18	FA-71500-35V16
FA-71325-35V13	FA-71375-35V12	FA-71425-30V16	FA-71475-30V09	FA-71500-35V17
FA-71325-35V15	FA-71375-35V13	FA-71425-30V18	FA-71475-30V10	FA-71500-35V18
FA-71350-25V09	FA-71375-35V15	FA-71425-30V19	FA-71475-30V11	FA-71500-35V19
FA-71350-25V10	FA-71375-35V16	FA-71425-35V09	FA-71475-30V12	FA-77250-10V10
FA-71350-25V13	FA-71375-35V18	FA-71425-35V11	FA-71475-30V13	FA-77250-10V11
FA-71350-25V15	FA-71400-25V09	FA-71425-35V12	FA-71475-30V15	FA-77250-10V12
FA-77250-10V13	FA-77275-14V11	FA-77300-16V19	FA-77325-20V11	FA-77350-20V16
FA-77250-10V19	FA-77275-14V12	FA-77300-18V10	FA-77325-20V12	FA-77350-20V17
FA-77250-12V11	FA-77275-14V13	FA-77300-18V11	FA-77325-20V13	FA-77350-20V18
FA-77250-12V12	FA-77275-14V17	FA-77300-18V13	FA-77325-20V15	FA-77375-10V11
FA-77250-12V13	FA-77275-14V18	FA-77300-18V16	FA-77325-20V16	FA-77375-10V12
FA-77250-12V18	FA-77275-14V19	FA-77300-18V18	FA-77325-20V18	FA-77375-10V13
FA-77250-14V12	FA-77275-16V12	FA-77300-18V19	FA-77350-10V11	FA-77375-10V16
FA-77250-14V13	FA-77275-16V13	FA-77300-20V09	FA-77350-10V12	FA-77375-10V18
FA-77250-14V15	FA-77275-16V16	FA-77300-20V10	FA-77350-10V13	FA-77375-10V19
FA-77250-14V18	FA-77275-16V18	FA-77300-20V11	FA-77350-10V18	FA-77375-12V11
FA-77250-14V19	FA-77275-16V19	FA-77300-20V12	FA-77350-10V19	FA-77375-12V12
FA-77250-16V11	FA-77275-18V10	FA-77300-20V13	FA-77350-12V11	FA-77375-12V13
FA-77250-16V12	FA-77275-18V12	FA-77300-20V15	FA-77350-12V12	FA-77375-12V17
FA-77250-16V13	FA-77275-18V16	FA-77325-10V11	FA-77350-12V13	FA-77375-12V19
FA-77250-16V15	FA-77275-18V18	FA-77325-10V12	FA-77350-12V19	FA-77375-14V11
FA-77250-16V16	FA-77275-18V19	FA-77325-10V18	FA-77350-14V10	FA-77375-14V12
FA-77250-16V19	FA-77275-20V10	FA-77325-10V19	FA-77350-14V11	FA-77375-14V13
FA-77250-18V09	FA-77275-20V13	FA-77325-12V11	FA-77350-14V12	FA-77375-14V18
FA-77250-18V10	FA-77275-20V15	FA-77325-12V12	FA-77350-14V13	FA-77375-14V19
FA-77250-18V11	FA-77275-20V16	FA-77325-12V13	FA-77350-14V16	FA-77375-16V10
FA-77250-18V12	FA-77275-20V18	FA-77325-12V18	FA-77350-14V18	FA-77375-16V11
FA-77250-18V13	FA-77300-10V11	FA-77325-12V19	FA-77350-14V19	FA-77375-16V12
FA-77250-18V16	FA-77300-10V12	FA-77325-14V11	FA-77350-16V10	FA-77375-16V13
FA-77250-18V18	FA-77300-10V13	FA-77325-14V12	FA-77350-16V11	FA-77375-16V16
FA-77250-18V19	FA-77300-10V16	FA-77325-14V13	FA-77350-16V12	FA-77375-16V17
FA-77250-20V10	FA-77300-10V19	FA-77325-14V18	FA-77350-16V13	FA-77375-16V18
FA-77250-20V11	FA-77300-12V11	FA-77325-14V19	FA-77350-16V17	FA-77375-16V19

Appendix 1 – List of Affected Lots (cont.):

Pipeline Classic Catalog Numbers

FA-77250-20V12	FA-77300-12V12	FA-77325-16V10	FA-77350-16V18	FA-77375-18V09
FA-77250-20V13	FA-77300-12V13	FA-77325-16V11	FA-77350-18V10	FA-77375-18V10
FA-77250-20V15	FA-77300-12V18	FA-77325-16V12	FA-77350-18V11	FA-77375-18V11
FA-77250-20V16	FA-77300-12V19	FA-77325-16V13	FA-77350-18V12	FA-77375-18V12
FA-77250-20V18	FA-77300-14V11	FA-77325-16V18	FA-77350-18V13	FA-77375-18V13
FA-77275-10V12	FA-77300-14V12	FA-77325-18V10	FA-77350-18V17	FA-77375-18V16
FA-77275-10V13	FA-77300-14V13	FA-77325-18V12	FA-77350-18V18	FA-77375-18V17
FA-77275-10V16	FA-77300-14V16	FA-77325-18V13	FA-77350-18V19	FA-77375-18V18
FA-77275-10V18	FA-77300-14V19	FA-77325-18V16	FA-77350-20V09	FA-77375-18V19
FA-77275-10V19	FA-77300-16V10	FA-77325-18V17	FA-77350-20V10	FA-77375-20V09
FA-77275-12V11	FA-77300-16V11	FA-77325-18V18	FA-77350-20V11	FA-77375-20V10
FA-77275-12V12	FA-77300-16V12	FA-77325-18V19	FA-77350-20V12	FA-77375-20V11
FA-77275-12V13	FA-77300-16V13	FA-77325-20V09	FA-77350-20V13	FA-77375-20V12
FA-77275-12V19	FA-77300-16V18	FA-77325-20V10	FA-77350-20V15	FA-77375-20V13
FA-77375-20V15	FA-77400-18V18	FA-77425-18V18	FA-77450-18V18	FA-77475-20V11
FA-77375-20V16	FA-77400-18V19	FA-77425-18V19	FA-77450-18V19	FA-77475-20V12
FA-77375-20V17	FA-77400-20V09	FA-77425-20V09	FA-77450-20V09	FA-77475-20V13
FA-77375-20V18	FA-77400-20V10	FA-77425-20V10	FA-77450-20V10	FA-77475-20V15
FA-77400-10V09	FA-77400-20V11	FA-77425-20V11	FA-77450-20V11	FA-77475-20V16
FA-77400-10V11	FA-77400-20V12	FA-77425-20V12	FA-77450-20V12	FA-77475-20V17
FA-77400-10V12	FA-77400-20V13	FA-77425-20V13	FA-77450-20V13	FA-77475-20V18
FA-77400-10V13	FA-77400-20V15	FA-77425-20V15	FA-77450-20V15	FA-77500-10V12
FA-77400-10V19	FA-77400-20V16	FA-77425-20V16	FA-77450-20V16	FA-77500-10V13
FA-77400-12V09	FA-77400-20V17	FA-77425-20V17	FA-77450-20V17	FA-77500-10V18
FA-77400-12V11	FA-77400-20V18	FA-77425-20V18	FA-77450-20V18	FA-77500-10V19
FA-77400-12V12	FA-77425-10V11	FA-77450-10V11	FA-77450-20V19	FA-77500-12V11
FA-77400-12V13	FA-77425-10V12	FA-77450-10V12	FA-77475-10V12	FA-77500-12V12
FA-77400-12V17	FA-77425-10V13	FA-77450-10V13	FA-77475-10V13	FA-77500-12V13
FA-77400-12V18	FA-77425-10V19	FA-77450-10V18	FA-77475-10V18	FA-77500-12V18
FA-77400-14V09	FA-77425-12V11	FA-77450-10V19	FA-77475-10V19	FA-77500-14V11
FA-77400-14V10	FA-77425-12V12	FA-77450-12V11	FA-77475-12V12	FA-77500-14V12
FA-77400-14V11	FA-77425-12V13	FA-77450-12V12	FA-77475-12V13	FA-77500-14V13
FA-77400-14V12	FA-77425-12V17	FA-77450-12V13	FA-77475-12V18	FA-77500-14V18
FA-77400-14V13	FA-77425-12V18	FA-77450-12V18	FA-77475-14V11	FA-77500-14V19

Appendix 1 – List of Affected Lots (cont.):

Pipeline Classic Catalog Numbers

FA-77400-14V15	FA-77425-14V10	FA-77450-14V11	FA-77475-14V12	FA-77500-16V09
FA-77400-14V18	FA-77425-14V11	FA-77450-14V12	FA-77475-14V18	FA-77500-16V10
FA-77400-14V19	FA-77425-14V12	FA-77450-14V13	FA-77475-14V19	FA-77500-16V11
FA-77400-16V09	FA-77425-14V13	FA-77450-14V19	FA-77475-16V10	FA-77500-16V12
FA-77400-16V10	FA-77425-14V19	FA-77450-16V09	FA-77475-16V11	FA-77500-16V13
FA-77400-16V11	FA-77425-16V09	FA-77450-16V10	FA-77475-16V12	FA-77500-16V15
FA-77400-16V12	FA-77425-16V10	FA-77450-16V11	FA-77475-16V13	FA-77500-16V17
FA-77400-16V13	FA-77425-16V11	FA-77450-16V12	FA-77475-16V16	FA-77500-18V09
FA-77400-16V15	FA-77425-16V12	FA-77450-16V13	FA-77475-16V17	FA-77500-18V10
FA-77400-16V16	FA-77425-16V13	FA-77450-16V15	FA-77475-16V18	FA-77500-18V11
FA-77400-16V17	FA-77425-16V15	FA-77450-16V16	FA-77475-18V09	FA-77500-18V12
FA-77400-16V18	FA-77425-16V16	FA-77450-16V17	FA-77475-18V10	FA-77500-18V13
FA-77400-16V19	FA-77425-16V17	FA-77450-16V18	FA-77475-18V11	FA-77500-18V15
FA-77400-18V09	FA-77425-16V18	FA-77450-18V09	FA-77475-18V12	FA-77500-18V19
FA-77400-18V10	FA-77425-18V09	FA-77450-18V10	FA-77475-18V15	FA-77500-20V09
FA-77400-18V11	FA-77425-18V10	FA-77450-18V11	FA-77475-18V16	FA-77500-20V10
FA-77400-18V12	FA-77425-18V11	FA-77450-18V12	FA-77475-18V17	FA-77500-20V11
FA-77400-18V13	FA-77425-18V12	FA-77450-18V13	FA-77475-18V18	FA-77500-20V12
FA-77400-18V15	FA-77425-18V13	FA-77450-18V15	FA-77475-18V19	FA-77500-20V13
FA-77400-18V16	FA-77425-18V15	FA-77450-18V16	FA-77475-20V09	FA-77500-20V15
FA-77400-18V17	FA-77425-18V16	FA-77450-18V17	FA-77475-20V10	FA-77500-20V16
				FA-77500-20V17
				FA-77500-20V18

Appendix 1 – List of Affected Lots (cont.):

Alligator Retrieval Device Catalog Numbers

FA-88810-20V04
FA-88810-20V06
FA-88810-30V04
FA-88810-30V06
FA-88810-40V04
FA-88810-40V06
FA-88810-50V04
FA-88810-50V06

Marathon Device Catalog Numbers

105-5055V02
105-5055V02
105-5055V05

UltraFlow Device Catalog Numbers

105-5065V02
105-5065V05
105-5066V02
105-5066V05

Xcelerator Guidewire Device Catalog Numbers

103-0601-300	103-0601-350V10	103-0602-350V06
103-0601-300V01	103-0602-300	103-0602-350V09
103-0601-300V02	103-0602-300V01	103-0602-350V10
103-0601-300V03	103-0602-300V02	
103-0601-300V06	103-0602-300v06	
103-0601-300V10	103-0602-300V09	
103-0601-300V11	103-0602-300V10	
103-0601-350	103-0602-350	
103-0601-350V01	103-0602-350V01	
103-0601-350V06	103-0602-350V02	

Appendix 1 – List of Affected Lots (cont.):

Pipeline Classic Lots

A010308	A020911	A026317	A032375	A036898	A043417	A046146	A053154
A010309	A020912	A026319	A032378	A036899	A043663	A047160	A053878
A010310	A020913	A026835	A032381	A036900	A043667	A047161	A053881
A010312	A020914	A026836	A032383	A036901	A043668	A047162	A053882
A010314	A020916	A026842	A033000	A036904	A043911	A047808	A053993
A010316	A021047	A026843	A033001	A036905	A043912	A047809	A053994
A013030	A021410	A026939	A033237	A036906	A043913	A047810	A054421
A013031	A021411	A026940	A033238	A037896	A043914	A047811	A054422
A013591	A021412	A026941	A033240	A037897	A043915	A048395	A054423
A013592	A021479	A026942	A033242	A038406	A044125	A048396	A054424
A013594	A021483	A027316	A033243	A038407	A044126	A048397	A054425
A014538	A021484	A027317	A033246	A038410	A044127	A048398	A054426
A014542	A021677	A027319	A033738	A038411	A044128	A049717	A054427
A014543	A021678	A028325	A033740	A039064	A044129	A049718	A054714
A015450	A021679	A028326	A033831	A039066	A044130	A049721	A054777
A015451	A021680	A028327	A033832	A039067	A044131	A049722	A054778
A015557	A021681	A028328	A033833	A039068	A044278	A049724	A054779
A016139	A021688	A028588	A033879	A039727	A044279	A049726	A054780
A016140	A021691	A028591	A033880	A040185	A044284	A049727	A055473
A016152	A022538	A028609	A033881	A040188	A044285	A049933	A055475
A016155	A022707	A029206	A033882	A040330	A044287	A049938	A055476
A016797	A022710	A029208	A034154	A040331	A044378	A049939	A055477
A016798	A023091	A029223	A034156	A040480	A044379	A049942	A056174
A016801	A023092	A029315	A034201	A041006	A044380	A049945	A056178
A016802	A023094	A029316	A034203	A041009	A044381	A049947	A056180
A017024	A023255	A029965	A034334	A041510	A044382	A049949	A056181
A017485	A023256	A029968	A034335	A041511	A044826	A050932	A056182
A017486	A023261	A030489	A034336	A041513	A044827	A050933	A056768
A017487	A023269	A030490	A034337	A041514	A044829	A050935	A056771
A017489	A024218	A030542	A036036	A041870	A044831	A050936	A056775
A017861	A024228	A030543	A036189	A041871	A044834	A050973	A056778
A017862	A024268	A030544	A036338	A042005	A044873	A050974	A056794
A018467	A024577	A030545	A036339	A042392	A044883	A050976	A057026
A018469	A024908	A031358	A036340	A042396	A045690	A051066	A057027
A018470	A024917	A031360	A036391	A042399	A045691	A051067	A058140
A019062	A024920	A031361	A036404	A042400	A045692	A051069	A058142
A019063	A024923	A031364	A036514	A042658	A045693	A052643	A058143
A019065	A025001	A031794	A036606	A043061	A045695	A052644	A058145
A019499	A025002	A031795	A036607	A043063	A045697	A052645	A058146
A019500	A025634	A031948	A036608	A043064	A045698	A052646	A058147
A019616	A025635	A031949	A036749	A043065	A046133	A052648	A058148
A020032	A025654	A031951	A036750	A043272	A046134	A052649	A058505
A020186	A025655	A031953	A036777	A043273	A046139	A053151	A058506
A020488	A026312	A032012	A036778	A043285	A046141	A053152	A058507
A020909	A026314	A032372	A036780	A043416	A046142	A053153	A058508

Appendix 1 – List of Affected Lots (cont.):

Pipeline Classic Lots

A058509	A069632	A081553	A098121	A113419	A125915	A134228	A139869
A058510	A069633	A081557	A098122	A113420	A125916	A134229	A139870
A059244	A070370	A081558	A098123	A113421	A125919	A134230	A139871
A059245	A070380	A086827	A098124	A113791	A125921	A134886	A139872
A059246	A070552	A086831	A100283	A113792	A125922	A134887	A140870
A059247	A070553	A086837	A100284	A113793	A125924	A134888	A140872
A061332	A070554	A086840	A100295	A113794	A125925	A134889	A140873
A061934	A070555	A087759	A100297	A113795	A125926	A134890	A140875
A061935	A073701	A087760	A100298	A113796	A125927	A134891	A140876
A062283	A073702	A087761	A100299	A115088	A127285	A135409	A140877
A062284	A073704	A088653	A100323	A115090	A127287	A135410	A140878
A062752	A073705	A088654	A100324	A115091	A127291	A135411	A142033
A062988	A073707	A088655	A100327	A115605	A127292	A135412	A142036
A062989	A073709	A088658	A100328	A115606	A128359	A135413	A142037
A063057	A073711	A088659	A106144	A115607	A128360	A135414	A142038
A063745	A074018	A088660	A106155	A117448	A128361	A135415	A142039
A063747	A074019	A090010	A106160	A117449	A128363	A135417	A142040
A064308	A074021	A090012	A106223	A117665	A128364	A135418	A142041
A064309	A074022	A090013	A106880	A117667	A129828	A135419	A142042
A064310	A074023	A090016	A107034	A117696	A129830	A137811	A142043
A064313	A074026	A090018	A107035	A117697	A129831	A137812	A143592
A064973	A074892	A090019	A107036	A117698	A131570	A137813	A143593
A064975	A074893	A091998	A107039	A117700	A131573	A137814	A143594
A065074	A074895	A092000	A107673	A117701	A131575	A137815	A143595
A065075	A074896	A092002	A107674	A117702	A131576	A137819	A148018
A065522	A076813	A092003	A107675	A119160	A131577	A137820	A148019
A065523	A076815	A092005	A107679	A119161	A131579	A137821	A148022
A065524	A076817	A092006	A107680	A119162	A131581	A137823	A148024
A065526	A076818	A092007	A107681	A119164	A131582	A138559	A148025
A065639	A076820	A093615	A107682	A120209	A132269	A138560	A148026
A065642	A076821	A093616	A108850	A120210	A132270	A138561	A148028
A065649	A076822	A093617	A108856	A120212	A132271	A138562	A148035
A065651	A076823	A093618	A108857	A120213	A132272	A138563	A148042
A065653	A078964	A093619	A108858	A120215	A132276	A139221	A148043
A066073	A078965	A094468	A108859	A122026	A133598	A139222	A148044
A066074	A080848	A094469	A108860	A122027	A133599	A139223	A148045
A066078	A080849	A094470	A110200	A122028	A133600	A139224	A150199
A066079	A080850	A094473	A110201	A122029	A133602	A139225	A150200
A066692	A080851	A096665	A110203	A122033	A133605	A139226	A150201
A066694	A080852	A096669	A110206	A122035	A133606	A139227	A150202
A066695	A080853	A096670	A110207	A123930	A133607	A139229	A150203
A066696	A080854	A096671	A110208	A123931	A134224	A139853	A150204
A069628	A080855	A096673	A110209	A123932	A134225	A139854	A150205
A069630	A081550	A096675	A110210	A123933	A134226	A139864	A150894
A069631	A081551	A098120	A113417	A123935	A134227	A139867	A150898

Appendix 1 – List of Affected Lots (cont.):

Pipeline Classic Lots

A150899	A162848	A185769	A225554	A261101	A297968	A317168	A015452
A150900	A162850	A185771	A225555	A261102	A297970	A317171	A017863
A150901	A162851	A185772	A226043	A261106	A297971	A317937	A019067
A150902	A162853	A185773	A226051	A261109	A297973	A317940	A019687
A151239	A162854	A188190	A227129	A263738	A300193	A318732	A019498
A151240	A162855	A188195	A229179	A265492	A300194	A319040	A019689
A151241	A164383	A190077	A232637	A266994	A300991	A319043	A019686
A151242	A164384	A190079	A232638	A266997	A300992	A319130	A019688
A151243	A164385	A190080	A232640	A266999	A300993	A320823	A020051
A151244	A165179	A190081	A234756	A267000	A302644	A320826	A020185
A151245	A165180	A190082	A238664	A267001	A302645	A320827	A020487
A151246	A166121	A193900	A238666	A269767	A302646	A322378	A020486
A151248	A166122	A193902	A238667	A269771	A304830	A322381	A020485
A152139	A166123	A193904	A239875	A269773	A304842	A322384	A020484
A152140	A167455	A193905	A239877	A269774	A304844	A327069	A020184
A152141	A167457	A193906	A241057	A269775	A304846	9961450	A020035
A152142	A169229	A197312	A242547	A271509	A304847	9961902	A020047
A152143	A169230	A197314	A242548	A272748	A306243	9961900	A020910
A153769	A169231	A199979	A242550	A272750	A306245	9961899	A021481
A153770	A171151	A199983	A242559	A272752	A306251	9962409	A021690
A153771	A171156	A204514	A242561	A274216	A306253	9962411	A021689
A153772	A171163	A204515	A245762	A274223	A306254	9962666	A022709
A153773	A171165	A204516	A245765	A276038	A306256	9962420	A022520
A154448	A171168	A211767	A245767	A276041	A307932	9962984	A022522
A154450	A171248	A211769	A245768	A278715	A307935	9963258	A022519
A154452	A174423	A214097	A245770	A278718	A307940	9963003	A024283
A154453	A174424	A214779	A245771	A281388	A307941	9963797	A024276
A154454	A174426	A214781	A249221	A281389	A309598	9963255	A023270
A156152	A174427	A214782	A249225	A281390	A309599	9963792	A024213
A156157	A174428	A214783	A250673	A283037	A309601	9962996	A021046
A157201	A174430	A214784	A250675	A283040	A309603	9964272	A029219
A157202	A174432	A215491	A252980	A283041	A311498	9963787	A029216
A157203	A177772	A215494	A252981	A283042	A311499	9963786	A033834
A157204	A177775	A216381	A252983	A287710	A311500	9964271	A044281
A157205	A180608	A217701	A253701	A287712	A311502	9964269	A044288
A157206	A180610	A219813	A253703	A289358	A311503	9965334	A044282
A157207	A180654	A219814	A256125	A289361	A313697	A005065	A044280
A157210	A180656	A220212	A256127	A291847	A313698	A005067	A043664
A157211	A183853	A220502	A256129	A294769	A313699	A008831	A043301
A160075	A183854	A220506	A256855	A294772	A313700	A008829	A044283
A160076	A183875	A221351	A258774	A294776	A313701	A008832	A043661
A160077	A183876	A221352	A258777	A294786	A316285	A010289	A047036
A160919	A183877	A221354	A258784	A296212	A316286	A010917	A047751
A160920	A183879	A223464	A258846	A296214	A316436	A013596	A049094
A160921	A185767	A223465	A258848	A296216	A316437	A014546	A047788

Appendix 1 – List of Affected Lots (cont.):

Pipeline Classic Lots

A047786	A157208	A230275	A239181	A033899	A036714	A043575	A059875
A049098	A162846	A230621	A250676	A033921	A036761	A043607	A059883
A047787	A167458	A300192	A263745	A034028	A036762	A043608	A059889
A047790	A169227	A300191	A294784	A034047	A037727	A043620	A059898
A049087	A167454	A256868	A307933	A034056	A037783	A043622	A059901
A047752	A167459	A296219	A024002	A034069	A037792	A043671	A059907
A047789	A167456	A252984	A025079	A034331	A037890	A044712	A059909
A049089	A177776	A253705	A025131	A034348	A037944	A044750	A061683
A049106	A177773	A256867	A025453	A034348	A037949	A047035	A064065
A049780	A177774	A249223	A025469	A034378	A037957	A047070	A064101
A049080	A177777	A268268	A025490	A034387	A037968	A047418	A064122
A043418	A180655	A287709	A025816	A034405	A037975	A047452	A064134
A049104	A185764	A241048	A027429	A034430	A037994	A047473	A064147
A049102	A185766	A249224	A027949	A034763	A038333	A047492	A064154
A049107	A185765	A261105	A027958	A034790	A038529	A047517	A065277
A049096	A185774	A271508	A027960	A034805	A038550	A048296	A065305
A049100	A198686	A287711	A027966	A034818	A034405	A048336	A065317
A050118	A188192	A230623	A028002	A034887	A039292	A048366	A065330
A049781	A188191	A263747	A028019	A035376	A039800	A050293	A065342
A049779	A188194	A249226	A028025	A035407	A039819	A050314	A065354
A049782	A193910	A278717	A028748	A035474	A039819	A050341	A071060
A049105	A199982	A289359	A028816	A035474	A039877	A050375	A071079
A050117	A193907	A289360	A028845	A035484	A039916	A050416	A071095
A086855	A215497	A297969	A029783	A035526	A040676	A050461	A071130
A096668	A204517	A263739	A029776	A035547	A040683	A051858	A071141
A131037	A214100	A263748	A030150	A035554	A040714	A051869	A071167
A131039	A215506	A234754	A030271	A035731	A040803	A051869	A071786
A131040	A215504	A256856	A030287	A035774	A040804	A052182	A071791
A132274	A216380	A271507	A031103	A036013	A040852	A052211	A072457
A132273	A218414	A239184	A031136	A036044	A040853	A052245	A074811
A133596	A216377	A241045	A031151	A036084	A040880	A052282	A074829
A133597	A216378	A250671	A031182	A036127	A040884	A052296	A074839
A137824	A218224	A263742	A031224	A036136	A040948	A052305	A074877
A148031	A217703	A238668	A031240	A036137	A040955	A053710	A074897
A151238	A218226	A239185	A031274	A036176	A040967	A053760	A067913
A152133	A216379	A250672	A032693	A036188	A040968	A056706	A076175
A152136	A215496	A263744	A033113	A036177	A041392	A056721	A076233
A152135	A218413	A230622	A033124	A036219	A041713	A056730	A076207
A152134	A220211	A250674	A033445	A036228	A041734	A056737	A076238
A152137	A227130	A261103	A033490	A036245	A041746	A057442	A076248
A156153	A226053	A271511	A033809	A036603	A041756	A057449	A076258
A156154	A226046	A230625	A033819	A036613	A041760	A057454	A078408
A156158	A229181	A241049	A033850	A036615	A043549	A057461	A078418
A156155	A228013	A253706	A033857	A036696	A043548	A057470	A078451
A157209	A228014	A265490	A033887	A036697	A043567	A057481	A079821

Appendix 1 – List of Affected Lots (cont.):

Pipeline Classic Lots

A080450	A089979	A090139	A095566	A097029	A099768	A104413	A110371
A083417	A089991	A090167	A095618	A097729	A100859	A105772	A110376
A084728	A090011	A090181	A095637	A099267	A103256	A112687	A110383
A084761	A090047	A090209	A095694	A099348	A103270	A112690	A110389
A084820	A090070	A090241	A095716	A099367	A104383	A112693	A110396
A084834	A090076	A092811	A095724	A099374	A104400	A105810	A112046
A087689	A090107	A095509	A097001	A099734	A104409	A110354	A112674
A089934	9947374	9956715	9961032	9969715	9972030	9974427	9977717
9936174	9947375	9956717	9965803	9969716	9972031	9974429	9977718
9938400	9947376	9956719	9966844	9969717	9972037	9974430	9979347
9939295	9947429	9957296	9966845	9969762	9972507	9974926	9979595
9939297	9948948	9958024	9966850	9971070	9973294	9974927	9979597
9940045	9948952	9958034	9966852	9971071	9973296	9976336	A000429
9940046	9952154	9958040	9966853	9971072	9973298	9976338	A003369
9941533	9952183	9958091	9967670	9971303	9973299	9976339	A003372
9941535	9953585	9958112	9967674	9971309	9974011	9976342	A010293
9941539	9954488	9958114	9967677	9971310	9974012	9977016	A010297
9941540	9954489	9959779	9969487	9971644	9974026	9977711	A010297
9944947	9955853	9960360	9969489	9971646	9974028	9977712	A010914
9944948	9956647	9960363	9969507	9971648	9974419	9977715	A011604
A011605	A011611	A012315					

Appendix 1 – List of Affected Lots (cont.):

Alligator Retrieval Device (ARD) Lots

9948077	9972887	A010848	A053541	A116064	A131036
9948396	9972889	A011322	A053542	A117476	A132044
9952512	9972890	A011359	A053543	A118368	A154323
9957329	9974418	A013554	A053763	A120498	A154324
9957336	9976051	A013556	A053764	A122265	A156580
9957342	9976060	A014832	A053766	A123213	A156581
9944965	9977018	A014868	A066445	A124445	A161317
9945280	9977019	A014953	A066446	A125283	A161318
9951876	9977396	A016804	A068660	A132937	A163895
9954273	9977397	A018473	A068664	A134257	A163896
9957330	9978556	A019068	A071295	A134900	A166199
9957339	9978557	A021485	A072578	A135945	A168759
9957343	9979590	A022023	A091991	A137794	A168760
9957347	9980261	A022527	A093697	A138489	A173721
9945281	9980272	A024541	A094487	A142716	A173722
9946234	A002055	A026133	A095510	A147975	A177071
9952511	A002875	A027036	A096506	A147976	A177073
9957331	A002976	A028730	A098665	A147978	A197692
9957340	A004046	A028731	A098667	A150526	A242153
9957344	A004049	A030919	A102229	A150527	A242931
9948387	A005156	A030920	A102232	A150529	A242936
9950866	A006168	A033160	A104111	A152100	A242937
9953649	A006169	A033161	A104113	A153748	A242939
9954974	A007336	A036879	A107473	A153749	A247199
9955940	A007955	A037969	A110323	A127224	A247612
9957333	A008188	A039135	A111254	A127844	A254814
9957341	A009329	A039136	A112093	A129002	A282350
9957346	A010342	A039142	A114200	A130573	A290956
A311410	A312892	A321463			

Appendix 1 – List of Affected Lots (cont.):

Marathon

9946730	9960999	9965968	9974493	A004327	A011297	A016502	A024591
9952156	9961016	9965970	9975117	A004329	A011515	A016504	A024594
9952519	9961297	9965980	9975118	A004435	A011522	A016507	A024893
9952892	9961308	9965983	9976018	A004436	A011731	A016510	A024894
9953171	9961313	9965984	9976024	A004437	A011735	A016632	A024896
9953177	9963245	9966022	9976049	A004535	A011818	A017149	A025130
9953699	9963247	9966200	9976052	A004537	A011845	A017150	A025857
9953892	9963249	9966203	9976055	A004542	A011972	A017154	A025862
9955091	9963378	9966205	9978598	A004543	A012088	A018239	A026331
9955406	9963380	9966550	A000215	A004544	A012341	A018301	A026336
9955408	9963382	9966562	A000219	A004614	A012591	A018307	A026337
9955530	9963386	9966563	A000220	A004616	A012691	A018310	A026375
9955561	9963622	9966674	A000369	A005866	A012738	A018721	A027449
9955726	9963627	9966674	A000453	A005870	A012743	A018722	A027805
9955727	9963633	9966676	A002304	A005980	A012755	A018741	A027862
9955845	9963851	9966934	A002305	A006045	A013058	A018885	A027922
9956028	9964145	9966945	A002307	A006222	A013211	A019293	A028081
9956396	9964148	9966946	A002404	A006225	A013685	A019401	A028468
9956460	9964157	9971999	A002406	A006229	A013686	A019539	A028566
9957242	9964161	9972000	A002408	A006237	A013688	A019550	A028760
9957512	9964259	9972002	A002595	A006829	A013689	A019567	A028777
9957684	9964494	9972610	A002606	A006831	A013796	A019573	A028780
9957685	9964496	9972613	A002616	A006834	A013969	A019834	A028787
9957927	9964498	9972719	A002898	A007148	A014018	A019912	A029143
9957947	9964756	9972720	A003029	A007150	A014061	A020026	A029145
9957949	9964810	9973291	A003038	A007304	A014063	A020056	A029146
9958187	9964826	9973322	A003039	A007344	A014351	A020244	A029318
9958345	9965055	9973329	A003044	A007345	A014352	A020313	A029460
9958346	9965062	9973567	A003164	A007349	A014359	A020315	A029508
9958393	9965067	9973570	A003291	A007354	A014450	A020457	A029610
9958740	9965070	9973597	A003296	A007462	A014753	A020464	A029613
9958741	9965483	9973779	A003311	A007463	A014768	A020564	A029617
9958750	9965488	9973880	A003316	A008617	A014795	A020569	A029794
9959038	9965489	9974008	A003603	A008767	A014840	A020699	A029862
9959490	9965768	9974013	A003692	A008886	A015227	A021233	A029863
9959502	9965771	9974068	A003703	A009203	A015228	A021384	A029865
9959520	9965805	9974070	A003706	A009221	A015332	A021386	A029866
9959884	9965946	9974071	A003710	A010851	A015362	A021416	A031039
9959885	9965947	9974165	A003805	A010947	A015541	A021885	A031102
9960141	9965950	9974226	A004047	A010948	A015559	A021887	A031106
9960147	9965956	9974230	A004085	A010956	A015563	A021913	A031385
9960150	9965958	9974247	A004116	A010990	A015570	A022134	A031387
9960720	9965962	9974250	A004259	A011059	A015957	A022404	A031391
9960761	9965963	9974426	A004318	A011267	A016217	A022406	A031649
9960776	9965965	9974433	A004321	A011295	A016477	A024199	A031663

Appendix 1 – List of Affected Lots (cont.):

Marathon

A031670	A040153	A049233	A058253	A067173	A075715	A081223	A091333
A031677	A040155	A049236	A059305	A067407	A075721	A081718	A091341
A031963	A040163	A049237	A059307	A067409	A076320	A081725	A091589
A032065	A040168	A050083	A059311	A067412	A076402	A081859	A091708
A032068	A040705	A050084	A059577	A068091	A076406	A081893	A091710
A032070	A040718	A050119	A059691	A068094	A076410	A081895	A091713
A032072	A040725	A050263	A059727	A068098	A076412	A082393	A093603
A032398	A041039	A050356	A059731	A068222	A077046	A083144	A093607
A032835	A041157	A050361	A059733	A068382	A077048	A083166	A093787
A032837	A041166	A050399	A061328	A068419	A077049	A083203	A093789
A032840	A041283	A050404	A061329	A068436	A077557	A083206	A093795
A033206	A041359	A050744	A062218	A068563	A077565	A083545	A093797
A033229	A041404	A050745	A062266	A068565	A077607	A083546	A094057
A033230	A041456	A050756	A062267	A068773	A077650	A085047	A094083
A033233	A041488	A050855	A062268	A068835	A077664	A085059	A094085
A034330	A041495	A050863	A062269	A068839	A078031	A085063	A094100
A034351	A041754	A051742	A062987	A068842	A078033	A085066	A094794
A034352	A041861	A051875	A062991	A068850	A078870	A085070	A094796
A034353	A041866	A051939	A063168	A068966	A078879	A085073	A094801
A034354	A041867	A051940	A063235	A068981	A078950	A085329	A095119
A034509	A042017	A051941	A063275	A072687	A078962	A086510	A095305
A034511	A043949	A051943	A063399	A072735	A078963	A086658	A095308
A034512	A043960	A051958	A063400	A072738	A079311	A087223	A095397
A034517	A043964	A052051	A063960	A072740	A079320	A087228	A096868
A034568	A046465	A052195	A063968	A072810	A079321	A087231	A096899
A034585	A046470	A052217	A063976	A073361	A079324	A087234	A097061
A034586	A046485	A052227	A064156	A073364	A079334	A088753	A097090
A034587	A046487	A052237	A064238	A073371	A079456	A088756	A097095
A037007	A046488	A052246	A064250	A073382	A079606	A089280	A097098
A037008	A046610	A052542	A064251	A073384	A080017	A089296	A097100
A037009	A046611	A052828	A064252	A074651	A080018	A089334	A097427
A037010	A046735	A052873	A064254	A075001	A080019	A089673	A097429
A037014	A046739	A052874	A064373	A075030	A080021	A089675	A097430
A037410	A046868	A052876	A064908	A075034	A080022	A089837	A097433
A037415	A046873	A053105	A064910	A075040	A080241	A090423	A097435
A037416	A047032	A053180	A064912	A075046	A080269	A090577	A097436
A037513	A048479	A053183	A064999	A075177	A080272	A090598	A097437
A038097	A048480	A053452	A065056	A075183	A080276	A090712	A098264
A038103	A048487	A053459	A065597	A075298	A080387	A091142	A100171
A038785	A048645	A058121	A065598	A075361	A080592	A091143	A100172
A038794	A048646	A058149	A065602	A075362	A080600	A091145	A100240
A038802	A048808	A058152	A066531	A075366	A080613	A091147	A100243
A039450	A048815	A058154	A066948	A075619	A080712	A091224	A100251
A039461	A048821	A058157	A066965	A075697	A081143	A091312	A100584
A040138	A048825	A058160	A066966	A075709	A081213	A091331	A100588

Appendix 1 – List of Affected Lots (cont.):

Marathon

A100593	A117561	A123429	A135805	A142934	A156901	A163572	A179306
A100895	A117826	A123434	A136017	A142937	A156904	A164739	A181904
A101575	A117907	A123805	A136086	A143067	A156909	A164745	A181919
A101628	A117969	A123807	A136448	A143074	A157053	A164825	A181921
A101634	A117971	A124180	A136451	A144556	A157135	A164903	A181952
A101638	A118192	A124276	A136456	A144585	A157322	A164981	A182269
A101641	A118425	A124277	A136580	A144591	A157325	A164989	A182274
A101941	A118687	A124464	A136583	A144597	A157367	A165291	A182284
A101944	A118689	A124476	A137622	A144736	A157895	A165294	A182289
A101950	A118708	A125090	A137641	A144873	A157896	A165297	A183791
A102234	A118842	A126669	A137643	A145017	A157898	A165563	A183801
A102251	A118850	A126670	A137645	A145023	A158504	A165564	A183802
A102287	A118854	A127819	A137757	A145031	A158512	A165565	A183928
A102289	A119044	A127937	A137803	A145039	A158557	A165932	A184118
A102293	A119046	A127949	A137873	A145165	A159131	A165935	A184221
A102646	A119144	A128031	A137879	A153798	A159132	A165946	A184222
A102685	A119346	A128212	A137885	A153890	A159133	A165985	A184357
A102686	A119879	A128333	A137886	A153892	A159134	A165986	A185869
A103326	A119891	A128942	A137888	A153894	A159337	A172223	A185871
A103395	A120131	A128953	A138208	A153896	A159338	A172229	A185872
A103398	A120204	A129090	A138219	A154486	A159813	A172232	A185875
A103401	A120240	A129091	A138220	A154498	A159830	A172236	A186101
A103406	A120283	A129327	A138225	A154505	A159844	A174414	A186135
A103410	A120400	A129335	A138231	A154514	A159848	A174416	A186140
A103412	A120413	A129898	A138490	A154522	A159851	A175520	A186333
A104092	A121145	A130244	A138906	A154633	A160096	A175556	A186481
A104114	A121151	A130945	A138908	A154634	A161031	A175695	A189181
A104120	A121155	A130985	A138921	A154813	A161039	A175710	A189190
A105499	A121156	A131044	A138923	A154816	A161042	A175715	A189192
A105549	A121772	A131047	A138926	A155028	A161051	A175935	A189194
A105553	A121779	A131939	A139599	A155131	A161142	A176084	A189815
A105757	A121785	A131943	A139607	A155133	A161279	A176157	A189816
A105880	A122401	A131966	A139610	A155140	A161288	A176158	A189818
A105883	A122468	A132855	A139615	A155147	A161289	A176216	A189825
A105887	A122621	A133983	A139617	A155157	A161292	A178717	A190048
A105895	A122651	A134515	A139993	A155434	A161463	A178718	A190117
A105913	A122655	A134519	A140160	A155668	A162951	A178726	A190118
A106016	A122656	A134630	A140235	A156090	A162954	A178883	A190283
A106026	A123032	A135227	A140561	A156286	A163149	A178912	A190446
A116756	A123037	A135229	A140671	A156380	A163157	A178918	A190453
A116759	A123039	A135232	A142410	A156402	A163168	A178921	A190455
A117473	A123202	A135380	A142419	A156409	A163353	A179058	A190456
A117485	A123205	A135786	A142422	A156415	A163499	A179077	A191757
A117494	A123211	A135789	A142552	A156419	A163513	A179079	A191759
A117516	A123427	A135800	A142556	A156421	A163526	A179304	A191767

Appendix 1 – List of Affected Lots (cont.):

Marathon

A191864	A213290	A228170	A235536	A246111	A253655	A264501	A275676
A200836	A214210	A228172	A235538	A246120	A253663	A264962	A275677
A200884	A214214	A228208	A236147	A246128	A253664	A264965	A275678
A201108	A214215	A228603	A236466	A246131	A255799	A264967	A275679
A201216	A214558	A228604	A237340	A246255	A255814	A264971	A276176
A201220	A214561	A228859	A237344	A246256	A255821	A265128	A276177
A201224	A214564	A228862	A237349	A246570	A255823	A265129	A276284
A201331	A214898	A230195	A237353	A246571	A255827	A265691	A276285
A201490	A214905	A230269	A237466	A246633	A255985	A265699	A277023
A201510	A214908	A230847	A237467	A246634	A255986	A265706	A277024
A201532	A214909	A230930	A237561	A247155	A256542	A265773	A277025
A201542	A214965	A230931	A237564	A247267	A256548	A265775	A277172
A201907	A215131	A230984	A237642	A247287	A256550	A265873	A278439
A203970	A215149	A230985	A237645	A248003	A256555	A267430	A278441
A203975	A215156	A230987	A237647	A248009	A256742	A267432	A278443
A203985	A215164	A230991	A238041	A248020	A256819	A267525	A278627
A203990	A215167	A231005	A238047	A248023	A256822	A267534	A279421
A204483	A215538	A231643	A239537	A248028	A256823	A267535	A279423
A204486	A215602	A231645	A239539	A248215	A256826	A267536	A279439
A204493	A215651	A231649	A239540	A248743	A256827	A267537	A279729
A204871	A215652	A231652	A239550	A248746	A257352	A268139	A279730
A204875	A222349	A231656	A239551	A248753	A257358	A268176	A280700
A204876	A222570	A232346	A240033	A248759	A257360	A268243	A280706
A205105	A222647	A232352	A240038	A248894	A257363	A268246	A280709
A209299	A222648	A232359	A240043	A248896	A257368	A268247	A280855
A209303	A223013	A232363	A240047	A249006	A258088	A268253	A280856
A209304	A223015	A232365	A240159	A249008	A258096	A268254	A280931
A209330	A223124	A232587	A240161	A249010	A258107	A268742	A281003
A209430	A223125	A232588	A243270	A249406	A258255	A272372	A281667
A209707	A223150	A232598	A243277	A249648	A260041	A272384	A282592
A209710	A224545	A233765	A243282	A249650	A260047	A272390	A282622
A209714	A224546	A233766	A243284	A249714	A260055	A272537	A282625
A209717	A224547	A233767	A243415	A249715	A260216	A272612	A282647
A210619	A225374	A233768	A243417	A249953	A260768	A272614	A282650
A210621	A225445	A234758	A243529	A250226	A260776	A272624	A282821
A210623	A225446	A234759	A243532	A250259	A260782	A274619	A282834
A210624	A225449	A234760	A243533	A250825	A260786	A274627	A282840
A210897	A225654	A234761	A243535	A250851	A260946	A274630	A282841
A210906	A227264	A235062	A243994	A250852	A262262	A274641	A282843
A210918	A227265	A235063	A243997	A250856	A262287	A274813	A283310
A212857	A227684	A235210	A244001	A251155	A262288	A275504	A283606
A212865	A227685	A235215	A244012	A253403	A262289	A275505	A283641
A213284	A227714	A235220	A244099	A253411	A262382	A275506	A283642
A213286	A227717	A235221	A244101	A253414	A262419	A275507	A283643
A213288	A227794	A235532	A244139	A253427	A264500	A275666	A283644

Appendix 1 – List of Affected Lots (cont.):

Marathon

A284103	A290773	A297537	A303119	A309234	A316750	A323580	A329235
A284108	A290791	A297612	A303136	A309394	A316751	A323582	A329356
A284109	A291365	A297616	A303137	A309522	A316752	A324273	A329358
A284113	A291405	A297618	A303138	A310667	A317330	A324274	A329359
A284119	A291406	A299056	A303140	A310668	A317394	A325147	A329761
A284907	A291407	A299057	A303142	A310803	A317412	A325150	A329765
A284920	A291567	A299059	A303388	A310809	A317437	A325360	A329891
A284927	A292912	A299329	A303393	A312420	A318015	A325361	A329892
A284934	A292918	A299330	A303411	A312421	A318019	A325772	A329935
A284939	A292921	A299331	A303412	A312461	A318187	A325773	A330004
A284944	A293023	A299593	A303413	A314037	A320445	A325949	A330008
A284948	A293329	A300051	A303688	A314041	A320448	A326557	A330010
A286589	A293332	A300060	A303840	A314042	A321077	A326564	A330011
A286591	A293333	A300237	A303853	A314043	A321079	A326565	A335114
A286649	A293540	A300319	A303859	A314044	A321202	A326723	A335138
A288843	A293541	A300569	A305247	A314477	A321291	A326787	A335146
A288845	A293856	A300570	A305283	A314479	A321294	A326805	A335151
A288854	A293858	A300595	A305285	A314492	A321295	A326806	A335171
A288946	A293861	A300982	A305287	A315198	A321297	A326808	A335537
A288960	A295763	A300994	A305291	A315204	A322497	A326809	A335543
A288961	A296040	A300995	A305950	A315207	A322498	A327070	A336022
A288965	A296042	A300996	A305960	A315210	A322504	A328301	A336359
A288966	A296043	A301001	A305962	A315214	A322669	A328310	A336361
A290567	A296134	A301556	A305966	A316607	A323235	A328319	A336692
A290575	A296135	A301558	A305969	A316610	A323237	A328442	A336700
A290582	A296701	A301560	A309144	A316612	A323241	A328474	A336702
A290599	A297361	A301563	A309230	A316746	A323579	A329194	A336704
9945407	9945399	9946910	9946737	9948666	9949668	9949448	9948915
9945725	9946271	9945991	9948660				

Appendix 1 – List of Affected Lots (cont.):

Ultraflow Lots

9963219	A006885	A028706	A046108	A092226	A125138	A197490	A217465
9972750	A006886	A032436	A046114	A092708	A125142	A197771	A218093
9972981	A017413	A032439	A049550	A092709	A125735	A198001	A218094
9972987	A017429	A033725	A049557	A092735	A125819	A198024	A218095
9973009	A017430	A033841	A060565	A092741	A125820	A198044	A221276
9976672	A017432	A033856	A060578	A092742	A125943	A198049	A221279
9976687	A020842	A033858	A072386	A119501	A126381	A205358	A221280
9976691	A020853	A033861	A072396	A119534	A126394	A205362	A221281
9977163	A020921	A034550	A081121	A119537	A126396	A211498	A229624
9977489	A021063	A036785	A081129	A119547	A126424	A211509	A229671
A001785	A021473	A036788	A081140	A119552	A127015	A211511	A265947
A001914	A025088	A036790	A091993	A119802	A127019	A217418	A265948
A001919	A025090	A044657	A091996	A120128	A127211	A217448	A266807
A001959	A025123	A044663	A091999	A120130	A153128	A217449	A266811
A002231	A025129	A044668	A092004	A124985	A197485	A217460	A280246
A003190	A028569						

Appendix 1 – List of Affected Lots (cont.):

Xcelerator Lots

9530849	9709603	9797735	A002136	A043280	A123449	A183571
9530851	9743324	9802698	A002460	A043281	A123450	A185838
9531882	9743328	9802700	A002610	A043969	A123452	A185841
9567732	9746399	9848318	A004560	A044251	A134837	A194324
9585211	9746403	9882269	A004561	A044252	A136164	A194327
9585247	9750614	9931088	A005008	A097470	A136165	A195347
9585249	9750615	9931089	A005011	A097472	A142394	A200386
9585282	9750618	9931093	A005477	A097473	A142395	A200387
9585283	9750619	9979485	A005478	A097474	A145004	A202560
9594629	9750621	9979495	A006080	A097475	A136161	A202563
9601922	9750624	9979497	A006171	A097476	A145006	A211246
9601924	9751548	9979498	A006176	A097477	A145149	A211248
9612762	9767048	9979500	A006182	A114384	A145150	A217138
9612763	9794407	201401398X	A006887	A114385	A145152	A217141
9628372	9794409	201601079E	A008199	A114389	A145153	A217142
9641711	9794795	201602502E	A008201	A114391	A150759	A221040
9643550	9794924	201603429E	A014137	A114393	A150762	A221041
9660868	9794925	201603430E	A014297	A114543	A150764	A221042
9666144	9795119	201605381E	A014819	A114545	A151577	A221043
9666146	9795120	201605382E	A014821	A114546	A152854	A237364
9666148	9795516	201605383E	A014825	A115726	A152856	A237366
9666835	9795517	201605384E	A019031	A116267	A159822	A237368
9667132	9796816	201606174E	A019337	A116270	A159826	A243986
9690683	9796817	201606178E	A035856	A116271	A165319	A243987
9690686	9796817	201606179E	A035858	A116272	A165320	A243988
9690687	9796818	201606180E	A037439	A116274	A165321	A247255
9701004	9797467	201606297E	A040174	A123439	A183549	A247256
9702147	9797469	201606298E	A040176	A123440	A183550	A247259
9706211	9797470	A002129	A040177	A123442	A183561	A247270
9709599	9797734	A002131	A040222	A123445	A183563	A258089

Medtronic	Field Action Plan	CF2120 Rev. A ECO-NV19829
------------------	--------------------------	--

Appendix 3: Customer Confirmation Certificate

Customer Confirmation Certificate Medtronic PTFE Delamination – October 2016

Account Number:
Account Name:
Address:
City, State, Zip
Sales Representative: *Rep Name*
Representative Phone: *Rep phone number*

By signing this form I confirm that I have read and understand the Urgent Medical Device Recall notice from Medtronic regarding the Medtronic PTFE delamination issue and that the information provided below is accurate.

Please mark the appropriate box that applies (select one only):

☐ I have quarantined all affected product in my inventory and will be returning it to Medtronic.

Contact Medtronic Customer Service at 1-800-716-6700 to initiate a return and credit of product referencing FCA-2016-1 for the following product from affected lots:

Pipeline™ embolization device
Alligator™ retrieval device
Marathon™ microcatheter
UltraFlow™ HPC microcatheter
Xcelerator™ hydrophilic guidewire

Product Return Address:

Medtronic
Attention: Field Returns
60 Middletown Ave
North Haven CT, 06473

☐ I have no affected product currently within my inventory.

Name (Print): _____

Title (Print): _____

Signature: _____ Date: _____

Telephone: _____

Please fax this completed form to 1-949-434-5020. For questions, contact your Medtronic Sales Representative.