



49 Changi South Avenue 2
Singapore 486056
www.medtronic.com

tel +65 6436 5000
fax +65 6776 6355

URGENT MEDICAL DEVICE RECALL
Covidien Endo GIA™ Articulating Loading Units

26 April 2018

Attention: Risk Management Director and O.R. Materials Management
CC: The Chairman Medical Board and relevant Head of Departments

Dear Valued Customer:

The purpose of this letter is to advise you that Medtronic is voluntarily recalling specific production lots of its

Covidien Endo GIA™ Articulating Loading Units.

Issue Description:

This voluntary recall is being conducted due to the potential for a device to be missing the sled component. The sled component is responsible for staple deployment. Without the sled, no staples would be deployed. Use of a product with this issue may result in failure to form a staple line when tissue is divided, leading to bleeding or leakage of luminal contents. Medtronic has received four reports of injury related to this issue.

This voluntary recall affects only the item codes and lots listed on Attachment A.

Medtronic requests that you immediately discontinue use of the affected products detailed on Attachment A and Attachment B (FOR USE IN EU ONLY). Please quarantine and return any unused products of the listed item codes and lots. Unused products from the affected item codes and lots should be returned as described in the Required Actions section below. If you have distributed Endo GIA™ devices listed on Attachment A and Attachment B (FOR USE IN EU ONLY) please promptly forward the information from this letter to those recipients. All unused products from the affected item codes and lots must be returned.

Required Actions:

1. Please immediately discontinue use, and quarantine, the affected item codes and lots listed on Attachment A and Attachment B (FOR USE IN EU ONLY).
2. Please return affected product as indicated below.
3. Complete the Customer confirmation form **even if you do not have inventory.**

- **CUSTOMERS WHO PURCHASED PRODUCT DIRECTLY FROM MEDTRONIC**

Please complete the Covidien Endo GIA™ Articulating Loading Units Customer Confirmation Form (attached) and return to your local Medtronic representative together with affected units if any.

- **CUSTOMERS WHO PURCHASED PRODUCT FROM A DISTRIBUTOR**

Please complete the Covidien Endo GIA™ Articulating Loading Units Customer Confirmation Form (attached) and contact your Distributor directly. All affected product and Customer Confirmation Form must be returned through the Distributor.

This field action is being issued or will be notified to relevant regulatory bodies according to applicable regulations. Please communicate this important information within your facility and or other facilities as required. We request that you contact Medtronic if you experienced quality problems or adverse events.

If you have any questions or concerns regarding this recall, please do not hesitate to contact your local Medtronic representative.

We appreciate your attention to this matter and apologize for any inconvenience this issue may have caused.

Sincerely,



Diana Teo
Quality System Manager
Medtronic

Attachment A

Item Code	Item Description	Affected Lots		
030450	Endo GIA™ Auto Suture™ Universal Articulating Loading Unit 30mm- 2.0mm	N7L0306KX	N8B0724KX	N8C1303KRX
		N7M0221KX	N8C0958KX	
030451	Endo GIA™ Auto Suture™ Universal Articulating Loading Unit 30mm - 2.5mm	N7L0071KX	N7M0464KX	N8A1014KX
		N7L0527KX	N7M0786KX	N8B0304KX
		N7L0792KX	N8A0154KX	N8C0830KX
		N7L0914KX	N8A0553KX	N8C0847KRX
		N7L1014KX	N8A0850KX	N8C0848KRX
		N7L1159KX		
030452	Endo GIA™ Auto Suture™ Universal Articulating Loading Unit 30mm - 3.5mm	N7L0191KX	N8B0234KX	N8C0147KX
		N7M0944KX		
030454	Endo GIA™ Auto Suture™ Universal Articulating Loading Unit 45mm - 2.5mm	N7K1146KX	N7M0509KX	N8A0434KX
		N7L0355KX	N7M0520KX	N8A0435KX
		N7L1119KX	N7M0521KX	N8B0487KX
		N7M0088KX	N8A0312KRX	N8C0026KX
		N7M0343KX	N8A0360KX	N8C0119KX
		N7M0477KRX	N8A0433KX	N8C0630KX
030455	Endo GIA™ Auto Suture™ Universal Articulating Loading Unit 45mm- 3.5mm	N7F0576KX	N7M0383KX	N8A0324KX
		N7L0032KX	N7M0584KX	N8A0477KRX
		N7L0239KX	N7M0696KX	N8A0551KX
		N7L0519KX	N8A0036KX	N8C0357KX
		N7M0104KX		
030457	Endo GIA™ Auto Suture™ Universal Articulating Loading Unit 60mm -2.5mm	N7D0492KRX	N7L0699KX	N8A0213KX
		N7F0518KX	N7L0732KX	N8B0002KX
		N7F0651KX	N7L1043KX	N8C0429KX
		N7L0270KX		
030458	Endo GIA™ Auto Suture™ Universal Articulating Loading Unit 60mm - 3.5mm	N7F0631KX	N7M0646KX	N8B0099KX
		N7L0117KX	N7M0740KX	N8C0008KX
		N7L0397KX	N7M0972KX	N8C0205KX
		N7L0465KX	N8A0172KX	N8C0532KX
		N7L0768KX	N8A0539KX	N8C0571KX
		N7L0911KX	N8A0679KX	N8C0758KX
		N7L1075KX	N8A0910KX	N8C0795KX
		N7M0067KX	N8A0961KX	N8C1284KRX
		N7M0314KX	N8B0054KX	
EGIA30CTAV	Endo GIA™ 30mm Gray Curved Tip Articulating Vascular Reload with Tri-Staple™ Technology	N7G0756KX		
EGIA30AV	Endo GIA™ 30mm Gray Articulating Vascular Reload with Tri-Staple™ Technology	N7G0757KX		
EGIA45CTAV	Endo GIA™ 45mm Gray Curved Tip Articulating Vascular Reload with Tri-Staple™ Technology	N7G0796KX		

Attachment B

For use in EU countries where Kits were distributed

Item Code	Item Description	Affected Lots
KIT1095V	BOX KIT1095V KIT OBESITAS	0214062622
01KBAR006116V PROC BYPASS	BOX 01KBAR006116V PROC BYPASS 01KBAR0061	0214710711
BOX1073076GASB0	BOX BOX1073076GASB0 BBR STVEIT GASTRIC B	0214043771
BOX1073076GASB0	BOX BOX1073076GASB0 BBR STVEIT GASTRIC B	0214232317
BOX1073076GASB0	BOX BOX1073076GASB0 BBR STVEIT GASTRIC B	0214375093
BOX1073076GASB0 BBR STVEIT GASTRIC B	BOX BOX1073076GASB0 BBR STVEIT GASTRIC B	0215113531
BOX1073076GASB0 BBR STVEIT GASTRIC B	BOX BOX1073076GASB0 BBR STVEIT GASTRIC B	0215132607
BOX1073076GASB1 BBR GASTRIC BYPASS S	BOX BOX1073076GASB1 BBR GASTRIC BYPASS S	0215227169
BOX1073180	BOX BOX1073180 APPE1 APPENDIX KIT ZAMS	0214050554
BOX1073180APPE1 APPENDIX KIT ZAMS	BOX BOX1073180APPE1 APPENDIX KIT ZAMS	0214907029
BOX1073180APPE1 APPENDIX KIT ZAMS	BOX BOX1073180APPE1 APPENDIX KIT ZAMS	0214991119
BOX1073180APPE1 APPENDIX KIT ZAMS	BOX BOX1073180APPE1 APPENDIX KIT ZAMS	0215130754
BOX1073180APPE1 APPENDIX KIT ZAMS	BOX BOX1073180APPE1 APPENDIX KIT ZAMS	0215158268
HEMIGX1 ZESTAW DO HEMICOLECTOMII	BOX HEMIGX1 ZESTAW DO HEMICOLECTOMII	0214909888
HEMIGX1 ZESTAW DO HEMICOLECTOMII	BOX HEMIGX1 ZESTAW DO HEMICOLECTOMII	0214916144
KIT00486R	KIT00486R OBESIDAD D. WORWARLD 2 X1	0213984473
KIT00486R	KIT00486R OBESIDAD D. WORWARLD 2 X1	0213583578
KIT00486R	KIT00486R OBESIDAD D. WORWARLD 2 X1	0213596057
KIT00486R	KIT00486R OBESIDAD D. WORWARLD 2 X1	0214068986
KIT00486R	KIT00486R OBESIDAD D. WORWARLD 2 X1	0213589227
KIT00713	BOX KIT00713 KIT BY PASS	0214364763
KIT00713	BOX KIT00713 KIT BY PASS	0214456590
KIT00713	BOX KIT00713 KIT BY PASS	0214435328
KIT00713	BOX KIT00713 KIT BY PASS	0214983139
KIT00713 KIT BY PASS	BOX KIT00713 KIT BY PASS	0215042012
KIT00713 KIT BY PASS	BOX KIT00713 KIT BY PASS	0215234651
KIT00713 KIT BY PASS	BOX KIT00713 KIT BY PASS	0215288296
KIT00714 KIT SLEEVE	BOX KIT00714 KIT SLEEVE	0214915642
KIT00714 KIT SLEEVE	BOX KIT00714 KIT SLEEVE	0214983138
KIT00714 KIT SLEEVE	BOX KIT00714 KIT SLEEVE	0215159225
KIT00714 KIT SLEEVE	BOX KIT00714 KIT SLEEVE	0215234652
KIT1171C	BOX KIT1171C KIT LAP CYSTECTOMY	0214090852
KIT1171C KIT LAP CYSTECTOMY	BOX KIT1171C KIT LAP CYSTECTOMY	0214522213
KIT1171C KIT LAP CYSTECTOMY	BOX KIT1171C KIT LAP CYSTECTOMY	0215009132
KIT1171C KIT LAP CYSTECTOMY	BOX KIT1171C KIT LAP CYSTECTOMY	0215153445
KIT4184 APEN LOBECTOMY KIT X0	KIT4184 APEN LOBECTOMY KIT X1	0214962218
KITBAR021- BARIATRIC KIT X0	KITBAR021- BARIATRIC KIT X1	0215065486
KIT-DE-0080	KIT-DE-0080 BYPASS KIT LOHNE X1	0214049646
KIT-DE-0080	KIT-DE-0080 BYPASS KIT LOHNE X1	0214370002
KIT-DE-0080	KIT-DE-0080 BYPASS KIT LOHNE X1	0214470664
KIT-DE-0080 BYPASS KIT LOHNE X0	KIT-DE-0080 BYPASS KIT LOHNE X1	0215208426
KIT-DE-0083	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0214798114
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0214618312
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0214795633
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0215002022
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0215010930
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0215010931

Attachment B

For use in EU countries where Kits were distributed

Item Code	Item Description	Affected Lots
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0215010932
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0215132611
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0215132612
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0215132613
KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	KIT-DE-0083 APPENEKTOMIE KIT NAUENX1	0215264173
KITDE0108	BOX KITDE0108 BYPASS KIT LUEBECK	0214894171
KITDE0143	BOX KITDE0143 BYPASS BOX	0214463634
KITDE0143 BYPASS BOX	BOX KITDE0143 BYPASS BOX	0215280952
KITDE0154	BOX KITDE0154 BYPASS KH NORDWEST	0214013887
KITDE0154	BOX KITDE0154 BYPASS KH NORDWEST	0214144177
KITDE0154	BOX KITDE0154 BYPASS KH NORDWEST	0214260194
KITDE0154	BOX KITDE0154 BYPASS KH NORDWEST	0214470767
KITDE0154	BOX KITDE0154 BYPASS KH NORDWEST	0214956044
KITDE0154 BYPASS KH NORDWEST	BOX KITDE0154 BYPASS KH NORDWEST	0214620382
KITDE0154 BYPASS KH NORDWEST	BOX KITDE0154 BYPASS KH NORDWEST	0214744957
KITDE0154 BYPASS KH NORDWEST	BOX KITDE0154 BYPASS KH NORDWEST	0215079103
KITDE0154 BYPASS KH NORDWEST	BOX KITDE0154 BYPASS KH NORDWEST	0215151206
KITDE0183	BOX KITDE0183 APPENDIX LIPPE DETMOLD	0214014173
KITDE0183 APPENDIX LIPPE DETMOLD	BOX KITDE0183 APPENDIX LIPPE DETMOLD	0214573291
KITDE0183 APPENDIX LIPPE DETMOLD	BOX KITDE0183 APPENDIX LIPPE DETMOLD	0214621781
KITDE0183 APPENDIX LIPPE DETMOLD	BOX KITDE0183 APPENDIX LIPPE DETMOLD	0214909631
KITDE0183 APPENDIX LIPPE DETMOLD	BOX KITDE0183 APPENDIX LIPPE DETMOLD	0215068744
KITDE0183 APPENDIX LIPPE DETMOLD	BOX KITDE0183 APPENDIX LIPPE DETMOLD	0215068743
KITDE0183 APPENDIX LIPPE DETMOLD	BOX KITDE0183 APPENDIX LIPPE DETMOLD	0215108744
KITDE0183 APPENDIX LIPPE DETMOLD	BOX KITDE0183 APPENDIX LIPPE DETMOLD	0215108743
KITDE0183 APPENDIX LIPPE DETMOLD	BOX KITDE0183 APPENDIX LIPPE DETMOLD	0215276369
LAPBARLECZNA ZESTAW BARIATRYCZNY LECZNA.	LAPBARLECZNA ZESTAW BARIATRYCZNY LECZNA.	0215274505
LAPGAR	BOX LAPGAR ZESTAW DO GB35	0214049704
LAPGAR	BOX LAPGAR ZESTAW DO GB35	0214529057
LAPGAR ZESTAW DO GB34	BOX LAPGAR ZESTAW DO GB35	0215176564
LAPIMBGP1	BOX LAPIMBGP1 ZESTAW DO MGBP	0214930629
LAPIMBGP1	BOX LAPIMBGP1 ZESTAW DO MGBP	0214993038
LAPIMBGP1	BOX LAPIMBGP1 ZESTAW DO MGBP	0215068848
LAPIMBGP1	BOX LAPIMBGP1 ZESTAW DO MGBP	0215124286
LAPIMBGP1	BOX LAPIMBGP1 ZESTAW DO MGBP	0215124285
LAPIMBGP1 ZESTAW DO MGBP	BOX LAPIMBGP1 ZESTAW DO MGBP	0214793818
LAPIMBGP1 ZESTAW DO MGBP	BOX LAPIMBGP1 ZESTAW DO MGBP	0215130759
LAPIMBGP1 ZESTAW DO MGBP	BOX LAPIMBGP1 ZESTAW DO MGBP	0215243237
LAPISG1	BOX LAPISG1 ZESTAW DO SG	0214993037
LAPISG1 ZESTAW DO SG	BOX LAPISG1 ZESTAW DO SG	0214793817
LAPISG1 ZESTAW DO SG	BOX LAPISG1 ZESTAW DO SG	0214930628
LAPISG1 ZESTAW DO SG	BOX LAPISG1 ZESTAW DO SG	0214930630
LAPISG1 ZESTAW DO SG	BOX LAPISG1 ZESTAW DO SG	0215068849
LAPISG1 ZESTAW DO SG	BOX LAPISG1 ZESTAW DO SG	0215124284
LAPISG1 ZESTAW DO SG	BOX LAPISG1 ZESTAW DO SG	0215243236
LAPISG1 ZESTAW DO SG	BOX LAPISG1 ZESTAW DO SG	0214925902
LAPISG1 ZESTAW DO SG	BOX LAPISG1 ZESTAW DO SG	0214925904

Attachment B

For use in EU countries where Kits were distributed


Item Code	Item Description	Affected Lots
LAPISG1 ZESTAW DO SG	BOX LAPISG1 ZESTAW DO SG	0214925903
LAPISGW ZESTAW DO SG	BOX LAPISGW ZESTAW DO SG	0215120700
LAPPH ZASTAW DO PRAWEJ HEMIKOLEKTOMI	BOX LAPPH ZASTAW DO PRAWEJ HEMIKOLEKTOMI	0215015272
LAPPH ZASTAW DO PRAWEJ HEMIKOLEKTOMI	BOX LAPPH ZASTAW DO PRAWEJ HEMIKOLEKTOMI	0215298205
LAPZD ZESTAW DO SG	BOX LAPZD ZESTAW DO SG	0215082707
LAPZD ZESTAW DO SG	BOX LAPZD ZESTAW DO SG	0215176563
MGBGLOGX1	BOX MGBGLOGX1 ZESTAW MINI GASTRIC BYPASS	0215124172
MGBGLOGX1 ZESTAW MINI GASTRIC BYPASS	BOX MGBGLOGX1 ZESTAW MINI GASTRIC BYPASS	0215253056
MGBX ZESTAW GB LS X1	MGBX ZESTAW GB LS X1	0215017845
MGBX ZESTAW GB LS X1	MGBX ZESTAW GB LS X1	0215002541
MGBX ZESTAW GB LS X1	MGBX ZESTAW GB LS X1	0215104507
PST01199	BOX PST01199 BARIATRIC TRISTAPLE	0214109631
PST01199 BARIATRIC TRISTAPLE	BOX PST01199 BARIATRIC TRISTAPLE	0214909770
PST03403	BOX PST03403 ALTERNATIVE APPENDIX KIT	0214975597
PST03403	BOX PST03533 KIT APP STAPLER	0214860061
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215036240
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215036239
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215036238
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215018854
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215018855
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215102248
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215102247
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215068747
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215068748
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215050427
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215102246
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215050429
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215050428
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215068746
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215117409
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215117411
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215117410
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215125963
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215125962
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215125964
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215275313
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215275315
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215275314
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215282505
PST03533 KIT APP STAPLER	BOX PST03533 KIT APP STAPLER	0215282506
PST03627 BYPASS AIWA CLINIC	BOX PST03627 BYPASS AIWA CLINIC	0214701832
PST03627 BYPASS AIWA CLINIC	BOX PST03627 BYPASS AIWA CLINIC	0214993796
PST03627 BYPASS AIWA CLINIC	BOX PST03627 BYPASS AIWA CLINIC	0214999435
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215050309
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215108007
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215050310
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215050308

Attachment B

For use in EU countries where Kits were distributed

Item Code	Item Description	Affected Lots
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215050311
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215050312
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215050313
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215037144
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215037143
PST03699 KIT APP MED ALTERNATIV STAP	BOX PST03699 KIT APP MED ALTERNATIV STAP	0215108008
PST04089 SLEEVE AIWA CLINIC	BOX PST04089 SLEEVE AIWA CLINIC	0215237526
PST04089 SLEEVE AIWA CLINIC	BOX PST04089 SLEEVE AIWA CLINIC	0214991351
PST04089 SLEEVE AIWA CLINIC	BOX PST04089 SLEEVE AIWA CLINIC	0214991350
PST04089 SLEEVE AIWA CLINIC	BOX PST04089 SLEEVE AIWA CLINIC	0215157998
PST04089 SLEEVE AIWA CLINIC	BOX PST04089 SLEEVE AIWA CLINIC	0215157997
PST04090 BYPASS AIWA CLINIC	BOX PST04090 BYPASS AIWA CLINIC	0215237527
PST04090 BYPASS AIWA CLINIC	BOX PST04090 BYPASS AIWA CLINIC	0215282382
PST04090 BYPASS AIWA CLINIC	BOX PST04090 BYPASS AIWA CLINIC	0214991349
PST04090 BYPASS AIWA CLINIC	BOX PST04090 BYPASS AIWA CLINIC	0214991348
PST04090 BYPASS AIWA CLINIC	BOX PST04090 BYPASS AIWA CLINIC	0214991347
SGGLOGX1 ZESTAW DO BARIATRII	BOX SGGLOGX1 ZESTAW DO BARIATRII	0215017846
SGGLOGX1 ZESTAW DO BARIATRII	BOX SGGLOGX1 ZESTAW DO BARIATRII	0215191012
SGSWP1 ZESTAW DO RESEKCJI REKAWOJEJ	BOX SGSWP1 ZESTAW DO RESEKCJI REKAWOJEJ	0215218687
SLEEVEGK1 ZESTAW DO ZABIEGOW SLEEVEX0	SLEEVEGK1 ZESTAW DO ZABIEGOW SLEEVEX1	0215244887
SLEEVEGK3	BOX SLEEVEGK3	0215050387
SLEEVEGK3	BOX SLEEVEGK3	0215193250
SLEEVEGK4 SLEEVE GASTRECTOMY	BOX SLEEVEGK4 SLEEVE GASTRECTOMY	0214933111
SLEEVEJAWI	BOX SLEEVEJAWI SLEEVE GASTRECTOMY BOX	0214030210
SLEEVEJAWI SLEEVE GASTRECTOMY BOX	BOX SLEEVEJAWI SLEEVE GASTRECTOMY BOX	0215041719
SLEEVEJAWI SLEEVE GASTRECTOMY BOX	BOX SLEEVEJAWI SLEEVE GASTRECTOMY BOX	0215161529
SLEEVEJAWI SLEEVE GASTRECTOMY BOX	BOX SLEEVEJAWI SLEEVE GASTRECTOMY BOX	0215272924
SSLEEVE ZEST DO SLEEVE GASTREKTOMIIX1	SSLEEVE ZEST DO SLEEVE GASTREKTOMIIX1	0214985253
SSLEEVE ZEST DO SLEEVE GASTREKTOMIIX1	SSLEEVE ZEST DO SLEEVE GASTREKTOMIIX1	0215104506
SSLEEVE ZEST DO SLEEVE GASTREKTOMIIX1	SSLEEVE ZEST DO SLEEVE GASTREKTOMIIX1	0215298203
TGGLOGOWX1 DO TOTAL GASTRECTOMIII	BOX TGGLOGOWX1 DO TOTAL GASTRECTOMIII	0214726876
VATS1SCDX	BOX VATS1SCDX DO ZABIEGOW VATS	0214874566
VATS1SCDX DO ZABIEGOW VATS	BOX VATS1SCDX DO ZABIEGOW VATS	0215001737
VATS1SCDX DO ZABIEGOW VATS	BOX VATS1SCDX DO ZABIEGOW VATS	0215151361
VATSLOB11	BOX VATSLOB11 ZESTAW DO VATS LOBECTOMII	0214866335
VATSLOB11 ZESTAW DO VATS LOBECTOMII	BOX VATSLOB11 ZESTAW DO VATS LOBECTOMII	0214993030
VATSLOB11 ZESTAW DO VATS LOBECTOMII	BOX VATSLOB11 ZESTAW DO VATS LOBECTOMII	0215178525
VMIN ZESTAW DO VATS PODSTAWOWY	BOX VMIN ZESTAW DO VATS PODSTAWOWY	0215122341
VMIN ZESTAW DO VATS PODSTAWOWY	BOX VMIN ZESTAW DO VATS PODSTAWOWY	0215226207
VWR ZESTAW DO VATS KLIN	BOX VWR ZESTAW DO VATS KLIN	0215122344

Attachment C





49 Changi South Avenue 2
Singapore 486056
www.medtronic.com

tel +65 6436 5000
fax +65 6776 6355

Customer Confirmation Form
Covidien Endo GIA™ Articulating Reloads

ALL CUSTOMERS PLEASE COMPLETE THE FORM IN ITS ENTIRETY

Customer Contact Details	Medtronic Contact Details
Hospital / HCP:	By E-mail:
Address:	By Post:
Telephone no:	
Fax no:	
E-mail:	

Please tick (either one) only:

- ☐ No, no affected inventory to be returned.
- ☐ Yes, affected inventory to be returned (Please fill up the below table if chosen)

Item code	Product Description	Lot number	QTY (Each)

I have read and understand the instructions provided and acknowledge receipt of the Covidien Endo GIA™ Articulating Reloads recall notification dated 26 April 2018 by signing below:

Name: _____ (print) Signature: _____ Stamp: _____ Date: _____