



Edwards

**MEDICAL DEVICE CORRECTION
ACTION REQUIRED**

Field Corrective Action #148

Product: HemoSphere Oximetry Cable

Model Number: HEMOXSC100

Serial Numbers: See Appendix A

(Note: No affected lots distributed into Singapore market)

Dec 12, 2019

<Customer #>

<Contact name or Dept.>

<Firm Name>

<Attention: RISK MANAGEMENT>

cc: Chairman Medical Board and relevant Head of Departments

<Address>

<City/state/zip>

Dear Valued Customer:

As part of our strong commitment to quality, Edwards Lifesciences continuously monitors our products throughout their life cycle to quickly identify and correct any potential issues. We became aware that some HemoSphere oximetry cables have not met customer expectations.

Description of the problem:


This issue involves HemoSphere oximetry cables that have operated at a temperature which may cause an Oximetry Cable Malfunction fault/error message and has the potential to stop producing Oximetry values. Although there have been complaints of the cable being warm to the touch, there have been no patient injuries noted to date. This non-conformance may cause a delay in the procedure while the product is exchanged.

Per investigation, the units were determined to be warmer than normal resulting in the cable shutting off due to internal safety checks, preventing the overheating from reaching a temperature that could cause injury/burn.

Edwards has designed a new hardware for the HemoSphere oximetry cable that is now available to address this issue.

The product with the previous hardware design is being voluntarily recalled. You may contact Edwards Technical Support at 65 9658 4570 for a Return Goods Authorization (RGA) number to return identified inventory for a credit. Technical Support can assist you with placing a new order to replace this inventory, or you may be contacted for an onsite visit to replace your inventory.

We request that you complete the attached Acknowledgement Form and return to Edwards Lifesciences per the instructions on the form.



Edwards

Description and Indication of product being recalled:

The HemoSphere oximetry cable is a multi-use device that connects to a compatible Edwards monitor on one end and any approved Edwards oximetry catheter on the other to continuously measure venous oxygen saturation by reflectance spectrophotometry. LEDs within the oximetry cable transmit light fiber optically to the distal end of the catheter. The amount of light absorbed, refracted, and reflected depends on the relative amounts of oxygenated and deoxygenated hemoglobin in the blood. This optical intensity data is gathered by the oximetry catheter, processed by the HemoSphere oximetry cable and displayed on a compatible monitoring platform.

Parameter output is mixed venous oxygen saturation (SvO₂) or central venous oxygen saturation (ScvO₂).

The HemoSphere advanced monitoring platform is intended to be used by qualified personnel or trained clinicians in a hospital setting. The HemoSphere advanced monitoring platform is intended for use with compatible Edwards oximetry and Swan-Ganz catheters.

If you have any questions, please contact Edwards Technical Support at 65 9658 4570.

Sincerely,

Sunita Das
Director | Quality Assurance
Edwards Lifesciences



Edwards

**MEDICAL DEVICE CORRECTION
ACTION REQUIRED**

Field Corrective Action #148

Product: HemoSphere Oximetry Cable

Model Number: HEMOXSC100

Serial Numbers: See Appendix A

(Note: No affected lots distributed into Singapore market)

CUSTOMER ACKNOWLEDGEMENT

<Customer #>

<Firm Name>

<Attention: RISK MANAGEMENT>

cc: Chairman Medical Board and relevant Head of Departments

<Address>

<City/state/zip>

Please follow all instructions below to complete the acknowledgement process.

Complete this acknowledgement form with the following information:

- Verify your inventory
- Complete all sections of the table below, indicate "0" if you do not have product
- If you have any product to return, call Edwards Technical Support at 65 9658 4570 to obtain a Returned Good Authorization (RGA) number.
- Complete and return the attached Acknowledgement Form within five (5) business days of receiving this notice to Technical Support via fax at 65 9658 4570.

Model	Serial Number	Ship to Date	Quantity Shipped From EW	Number of units to be returned	RGA Number
HEMOXSC100					

Name (Print): _____

Title/Dept. _____

Telephone Number: _____

Signature: _____

Date: _____

Item Number	Serial Number
HEMOXSC100	13703587
HEMOXSC100	13700534
HEMOXSC100	13700244
HEMOXSC100	13700371
HEMOXSC100	13700967
HEMOXSC100	13701335
HEMOXSC100	13701455
HEMOXSC100	13703371
HEMOXSC100	13703381
HEMOXSC100	13703396
HEMOXSC100	13704342
HEMOXSC100	13704472
HEMOXSC100	13703175
HEMOXSC100	13701375
HEMOXSC100	13700278
HEMOXSC100	13700234
HEMOXSC100	13700973
HEMOXSC100	13701185
HEMOXSC100	13700970
HEMOXSC100	13703161
HEMOXSC100	13703595
HEMOXSC100	13702094
HEMOXSC100	13702144
HEMOXSC100	13703129
HEMOXSC100	13701054
HEMOXSC100	13701208
HEMOXSC100	13700256
HEMOXSC100	13701216
HEMOXSC100	13700054
HEMOXSC100	13700719
HEMOXSC100	13700405
HEMOXSC100	13700969
HEMOXSC100	13700975
HEMOXSC100	13700987
HEMOXSC100	13701009
HEMOXSC100	13701015
HEMOXSC100	13701459
HEMOXSC100	13702761
HEMOXSC100	13704495
HEMOXSC100	13700334
HEMOXSC100	13700321
HEMOXSC100	13703584
HEMOXSC100	13703585
HEMOXSC100	13703594
HEMOXSC100	13700152
HEMOXSC100	13700341
HEMOXSC100	13700453
HEMOXSC100	13701494

HEMOXSC100	13704488
HEMOXSC100	13700486
HEMOXSC100	13700991
HEMOXSC100	13700281
HEMOXSC100	13700172
HEMOXSC100	13700057
HEMOXSC100	13700055
HEMOXSC100	13701512
HEMOXSC100	13701542
HEMOXSC100	13701545
HEMOXSC100	13705280
HEMOXSC100	13705283
HEMOXSC100	13705270
HEMOXSC100	13705274
HEMOXSC100	13705279
HEMOXSC100	13705281
HEMOXSC100	13705287
HEMOXSC100	13705290
HEMOXSC100	13705291
HEMOXSC100	13705292
HEMOXSC100	13705295
HEMOXSC100	13705307
HEMOXSC100	13704967
HEMOXSC100	13704970
HEMOXSC100	13704978
HEMOXSC100	13704980
HEMOXSC100	13704981
HEMOXSC100	13704986
HEMOXSC100	13701983
HEMOXSC100	13702014
HEMOXSC100	13705203
HEMOXSC100	13703357
HEMOXSC100	13703358
HEMOXSC100	13703361
HEMOXSC100	13703363
HEMOXSC100	13705120
HEMOXSC100	13710012
HEMOXSC100	13710291
HEMOXSC100	13710292
HEMOXSC100	13705194
HEMOXSC100	13705195
HEMOXSC100	13705193
HEMOXSC100	13701359
HEMOXSC100	13701594
HEMOXSC100	13702847
HEMOXSC100	13702851
HEMOXSC100	13702855
HEMOXSC100	13704791
HEMOXSC100	13704806
HEMOXSC100	13704809

HEMOXSC100	13704825
HEMOXSC100	13704829
HEMOXSC100	13704786
HEMOXSC100	13704816
HEMOXSC100	13704827
HEMOXSC100	13705211
HEMOXSC100	13705320
HEMOXSC100	13705321
HEMOXSC100	13705319
HEMOXSC100	13705323
HEMOXSC100	13703377
HEMOXSC100	13703370
HEMOXSC100	13703372
HEMOXSC100	13710192
HEMOXSC100	13710205
HEMOXSC100	13710209
HEMOXSC100	13710217
HEMOXSC100	13710193
HEMOXSC100	13710195
HEMOXSC100	13710207
HEMOXSC100	13700768
HEMOXSC100	13700774
HEMOXSC100	13701623
HEMOXSC100	13701626
HEMOXSC100	13700670
HEMOXSC100	13701202
HEMOXSC100	13701338
HEMOXSC100	13701457
HEMOXSC100	13710057
HEMOXSC100	13710058
HEMOXSC100	13704460
HEMOXSC100	13705271
HEMOXSC100	13702313
HEMOXSC100	13702317
HEMOXSC100	13702320
HEMOXSC100	13703571
HEMOXSC100	13704310
HEMOXSC100	13701426
HEMOXSC100	13701429
HEMOXSC100	13701431
HEMOXSC100	13702265
HEMOXSC100	13702614
HEMOXSC100	13702616
HEMOXSC100	13704026
HEMOXSC100	13705007
HEMOXSC100	13705011
HEMOXSC100	13705012
HEMOXSC100	13705020
HEMOXSC100	13705021
HEMOXSC100	13705022

HEMOXSC100	13705024
HEMOXSC100	13705026
HEMOXSC100	13710287
HEMOXSC100	13710289
HEMOXSC100	13704798
HEMOXSC100	13704797
HEMOXSC100	13700459
HEMOXSC100	13701561
HEMOXSC100	13700210
HEMOXSC100	13703214
HEMOXSC100	13703216
HEMOXSC100	13703217
HEMOXSC100	13703218
HEMOXSC100	13703219
HEMOXSC100	13703221
HEMOXSC100	13703222
HEMOXSC100	13701378
HEMOXSC100	13701379
HEMOXSC100	13701885
HEMOXSC100	13704728
HEMOXSC100	13703287
HEMOXSC100	13703291
HEMOXSC100	13703295
HEMOXSC100	13703297
HEMOXSC100	13703281
HEMOXSC100	13703283
HEMOXSC100	13703284
HEMOXSC100	13703286
HEMOXSC100	13703288
HEMOXSC100	13703290
HEMOXSC100	13703293
HEMOXSC100	13703294
HEMOXSC100	13705216
HEMOXSC100	13705255
HEMOXSC100	13705259
HEMOXSC100	13705213
HEMOXSC100	13705221
HEMOXSC100	13705230
HEMOXSC100	13705231
HEMOXSC100	13705232
HEMOXSC100	13705240
HEMOXSC100	13705248
HEMOXSC100	13705249
HEMOXSC100	13705250
HEMOXSC100	13705251
HEMOXSC100	13705252
HEMOXSC100	13705253
HEMOXSC100	13710335
HEMOXSC100	13710340
HEMOXSC100	13701675

HEMOXSC100	13701671
HEMOXSC100	13701674
HEMOXSC100	13701656
HEMOXSC100	13701496
HEMOXSC100	13701411
HEMOXSC100	13701412
HEMOXSC100	13701413
HEMOXSC100	13701414
HEMOXSC100	13701415
HEMOXSC100	13701416
HEMOXSC100	13703163
HEMOXSC100	13703183
HEMOXSC100	13703203
HEMOXSC100	13703205
HEMOXSC100	13703355
HEMOXSC100	13703356
HEMOXSC100	13705303
HEMOXSC100	13705304
HEMOXSC100	13703401
HEMOXSC100	13703392
HEMOXSC100	13703394
HEMOXSC100	13703395
HEMOXSC100	13703397
HEMOXSC100	13703398
HEMOXSC100	13703400
HEMOXSC100	13703402
HEMOXSC100	13705165
HEMOXSC100	13705163
HEMOXSC100	13705166
HEMOXSC100	13705167
HEMOXSC100	13701844
HEMOXSC100	13702396
HEMOXSC100	13702332
HEMOXSC100	13702381
HEMOXSC100	13702844
HEMOXSC100	13704672
HEMOXSC100	13704666
HEMOXSC100	13704665
HEMOXSC100	13705261
HEMOXSC100	13702558
HEMOXSC100	13702562
HEMOXSC100	13702563
HEMOXSC100	13702566
HEMOXSC100	13702567
HEMOXSC100	13702568
HEMOXSC100	13702443
HEMOXSC100	13702470
HEMOXSC100	13702517
HEMOXSC100	13702519
HEMOXSC100	13702582

HEMOXSC100	13704173
HEMOXSC100	13704174
HEMOXSC100	13704175
HEMOXSC100	13704176
HEMOXSC100	13704177
HEMOXSC100	13704182
HEMOXSC100	13703788
HEMOXSC100	13703823
HEMOXSC100	13703825
HEMOXSC100	13703914
HEMOXSC100	13703918
HEMOXSC100	13703922
HEMOXSC100	13703929
HEMOXSC100	13703945
HEMOXSC100	13703961
HEMOXSC100	13704278
HEMOXSC100	13704281
HEMOXSC100	13704283
HEMOXSC100	13704284
HEMOXSC100	13704266
HEMOXSC100	13703643
HEMOXSC100	13703656
HEMOXSC100	13703657
HEMOXSC100	13703660
HEMOXSC100	13703669
HEMOXSC100	13703637
HEMOXSC100	13703638
HEMOXSC100	13703661
HEMOXSC100	13703689
HEMOXSC100	13703647
HEMOXSC100	13703650
HEMOXSC100	13703688
HEMOXSC100	13703701
HEMOXSC100	13701649
HEMOXSC100	13701868
HEMOXSC100	13710458
HEMOXSC100	13710459
HEMOXSC100	13710461
HEMOXSC100	13710463
HEMOXSC100	13710466
HEMOXSC100	13710467
HEMOXSC100	13710469
HEMOXSC100	13710471
HEMOXSC100	13710472
HEMOXSC100	13710473
HEMOXSC100	13710474
HEMOXSC100	13710475
HEMOXSC100	13710476
HEMOXSC100	13710477
HEMOXSC100	13710480

HEMOXSC100	13710481
HEMOXSC100	13710483
HEMOXSC100	13710486
HEMOXSC100	13710487
HEMOXSC100	13710488
HEMOXSC100	13710489
HEMOXSC100	13710490
HEMOXSC100	13710491
HEMOXSC100	13710495
HEMOXSC100	13710389
HEMOXSC100	13710417
HEMOXSC100	13710453
HEMOXSC100	13710454
HEMOXSC100	13710456
HEMOXSC100	13710457
HEMOXSC100	13710300
HEMOXSC100	13710304
HEMOXSC100	13702461
HEMOXSC100	13702466
HEMOXSC100	13702564
HEMOXSC100	13701204
HEMOXSC100	13701923
HEMOXSC100	13702624
HEMOXSC100	13702626
HEMOXSC100	13702645
HEMOXSC100	13702646
HEMOXSC100	13702647
HEMOXSC100	13710167
HEMOXSC100	13710175
HEMOXSC100	13705199
HEMOXSC100	13705205
HEMOXSC100	13705207
HEMOXSC100	13705209
HEMOXSC100	13703837
HEMOXSC100	13703068
HEMOXSC100	13703074
HEMOXSC100	13703075
HEMOXSC100	13703054
HEMOXSC100	13703064
HEMOXSC100	13703065
HEMOXSC100	13703057
HEMOXSC100	13703062
HEMOXSC100	13703066
HEMOXSC100	13703067
HEMOXSC100	13703073
HEMOXSC100	13704789
HEMOXSC100	13702080
HEMOXSC100	13705305
HEMOXSC100	13705309
HEMOXSC100	13701583

HEMOXSC100	13701584
HEMOXSC100	13700139
HEMOXSC100	13700140
HEMOXSC100	13700141
HEMOXSC100	13700143
HEMOXSC100	13700144
HEMOXSC100	13700145
HEMOXSC100	13700147
HEMOXSC100	13700148
HEMOXSC100	13700309
HEMOXSC100	13700375
HEMOXSC100	13700379
HEMOXSC100	13700404
HEMOXSC100	13700400
HEMOXSC100	13700799
HEMOXSC100	13700802
HEMOXSC100	13700557
HEMOXSC100	13701563
HEMOXSC100	13710371
HEMOXSC100	13704659
HEMOXSC100	13704663
HEMOXSC100	13702386
HEMOXSC100	13702408
HEMOXSC100	13700463
HEMOXSC100	13700423
HEMOXSC100	13704903
HEMOXSC100	13700981
HEMOXSC100	13700983
HEMOXSC100	13700986
HEMOXSC100	13700998
HEMOXSC100	13701014
HEMOXSC100	13701027
HEMOXSC100	13701053
HEMOXSC100	13700251
HEMOXSC100	13700268
HEMOXSC100	13702739
HEMOXSC100	13703557
HEMOXSC100	13703519
HEMOXSC100	13710027
HEMOXSC100	13710028
HEMOXSC100	13710029
HEMOXSC100	13710030
HEMOXSC100	13710031
HEMOXSC100	13710041
HEMOXSC100	13700635
HEMOXSC100	13700638
HEMOXSC100	13700644
HEMOXSC100	13700651
HEMOXSC100	13700658
HEMOXSC100	13700632

HEMOXSC100	13700633
HEMOXSC100	13700636
HEMOXSC100	13700639
HEMOXSC100	13703271
HEMOXSC100	13702858
HEMOXSC100	13702861
HEMOXSC100	13704021
HEMOXSC100	13704024
HEMOXSC100	13704025
HEMOXSC100	13704028
HEMOXSC100	13704037
HEMOXSC100	13704038
HEMOXSC100	13704042
HEMOXSC100	13704043
HEMOXSC100	13704045
HEMOXSC100	13704054
HEMOXSC100	13704055
HEMOXSC100	13704056
HEMOXSC100	13704057
HEMOXSC100	13704058
HEMOXSC100	13704059
HEMOXSC100	13704061
HEMOXSC100	13704062
HEMOXSC100	13704064
HEMOXSC100	13704671
HEMOXSC100	13703704
HEMOXSC100	13703987
HEMOXSC100	13704027
HEMOXSC100	13704418
HEMOXSC100	13704454
HEMOXSC100	13704471
HEMOXSC100	13704620
HEMOXSC100	13704640
HEMOXSC100	13704550
HEMOXSC100	13704565
HEMOXSC100	13704570
HEMOXSC100	13704578
HEMOXSC100	13704580
HEMOXSC100	13704581
HEMOXSC100	13704585
HEMOXSC100	13704587
HEMOXSC100	13704592
HEMOXSC100	13704593
HEMOXSC100	13704595
HEMOXSC100	13704597
HEMOXSC100	13704598
HEMOXSC100	13704599
HEMOXSC100	13704600
HEMOXSC100	13704601
HEMOXSC100	13704604

HEMOXSC100	13704607
HEMOXSC100	13704608
HEMOXSC100	13704609
HEMOXSC100	13704610
HEMOXSC100	13704613
HEMOXSC100	13704617
HEMOXSC100	13704618
HEMOXSC100	13704622
HEMOXSC100	13701322
HEMOXSC100	13701358
HEMOXSC100	13701423
HEMOXSC100	13701317
HEMOXSC100	13700540
HEMOXSC100	13701838
HEMOXSC100	13703447
HEMOXSC100	13703450
HEMOXSC100	13703451
HEMOXSC100	13703718
HEMOXSC100	13700399
HEMOXSC100	13702472
HEMOXSC100	13702477
HEMOXSC100	13702772
HEMOXSC100	13703299
HEMOXSC100	13703456
HEMOXSC100	13703458
HEMOXSC100	13703459
HEMOXSC100	13704318
HEMOXSC100	13704319
HEMOXSC100	13704688
HEMOXSC100	13704690
HEMOXSC100	13702326
HEMOXSC100	13702391
HEMOXSC100	13701971
HEMOXSC100	13702058
HEMOXSC100	13702199
HEMOXSC100	13702207
HEMOXSC100	13702218
HEMOXSC100	13702250
HEMOXSC100	13702323
HEMOXSC100	13702370
HEMOXSC100	13701795
HEMOXSC100	13701824
HEMOXSC100	13701860
HEMOXSC100	13701962
HEMOXSC100	13702304
HEMOXSC100	13702436
HEMOXSC100	13702437
HEMOXSC100	13702450
HEMOXSC100	13702453
HEMOXSC100	13702471

HEMOXSC100	13702526
HEMOXSC100	13702551
HEMOXSC100	13702569
HEMOXSC100	13702572
HEMOXSC100	13702603
HEMOXSC100	13702604
HEMOXSC100	13702608
HEMOXSC100	13702609
HEMOXSC100	13702611
HEMOXSC100	13701833
HEMOXSC100	13702630
HEMOXSC100	13702631
HEMOXSC100	13702633
HEMOXSC100	13702639
HEMOXSC100	13702654
HEMOXSC100	13702693
HEMOXSC100	13702741
HEMOXSC100	13702746
HEMOXSC100	13702710
HEMOXSC100	13702748
HEMOXSC100	13702816
HEMOXSC100	13702836
HEMOXSC100	13702762
HEMOXSC100	13702776
HEMOXSC100	13702789
HEMOXSC100	13702807
HEMOXSC100	13702812
HEMOXSC100	13702819
HEMOXSC100	13702797
HEMOXSC100	13702800
HEMOXSC100	13704295
HEMOXSC100	13704302
HEMOXSC100	13704299
HEMOXSC100	13704301
HEMOXSC100	13704317
HEMOXSC100	13704320
HEMOXSC100	13704297
HEMOXSC100	13704312
HEMOXSC100	13704422
HEMOXSC100	13704325
HEMOXSC100	13702281
HEMOXSC100	13702266
HEMOXSC100	13702272
HEMOXSC100	13702306
HEMOXSC100	13702318
HEMOXSC100	13702246
HEMOXSC100	13702267
HEMOXSC100	13702268
HEMOXSC100	13702284
HEMOXSC100	13702342

HEMOXSC100	13702401
HEMOXSC100	13701935
HEMOXSC100	13702079
HEMOXSC100	13702086
HEMOXSC100	13702201
HEMOXSC100	13702211
HEMOXSC100	13702216
HEMOXSC100	13702217
HEMOXSC100	13702219
HEMOXSC100	13702678
HEMOXSC100	13702714
HEMOXSC100	13704576
HEMOXSC100	13704647
HEMOXSC100	13700303
HEMOXSC100	13704803
HEMOXSC100	13704801
HEMOXSC100	13700326
HEMOXSC100	13700665
HEMOXSC100	13700666
HEMOXSC100	13700680
HEMOXSC100	13700682
HEMOXSC100	13700686
HEMOXSC100	13700723
HEMOXSC100	13700727
HEMOXSC100	13700745
HEMOXSC100	13700747
HEMOXSC100	13700750
HEMOXSC100	13700755
HEMOXSC100	13700756
HEMOXSC100	13700758
HEMOXSC100	13700412
HEMOXSC100	13700415
HEMOXSC100	13700422
HEMOXSC100	13700203
HEMOXSC100	13700227
HEMOXSC100	13700253
HEMOXSC100	13700255
HEMOXSC100	13700302
HEMOXSC100	13700312
HEMOXSC100	13700318
HEMOXSC100	13700324
HEMOXSC100	13700350
HEMOXSC100	13700359
HEMOXSC100	13700377
HEMOXSC100	13700388
HEMOXSC100	13700622
HEMOXSC100	13702171
HEMOXSC100	13705224
HEMOXSC100	13705239
HEMOXSC100	13705223

HEMOXSC100	13705225
HEMOXSC100	13705226
HEMOXSC100	13705228
HEMOXSC100	13705234
HEMOXSC100	13705238
HEMOXSC100	13702084
HEMOXSC100	13702210
HEMOXSC100	13701889
HEMOXSC100	13701890
HEMOXSC100	13701895
HEMOXSC100	13701897
HEMOXSC100	13701097
HEMOXSC100	13702101
HEMOXSC100	13702102
HEMOXSC100	13700815
HEMOXSC100	13702651
HEMOXSC100	13702930
HEMOXSC100	13703809
HEMOXSC100	13705289
HEMOXSC100	13705294
HEMOXSC100	13705297
HEMOXSC100	13705288
HEMOXSC100	13710153
HEMOXSC100	13710154
HEMOXSC100	13710155
HEMOXSC100	13710156
HEMOXSC100	13704861
HEMOXSC100	13704862
HEMOXSC100	13700982
HEMOXSC100	13702234
HEMOXSC100	13702236
HEMOXSC100	13704132
HEMOXSC100	13700742
HEMOXSC100	13700776
HEMOXSC100	13700779
HEMOXSC100	13701565
HEMOXSC100	13701307
HEMOXSC100	13701303
HEMOXSC100	13701316
HEMOXSC100	13702514
HEMOXSC100	13702515
HEMOXSC100	13702520
HEMOXSC100	13702523
HEMOXSC100	13702524
HEMOXSC100	13702559
HEMOXSC100	13701664
HEMOXSC100	13701663
HEMOXSC100	13701873
HEMOXSC100	13702325
HEMOXSC100	13702165

HEMOXSC100	13704287
HEMOXSC100	13704285
HEMOXSC100	13700179
HEMOXSC100	13700199
HEMOXSC100	13701514
HEMOXSC100	13701518
HEMOXSC100	13701519
HEMOXSC100	13702508
HEMOXSC100	13704871
HEMOXSC100	13704869
HEMOXSC100	13704873
HEMOXSC100	13704883
HEMOXSC100	13704872
HEMOXSC100	13704879
HEMOXSC100	13704880
HEMOXSC100	13704881
HEMOXSC100	13710513
HEMOXSC100	13710514
HEMOXSC100	13700703
HEMOXSC100	13700705
HEMOXSC100	13700706
HEMOXSC100	13700895
HEMOXSC100	13700898
HEMOXSC100	13700900
HEMOXSC100	13700908
HEMOXSC100	13700911
HEMOXSC100	13700914
HEMOXSC100	13700896
HEMOXSC100	13700904
HEMOXSC100	13700913
HEMOXSC100	13702576
HEMOXSC100	13701955
HEMOXSC100	13701956
HEMOXSC100	13704303
HEMOXSC100	13704311
HEMOXSC100	13704304
HEMOXSC100	13704309
HEMOXSC100	13704702
HEMOXSC100	13704727
HEMOXSC100	13705227
HEMOXSC100	13705229
HEMOXSC100	13705237
HEMOXSC100	13710407
HEMOXSC100	13710325
HEMOXSC100	13710330
HEMOXSC100	13700759
HEMOXSC100	13701842
HEMOXSC100	13702221
HEMOXSC100	13702560
HEMOXSC100	13710329

HEMOXSC100	13710334
HEMOXSC100	13710310
HEMOXSC100	13710374
HEMOXSC100	13710375
HEMOXSC100	13710377
HEMOXSC100	13710384
HEMOXSC100	13710045
HEMOXSC100	13710051
HEMOXSC100	13710056
HEMOXSC100	13710492
HEMOXSC100	13702857
HEMOXSC100	13703836
HEMOXSC100	13704631
HEMOXSC100	13704094
HEMOXSC100	13704129
HEMOXSC100	13710176
HEMOXSC100	13702574
HEMOXSC100	13703469
HEMOXSC100	13703470
HEMOXSC100	13703471
HEMOXSC100	13703475
HEMOXSC100	13703477
HEMOXSC100	13703480
HEMOXSC100	13703409
HEMOXSC100	13703415
HEMOXSC100	13703417
HEMOXSC100	13703423
HEMOXSC100	13703433
HEMOXSC100	13703406
HEMOXSC100	13703408
HEMOXSC100	13703410
HEMOXSC100	13703412
HEMOXSC100	13703413
HEMOXSC100	13703424
HEMOXSC100	13703426
HEMOXSC100	13703429
HEMOXSC100	13703431
HEMOXSC100	13703435
HEMOXSC100	13703436
HEMOXSC100	13703100
HEMOXSC100	13703527
HEMOXSC100	13703474
HEMOXSC100	13704161
HEMOXSC100	13701691
HEMOXSC100	13703168
HEMOXSC100	13703573
HEMOXSC100	13704330
HEMOXSC100	13702418
HEMOXSC100	13702420
HEMOXSC100	13700497

HEMOXSC100	13700888
HEMOXSC100	13700890
HEMOXSC100	13700899
HEMOXSC100	13700909
HEMOXSC100	13700910
HEMOXSC100	13700935
HEMOXSC100	13700942
HEMOXSC100	13700952
HEMOXSC100	13700965
HEMOXSC100	13700976
HEMOXSC100	13700984
HEMOXSC100	13700996
HEMOXSC100	13701001
HEMOXSC100	13701004
HEMOXSC100	13701005
HEMOXSC100	13701006
HEMOXSC100	13701007
HEMOXSC100	13701011
HEMOXSC100	13701012
HEMOXSC100	13701022
HEMOXSC100	13701030
HEMOXSC100	13701033
HEMOXSC100	13701109
HEMOXSC100	13701115
HEMOXSC100	13701119
HEMOXSC100	13701121
HEMOXSC100	13701125
HEMOXSC100	13701114
HEMOXSC100	13700868
HEMOXSC100	13700869
HEMOXSC100	13701106
HEMOXSC100	13701110
HEMOXSC100	13701116
HEMOXSC100	13701123
HEMOXSC100	13701124
HEMOXSC100	13700181
HEMOXSC100	13700207
HEMOXSC100	13700219
HEMOXSC100	13700260
HEMOXSC100	13700535
HEMOXSC100	13700209
HEMOXSC100	13700220
HEMOXSC100	13700226
HEMOXSC100	13700263
HEMOXSC100	13700322
HEMOXSC100	13700349
HEMOXSC100	13700368
HEMOXSC100	13700478
HEMOXSC100	13700626
HEMOXSC100	13700669

HEMOXSC100	13700204
HEMOXSC100	13700221
HEMOXSC100	13700238
HEMOXSC100	13700292
HEMOXSC100	13700374
HEMOXSC100	13700396
HEMOXSC100	13700481
HEMOXSC100	13700487
HEMOXSC100	13700571
HEMOXSC100	13700576
HEMOXSC100	13700595
HEMOXSC100	13700677
HEMOXSC100	13701118
HEMOXSC100	13701107
HEMOXSC100	13701108
HEMOXSC100	13701112
HEMOXSC100	13701464
HEMOXSC100	13701484
HEMOXSC100	13701485
HEMOXSC100	13701608
HEMOXSC100	13701894
HEMOXSC100	13701901
HEMOXSC100	13701798
HEMOXSC100	13701810
HEMOXSC100	13701817
HEMOXSC100	13701828
HEMOXSC100	13701862
HEMOXSC100	13701863
HEMOXSC100	13701865
HEMOXSC100	13701869
HEMOXSC100	13701883
HEMOXSC100	13701899
HEMOXSC100	13701900
HEMOXSC100	13701905
HEMOXSC100	13701906
HEMOXSC100	13701910
HEMOXSC100	13701981
HEMOXSC100	13702055
HEMOXSC100	13702118
HEMOXSC100	13702279
HEMOXSC100	13702557
HEMOXSC100	13703274
HEMOXSC100	13703315
HEMOXSC100	13703254
HEMOXSC100	13703302
HEMOXSC100	13703303
HEMOXSC100	13703565
HEMOXSC100	13704000
HEMOXSC100	13704116
HEMOXSC100	13704117

HEMOXSC100	13702425
HEMOXSC100	13702431
HEMOXSC100	13700676
HEMOXSC100	13700687
HEMOXSC100	13700674
HEMOXSC100	13701977
HEMOXSC100	13702197
HEMOXSC100	13702200
HEMOXSC100	13710305
HEMOXSC100	13710306
HEMOXSC100	13710308
HEMOXSC100	13700171
HEMOXSC100	13700173
HEMOXSC100	13700176
HEMOXSC100	13700178
HEMOXSC100	13700200
HEMOXSC100	13700232
HEMOXSC100	13703298
HEMOXSC100	13703171
HEMOXSC100	13703182
HEMOXSC100	13703184
HEMOXSC100	13703188
HEMOXSC100	13703189
HEMOXSC100	13703190
HEMOXSC100	13703191
HEMOXSC100	13703192
HEMOXSC100	13703162
HEMOXSC100	13703164
HEMOXSC100	13703165
HEMOXSC100	13703166
HEMOXSC100	13703170
HEMOXSC100	13703173
HEMOXSC100	13703174
HEMOXSC100	13703176
HEMOXSC100	13703185
HEMOXSC100	13703197
HEMOXSC100	13703199
HEMOXSC100	13703201
HEMOXSC100	13703206
HEMOXSC100	13703208
HEMOXSC100	13703209
HEMOXSC100	13703210
HEMOXSC100	13703212
HEMOXSC100	13704792
HEMOXSC100	13700818
HEMOXSC100	13700819
HEMOXSC100	13700822
HEMOXSC100	13700827
HEMOXSC100	13700839
HEMOXSC100	13700841

HEMOXSC100	13700816
HEMOXSC100	13700817
HEMOXSC100	13700820
HEMOXSC100	13700838
HEMOXSC100	13700842
HEMOXSC100	13700825
HEMOXSC100	13700823
HEMOXSC100	13700834
HEMOXSC100	13700840
HEMOXSC100	13700824
HEMOXSC100	13700826
HEMOXSC100	13703326
HEMOXSC100	13703328
HEMOXSC100	13703321
HEMOXSC100	13703323
HEMOXSC100	13703324
HEMOXSC100	13703325
HEMOXSC100	13703327
HEMOXSC100	13700300
HEMOXSC100	13701546
HEMOXSC100	13703403
HEMOXSC100	13703562
HEMOXSC100	13703569
HEMOXSC100	13704648
HEMOXSC100	13704654
HEMOXSC100	13704655
HEMOXSC100	13704656
HEMOXSC100	13704657
HEMOXSC100	13704822
HEMOXSC100	13704815
HEMOXSC100	13704817
HEMOXSC100	13704818
HEMOXSC100	13704819
HEMOXSC100	13704820
HEMOXSC100	13704828
HEMOXSC100	13704823
HEMOXSC100	13704821
HEMOXSC100	13704741
HEMOXSC100	13704866
HEMOXSC100	13704875
HEMOXSC100	13710119
HEMOXSC100	13710122
HEMOXSC100	13710123
HEMOXSC100	13710127
HEMOXSC100	13710134
HEMOXSC100	13710135
HEMOXSC100	13710143
HEMOXSC100	13710144
HEMOXSC100	13710114
HEMOXSC100	13710125

HEMOXSC100	13710141
HEMOXSC100	13710145
HEMOXSC100	13700290
HEMOXSC100	13702227
HEMOXSC100	13702238
HEMOXSC100	13702256
HEMOXSC100	13702259
HEMOXSC100	13702282
HEMOXSC100	13702296
HEMOXSC100	13702511
HEMOXSC100	13710341
HEMOXSC100	13710549
HEMOXSC100	13710574
HEMOXSC100	13703448
HEMOXSC100	13703449
HEMOXSC100	13703453
HEMOXSC100	13703457
HEMOXSC100	13703460
HEMOXSC100	13703461
HEMOXSC100	13703482
HEMOXSC100	13703487
HEMOXSC100	13703493
HEMOXSC100	13703494
HEMOXSC100	13703495
HEMOXSC100	13703498
HEMOXSC100	13703500
HEMOXSC100	13703501
HEMOXSC100	13703486
HEMOXSC100	13703784
HEMOXSC100	13710321
HEMOXSC100	13700714
HEMOXSC100	13700737
HEMOXSC100	13700709
HEMOXSC100	13700720
HEMOXSC100	13700748
HEMOXSC100	13700772
HEMOXSC100	13700780
HEMOXSC100	13700786
HEMOXSC100	13700796
HEMOXSC100	13700800
HEMOXSC100	13701703
HEMOXSC100	13701704
HEMOXSC100	13701066
HEMOXSC100	13701196
HEMOXSC100	13701374
HEMOXSC100	13701406
HEMOXSC100	13701203
HEMOXSC100	13701463
HEMOXSC100	13701465
HEMOXSC100	13701500

HEMOXSC100	13701391
HEMOXSC100	13701257
HEMOXSC100	13701289
HEMOXSC100	13701258
HEMOXSC100	13701267
HEMOXSC100	13701282
HEMOXSC100	13701286
HEMOXSC100	13701290
HEMOXSC100	13701291
HEMOXSC100	13703916
HEMOXSC100	13703928
HEMOXSC100	13703919
HEMOXSC100	13703925
HEMOXSC100	13703930
HEMOXSC100	13703936
HEMOXSC100	13703442
HEMOXSC100	13703506
HEMOXSC100	13703509
HEMOXSC100	13703512
HEMOXSC100	13703515
HEMOXSC100	13703517
HEMOXSC100	13703464
HEMOXSC100	13703503
HEMOXSC100	13703505
HEMOXSC100	13703180
HEMOXSC100	13710428
HEMOXSC100	13710429
HEMOXSC100	13710430
HEMOXSC100	13704974
HEMOXSC100	13701036
HEMOXSC100	13701039
HEMOXSC100	13701037
HEMOXSC100	13701038
HEMOXSC100	13702581
HEMOXSC100	13702586
HEMOXSC100	13704226
HEMOXSC100	13704231
HEMOXSC100	13704236
HEMOXSC100	13704239
HEMOXSC100	13704240
HEMOXSC100	13704241
HEMOXSC100	13704242
HEMOXSC100	13704244
HEMOXSC100	13704248
HEMOXSC100	13704249
HEMOXSC100	13704252
HEMOXSC100	13704253
HEMOXSC100	13704254
HEMOXSC100	13704256
HEMOXSC100	13704258

HEMOXSC100	13704263
HEMOXSC100	13704265
HEMOXSC100	13704228
HEMOXSC100	13704232
HEMOXSC100	13704233
HEMOXSC100	13704237
HEMOXSC100	13704238
HEMOXSC100	13704243
HEMOXSC100	13704251
HEMOXSC100	13704260
HEMOXSC100	13704987
HEMOXSC100	13704989
HEMOXSC100	13704992
HEMOXSC100	13704994
HEMOXSC100	13705009
HEMOXSC100	13705015
HEMOXSC100	13700136
HEMOXSC100	13700137
HEMOXSC100	13701057
HEMOXSC100	13701425
HEMOXSC100	13700767
HEMOXSC100	13700753
HEMOXSC100	13700760
HEMOXSC100	13700762
HEMOXSC100	13700763
HEMOXSC100	13701120
HEMOXSC100	13701122
HEMOXSC100	13701127
HEMOXSC100	13701136
HEMOXSC100	13701158
HEMOXSC100	13705260
HEMOXSC100	13705236
HEMOXSC100	13705244
HEMOXSC100	13705266
HEMOXSC100	13705267
HEMOXSC100	13705269
HEMOXSC100	13702696
HEMOXSC100	13702698
HEMOXSC100	13704163
HEMOXSC100	13704168
HEMOXSC100	13704134
HEMOXSC100	13704178
HEMOXSC100	13704181
HEMOXSC100	13704187
HEMOXSC100	13704189
HEMOXSC100	13704203
HEMOXSC100	13704204
HEMOXSC100	13704205
HEMOXSC100	13710190
HEMOXSC100	13710191

HEMOXSC100	13710198
HEMOXSC100	13710199
HEMOXSC100	13710201
HEMOXSC100	13703149
HEMOXSC100	13705045
HEMOXSC100	13705052
HEMOXSC100	13705053
HEMOXSC100	13705057
HEMOXSC100	13700696
HEMOXSC100	13700698
HEMOXSC100	13700699
HEMOXSC100	13700701
HEMOXSC100	13700691
HEMOXSC100	13700689
HEMOXSC100	13700694
HEMOXSC100	13705072
HEMOXSC100	13705075
HEMOXSC100	13705076
HEMOXSC100	13701468
HEMOXSC100	13705296
HEMOXSC100	13705298
HEMOXSC100	13705246
HEMOXSC100	13705243
HEMOXSC100	13705247
HEMOXSC100	13705254
HEMOXSC100	13705257
HEMOXSC100	13705258
HEMOXSC100	13700283
HEMOXSC100	13701321
HEMOXSC100	13701326
HEMOXSC100	13701341
HEMOXSC100	13703786
HEMOXSC100	13703787
HEMOXSC100	13703791
HEMOXSC100	13703813
HEMOXSC100	13703815
HEMOXSC100	13703818
HEMOXSC100	13703910
HEMOXSC100	13703128
HEMOXSC100	13703374
HEMOXSC100	13703432
HEMOXSC100	13703602
HEMOXSC100	13703648
HEMOXSC100	13703654
HEMOXSC100	13703663
HEMOXSC100	13703759
HEMOXSC100	13703762
HEMOXSC100	13703764
HEMOXSC100	13703766
HEMOXSC100	13703772

HEMOXSC100	13703777
HEMOXSC100	13703792
HEMOXSC100	13703797
HEMOXSC100	13703798
HEMOXSC100	13703800
HEMOXSC100	13703826
HEMOXSC100	13703830
HEMOXSC100	13703843
HEMOXSC100	13703846
HEMOXSC100	13703850
HEMOXSC100	13703857
HEMOXSC100	13703869
HEMOXSC100	13703876
HEMOXSC100	13703877
HEMOXSC100	13703880
HEMOXSC100	13703883
HEMOXSC100	13703884
HEMOXSC100	13703886
HEMOXSC100	13703738
HEMOXSC100	13703768
HEMOXSC100	13703778
HEMOXSC100	13703780
HEMOXSC100	13703799
HEMOXSC100	13703817
HEMOXSC100	13703820
HEMOXSC100	13703824
HEMOXSC100	13703832
HEMOXSC100	13703844
HEMOXSC100	13703847
HEMOXSC100	13703854
HEMOXSC100	13703858
HEMOXSC100	13703861
HEMOXSC100	13703863
HEMOXSC100	13703864
HEMOXSC100	13703871
HEMOXSC100	13703873
HEMOXSC100	13703882
HEMOXSC100	13703885
HEMOXSC100	13703603
HEMOXSC100	13703607
HEMOXSC100	13703621
HEMOXSC100	13703630
HEMOXSC100	13703838
HEMOXSC100	13703839
HEMOXSC100	13703840
HEMOXSC100	13704590
HEMOXSC100	13710399
HEMOXSC100	13701483
HEMOXSC100	13702152
HEMOXSC100	13702157

HEMOXSC100	13702167
HEMOXSC100	13702170
HEMOXSC100	13702181
HEMOXSC100	13702164
HEMOXSC100	13702173
HEMOXSC100	13702175
HEMOXSC100	13702182
HEMOXSC100	13702150
HEMOXSC100	13702154
HEMOXSC100	13702155
HEMOXSC100	13702177
HEMOXSC100	13702178
HEMOXSC100	13702283
HEMOXSC100	13701930
HEMOXSC100	13702928
HEMOXSC100	13703156
HEMOXSC100	13703896
HEMOXSC100	13703892
HEMOXSC100	13704034
HEMOXSC100	13703796
HEMOXSC100	13703805
HEMOXSC100	13705151
HEMOXSC100	13710256
HEMOXSC100	13710339
HEMOXSC100	13704429
HEMOXSC100	13704430
HEMOXSC100	13704433
HEMOXSC100	13704441
HEMOXSC100	13701677
HEMOXSC100	13701676
HEMOXSC100	13701680
HEMOXSC100	13701681
HEMOXSC100	13703841
HEMOXSC100	13704087
HEMOXSC100	13704089
HEMOXSC100	13704090
HEMOXSC100	13704091
HEMOXSC100	13700191
HEMOXSC100	13700192
HEMOXSC100	13700193
HEMOXSC100	13700194
HEMOXSC100	13702341
HEMOXSC100	13702346
HEMOXSC100	13702344
HEMOXSC100	13702525
HEMOXSC100	13702528
HEMOXSC100	13702532
HEMOXSC100	13702534
HEMOXSC100	13702535
HEMOXSC100	13702538

HEMOXSC100	13702540
HEMOXSC100	13702948
HEMOXSC100	13703224
HEMOXSC100	13704912
HEMOXSC100	13705028
HEMOXSC100	13705029
HEMOXSC100	13705030
HEMOXSC100	13705034
HEMOXSC100	13705036
HEMOXSC100	13705043
HEMOXSC100	13705046
HEMOXSC100	13705032
HEMOXSC100	13705037
HEMOXSC100	13705038
HEMOXSC100	13705040
HEMOXSC100	13705276
HEMOXSC100	13705215
HEMOXSC100	13705214
HEMOXSC100	13705217
HEMOXSC100	13701420
HEMOXSC100	13701434
HEMOXSC100	13701438
HEMOXSC100	13702410
HEMOXSC100	13702117
HEMOXSC100	13702409
HEMOXSC100	13702413
HEMOXSC100	13701940
HEMOXSC100	13701946
HEMOXSC100	13701951
HEMOXSC100	13701988
HEMOXSC100	13701997
HEMOXSC100	13702026
HEMOXSC100	13702034
HEMOXSC100	13702037
HEMOXSC100	13702038
HEMOXSC100	13702039
HEMOXSC100	13702041
HEMOXSC100	13702044
HEMOXSC100	13702046
HEMOXSC100	13702047
HEMOXSC100	13702049
HEMOXSC100	13702060
HEMOXSC100	13702933
HEMOXSC100	13702934
HEMOXSC100	13702946
HEMOXSC100	13703085
HEMOXSC100	13703155
HEMOXSC100	13703390
HEMOXSC100	13703537
HEMOXSC100	13703538

HEMOXSC100	13703819
HEMOXSC100	13704591
HEMOXSC100	13704612
HEMOXSC100	13704790
HEMOXSC100	13710194
HEMOXSC100	13701967
HEMOXSC100	13702176
HEMOXSC100	13702179
HEMOXSC100	13702196
HEMOXSC100	13710219
HEMOXSC100	13702287
HEMOXSC100	13702480
HEMOXSC100	13702481
HEMOXSC100	13702482
HEMOXSC100	13702483
HEMOXSC100	13702484
HEMOXSC100	13702485
HEMOXSC100	13702486
HEMOXSC100	13702487
HEMOXSC100	13702488
HEMOXSC100	13702489
HEMOXSC100	13704639
HEMOXSC100	13704644
HEMOXSC100	13710221
HEMOXSC100	13710231
HEMOXSC100	13710240
HEMOXSC100	13710245
HEMOXSC100	13710267
HEMOXSC100	13710270
HEMOXSC100	13710274
HEMOXSC100	13710277
HEMOXSC100	13710284
HEMOXSC100	13710286
HEMOXSC100	13710312
HEMOXSC100	13710316
HEMOXSC100	13710047
HEMOXSC100	13702667
HEMOXSC100	13702679
HEMOXSC100	13702687
HEMOXSC100	13702697
HEMOXSC100	13702701
HEMOXSC100	13702432
HEMOXSC100	13702433
HEMOXSC100	13702441
HEMOXSC100	13702442
HEMOXSC100	13702445
HEMOXSC100	13702454
HEMOXSC100	13710257
HEMOXSC100	13710285
HEMOXSC100	13701989

HEMOXSC100	13701990
HEMOXSC100	13701995
HEMOXSC100	13701996
HEMOXSC100	13700276
HEMOXSC100	13702056
HEMOXSC100	13702062
HEMOXSC100	13702065
HEMOXSC100	13702125
HEMOXSC100	13702592
HEMOXSC100	13702593
HEMOXSC100	13702594
HEMOXSC100	13701792
HEMOXSC100	13701800
HEMOXSC100	13704630
HEMOXSC100	13704635
HEMOXSC100	13704646
HEMOXSC100	13704214
HEMOXSC100	13704229
HEMOXSC100	13703632
HEMOXSC100	13703670
HEMOXSC100	13703895
HEMOXSC100	13703997
HEMOXSC100	13704040
HEMOXSC100	13704041
HEMOXSC100	13704092
HEMOXSC100	13704120
HEMOXSC100	13704125
HEMOXSC100	13704730
HEMOXSC100	13700213
HEMOXSC100	13700214
HEMOXSC100	13700224
HEMOXSC100	13700231
HEMOXSC100	13704652
HEMOXSC100	13704020
HEMOXSC100	13704577
HEMOXSC100	13704574
HEMOXSC100	13704564
HEMOXSC100	13704626
HEMOXSC100	13703870
HEMOXSC100	13703897
HEMOXSC100	13704298
HEMOXSC100	13704347
HEMOXSC100	13704224
HEMOXSC100	13703879
HEMOXSC100	13704461
HEMOXSC100	13704468
HEMOXSC100	13704443
HEMOXSC100	13704446
HEMOXSC100	13704450
HEMOXSC100	13704959

HEMOXSC100	13703516
HEMOXSC100	13703523
HEMOXSC100	13703526
HEMOXSC100	13703507
HEMOXSC100	13703524
HEMOXSC100	13703530
HEMOXSC100	13703532
HEMOXSC100	13703540
HEMOXSC100	13702105
HEMOXSC100	13704938
HEMOXSC100	13704943
HEMOXSC100	13705208
HEMOXSC100	13702126
HEMOXSC100	13703362
HEMOXSC100	13703349
HEMOXSC100	13703352
HEMOXSC100	13703354
HEMOXSC100	13703359
HEMOXSC100	13703360
HEMOXSC100	13703365
HEMOXSC100	13700301
HEMOXSC100	13700813
HEMOXSC100	13700892
HEMOXSC100	13701559
HEMOXSC100	13703803
HEMOXSC100	13703785
HEMOXSC100	13703795
HEMOXSC100	13702860
HEMOXSC100	13703698
HEMOXSC100	13704140
HEMOXSC100	13704144
HEMOXSC100	13700509
HEMOXSC100	13700920
HEMOXSC100	13700923
HEMOXSC100	13701392
HEMOXSC100	13703626
HEMOXSC100	13703642
HEMOXSC100	13703874
HEMOXSC100	13703020
HEMOXSC100	13703021
HEMOXSC100	13703026
HEMOXSC100	13703028
HEMOXSC100	13703029
HEMOXSC100	13703030
HEMOXSC100	13700724
HEMOXSC100	13704718
HEMOXSC100	13704719
HEMOXSC100	13704721
HEMOXSC100	13703472
HEMOXSC100	13704778

HEMOXSC100	13704761
HEMOXSC100	13704743
HEMOXSC100	13704782
HEMOXSC100	13704788
HEMOXSC100	13704775
HEMOXSC100	13704794
HEMOXSC100	13702113
HEMOXSC100	13702198
HEMOXSC100	13702427
HEMOXSC100	13702438
HEMOXSC100	13705003
HEMOXSC100	13705004
HEMOXSC100	13705005
HEMOXSC100	13700743
HEMOXSC100	13702407
HEMOXSC100	13702412
HEMOXSC100	13703096
HEMOXSC100	13703087
HEMOXSC100	13703105
HEMOXSC100	13703097
HEMOXSC100	13703101
HEMOXSC100	13703106
HEMOXSC100	13703091
HEMOXSC100	13703855
HEMOXSC100	13703806
HEMOXSC100	13701403
HEMOXSC100	13701439
HEMOXSC100	13701443
HEMOXSC100	13701445
HEMOXSC100	13701447
HEMOXSC100	13701448
HEMOXSC100	13702981
HEMOXSC100	13703202
HEMOXSC100	13703207
HEMOXSC100	13703211
HEMOXSC100	13703198
HEMOXSC100	13703200
HEMOXSC100	13703204
HEMOXSC100	13703220
HEMOXSC100	13701031
HEMOXSC100	13702273
HEMOXSC100	13702242
HEMOXSC100	13702269
HEMOXSC100	13703069
HEMOXSC100	13703116
HEMOXSC100	13703257
HEMOXSC100	13703993
HEMOXSC100	13704162
HEMOXSC100	13705204
HEMOXSC100	13703261

HEMOXSC100	13703258
HEMOXSC100	13703259
HEMOXSC100	13703262
HEMOXSC100	13703263
HEMOXSC100	13703264
HEMOXSC100	13710560
HEMOXSC100	13702255
HEMOXSC100	13702297
HEMOXSC100	13702305
HEMOXSC100	13702351
HEMOXSC100	13702991
HEMOXSC100	13702992
HEMOXSC100	13702985
HEMOXSC100	13703114
HEMOXSC100	13700647
HEMOXSC100	13700642
HEMOXSC100	13700643
HEMOXSC100	13710311
HEMOXSC100	13710386
HEMOXSC100	13710388
HEMOXSC100	13710391
HEMOXSC100	13710392
HEMOXSC100	13710393
HEMOXSC100	13710395
HEMOXSC100	13710396
HEMOXSC100	13710397
HEMOXSC100	13710400
HEMOXSC100	13710401
HEMOXSC100	13710402
HEMOXSC100	13710403
HEMOXSC100	13710405
HEMOXSC100	13710406
HEMOXSC100	13700879
HEMOXSC100	13701507
HEMOXSC100	13702339
HEMOXSC100	13705062
HEMOXSC100	13705067
HEMOXSC100	13700514
HEMOXSC100	13700523
HEMOXSC100	13700351
HEMOXSC100	13700442
HEMOXSC100	13700474
HEMOXSC100	13700488
HEMOXSC100	13700490
HEMOXSC100	13705212
HEMOXSC100	13700867
HEMOXSC100	13703427
HEMOXSC100	13703428
HEMOXSC100	13704583
HEMOXSC100	13705302

HEMOXSC100	13705300
HEMOXSC100	13704245
HEMOXSC100	13702521
HEMOXSC100	13702537
HEMOXSC100	13702543
HEMOXSC100	13703613
HEMOXSC100	13701796
HEMOXSC100	13701799
HEMOXSC100	13701645
HEMOXSC100	13701650
HEMOXSC100	13701652
HEMOXSC100	13701653
HEMOXSC100	13701672
HEMOXSC100	13701673
HEMOXSC100	13701678
HEMOXSC100	13701679
HEMOXSC100	13700249
HEMOXSC100	13700520
HEMOXSC100	13704268
HEMOXSC100	13700556
HEMOXSC100	13700627
HEMOXSC100	13703336
HEMOXSC100	13703337
HEMOXSC100	13704824
HEMOXSC100	13704865
HEMOXSC100	13704867
HEMOXSC100	13704868
HEMOXSC100	13703723
HEMOXSC100	13703724
HEMOXSC100	13703875
HEMOXSC100	13704050
HEMOXSC100	13704195
HEMOXSC100	13704202
HEMOXSC100	13704222
HEMOXSC100	13700925
HEMOXSC100	13700926
HEMOXSC100	13700928
HEMOXSC100	13700929
HEMOXSC100	13700932
HEMOXSC100	13700933
HEMOXSC100	13700934
HEMOXSC100	13700937
HEMOXSC100	13701047
HEMOXSC100	13701049
HEMOXSC100	13701050
HEMOXSC100	13701056
HEMOXSC100	13702075
HEMOXSC100	13704722
HEMOXSC100	13704725
HEMOXSC100	13704726

HEMOXSC100	13703667
HEMOXSC100	13703674
HEMOXSC100	13703675
HEMOXSC100	13703681
HEMOXSC100	13703684
HEMOXSC100	13705326
HEMOXSC100	13710013
HEMOXSC100	13710014
HEMOXSC100	13710015
HEMOXSC100	13710016
HEMOXSC100	13710021
HEMOXSC100	13710023
HEMOXSC100	13710034
HEMOXSC100	13700765
HEMOXSC100	13700766
HEMOXSC100	13700471
HEMOXSC100	13701393
HEMOXSC100	13704925
HEMOXSC100	13704926
HEMOXSC100	13710025
HEMOXSC100	13704563
HEMOXSC100	13704596
HEMOXSC100	13704603
HEMOXSC100	13704605
HEMOXSC100	13704606
HEMOXSC100	13700183
HEMOXSC100	13700187
HEMOXSC100	13700583
HEMOXSC100	13700590
HEMOXSC100	13710110
HEMOXSC100	13710111
HEMOXSC100	13710065
HEMOXSC100	13710075
HEMOXSC100	13710121
HEMOXSC100	13710124
HEMOXSC100	13710133
HEMOXSC100	13700664
HEMOXSC100	13700810
HEMOXSC100	13700858
HEMOXSC100	13700894
HEMOXSC100	13700974
HEMOXSC100	13701046
HEMOXSC100	13701436
HEMOXSC100	13701437
HEMOXSC100	13703304
HEMOXSC100	13705301
HEMOXSC100	13705313
HEMOXSC100	13705299
HEMOXSC100	13705306
HEMOXSC100	13705311

HEMOXSC100	13704845
HEMOXSC100	13705059
HEMOXSC100	13705178
HEMOXSC100	13705206
HEMOXSC100	13705233
HEMOXSC100	13705235
HEMOXSC100	13710213
HEMOXSC100	13701870
HEMOXSC100	13700243
HEMOXSC100	13703306
HEMOXSC100	13703612
HEMOXSC100	13703614
HEMOXSC100	13703615
HEMOXSC100	13703619
HEMOXSC100	13703351
HEMOXSC100	13703383
HEMOXSC100	13703387
HEMOXSC100	13703393
HEMOXSC100	13703399
HEMOXSC100	13703411
HEMOXSC100	13703414
HEMOXSC100	13703416
HEMOXSC100	13703418
HEMOXSC100	13703419
HEMOXSC100	13703420
HEMOXSC100	13703421
HEMOXSC100	13703422
HEMOXSC100	13710177
HEMOXSC100	13710180
HEMOXSC100	13710183
HEMOXSC100	13700903
HEMOXSC100	13700711
HEMOXSC100	13700708
HEMOXSC100	13700710
HEMOXSC100	13700919
HEMOXSC100	13701959
HEMOXSC100	13702220
HEMOXSC100	13703151
HEMOXSC100	13703234
HEMOXSC100	13703368
HEMOXSC100	13703369
HEMOXSC100	13705060
HEMOXSC100	13700280
HEMOXSC100	13700282
HEMOXSC100	13700291
HEMOXSC100	13703059
HEMOXSC100	13703104
HEMOXSC100	13703109
HEMOXSC100	13710452
HEMOXSC100	13704270

HEMOXSC100	13705245
HEMOXSC100	13705256
HEMOXSC100	13705293
HEMOXSC100	13705315
HEMOXSC100	13705322
HEMOXSC100	13705324
HEMOXSC100	13705325
HEMOXSC100	13710318
HEMOXSC100	13702121
HEMOXSC100	13702124
HEMOXSC100	13702127
HEMOXSC100	13703117
HEMOXSC100	13703118
HEMOXSC100	13703119
HEMOXSC100	13703004
HEMOXSC100	13703024
HEMOXSC100	13703025
HEMOXSC100	13703027
HEMOXSC100	13703031
HEMOXSC100	13704076
HEMOXSC100	13704077
HEMOXSC100	13704078
HEMOXSC100	13704079
HEMOXSC100	13704081
HEMOXSC100	13704082
HEMOXSC100	13704083
HEMOXSC100	13704084
HEMOXSC100	13704085
HEMOXSC100	13704086
HEMOXSC100	13703536
HEMOXSC100	13704003
HEMOXSC100	13704101
HEMOXSC100	13704103
HEMOXSC100	13704104
HEMOXSC100	13704105
HEMOXSC100	13704107
HEMOXSC100	13704118
HEMOXSC100	13704119
HEMOXSC100	13704121
HEMOXSC100	13704122
HEMOXSC100	13704124
HEMOXSC100	13704126
HEMOXSC100	13704130
HEMOXSC100	13704131
HEMOXSC100	13704133
HEMOXSC100	13704135
HEMOXSC100	13704137
HEMOXSC100	13704138
HEMOXSC100	13704141
HEMOXSC100	13704142

HEMOXSC100	13704143
HEMOXSC100	13704145
HEMOXSC100	13704345
HEMOXSC100	13704353
HEMOXSC100	13704352
HEMOXSC100	13703032
HEMOXSC100	13702675
HEMOXSC100	13702680
HEMOXSC100	13702681
HEMOXSC100	13701170
HEMOXSC100	13702956
HEMOXSC100	13702957
HEMOXSC100	13702960
HEMOXSC100	13702962
HEMOXSC100	13703473
HEMOXSC100	13701685
HEMOXSC100	13710214
HEMOXSC100	13704048
HEMOXSC100	13704235
HEMOXSC100	13704416
HEMOXSC100	13701801
HEMOXSC100	13702012
HEMOXSC100	13702019
HEMOXSC100	13704707
HEMOXSC100	13704776
HEMOXSC100	13704961
HEMOXSC100	13704962
HEMOXSC100	13705210
HEMOXSC100	13705083
HEMOXSC100	13705102
HEMOXSC100	13705111
HEMOXSC100	13705119
HEMOXSC100	13704757
HEMOXSC100	13704758
HEMOXSC100	13701551
HEMOXSC100	13701557
HEMOXSC100	13700921
HEMOXSC100	13700924
HEMOXSC100	13700936
HEMOXSC100	13700939
HEMOXSC100	13700940
HEMOXSC100	13700941
HEMOXSC100	13700944
HEMOXSC100	13700946
HEMOXSC100	13700949
HEMOXSC100	13700951
HEMOXSC100	13700955
HEMOXSC100	13700957
HEMOXSC100	13700963
HEMOXSC100	13700971

HEMOXSC100	13700977
HEMOXSC100	13700978
HEMOXSC100	13701692
HEMOXSC100	13701693
HEMOXSC100	13701694
HEMOXSC100	13701695
HEMOXSC100	13701696
HEMOXSC100	13701698
HEMOXSC100	13701697
HEMOXSC100	13701699
HEMOXSC100	13702233
HEMOXSC100	13702254
HEMOXSC100	13702403
HEMOXSC100	13702402
HEMOXSC100	13710373
HEMOXSC100	13710378
HEMOXSC100	13710379
HEMOXSC100	13710380
HEMOXSC100	13710381
HEMOXSC100	13710382
HEMOXSC100	13710383
HEMOXSC100	13710385
HEMOXSC100	13710387
HEMOXSC100	13710390
HEMOXSC100	13710394
HEMOXSC100	13710398
HEMOXSC100	13710404
HEMOXSC100	13710408
HEMOXSC100	13710410
HEMOXSC100	13710411
HEMOXSC100	13710412
HEMOXSC100	13710413
HEMOXSC100	13710414
HEMOXSC100	13710415
HEMOXSC100	13710422
HEMOXSC100	13710423
HEMOXSC100	13710425
HEMOXSC100	13710427
HEMOXSC100	13710499
HEMOXSC100	13703120
HEMOXSC100	13703167
HEMOXSC100	13703478
HEMOXSC100	13703525
HEMOXSC100	13704272
HEMOXSC100	13704273
HEMOXSC100	13710042
HEMOXSC100	13710170
HEMOXSC100	13703023
HEMOXSC100	13703047
HEMOXSC100	13703022

HEMOXSC100	13703035
HEMOXSC100	13703048
HEMOXSC100	13703455
HEMOXSC100	13703467
HEMOXSC100	13700744
HEMOXSC100	13701080
HEMOXSC100	13701085
HEMOXSC100	13701087
HEMOXSC100	13701089
HEMOXSC100	13701093
HEMOXSC100	13701095
HEMOXSC100	13701101
HEMOXSC100	13701086
HEMOXSC100	13701090
HEMOXSC100	13701100
HEMOXSC100	13701088
HEMOXSC100	13701094
HEMOXSC100	13701102
HEMOXSC100	13701099
HEMOXSC100	13701558
HEMOXSC100	13701395
HEMOXSC100	13703822
HEMOXSC100	13710043
HEMOXSC100	13710523
HEMOXSC100	13701603
HEMOXSC100	13701604
HEMOXSC100	13701606
HEMOXSC100	13701612
HEMOXSC100	13701619
HEMOXSC100	13700382
HEMOXSC100	13700475
HEMOXSC100	13700582
HEMOXSC100	13700695
HEMOXSC100	13700773
HEMOXSC100	13700778
HEMOXSC100	13700832
HEMOXSC100	13701131
HEMOXSC100	13701152
HEMOXSC100	13701189
HEMOXSC100	13701193
HEMOXSC100	13701195
HEMOXSC100	13701198
HEMOXSC100	13701310
HEMOXSC100	13701362
HEMOXSC100	13701364
HEMOXSC100	13701367
HEMOXSC100	13701369
HEMOXSC100	13701372
HEMOXSC100	13701376
HEMOXSC100	13702974

HEMOXSC100	13702977
HEMOXSC100	13702978
HEMOXSC100	13702979
HEMOXSC100	13702360
HEMOXSC100	13702411
HEMOXSC100	13702406
HEMOXSC100	13702048
HEMOXSC100	13702081
HEMOXSC100	13701470
HEMOXSC100	13701501
HEMOXSC100	13701806
HEMOXSC100	13701827
HEMOXSC100	13701846
HEMOXSC100	13702843
HEMOXSC100	13702838
HEMOXSC100	13702842
HEMOXSC100	13702859
HEMOXSC100	13702862
HEMOXSC100	13702863
HEMOXSC100	13702864
HEMOXSC100	13702136
HEMOXSC100	13702139
HEMOXSC100	13705082
HEMOXSC100	13705087
HEMOXSC100	13705084
HEMOXSC100	13705086
HEMOXSC100	13702202
HEMOXSC100	13702206
HEMOXSC100	13702290
HEMOXSC100	13702302
HEMOXSC100	13702641
HEMOXSC100	13702650
HEMOXSC100	13703894
HEMOXSC100	13703901
HEMOXSC100	13703902
HEMOXSC100	13703903
HEMOXSC100	13703905
HEMOXSC100	13703907
HEMOXSC100	13703908
HEMOXSC100	13703891
HEMOXSC100	13703912
HEMOXSC100	13704611
HEMOXSC100	13704621
HEMOXSC100	13704624
HEMOXSC100	13704625
HEMOXSC100	13700883
HEMOXSC100	13701117
HEMOXSC100	13701259
HEMOXSC100	13701272
HEMOXSC100	13701280

HEMOXSC100	13701304
HEMOXSC100	13701540
HEMOXSC100	13701794
HEMOXSC100	13702073
HEMOXSC100	13703444
HEMOXSC100	13700592
HEMOXSC100	13700597
HEMOXSC100	13700631
HEMOXSC100	13700653
HEMOXSC100	13700267
HEMOXSC100	13702303
HEMOXSC100	13704314
HEMOXSC100	13710077
HEMOXSC100	13710078
HEMOXSC100	13710079
HEMOXSC100	13710080
HEMOXSC100	13710081
HEMOXSC100	13710082
HEMOXSC100	13710083
HEMOXSC100	13710087
HEMOXSC100	13710088
HEMOXSC100	13710090
HEMOXSC100	13710091
HEMOXSC100	13710095
HEMOXSC100	13710096
HEMOXSC100	13710097
HEMOXSC100	13710104
HEMOXSC100	13710019
HEMOXSC100	13702649
HEMOXSC100	13702655
HEMOXSC100	13702240
HEMOXSC100	13702277
HEMOXSC100	13702288
HEMOXSC100	13704165
HEMOXSC100	13701526
HEMOXSC100	13703810
HEMOXSC100	13710347
HEMOXSC100	13710357
HEMOXSC100	13710361
HEMOXSC100	13710366
HEMOXSC100	13710432
HEMOXSC100	13710350
HEMOXSC100	13710369
HEMOXSC100	13701929
HEMOXSC100	13701809
HEMOXSC100	13701830
HEMOXSC100	13701832
HEMOXSC100	13701835
HEMOXSC100	13701841
HEMOXSC100	13704960

HEMOXSC100	13704979
HEMOXSC100	13704983
HEMOXSC100	13704988
HEMOXSC100	13704990
HEMOXSC100	13704954
HEMOXSC100	13704947
HEMOXSC100	13700572
HEMOXSC100	13701547
HEMOXSC100	13702213
HEMOXSC100	13702225
HEMOXSC100	13702228
HEMOXSC100	13702230
HEMOXSC100	13702237
HEMOXSC100	13702239
HEMOXSC100	13702299
HEMOXSC100	13702300
HEMOXSC100	13702340
HEMOXSC100	13702478
HEMOXSC100	13702504
HEMOXSC100	13702505
HEMOXSC100	13702506
HEMOXSC100	13702507
HEMOXSC100	13702588
HEMOXSC100	13704172
HEMOXSC100	13710225
HEMOXSC100	13710226
HEMOXSC100	13710228
HEMOXSC100	13710229
HEMOXSC100	13710232
HEMOXSC100	13710233
HEMOXSC100	13710234
HEMOXSC100	13710236
HEMOXSC100	13710237
HEMOXSC100	13710238
HEMOXSC100	13710239
HEMOXSC100	13710241
HEMOXSC100	13710242
HEMOXSC100	13710243
HEMOXSC100	13710244
HEMOXSC100	13710246
HEMOXSC100	13710247
HEMOXSC100	13710248
HEMOXSC100	13710252
HEMOXSC100	13710255
HEMOXSC100	13710032
HEMOXSC100	13701350
HEMOXSC100	13701352
HEMOXSC100	13700225
HEMOXSC100	13700367
HEMOXSC100	13700398

HEMOXSC100	13700550
HEMOXSC100	13702352
HEMOXSC100	13702353
HEMOXSC100	13702355
HEMOXSC100	13702356
HEMOXSC100	13700721
HEMOXSC100	13700891
HEMOXSC100	13701442
HEMOXSC100	13701462
HEMOXSC100	13702849
HEMOXSC100	13705136
HEMOXSC100	13705138
HEMOXSC100	13705135
HEMOXSC100	13705143
HEMOXSC100	13705140
HEMOXSC100	13705146
HEMOXSC100	13705147
HEMOXSC100	13703046
HEMOXSC100	13703049
HEMOXSC100	13703050
HEMOXSC100	13703053
HEMOXSC100	13704036
HEMOXSC100	13704046
HEMOXSC100	13704717
HEMOXSC100	13704888
HEMOXSC100	13704886
HEMOXSC100	13704889
HEMOXSC100	13704885
HEMOXSC100	13704890
HEMOXSC100	13700906
HEMOXSC100	13701700
HEMOXSC100	13704700
HEMOXSC100	13704704
HEMOXSC100	13705241
HEMOXSC100	13705242
HEMOXSC100	13701701
HEMOXSC100	13701849
HEMOXSC100	13701853
HEMOXSC100	13701854
HEMOXSC100	13701855
HEMOXSC100	13701857
HEMOXSC100	13701864
HEMOXSC100	13702927
HEMOXSC100	13704428
HEMOXSC100	13704348
HEMOXSC100	13704417
HEMOXSC100	13704420
HEMOXSC100	13704435
HEMOXSC100	13704436
HEMOXSC100	13704439

HEMOXSC100	13704447
HEMOXSC100	13704465
HEMOXSC100	13704466
HEMOXSC100	13704473
HEMOXSC100	13704492
HEMOXSC100	13704498
HEMOXSC100	13704553
HEMOXSC100	13704555
HEMOXSC100	13704557
HEMOXSC100	13704560
HEMOXSC100	13704571
HEMOXSC100	13704573
HEMOXSC100	13704575
HEMOXSC100	13704997
HEMOXSC100	13704998
HEMOXSC100	13704999
HEMOXSC100	13705000
HEMOXSC100	13704993
HEMOXSC100	13704995
HEMOXSC100	13710018
HEMOXSC100	13710000
HEMOXSC100	13710003
HEMOXSC100	13710005
HEMOXSC100	13710007
HEMOXSC100	13710008
HEMOXSC100	13710009
HEMOXSC100	13710011
HEMOXSC100	13701916
HEMOXSC100	13701918
HEMOXSC100	13701920
HEMOXSC100	13703279
HEMOXSC100	13701384
HEMOXSC100	13701385
HEMOXSC100	13701387
HEMOXSC100	13701389
HEMOXSC100	13701597
HEMOXSC100	13701601
HEMOXSC100	13701600
HEMOXSC100	13701598
HEMOXSC100	13704315
HEMOXSC100	13702853
HEMOXSC100	13702820
HEMOXSC100	13702834
HEMOXSC100	13702835
HEMOXSC100	13702837
HEMOXSC100	13702839
HEMOXSC100	13702841
HEMOXSC100	13702845
HEMOXSC100	13702846
HEMOXSC100	13704068

HEMOXSC100	13704070
HEMOXSC100	13704071
HEMOXSC100	13704072
HEMOXSC100	13704074
HEMOXSC100	13704075
HEMOXSC100	13700195
HEMOXSC100	13700196
HEMOXSC100	13700197
HEMOXSC100	13700198
HEMOXSC100	13705035
HEMOXSC100	13705044
HEMOXSC100	13705033
HEMOXSC100	13704060
HEMOXSC100	13702241
HEMOXSC100	13702232
HEMOXSC100	13703463
HEMOXSC100	13704667
HEMOXSC100	13704660
HEMOXSC100	13704664
HEMOXSC100	13703055
HEMOXSC100	13703821
HEMOXSC100	13705156
HEMOXSC100	13705158
HEMOXSC100	13704841
HEMOXSC100	13702372
HEMOXSC100	13702388
HEMOXSC100	13701199
HEMOXSC100	13701811
HEMOXSC100	13702007
HEMOXSC100	13702120
HEMOXSC100	13702275
HEMOXSC100	13702293
HEMOXSC100	13702335
HEMOXSC100	13702336
HEMOXSC100	13702430
HEMOXSC100	13702456
HEMOXSC100	13702462
HEMOXSC100	13702465
HEMOXSC100	13702467
HEMOXSC100	13702468
HEMOXSC100	13702479
HEMOXSC100	13702518
HEMOXSC100	13702545
HEMOXSC100	13702573
HEMOXSC100	13702578
HEMOXSC100	13702584
HEMOXSC100	13702585
HEMOXSC100	13702587
HEMOXSC100	13701055
HEMOXSC100	13701502

HEMOXSC100	13701848
HEMOXSC100	13701858
HEMOXSC100	13702122
HEMOXSC100	13710317
HEMOXSC100	13700285
HEMOXSC100	13700306
HEMOXSC100	13700348
HEMOXSC100	13700362
HEMOXSC100	13700387
HEMOXSC100	13700392
HEMOXSC100	13700397
HEMOXSC100	13700401
HEMOXSC100	13700483
HEMOXSC100	13700485
HEMOXSC100	13700489
HEMOXSC100	13700491
HEMOXSC100	13700492
HEMOXSC100	13700498
HEMOXSC100	13700502
HEMOXSC100	13700503
HEMOXSC100	13700620
HEMOXSC100	13702286
HEMOXSC100	13702204
HEMOXSC100	13703292
HEMOXSC100	13702189
HEMOXSC100	13702223
HEMOXSC100	13701175
HEMOXSC100	13700688
HEMOXSC100	13702032
HEMOXSC100	13710460
HEMOXSC100	13710506
HEMOXSC100	13710507
HEMOXSC100	13710521
HEMOXSC100	13710498
HEMOXSC100	13700424
HEMOXSC100	13701327
HEMOXSC100	13701336
HEMOXSC100	13701164
HEMOXSC100	13701278
HEMOXSC100	13710493
HEMOXSC100	13710494
HEMOXSC100	13710033
HEMOXSC100	13702828
HEMOXSC100	13704444
HEMOXSC100	13704448
HEMOXSC100	13704486
HEMOXSC100	13704412
HEMOXSC100	13704321
HEMOXSC100	13703705
HEMOXSC100	13703709

HEMOXSC100	13703711
HEMOXSC100	13703712
HEMOXSC100	13703714
HEMOXSC100	13703730
HEMOXSC100	13703735
HEMOXSC100	13703741
HEMOXSC100	13703706
HEMOXSC100	13703720
HEMOXSC100	13703722
HEMOXSC100	13703732
HEMOXSC100	13703733
HEMOXSC100	13703736
HEMOXSC100	13703739
HEMOXSC100	13703740
HEMOXSC100	13705079
HEMOXSC100	13705080
HEMOXSC100	13710462
HEMOXSC100	13710482
HEMOXSC100	13701300
HEMOXSC100	13701595
HEMOXSC100	13701599
HEMOXSC100	13701042
HEMOXSC100	13701045
HEMOXSC100	13702071
HEMOXSC100	13703305
HEMOXSC100	13710530
HEMOXSC100	13710531
HEMOXSC100	13710533
HEMOXSC100	13710550
HEMOXSC100	13710554
HEMOXSC100	13710555
HEMOXSC100	13701625
HEMOXSC100	13701684
HEMOXSC100	13703744
HEMOXSC100	13703745
HEMOXSC100	13704710
HEMOXSC100	13704732
HEMOXSC100	13704720
HEMOXSC100	13704733
HEMOXSC100	13704734
HEMOXSC100	13700541
HEMOXSC100	13700544
HEMOXSC100	13700558
HEMOXSC100	13700589
HEMOXSC100	13700591
HEMOXSC100	13700553
HEMOXSC100	13700584
HEMOXSC100	13700681
HEMOXSC100	13700618
HEMOXSC100	13700634

HEMOXSC100	13700637
HEMOXSC100	13700654
HEMOXSC100	13700655
HEMOXSC100	13700679
HEMOXSC100	13700522
HEMOXSC100	13700649
HEMOXSC100	13700623
HEMOXSC100	13700625
HEMOXSC100	13700640
HEMOXSC100	13700645
HEMOXSC100	13700646
HEMOXSC100	13700648
HEMOXSC100	13700652
HEMOXSC100	13700656
HEMOXSC100	13700661
HEMOXSC100	13700663
HEMOXSC100	13701184
HEMOXSC100	13703445
HEMOXSC100	13710035
HEMOXSC100	13710037
HEMOXSC100	13710038
HEMOXSC100	13710046
HEMOXSC100	13710053
HEMOXSC100	13710054
HEMOXSC100	13710055
HEMOXSC100	13710059
HEMOXSC100	13710060
HEMOXSC100	13710067
HEMOXSC100	13710072
HEMOXSC100	13710073
HEMOXSC100	13704632
HEMOXSC100	13704642
HEMOXSC100	13704643
HEMOXSC100	13704645
HEMOXSC100	13704649
HEMOXSC100	13704650
HEMOXSC100	13704673
HEMOXSC100	13704661
HEMOXSC100	13704670
HEMOXSC100	13704030
HEMOXSC100	13705063
HEMOXSC100	13705065
HEMOXSC100	13705069
HEMOXSC100	13705078
HEMOXSC100	13705061
HEMOXSC100	13705073
HEMOXSC100	13705074
HEMOXSC100	13705077
HEMOXSC100	13705081
HEMOXSC100	13705064

HEMOXSC100	13705066
HEMOXSC100	13710497
HEMOXSC100	13702404
HEMOXSC100	13704146
HEMOXSC100	13704147
HEMOXSC100	13704148
HEMOXSC100	13704149
HEMOXSC100	13704150
HEMOXSC100	13704152
HEMOXSC100	13704153
HEMOXSC100	13704155
HEMOXSC100	13704156
HEMOXSC100	13704157
HEMOXSC100	13704209
HEMOXSC100	13704211
HEMOXSC100	13704216
HEMOXSC100	13704246
HEMOXSC100	13704255
HEMOXSC100	13704257
HEMOXSC100	13704261
HEMOXSC100	13704264
HEMOXSC100	13704267
HEMOXSC100	13704207
HEMOXSC100	13704218
HEMOXSC100	13704221
HEMOXSC100	13704199
HEMOXSC100	13704200
HEMOXSC100	13704210
HEMOXSC100	13704212
HEMOXSC100	13704215
HEMOXSC100	13704220
HEMOXSC100	13704185
HEMOXSC100	13704188
HEMOXSC100	13704192
HEMOXSC100	13704198
HEMOXSC100	13704213
HEMOXSC100	13704223
HEMOXSC100	13704225
HEMOXSC100	13703636
HEMOXSC100	13703653
HEMOXSC100	13703783
HEMOXSC100	13703794
HEMOXSC100	13703804
HEMOXSC100	13703808
HEMOXSC100	13703890
HEMOXSC100	13703465
HEMOXSC100	13703468
HEMOXSC100	13703466
HEMOXSC100	13703715
HEMOXSC100	13703801

HEMOXSC100	13704456
HEMOXSC100	13710446
HEMOXSC100	13700383
HEMOXSC100	13700421
HEMOXSC100	13700546
HEMOXSC100	13700777
HEMOXSC100	13700783
HEMOXSC100	13700785
HEMOXSC100	13700787
HEMOXSC100	13700788
HEMOXSC100	13700790
HEMOXSC100	13700793
HEMOXSC100	13700794
HEMOXSC100	13702224
HEMOXSC100	13702231
HEMOXSC100	13702327
HEMOXSC100	13702063
HEMOXSC100	13702066
HEMOXSC100	13702077
HEMOXSC100	13702078
HEMOXSC100	13702949
HEMOXSC100	13702952
HEMOXSC100	13703014
HEMOXSC100	13703015
HEMOXSC100	13703016
HEMOXSC100	13703017
HEMOXSC100	13703019
HEMOXSC100	13703236
HEMOXSC100	13703348
HEMOXSC100	13703616
HEMOXSC100	13710039
HEMOXSC100	13710189
HEMOXSC100	13700812
HEMOXSC100	13701504
HEMOXSC100	13701498
HEMOXSC100	13703037
HEMOXSC100	13703038
HEMOXSC100	13703039
HEMOXSC100	13703040
HEMOXSC100	13703041
HEMOXSC100	13703042
HEMOXSC100	13703043
HEMOXSC100	13703044
HEMOXSC100	13703045
HEMOXSC100	13703434
HEMOXSC100	13705089
HEMOXSC100	13705090
HEMOXSC100	13705088
HEMOXSC100	13700801
HEMOXSC100	13703058

HEMOXSC100	13703835
HEMOXSC100	13704458
HEMOXSC100	13704462
HEMOXSC100	13701665
HEMOXSC100	13701666
HEMOXSC100	13703093
HEMOXSC100	13703094
HEMOXSC100	13703098
HEMOXSC100	13703099
HEMOXSC100	13703102
HEMOXSC100	13701396
HEMOXSC100	13701397
HEMOXSC100	13701398
HEMOXSC100	13701400
HEMOXSC100	13701401
HEMOXSC100	13701402
HEMOXSC100	13701405
HEMOXSC100	13701407
HEMOXSC100	13701410
HEMOXSC100	13702361
HEMOXSC100	13702405
HEMOXSC100	13702490
HEMOXSC100	13702493
HEMOXSC100	13702494
HEMOXSC100	13702495
HEMOXSC100	13702496
HEMOXSC100	13702497
HEMOXSC100	13702498
HEMOXSC100	13702501
HEMOXSC100	13702502
HEMOXSC100	13702503
HEMOXSC100	13703334
HEMOXSC100	13701576
HEMOXSC100	13703335
HEMOXSC100	13700188
HEMOXSC100	13700189
HEMOXSC100	13700190
HEMOXSC100	13704807
HEMOXSC100	13704811
HEMOXSC100	13704812
HEMOXSC100	13704814
HEMOXSC100	13704805
HEMOXSC100	13700713
HEMOXSC100	13702035
HEMOXSC100	13702138
HEMOXSC100	13702180
HEMOXSC100	13702301
HEMOXSC100	13702419
HEMOXSC100	13702421
HEMOXSC100	13702423

HEMOXSC100	13700233
HEMOXSC100	13700235
HEMOXSC100	13703033
HEMOXSC100	13703034
HEMOXSC100	13703036
HEMOXSC100	13703060
HEMOXSC100	13703063
HEMOXSC100	13700675
HEMOXSC100	13701581
HEMOXSC100	13701582
HEMOXSC100	13701041
HEMOXSC100	13701182
HEMOXSC100	13702531
HEMOXSC100	13703052
HEMOXSC100	13710258
HEMOXSC100	13710263
HEMOXSC100	13710269
HEMOXSC100	13710271
HEMOXSC100	13710272
HEMOXSC100	13703310
HEMOXSC100	13703311
HEMOXSC100	13703314
HEMOXSC100	13703320
HEMOXSC100	13703307
HEMOXSC100	13703308
HEMOXSC100	13703312
HEMOXSC100	13703317
HEMOXSC100	13703319
HEMOXSC100	13703350
HEMOXSC100	13701845
HEMOXSC100	13703856
HEMOXSC100	13703845
HEMOXSC100	13703849
HEMOXSC100	13703851
HEMOXSC100	13703852
HEMOXSC100	13703853
HEMOXSC100	13703860
HEMOXSC100	13703868
HEMOXSC100	13703878
HEMOXSC100	13703887
HEMOXSC100	13703904
HEMOXSC100	13703906
HEMOXSC100	13703913
HEMOXSC100	13703917
HEMOXSC100	13703848
HEMOXSC100	13703865
HEMOXSC100	13703866
HEMOXSC100	13703867
HEMOXSC100	13703881
HEMOXSC100	13703893

HEMOXSC100	13703898
HEMOXSC100	13703909
HEMOXSC100	13703915
HEMOXSC100	13704316
HEMOXSC100	13702832
HEMOXSC100	13702840
HEMOXSC100	13703082
HEMOXSC100	13703076
HEMOXSC100	13703083
HEMOXSC100	13710313
HEMOXSC100	13710314
HEMOXSC100	13710322
HEMOXSC100	13710323
HEMOXSC100	13710324
HEMOXSC100	13710328
HEMOXSC100	13700835
HEMOXSC100	13700837
HEMOXSC100	13700885
HEMOXSC100	13703366
HEMOXSC100	13704740
HEMOXSC100	13704908
HEMOXSC100	13704909
HEMOXSC100	13704910
HEMOXSC100	13710525
HEMOXSC100	13710419
HEMOXSC100	13710420
HEMOXSC100	13710433
HEMOXSC100	13704739
HEMOXSC100	13704749
HEMOXSC100	13704753
HEMOXSC100	13704756
HEMOXSC100	13705327
HEMOXSC100	13705328
HEMOXSC100	13705329
HEMOXSC100	13710001
HEMOXSC100	13710002
HEMOXSC100	13710004
HEMOXSC100	13710006
HEMOXSC100	13710010
HEMOXSC100	13710017
HEMOXSC100	13702619
HEMOXSC100	13702620
HEMOXSC100	13702622
HEMOXSC100	13702625
HEMOXSC100	13702627
HEMOXSC100	13702629
HEMOXSC100	13702632
HEMOXSC100	13703380
HEMOXSC100	13703388
HEMOXSC100	13703378

HEMOXSC100	13703379
HEMOXSC100	13703382
HEMOXSC100	13703384
HEMOXSC100	13703385
HEMOXSC100	13703386
HEMOXSC100	13704561
HEMOXSC100	13704562
HEMOXSC100	13704194
HEMOXSC100	13704306
HEMOXSC100	13704351
HEMOXSC100	13704186
HEMOXSC100	13704459
HEMOXSC100	13704554
HEMOXSC100	13702848
HEMOXSC100	13700700
HEMOXSC100	13700704
HEMOXSC100	13704637
HEMOXSC100	13704638
HEMOXSC100	13705159
HEMOXSC100	13701132
HEMOXSC100	13702784
HEMOXSC100	13702809
HEMOXSC100	13702814
HEMOXSC100	13702822
HEMOXSC100	13702833
HEMOXSC100	13704920
HEMOXSC100	13704921
HEMOXSC100	13704922
HEMOXSC100	13704924
HEMOXSC100	13702375
HEMOXSC100	13701941
HEMOXSC100	13701994
HEMOXSC100	13702229
HEMOXSC100	13702338
HEMOXSC100	13702357
HEMOXSC100	13702377
HEMOXSC100	13701802
HEMOXSC100	13700725
HEMOXSC100	13700728
HEMOXSC100	13700730
HEMOXSC100	13700733
HEMOXSC100	13700736
HEMOXSC100	13700738
HEMOXSC100	13700729
HEMOXSC100	13700735
HEMOXSC100	13701357
HEMOXSC100	13701363
HEMOXSC100	13701365
HEMOXSC100	13701370
HEMOXSC100	13701373

HEMOXSC100	13701383
HEMOXSC100	13701418
HEMOXSC100	13701419
HEMOXSC100	13701424
HEMOXSC100	13701330
HEMOXSC100	13701254
HEMOXSC100	13701016
HEMOXSC100	13701058
HEMOXSC100	13701111
HEMOXSC100	13701141
HEMOXSC100	13701417
HEMOXSC100	13701430
HEMOXSC100	13701577
HEMOXSC100	13701586
HEMOXSC100	13701587
HEMOXSC100	13701588
HEMOXSC100	13701589
HEMOXSC100	13701611
HEMOXSC100	13701580
HEMOXSC100	13701590
HEMOXSC100	13701574
HEMOXSC100	13701933
HEMOXSC100	13701936
HEMOXSC100	13701937
HEMOXSC100	13701939
HEMOXSC100	13703340
HEMOXSC100	13703341
HEMOXSC100	13703342
HEMOXSC100	13703343
HEMOXSC100	13703130
HEMOXSC100	13703133
HEMOXSC100	13703139
HEMOXSC100	13703152
HEMOXSC100	13703123
HEMOXSC100	13703124
HEMOXSC100	13703634
HEMOXSC100	13703814
HEMOXSC100	13704114
HEMOXSC100	13704425
HEMOXSC100	13705058
HEMOXSC100	13710113
HEMOXSC100	13710115
HEMOXSC100	13710129
HEMOXSC100	13710130
HEMOXSC100	13710131
HEMOXSC100	13702830
HEMOXSC100	13703623
HEMOXSC100	13703624
HEMOXSC100	13703633
HEMOXSC100	13703598

HEMOXSC100	13703622
HEMOXSC100	13703625
HEMOXSC100	13703627
HEMOXSC100	13703628
HEMOXSC100	13703629
HEMOXSC100	13704276
HEMOXSC100	13703685
HEMOXSC100	13705308
HEMOXSC100	13705312
HEMOXSC100	13705164
HEMOXSC100	13705171
HEMOXSC100	13704774
HEMOXSC100	13704633
HEMOXSC100	13703976
HEMOXSC100	13702993
HEMOXSC100	13702997
HEMOXSC100	13704895
HEMOXSC100	13704696
HEMOXSC100	13705001
HEMOXSC100	13705002
HEMOXSC100	13703484
HEMOXSC100	13703485
HEMOXSC100	13703488
HEMOXSC100	13703490
HEMOXSC100	13703491
HEMOXSC100	13703492
HEMOXSC100	13703496
HEMOXSC100	13703497
HEMOXSC100	13703499
HEMOXSC100	13703502
HEMOXSC100	13703504
HEMOXSC100	13703511
HEMOXSC100	13703513
HEMOXSC100	13703514
HEMOXSC100	13700791
HEMOXSC100	13700792
HEMOXSC100	13700781
HEMOXSC100	13700789
HEMOXSC100	13701487
HEMOXSC100	13701490
HEMOXSC100	13701491
HEMOXSC100	13703333
HEMOXSC100	13702935
HEMOXSC100	13702937
HEMOXSC100	13702315
HEMOXSC100	13703255
HEMOXSC100	13703256
HEMOXSC100	13703781
HEMOXSC100	13703782
HEMOXSC100	13702541

HEMOXSC100	13702542
HEMOXSC100	13702546
HEMOXSC100	13702548
HEMOXSC100	13705219
HEMOXSC100	13705222
HEMOXSC100	13705218
HEMOXSC100	13705220
HEMOXSC100	13702565
HEMOXSC100	13704002
HEMOXSC100	13704004
HEMOXSC100	13704005
HEMOXSC100	13704006
HEMOXSC100	13704007
HEMOXSC100	13704008
HEMOXSC100	13704009
HEMOXSC100	13704012
HEMOXSC100	13704013
HEMOXSC100	13704014
HEMOXSC100	13704015
HEMOXSC100	13704016
HEMOXSC100	13704018
HEMOXSC100	13704019
HEMOXSC100	13704169
HEMOXSC100	13704171
HEMOXSC100	13704180
HEMOXSC100	13704184
HEMOXSC100	13710273
HEMOXSC100	13710275
HEMOXSC100	13710281
HEMOXSC100	13710283
HEMOXSC100	13705268
HEMOXSC100	13704793
HEMOXSC100	13704799
HEMOXSC100	13704800
HEMOXSC100	13704802
HEMOXSC100	13704808
HEMOXSC100	13704810
HEMOXSC100	13703920
HEMOXSC100	13703923
HEMOXSC100	13703927
HEMOXSC100	13703939
HEMOXSC100	13703940
HEMOXSC100	13703957
HEMOXSC100	13703958
HEMOXSC100	13703966
HEMOXSC100	13703921
HEMOXSC100	13703924
HEMOXSC100	13703934
HEMOXSC100	13703935
HEMOXSC100	13703937

HEMOXSC100	13703941
HEMOXSC100	13703965
HEMOXSC100	13703811
HEMOXSC100	13703827
HEMOXSC100	13703829
HEMOXSC100	13703831
HEMOXSC100	13703834
HEMOXSC100	13703944
HEMOXSC100	13703946
HEMOXSC100	13703947
HEMOXSC100	13703948
HEMOXSC100	13703949
HEMOXSC100	13703950
HEMOXSC100	13703951
HEMOXSC100	13703953
HEMOXSC100	13703954
HEMOXSC100	13703964
HEMOXSC100	13703641
HEMOXSC100	13703651
HEMOXSC100	13703655
HEMOXSC100	13703666
HEMOXSC100	13703679
HEMOXSC100	13703700
HEMOXSC100	13703728
HEMOXSC100	13703597
HEMOXSC100	13703713
HEMOXSC100	13703717
HEMOXSC100	13703731
HEMOXSC100	13703816
HEMOXSC100	13703579
HEMOXSC100	13703726
HEMOXSC100	13703729
HEMOXSC100	13703773
HEMOXSC100	13703981
HEMOXSC100	13703996
HEMOXSC100	13704011
HEMOXSC100	13704073
HEMOXSC100	13704098
HEMOXSC100	13703533
HEMOXSC100	13703581
HEMOXSC100	13703610
HEMOXSC100	13703644
HEMOXSC100	13703696
HEMOXSC100	13703776
HEMOXSC100	13703790
HEMOXSC100	13703793
HEMOXSC100	13703828
HEMOXSC100	13703859
HEMOXSC100	13703862
HEMOXSC100	13703900

HEMOXSC100	13703926
HEMOXSC100	13703931
HEMOXSC100	13703932
HEMOXSC100	13703933
HEMOXSC100	13703938
HEMOXSC100	13703943
HEMOXSC100	13703952
HEMOXSC100	13703960
HEMOXSC100	13703963
HEMOXSC100	13703967
HEMOXSC100	13703968
HEMOXSC100	13703969
HEMOXSC100	13703970
HEMOXSC100	13703971
HEMOXSC100	13703973
HEMOXSC100	13703974
HEMOXSC100	13703975
HEMOXSC100	13703977
HEMOXSC100	13703978
HEMOXSC100	13703979
HEMOXSC100	13710505
HEMOXSC100	13710509
HEMOXSC100	13710510
HEMOXSC100	13710511
HEMOXSC100	13701812
HEMOXSC100	13701856
HEMOXSC100	13702010
HEMOXSC100	13702097
HEMOXSC100	13703888
HEMOXSC100	13703889
HEMOXSC100	13702661
HEMOXSC100	13702662
HEMOXSC100	13702663
HEMOXSC100	13702666
HEMOXSC100	13702668
HEMOXSC100	13702669
HEMOXSC100	13702670
HEMOXSC100	13702674
HEMOXSC100	13702245
HEMOXSC100	13702262
HEMOXSC100	13702263
HEMOXSC100	13702264
HEMOXSC100	13702294
HEMOXSC100	13702295
HEMOXSC100	13702308
HEMOXSC100	13702309
HEMOXSC100	13702310
HEMOXSC100	13702311
HEMOXSC100	13702276
HEMOXSC100	13702278

HEMOXSC100	13702280
HEMOXSC100	13702285
HEMOXSC100	13702312
HEMOXSC100	13703956
HEMOXSC100	13703955
HEMOXSC100	13703446
HEMOXSC100	13701615
HEMOXSC100	13701622
HEMOXSC100	13701621
HEMOXSC100	13700473
HEMOXSC100	13701957
HEMOXSC100	13702104
HEMOXSC100	13701851
HEMOXSC100	13701911
HEMOXSC100	13701931
HEMOXSC100	13701969
HEMOXSC100	13701987
HEMOXSC100	13702108
HEMOXSC100	13702112
HEMOXSC100	13702119
HEMOXSC100	13702131
HEMOXSC100	13702133
HEMOXSC100	13701805
HEMOXSC100	13701814
HEMOXSC100	13701850
HEMOXSC100	13703275
HEMOXSC100	13703277
HEMOXSC100	13703278
HEMOXSC100	13703404
HEMOXSC100	13703405
HEMOXSC100	13703564
HEMOXSC100	13703567
HEMOXSC100	13703443
HEMOXSC100	13704201
HEMOXSC100	13703707
HEMOXSC100	13703774
HEMOXSC100	13701295
HEMOXSC100	13701297
HEMOXSC100	13705263
HEMOXSC100	13705264
HEMOXSC100	13700915
HEMOXSC100	13700566
HEMOXSC100	13703983
HEMOXSC100	13703984
HEMOXSC100	13703985
HEMOXSC100	13703986
HEMOXSC100	13703989
HEMOXSC100	13703991
HEMOXSC100	13703000
HEMOXSC100	13703006

HEMOXSC100	13703011
HEMOXSC100	13702996
HEMOXSC100	13703009
HEMOXSC100	13703001
HEMOXSC100	13703018
HEMOXSC100	13704065
HEMOXSC100	13704044
HEMOXSC100	13704127
HEMOXSC100	13703140
HEMOXSC100	13703148
HEMOXSC100	13704029
HEMOXSC100	13704039
HEMOXSC100	13704049
HEMOXSC100	13704289
HEMOXSC100	13700662
HEMOXSC100	13701917
HEMOXSC100	13700726
HEMOXSC100	13700741
HEMOXSC100	13700784
HEMOXSC100	13700722
HEMOXSC100	13700734
HEMOXSC100	13700746
HEMOXSC100	13700479
HEMOXSC100	13700482
HEMOXSC100	13700499
HEMOXSC100	13700505
HEMOXSC100	13700406
HEMOXSC100	13700495
HEMOXSC100	13700515
HEMOXSC100	13700542
HEMOXSC100	13700391
HEMOXSC100	13700403
HEMOXSC100	13700407
HEMOXSC100	13700434
HEMOXSC100	13700444
HEMOXSC100	13700447
HEMOXSC100	13700449
HEMOXSC100	13700469
HEMOXSC100	13700524
HEMOXSC100	13700529
HEMOXSC100	13700538
HEMOXSC100	13700545
HEMOXSC100	13700548
HEMOXSC100	13700549
HEMOXSC100	13700552
HEMOXSC100	13700554
HEMOXSC100	13700561
HEMOXSC100	13700562
HEMOXSC100	13700574
HEMOXSC100	13700575

HEMOXSC100	13700578
HEMOXSC100	13701875
HEMOXSC100	13701877
HEMOXSC100	13701881
HEMOXSC100	13701882
HEMOXSC100	13701874
HEMOXSC100	13704785
HEMOXSC100	13702963
HEMOXSC100	13702967
HEMOXSC100	13702973
HEMOXSC100	13700241
HEMOXSC100	13700246
HEMOXSC100	13700247
HEMOXSC100	13700248
HEMOXSC100	13700804
HEMOXSC100	13702831
HEMOXSC100	13703213
HEMOXSC100	13703225
HEMOXSC100	13703235
HEMOXSC100	13703250
HEMOXSC100	13703187
HEMOXSC100	13703193
HEMOXSC100	13703194
HEMOXSC100	13703215
HEMOXSC100	13703226
HEMOXSC100	13703227
HEMOXSC100	13703228
HEMOXSC100	13703229
HEMOXSC100	13703230
HEMOXSC100	13703231
HEMOXSC100	13703232
HEMOXSC100	13703233
HEMOXSC100	13703237
HEMOXSC100	13703238
HEMOXSC100	13703239
HEMOXSC100	13703244
HEMOXSC100	13703246
HEMOXSC100	13703248
HEMOXSC100	13703251
HEMOXSC100	13703252
HEMOXSC100	13703253
HEMOXSC100	13703241
HEMOXSC100	13703245
HEMOXSC100	13703265
HEMOXSC100	13703270
HEMOXSC100	13703276
HEMOXSC100	13703280
HEMOXSC100	13703282
HEMOXSC100	13703285
HEMOXSC100	13703296

HEMOXSC100	13703300
HEMOXSC100	13703301
HEMOXSC100	13703223
HEMOXSC100	13710184
HEMOXSC100	13710185
HEMOXSC100	13703668
HEMOXSC100	13702134
HEMOXSC100	13702137
HEMOXSC100	13700413
HEMOXSC100	13703367
HEMOXSC100	13705054
HEMOXSC100	13704906
HEMOXSC100	13702088
HEMOXSC100	13704490
HEMOXSC100	13704489
HEMOXSC100	13704911
HEMOXSC100	13710512
HEMOXSC100	13710520
HEMOXSC100	13710529
HEMOXSC100	13710544
HEMOXSC100	13701032
HEMOXSC100	13703103
HEMOXSC100	13704329
HEMOXSC100	13704440
HEMOXSC100	13704455
HEMOXSC100	13704432
HEMOXSC100	13704431
HEMOXSC100	13704442
HEMOXSC100	13704452
HEMOXSC100	13704421
HEMOXSC100	13704423
HEMOXSC100	13704438
HEMOXSC100	13704445
HEMOXSC100	13704247
HEMOXSC100	13700650
HEMOXSC100	13705031
HEMOXSC100	13705017
HEMOXSC100	13704918
HEMOXSC100	13704919
HEMOXSC100	13702020
HEMOXSC100	13702022
HEMOXSC100	13702023
HEMOXSC100	13700138
HEMOXSC100	13700304
HEMOXSC100	13700395
HEMOXSC100	13700258
HEMOXSC100	13700431
HEMOXSC100	13700692
HEMOXSC100	13700366
HEMOXSC100	13700902

HEMOXSC100	13700875
HEMOXSC100	13700884
HEMOXSC100	13701351
HEMOXSC100	13701647
HEMOXSC100	13703347
HEMOXSC100	13700678
HEMOXSC100	13700702
HEMOXSC100	13700690
HEMOXSC100	13705071
HEMOXSC100	13701427
HEMOXSC100	13701428
HEMOXSC100	13704340
HEMOXSC100	13704344
HEMOXSC100	13704305
HEMOXSC100	13701222
HEMOXSC100	13701161
HEMOXSC100	13701162
HEMOXSC100	13701156
HEMOXSC100	13701163
HEMOXSC100	13701166
HEMOXSC100	13701155
HEMOXSC100	13701157
HEMOXSC100	13701218
HEMOXSC100	13702161
HEMOXSC100	13702980
HEMOXSC100	13703980
HEMOXSC100	13704765
HEMOXSC100	13703181
HEMOXSC100	13701480
HEMOXSC100	13701488
HEMOXSC100	13705047
HEMOXSC100	13702089
HEMOXSC100	13702091
HEMOXSC100	13702092
HEMOXSC100	13705027
HEMOXSC100	13701454
HEMOXSC100	13704277
HEMOXSC100	13710208
HEMOXSC100	13710211
HEMOXSC100	13710212
HEMOXSC100	13704891
HEMOXSC100	13704893
HEMOXSC100	13700066
HEMOXSC100	13700507
HEMOXSC100	13700508
HEMOXSC100	13700211
HEMOXSC100	13703179
HEMOXSC100	13701636
HEMOXSC100	13700376
HEMOXSC100	13702331

HEMOXSC100	13702512
HEMOXSC100	13702513
HEMOXSC100	13702516
HEMOXSC100	13700286
HEMOXSC100	13700218
HEMOXSC100	13700202
HEMOXSC100	13700521
HEMOXSC100	13700496
HEMOXSC100	13700624
HEMOXSC100	13700208
HEMOXSC100	13700671
HEMOXSC100	13700462
HEMOXSC100	13700223
HEMOXSC100	13703911
HEMOXSC100	13703998
HEMOXSC100	13700716
HEMOXSC100	13700715
HEMOXSC100	13701932
HEMOXSC100	13701934
HEMOXSC100	13701938
HEMOXSC100	13701942
HEMOXSC100	13701943
HEMOXSC100	13701944
HEMOXSC100	13701945
HEMOXSC100	13701947
HEMOXSC100	13701949
HEMOXSC100	13701950
HEMOXSC100	13701952
HEMOXSC100	13705273
HEMOXSC100	13705277
HEMOXSC100	13705278
HEMOXSC100	13705282
HEMOXSC100	13705284
HEMOXSC100	13705272
HEMOXSC100	13705285
HEMOXSC100	13702938
HEMOXSC100	13702491
HEMOXSC100	13702500
HEMOXSC100	13700295
HEMOXSC100	13700752
HEMOXSC100	13704763
HEMOXSC100	13704735
HEMOXSC100	13700455
HEMOXSC100	13700452
HEMOXSC100	13700448
HEMOXSC100	13701610
HEMOXSC100	13702549
HEMOXSC100	13701573
HEMOXSC100	13701607
HEMOXSC100	13701624

HEMOXSC100	13701807
HEMOXSC100	13701560
HEMOXSC100	13701567
HEMOXSC100	13701572
HEMOXSC100	13701605
HEMOXSC100	13701609
HEMOXSC100	13701620
HEMOXSC100	13701640
HEMOXSC100	13701186
HEMOXSC100	13700673
HEMOXSC100	13701515
HEMOXSC100	13701509
HEMOXSC100	13700707
HEMOXSC100	13704426
HEMOXSC100	13704415
HEMOXSC100	13700364
HEMOXSC100	13700228
HEMOXSC100	13700254
HEMOXSC100	13701181
HEMOXSC100	13701349
HEMOXSC100	13702070
HEMOXSC100	13702334
HEMOXSC100	13702337
HEMOXSC100	13702321
HEMOXSC100	13702330
HEMOXSC100	13702333
HEMOXSC100	13703107
HEMOXSC100	13710120
HEMOXSC100	13710504
HEMOXSC100	13700239
HEMOXSC100	13703150
HEMOXSC100	13703153
HEMOXSC100	13703154
HEMOXSC100	13703157
HEMOXSC100	13703158
HEMOXSC100	13703160
HEMOXSC100	13705051
HEMOXSC100	13700381
HEMOXSC100	13700385
HEMOXSC100	13700994
HEMOXSC100	13701003
HEMOXSC100	13702093
HEMOXSC100	13701912
HEMOXSC100	13703172
HEMOXSC100	13704478
HEMOXSC100	13704484
HEMOXSC100	13700361
HEMOXSC100	13700764
HEMOXSC100	13700438
HEMOXSC100	13700551

HEMOXSC100	13700560
HEMOXSC100	13700568
HEMOXSC100	13701178
HEMOXSC100	13701635
HEMOXSC100	13701634
HEMOXSC100	13701636
HEMOXSC100	13701638
HEMOXSC100	13701637
HEMOXSC100	13701629
HEMOXSC100	13703454
HEMOXSC100	13703452
HEMOXSC100	13704097
HEMOXSC100	13702777
HEMOXSC100	13702808
HEMOXSC100	13704414
HEMOXSC100	13704424
HEMOXSC100	13704836
HEMOXSC100	13704838
HEMOXSC100	13701450
HEMOXSC100	13701644
HEMOXSC100	13701466
HEMOXSC100	13703195
HEMOXSC100	13703196
HEMOXSC100	13703548
HEMOXSC100	13702384
HEMOXSC100	13702387
HEMOXSC100	13702389
HEMOXSC100	13702393
HEMOXSC100	13702394
HEMOXSC100	13702376
HEMOXSC100	13702392
HEMOXSC100	13702395
HEMOXSC100	13703249
HEMOXSC100	13703159
HEMOXSC100	13704752
HEMOXSC100	13701569
HEMOXSC100	13704857
HEMOXSC100	13704860
HEMOXSC100	13704854
HEMOXSC100	13704859
HEMOXSC100	13701432
HEMOXSC100	13701433
HEMOXSC100	13703692
HEMOXSC100	13700072
HEMOXSC100	13700355
HEMOXSC100	13700369
HEMOXSC100	13704300
HEMOXSC100	13704766
HEMOXSC100	13700464
HEMOXSC100	13702422

HEMOXSC100	13700180
HEMOXSC100	13700185
HEMOXSC100	13700186
HEMOXSC100	13700821
HEMOXSC100	13700175
HEMOXSC100	13700174
HEMOXSC100	13700177
HEMOXSC100	13704842
HEMOXSC100	13704847
HEMOXSC100	13704849
HEMOXSC100	13704852
HEMOXSC100	13704856
HEMOXSC100	13704864
HEMOXSC100	13704835
HEMOXSC100	13704855
HEMOXSC100	13704884
HEMOXSC100	13704892
HEMOXSC100	13704900
HEMOXSC100	13704915
HEMOXSC100	13704928
HEMOXSC100	13704936
HEMOXSC100	13700201
HEMOXSC100	13700205
HEMOXSC100	13700236
HEMOXSC100	13700265
HEMOXSC100	13700270
HEMOXSC100	13700271
HEMOXSC100	13700240
HEMOXSC100	13700215
HEMOXSC100	13703125
HEMOXSC100	13702763
HEMOXSC100	13704338
HEMOXSC100	13703558
HEMOXSC100	13701568
HEMOXSC100	13701575
HEMOXSC100	13700712
HEMOXSC100	13700717
HEMOXSC100	13701521
HEMOXSC100	13701513
HEMOXSC100	13701063
HEMOXSC100	13702552
HEMOXSC100	13702553
HEMOXSC100	13704190
HEMOXSC100	13701211
HEMOXSC100	13701052
HEMOXSC100	13700356
HEMOXSC100	13700384
HEMOXSC100	13700619
HEMOXSC100	13700621
HEMOXSC100	13702397

HEMOXSC100	13702399
HEMOXSC100	13702400
HEMOXSC100	13702398
HEMOXSC100	13700372
HEMOXSC100	13704843
HEMOXSC100	13704846
HEMOXSC100	13704851
HEMOXSC100	13704844
HEMOXSC100	13704848
HEMOXSC100	13704850
HEMOXSC100	13700985
HEMOXSC100	13701017
HEMOXSC100	13703725
HEMOXSC100	13700277
HEMOXSC100	13701474
HEMOXSC100	13701475
HEMOXSC100	13701476
HEMOXSC100	13701477
HEMOXSC100	13704556
HEMOXSC100	13704558
HEMOXSC100	13704567
HEMOXSC100	13704568
HEMOXSC100	13704569
HEMOXSC100	13704572
HEMOXSC100	13703646
HEMOXSC100	13703672
HEMOXSC100	13703686
HEMOXSC100	13703693
HEMOXSC100	13704691
HEMOXSC100	13704692
HEMOXSC100	13704695
HEMOXSC100	13704697
HEMOXSC100	13704699
HEMOXSC100	13704701
HEMOXSC100	13704703
HEMOXSC100	13704705
HEMOXSC100	13704706
HEMOXSC100	13704708
HEMOXSC100	13704709
HEMOXSC100	13701168
HEMOXSC100	13701176
HEMOXSC100	13701179
HEMOXSC100	13701167
HEMOXSC100	13701173
HEMOXSC100	13701174
HEMOXSC100	13701183
HEMOXSC100	13701171
HEMOXSC100	13701382
HEMOXSC100	13701019
HEMOXSC100	13703742

HEMOXSC100	13702769
HEMOXSC100	13702770
HEMOXSC100	13704951
HEMOXSC100	13704952
HEMOXSC100	13701371
HEMOXSC100	13704948
HEMOXSC100	13700685
HEMOXSC100	13704834
HEMOXSC100	13704853
HEMOXSC100	13701021
HEMOXSC100	13701630
HEMOXSC100	13701219
HEMOXSC100	13701223
HEMOXSC100	13701999
HEMOXSC100	13702001
HEMOXSC100	13705201
HEMOXSC100	13705202
HEMOXSC100	13703056
HEMOXSC100	13703699
HEMOXSC100	13703703
HEMOXSC100	13703676
HEMOXSC100	13703146
HEMOXSC100	13702190
HEMOXSC100	13702191
HEMOXSC100	13702192
HEMOXSC100	13702193
HEMOXSC100	13702195
HEMOXSC100	13702187
HEMOXSC100	13702188
HEMOXSC100	13702208
HEMOXSC100	13702209
HEMOXSC100	13702215
HEMOXSC100	13702214
HEMOXSC100	13701614
HEMOXSC100	13701618
HEMOXSC100	13704940
HEMOXSC100	13704941
HEMOXSC100	13704944
HEMOXSC100	13704470
HEMOXSC100	13703687
HEMOXSC100	13703737
HEMOXSC100	13703313
HEMOXSC100	13703664
HEMOXSC100	13703990
HEMOXSC100	13704093
HEMOXSC100	13704227
HEMOXSC100	13704346
HEMOXSC100	13704350
HEMOXSC100	13704360
HEMOXSC100	13700990

HEMOXSC100	13710040
HEMOXSC100	13703092
HEMOXSC100	13703070
HEMOXSC100	13703071
HEMOXSC100	13703077
HEMOXSC100	13703078
HEMOXSC100	13703081
HEMOXSC100	13703088
HEMOXSC100	13710346
HEMOXSC100	13704196
HEMOXSC100	13700441
HEMOXSC100	13703589
HEMOXSC100	13703590
HEMOXSC100	13703591
HEMOXSC100	13703593
HEMOXSC100	13703599
HEMOXSC100	13700539
HEMOXSC100	13704128
HEMOXSC100	13702825
HEMOXSC100	13704731
HEMOXSC100	13701467
HEMOXSC100	13701532
HEMOXSC100	13701556
HEMOXSC100	13701564
HEMOXSC100	13701578
HEMOXSC100	13701469
HEMOXSC100	13701473
HEMOXSC100	13701516
HEMOXSC100	13701520
HEMOXSC100	13701544
HEMOXSC100	13701550
HEMOXSC100	13701570
HEMOXSC100	13701591
HEMOXSC100	13701596
HEMOXSC100	13701008
HEMOXSC100	13704796
HEMOXSC100	13704795
HEMOXSC100	13704032
HEMOXSC100	13704863
HEMOXSC100	13701064
HEMOXSC100	13704956
HEMOXSC100	13704957
HEMOXSC100	13704965
HEMOXSC100	13704963
HEMOXSC100	13705157
HEMOXSC100	13704826
HEMOXSC100	13703005
HEMOXSC100	13703010
HEMOXSC100	13703079
HEMOXSC100	13704839

HEMOXSC100	13704877
HEMOXSC100	13704932
HEMOXSC100	13704779
HEMOXSC100	13704966
HEMOXSC100	13705265
HEMOXSC100	13705262
HEMOXSC100	13705275
HEMOXSC100	13705286
HEMOXSC100	13700992
HEMOXSC100	13710445
HEMOXSC100	13704942
HEMOXSC100	13704982
HEMOXSC100	13705016
HEMOXSC100	13702520
HEMOXSC100	13702570
HEMOXSC100	13702582
HEMOXSC100	13702626
HEMOXSC100	13702628
HEMOXSC100	13702644
HEMOXSC100	13702646
HEMOXSC100	13702627
HEMOXSC100	13705176
HEMOXSC100	13703752
HEMOXSC100	13700007
HEMOXSC100	13700018
HEMOXSC100	13702999
HEMOXSC100	13702950
HEMOXSC100	13701200
HEMOXSC100	13702280
HEMOXSC100	13701020
HEMOXSC100	13702650
HEMOXSC100	13702660
HEMOXSC100	13702674
HEMOXSC100	13703115
HEMOXSC100	13703121
HEMOXSC100	13703131
HEMOXSC100	13702426
HEMOXSC100	13703437
HEMOXSC100	13703438
HEMOXSC100	13701974
HEMOXSC100	13705170
HEMOXSC100	13703551
HEMOXSC100	13701356
HEMOXSC100	13700547
HEMOXSC100	13700563
HEMOXSC100	13700565
HEMOXSC100	13700567
HEMOXSC100	13700573
HEMOXSC100	13700577
HEMOXSC100	13700580
HEMOXSC100	13700581
HEMOXSC100	13700586
HEMOXSC100	13700587
HEMOXSC100	13700593
HEMOXSC100	13700596
HEMOXSC100	13705184
HEMOXSC100	13704063
HEMOXSC100	13701654
HEMOXSC100	13701655
HEMOXSC100	13701658

HEMOXSC100	13701660
HEMOXSC100	13701667
HEMOXSC100	13701012
HEMOXSC100	13700050
HEMOXSC100	13701000
HEMOXSC100	13703770
HEMOXSC100	13700512
HEMOXSC100	13705180
HEMOXSC100	13705181
HEMOXSC100	13703322
HEMOXSC100	13700446
HEMOXSC100	13700525
HEMOXSC100	13700458
HEMOXSC100	13700460
HEMOXSC100	13701388
HEMOXSC100	13702705
HEMOXSC100	13702643
HEMOXSC100	13704108
HEMOXSC100	13700226
HEMOXSC100	13710211
HEMOXSC100	13710215
HEMOXSC100	13710210
HEMOXSC100	13710252
HEMOXSC100	13703604
HEMOXSC100	13703606
HEMOXSC100	13703611
HEMOXSC100	13703747
HEMOXSC100	13703751
HEMOXSC100	13703746
HEMOXSC100	13703749
HEMOXSC100	13703750
HEMOXSC100	13703751
HEMOXSC100	13703779
HEMOXSC100	13704651
HEMOXSC100	13700504
HEMOXSC100	13704668
HEMOXSC100	13701861
HEMOXSC100	13700877
HEMOXSC100	13702108
HEMOXSC100	13705101
HEMOXSC100	13700854
HEMOXSC100	13710159
HEMOXSC100	13710160
HEMOXSC100	13710162
HEMOXSC100	13702976
HEMOXSC100	13702983
HEMOXSC100	13702984
HEMOXSC100	13700938
HEMOXSC100	13700855
HEMOXSC100	13702575
HEMOXSC100	13702061
HEMOXSC100	13702723
HEMOXSC100	13703012
HEMOXSC100	13710168
HEMOXSC100	13701651
HEMOXSC100	13705117
HEMOXSC100	13705172
HEMOXSC100	13705171
HEMOXSC100	13701260
HEMOXSC100	13701660
HEMOXSC100	13701110

HEMOXSC100	13700460
HEMOXSC100	13701451
HEMOXSC100	13705110
HEMOXSC100	13705110
HEMOXSC100	13700893
HEMOXSC100	13703755
HEMOXSC100	13703757
HEMOXSC100	13703758
HEMOXSC100	13703727
HEMOXSC100	13700304
HEMOXSC100	13703550
HEMOXSC100	13700005
HEMOXSC100	13700000
HEMOXSC100	13700000
HEMOXSC100	13700000
HEMOXSC100	13700000
HEMOXSC100	13700443
HEMOXSC100	13700570
HEMOXSC100	13705124
HEMOXSC100	13700524
HEMOXSC100	13700871
HEMOXSC100	13701148
HEMOXSC100	13700537
HEMOXSC100	13700555
HEMOXSC100	13700530
HEMOXSC100	13703243
HEMOXSC100	13703289
HEMOXSC100	13705000
HEMOXSC100	13705000
HEMOXSC100	13703719
HEMOXSC100	13700316
HEMOXSC100	13700300
HEMOXSC100	13700352
HEMOXSC100	13701909
HEMOXSC100	13701878
HEMOXSC100	13701927
HEMOXSC100	13700451
HEMOXSC100	13700457
HEMOXSC100	13703309
HEMOXSC100	13703316
HEMOXSC100	13703318
HEMOXSC100	13703645
HEMOXSC100	13700352
HEMOXSC100	13700600
HEMOXSC100	13700307
HEMOXSC100	13702986
HEMOXSC100	13702987
HEMOXSC100	13702989
HEMOXSC100	13702990
HEMOXSC100	13701115
HEMOXSC100	13700960
HEMOXSC100	13701610
HEMOXSC100	13702347
HEMOXSC100	13701241
HEMOXSC100	13701804
HEMOXSC100	13702448
HEMOXSC100	13702458
HEMOXSC100	13701642
HEMOXSC100	13701839
HEMOXSC100	13701192
HEMOXSC100	13703535

HEMOXSC100	13703539
HEMOXSC100	13703541
HEMOXSC100	13703765
HEMOXSC100	13703767
HEMOXSC100	13704159
HEMOXSC100	13704160
HEMOXSC100	13702635
HEMOXSC100	13702672
HEMOXSC100	13703528
HEMOXSC100	13703441
HEMOXSC100	13700052
HEMOXSC100	13705122
HEMOXSC100	13700510
HEMOXSC100	13700870
HEMOXSC100	13700871
HEMOXSC100	13702008
HEMOXSC100	13701815
HEMOXSC100	13701826
HEMOXSC100	13703547
HEMOXSC100	13703439
HEMOXSC100	13703440
HEMOXSC100	13701150
HEMOXSC100	13703110
HEMOXSC100	13705185
HEMOXSC100	13705186
HEMOXSC100	13700852
HEMOXSC100	13700205
HEMOXSC100	13700215
HEMOXSC100	13700227
HEMOXSC100	13700851
HEMOXSC100	13704294
HEMOXSC100	13704658
HEMOXSC100	13705096
HEMOXSC100	13701843
HEMOXSC100	13700257
HEMOXSC100	13700272
HEMOXSC100	13700280
HEMOXSC100	13701661
HEMOXSC100	13702760
HEMOXSC100	13705182
HEMOXSC100	13702658
HEMOXSC100	13702007
HEMOXSC100	13701108
HEMOXSC100	13701000
HEMOXSC100	13701100
HEMOXSC100	13701106
HEMOXSC100	13701206
HEMOXSC100	13702686
HEMOXSC100	13700845
HEMOXSC100	13705177
HEMOXSC100	13702291
HEMOXSC100	13710251
HEMOXSC100	13710255
HEMOXSC100	13710250
HEMOXSC100	13710261
HEMOXSC100	13710265
HEMOXSC100	13703003
HEMOXSC100	13701802
HEMOXSC100	13701100
HEMOXSC100	13701110
HEMOXSC100	13701112
HEMOXSC100	13705100
HEMOXSC100	13705108
HEMOXSC100	13705200
HEMOXSC100	13705196
HEMOXSC100	13705187
HEMOXSC100	13705189
HEMOXSC100	13705191
HEMOXSC100	13705192

HEMOXSC100	13705188
HEMOXSC100	13705197
HEMOXSC100	13710020
HEMOXSC100	13710036
HEMOXSC100	13710085
HEMOXSC100	13710118
HEMOXSC100	13710132
HEMOXSC100	13710136
HEMOXSC100	13710140
HEMOXSC100	13710146
HEMOXSC100	13700135
HEMOXSC100	13705093
HEMOXSC100	13702764
HEMOXSC100	13702415
HEMOXSC100	13710169
HEMOXSC100	13710170
HEMOXSC100	13710182
HEMOXSC100	13702502
HEMOXSC100	13701346
HEMOXSC100	13702207
HEMOXSC100	13702250
HEMOXSC100	13702251
HEMOXSC100	13702207
HEMOXSC100	13702057
HEMOXSC100	13702271
HEMOXSC100	13702261
HEMOXSC100	13702120
HEMOXSC100	13702117
HEMOXSC100	13702151
HEMOXSC100	13702556
HEMOXSC100	13702580
HEMOXSC100	13702500
HEMOXSC100	13702505
HEMOXSC100	13702506
HEMOXSC100	13702600
HEMOXSC100	13702602
HEMOXSC100	13702720
HEMOXSC100	13702715
HEMOXSC100	13702721
HEMOXSC100	13702720
HEMOXSC100	13702810
HEMOXSC100	13702821
HEMOXSC100	13702815
HEMOXSC100	13702850
HEMOXSC100	13702610
HEMOXSC100	13702662
HEMOXSC100	13702601
HEMOXSC100	13702710
HEMOXSC100	13702716
HEMOXSC100	13702742
HEMOXSC100	13703763
HEMOXSC100	13702464
HEMOXSC100	13703137
HEMOXSC100	13703242
HEMOXSC100	13703134
HEMOXSC100	13703136
HEMOXSC100	13700927
HEMOXSC100	13705182
HEMOXSC100	13701616
HEMOXSC100	13702673
HEMOXSC100	13700228
HEMOXSC100	13700220
HEMOXSC100	13700220
HEMOXSC100	13700252
HEMOXSC100	13700251
HEMOXSC100	13700270
HEMOXSC100	13700280
HEMOXSC100	13700211
HEMOXSC100	13710171
HEMOXSC100	13710173
HEMOXSC100	13702773
HEMOXSC100	13700101
HEMOXSC100	13700170

HEMOXSC100	13700480
HEMOXSC100	13700527
HEMOXSC100	13703113
HEMOXSC100	13703112
HEMOXSC100	13703529
HEMOXSC100	13703004
HEMOXSC100	13701144
HEMOXSC100	13700082
HEMOXSC100	13701129
HEMOXSC100	13701884
HEMOXSC100	13701125
HEMOXSC100	13701127
HEMOXSC100	13702689
HEMOXSC100	13703761
HEMOXSC100	13703760
HEMOXSC100	13702226
HEMOXSC100	13702711
HEMOXSC100	13702700
HEMOXSC100	13701018
HEMOXSC100	13702217
HEMOXSC100	13702264
HEMOXSC100	13705175
HEMOXSC100	13700169
HEMOXSC100	13701366
HEMOXSC100	13701368
HEMOXSC100	13702715
HEMOXSC100	13702215
HEMOXSC100	13702216
HEMOXSC100	13710151
HEMOXSC100	13705121
HEMOXSC100	13704021
HEMOXSC100	13705120
HEMOXSC100	13710186
HEMOXSC100	13700163
HEMOXSC100	13700164
HEMOXSC100	13700167
HEMOXSC100	13710149
HEMOXSC100	13705115
HEMOXSC100	13702657
HEMOXSC100	13710280
HEMOXSC100	13700847
HEMOXSC100	13700420
HEMOXSC100	13700849
HEMOXSC100	13700850
HEMOXSC100	13701976
HEMOXSC100	13702782
HEMOXSC100	13705170
HEMOXSC100	13701116
HEMOXSC100	13701117
HEMOXSC100	13700518
HEMOXSC100	13700517
HEMOXSC100	13702577
HEMOXSC100	13702580
HEMOXSC100	13702582
HEMOXSC100	13702712
HEMOXSC100	13702522
HEMOXSC100	13701876
HEMOXSC100	13702712
HEMOXSC100	13700426
HEMOXSC100	13702708
HEMOXSC100	13702682
HEMOXSC100	13702685
HEMOXSC100	13701893
HEMOXSC100	13705122

HEMOXSC100	13705125
HEMOXSC100	13705128
HEMOXSC100	13710367
HEMOXSC100	13702271
HEMOXSC100	13700439
HEMOXSC100	13700432
HEMOXSC100	13700436
HEMOXSC100	13700274
HEMOXSC100	13700317
HEMOXSC100	13700042
HEMOXSC100	13700045
HEMOXSC100	13700047
HEMOXSC100	13702756
HEMOXSC100	13704487
HEMOXSC100	13703143
HEMOXSC100	13704744
HEMOXSC100	13704745
HEMOXSC100	13704746
HEMOXSC100	13704747
HEMOXSC100	13704750
HEMOXSC100	13704751
HEMOXSC100	13704751
HEMOXSC100	13704760
HEMOXSC100	13702612
HEMOXSC100	13704724
HEMOXSC100	13704728
HEMOXSC100	13704767
HEMOXSC100	13700165
HEMOXSC100	13703132
HEMOXSC100	13704482
HEMOXSC100	13704483
HEMOXSC100	13704485
HEMOXSC100	13704477
HEMOXSC100	13704479
HEMOXSC100	13704496
HEMOXSC100	13704332
HEMOXSC100	13700275
HEMOXSC100	13701227
HEMOXSC100	13701228
HEMOXSC100	13700149
HEMOXSC100	13700150
HEMOXSC100	13700151
HEMOXSC100	13700153
HEMOXSC100	13700154
HEMOXSC100	13700155
HEMOXSC100	13700156
HEMOXSC100	13700157
HEMOXSC100	13700158
HEMOXSC100	13700159
HEMOXSC100	13700160
HEMOXSC100	13700161
HEMOXSC100	13700162
HEMOXSC100	13702759
HEMOXSC100	13703177
HEMOXSC100	13705013
HEMOXSC100	13701579
HEMOXSC100	13701616

HEMOXSC100	13701260
HEMOXSC100	13701264
HEMOXSC100	13701266
HEMOXSC100	13701269
HEMOXSC100	13701271
HEMOXSC100	13701273
HEMOXSC100	13703596
HEMOXSC100	13703658
HEMOXSC100	13703659
HEMOXSC100	13703665
HEMOXSC100	13704913
HEMOXSC100	13704736
HEMOXSC100	13700074
HEMOXSC100	13700294
HEMOXSC100	13700296
HEMOXSC100	13710061
HEMOXSC100	13710052
HEMOXSC100	13704035
HEMOXSC100	13704066
HEMOXSC100	13704167
HEMOXSC100	13704219
HEMOXSC100	13704259
HEMOXSC100	13704419
HEMOXSC100	13704437
HEMOXSC100	13704566
HEMOXSC100	13703807
HEMOXSC100	13703872
HEMOXSC100	13701237
HEMOXSC100	13701238
HEMOXSC100	13704742
HEMOXSC100	13704551
HEMOXSC100	13704552
HEMOXSC100	13700456
HEMOXSC100	13700930
HEMOXSC100	13700997
HEMOXSC100	13701002
HEMOXSC100	13701026
HEMOXSC100	13704715
HEMOXSC100	13700993
HEMOXSC100	13704493
HEMOXSC100	13700262
HEMOXSC100	13700811
HEMOXSC100	13700809
HEMOXSC100	13703600
HEMOXSC100	13703683
HEMOXSC100	13703697
HEMOXSC100	13701226
HEMOXSC100	13704882
HEMOXSC100	13704497
HEMOXSC100	13704494

HEMOXSC100	13700086
HEMOXSC100	13701298
HEMOXSC100	13704475
HEMOXSC100	13704476
HEMOXSC100	13702440
HEMOXSC100	13704737
HEMOXSC100	13700098
HEMOXSC100	13700526
HEMOXSC100	13700528
HEMOXSC100	13700531
HEMOXSC100	13700532
HEMOXSC100	13704111
HEMOXSC100	13704123
HEMOXSC100	13704687
HEMOXSC100	13704689
HEMOXSC100	13701190
HEMOXSC100	13701187
HEMOXSC100	13705055
HEMOXSC100	13705056
HEMOXSC100	13705144
HEMOXSC100	13701220
HEMOXSC100	13704713
HEMOXSC100	13702852
HEMOXSC100	13703169
HEMOXSC100	13703609
HEMOXSC100	13703518
HEMOXSC100	13703522
HEMOXSC100	13703543
HEMOXSC100	13703553
HEMOXSC100	13703578
HEMOXSC100	13703462
HEMOXSC100	13703601
HEMOXSC100	13703608
HEMOXSC100	13703671
HEMOXSC100	13710100
HEMOXSC100	13701461
HEMOXSC100	13701472
HEMOXSC100	13704463
HEMOXSC100	13703721
HEMOXSC100	13700533
HEMOXSC100	13703566
HEMOXSC100	13703568
HEMOXSC100	13703588
HEMOXSC100	13700343
HEMOXSC100	13700358
HEMOXSC100	13703178
HEMOXSC100	13700319
HEMOXSC100	13700216
HEMOXSC100	13704929
HEMOXSC100	13704937

HEMOXSC100	13701240
HEMOXSC100	13701239
HEMOXSC100	13704193
HEMOXSC100	13700418
HEMOXSC100	13704711
HEMOXSC100	13701613
HEMOXSC100	13703186
HEMOXSC100	13704474
HEMOXSC100	13705094
HEMOXSC100	13704464
HEMOXSC100	13704457
HEMOXSC100	13704840
HEMOXSC100	13701456
HEMOXSC100	13703364
HEMOXSC100	13703375
HEMOXSC100	13703138
HEMOXSC100	13703141
HEMOXSC100	13703142
HEMOXSC100	13703144
HEMOXSC100	13700311
HEMOXSC100	13700279
HEMOXSC100	13700310
HEMOXSC100	13701209
HEMOXSC100	13704971
HEMOXSC100	13704927
HEMOXSC100	13704969
HEMOXSC100	13704972
HEMOXSC100	13704976
HEMOXSC100	13704469
HEMOXSC100	13700284
HEMOXSC100	13701212
HEMOXSC100	13701538
HEMOXSC100	13701541
HEMOXSC100	13704964
HEMOXSC100	13704968
HEMOXSC100	13705123
HEMOXSC100	13705095
HEMOXSC100	13704491
HEMOXSC100	13700337
HEMOXSC100	13701460
HEMOXSC100	13704354
HEMOXSC100	13704355
HEMOXSC100	13703549
HEMOXSC100	13701231
HEMOXSC100	13700988
HEMOXSC100	13700989
HEMOXSC100	13703542
HEMOXSC100	13703544
HEMOXSC100	13703545
HEMOXSC100	13703546

HEMOXSC100	13704714
HEMOXSC100	13704723
HEMOXSC100	13704759
HEMOXSC100	13704946
HEMOXSC100	13703135
HEMOXSC100	13702612
HEMOXSC100	13702533
HEMOXSC100	13702575
HEMOXSC100	13701662
HEMOXSC100	13703556
HEMOXSC100	13703563
HEMOXSC100	13703570
HEMOXSC100	13703576
HEMOXSC100	13703582
HEMOXSC100	13701230
HEMOXSC100	13701035
HEMOXSC100	13701217
HEMOXSC100	13700501
HEMOXSC100	13703145
HEMOXSC100	13705070
HEMOXSC100	13705068
HEMOXSC100	13703962
HEMOXSC100	13700807
HEMOXSC100	13700808
HEMOXSC100	13700806
HEMOXSC100	13700805
HEMOXSC100	13702348
HEMOXSC100	13702383
HEMOXSC100	13702109
HEMOXSC100	13702114
HEMOXSC100	13702151
HEMOXSC100	13702428
HEMOXSC100	13702435
HEMOXSC100	13702446
HEMOXSC100	13702449
HEMOXSC100	13702460
HEMOXSC100	13702509
HEMOXSC100	13702547
HEMOXSC100	13702571
HEMOXSC100	13702580
HEMOXSC100	13702591
HEMOXSC100	13702597
HEMOXSC100	13702598
HEMOXSC100	13702601
HEMOXSC100	13702606
HEMOXSC100	13702607
HEMOXSC100	13702610
HEMOXSC100	13702615
HEMOXSC100	13702617
HEMOXSC100	13702618

HEMOXSC100	13701479
HEMOXSC100	13701275
HEMOXSC100	13702621
HEMOXSC100	13702743
HEMOXSC100	13702676
HEMOXSC100	13702691
HEMOXSC100	13702694
HEMOXSC100	13702813
HEMOXSC100	13702821
HEMOXSC100	13702854
HEMOXSC100	13702718
HEMOXSC100	13702818
HEMOXSC100	13702865
HEMOXSC100	13702642
HEMOXSC100	13704331
HEMOXSC100	13704293
HEMOXSC100	13704768
HEMOXSC100	13704770
HEMOXSC100	13704771
HEMOXSC100	13704773
HEMOXSC100	13704777
HEMOXSC100	13704833
HEMOXSC100	13704904
HEMOXSC100	13704905
HEMOXSC100	13704916
HEMOXSC100	13704917
HEMOXSC100	13704923
HEMOXSC100	13704945
HEMOXSC100	13705006
HEMOXSC100	13705008
HEMOXSC100	13705010
HEMOXSC100	13705023
HEMOXSC100	13705039
HEMOXSC100	13704907
HEMOXSC100	13704755
HEMOXSC100	13704762
HEMOXSC100	13704764
HEMOXSC100	13704772
HEMOXSC100	13704780
HEMOXSC100	13704781
HEMOXSC100	13704783
HEMOXSC100	13704784
HEMOXSC100	13704896
HEMOXSC100	13704933
HEMOXSC100	13704939
HEMOXSC100	13704935
HEMOXSC100	13704887
HEMOXSC100	13704934
HEMOXSC100	13710062
HEMOXSC100	13710063

HEMOXSC100	13710064
HEMOXSC100	13710066
HEMOXSC100	13710069
HEMOXSC100	13710070
HEMOXSC100	13710084
HEMOXSC100	13710086
HEMOXSC100	13710089
HEMOXSC100	13710092
HEMOXSC100	13710093
HEMOXSC100	13710094
HEMOXSC100	13710098
HEMOXSC100	13710099
HEMOXSC100	13701191
HEMOXSC100	13701486
HEMOXSC100	13700340
HEMOXSC100	13701194
HEMOXSC100	13703678
HEMOXSC100	13704313
HEMOXSC100	13704307
HEMOXSC100	13704308
HEMOXSC100	13701210
HEMOXSC100	13704914
HEMOXSC100	13700536
HEMOXSC100	13701232
HEMOXSC100	13700718
HEMOXSC100	13703680
HEMOXSC100	13705050
HEMOXSC100	13701235
HEMOXSC100	13703959
HEMOXSC100	13704878
HEMOXSC100	13704876
HEMOXSC100	13702095
HEMOXSC100	13703122
HEMOXSC100	13703126
HEMOXSC100	13703127
HEMOXSC100	13703982
HEMOXSC100	13704001
HEMOXSC100	13701523
HEMOXSC100	13701234
HEMOXSC100	13704712
HEMOXSC100	13702744
HEMOXSC100	13702087
HEMOXSC100	13702826