

9<sup>th</sup> October 2015

**FIELD SAFETY NOTICE**

**Thermo Scientific™ Oxoid™ M.I.C.Evaluator™ (M.I.C.E.™)**

**MA0101D/F - Penicillin 256µg/ml, lot 1708082**

Customers are to be advised of the following:

A technical investigation has confirmed that MA0101D/F Penicillin 256µg/ml, lot 1708082 contains insufficient concentration of penicillin.

Continued use of this product could result in the incorrect reporting of Minimum Inhibitory Concentration against penicillin.

Customers should review their internal quality control and test results reported using this lot and consider retesting and/or seek appropriate expert advice.

Accordingly, in keeping with our Quality Policy, we request that you destroy any remaining inventory of the lot listed above and contact Customer Services or your local distributor who will issue replacements.

This notice needs to be passed on to all who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

The Medicines and Healthcare products Regulatory Agency (MHRA) have been informed of this Field Safety Corrective Action.

If you experience any problems complying with this notice or you require any further information, please contact our Technical Support Department on +44 (0)1256 694238, or at [microbiology.techsupport.uk@thermofisher.com](mailto:microbiology.techsupport.uk@thermofisher.com)

Yours faithfully,



**Jim Filer**  
**Vice President, Quality and Regulatory, MBD**