

16 May 2018

**RE: Field Safety Notice –
Thermo Scientific™ Remel™ Salmonella paratyphi A-H Stained Suspension**

Dear Customer,
(cc Chairman Medical Board and Relevant Head of Departments)

Please be informed of the product notification for the batches below.

1. Item(s) affected

| Product Code | Bottle Lot No. | Carton Lot No, |
|--------------|----------------|--|
| R30855101 | 2169556 | 2184561, 2308897, 2228271, 2213610 & 2258333 |

2. Reason for notice

An internal technical investigation has determined that Thermo Scientific™ Remel™ Salmonella paratyphi A-H Stained Suspension (R30855101) is showing granularity, upon dispensing, which may be misinterpreted as a positive result when used according to the Instructions For Use (IFU).

Please refer to Field Safety Notice for full details.

3. Actions to be taken

- Review any remaining inventory of the affected lots
- Kindly fill up the acknowledgement section below and send back to us.

Please refer to attached memo for more information and feel free to contact us should you have any queries regarding this notice.

Thank you for your prompt attention.

Yours Sincerely,



Ji Wenxin
Tel: 6499 9992
Email: wenxin.ji@thermofisher.com

Acknowledgement (To be completed by end-user customer)

I confirm receipt of the F.S.N. for ThermoScientific™ Remel™ Salmonella paratyphi A-H Stained Suspension,
Dated 11th May 2018.

Signed.....Name (please print).....

Position.....

Date.....Email Address/Telephone Number.....

Company Name & Stamp: