

---

## Provisional Authorisation for COVID-19 Tests

|   |  |
|---|--|
| <i>Applicant</i>  | AITbiotech Pte Ltd   |
| <i>Name of test</i>   | abTES™ COVID-19 qPCR I Kit   |
| <i>Intended purpose (As per manufacturer's information for use)</i> | <i>The abTES™ COVID-19 qPCR I Kit is a qualitative real-time polymerase chain reaction (qPCR) kit which enables simultaneous detection of two COVID-19- specific signature regions from its non-structure polypeptide in a single reaction. It also includes detection of human housekeeping gene, GAPDH, as an Internal Control (IC) to identify possible PCR inhibitions from sample processing. This kit has been validated on samples extracted from sputum, nasopharyngeal and throat swabs</i> |
| <i>Date of Provisional Authorisation</i>                            | 5 March 2020   |

---