
Provisional Authorisation for COVID-19 Tests

| | |
|---|--|
| <i>Applicant</i> | Alcotech Pte Ltd |
| <i>Name of test</i> | 2019-nCoV IgG/IgM Rapid Test Cassette |
| <i>Intended purpose (As per manufacturer's information for use)</i> | <p><i>The 2019-nCoV IgG/IgM Rapid Test Cassette is a lateral flow chromatographic immunoassay for the qualitative detection of IgG and IgM antibodies to 2019-nCoV in human whole blood, serum or plasma specimen.</i></p> <p><i>The result from this test is not to be used for confirmatory testing or as sole basis for diagnosis. The results will have to be interpreted together with clinical presentation and are to be confirmed with supplemental testing (e.g. RT-PCR).</i></p> |
| <i>Date of Provisional Authorisation</i> | 5 June 2020 |
