

---

## Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	BioMérieux Singapore Pte Ltd
<i>Name of test</i>	BioFire® Respiratory Panel 2.1 (RP2.1)
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>The BioFire Respiratory Panel 2.1 (RP2.1) is a multiplexed nucleic acid test intended for the simultaneous qualitative detection and differentiation of nucleic acids from multiple viral and bacterial respiratory organisms, including nucleic acid from Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2), in nasopharyngeal swabs (NPS) obtained from individuals suspected of COVID-19 by their healthcare provider.</i>
<i>Date of Provisional Authorisation</i>	14 May 2020

---