

---

## Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Biomed Diagnostics Pte Ltd
<i>Name of test</i>	InnoCare SARS-CoV-2 RT-PCR Kit
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>The InnoCare SARS-CoV-2 RT-PCR Kit is a qualitative reverse transcription polymerase chain reaction (RT-PCR) kit which enables simultaneous detection of two SARS-CoV-2 specific signature regions. The kit contains all the necessary RT-PCR reagents for rapid, sensitive and specific detection using target-specific primers and double-labeled probes. The RNA should be extracted from respiratory specimen such as nasopharyngeal swabs. The InnoCare SARS-CoV-2 RT-PCR Kit enables detection of human housekeeping gene, as an Internal Control (IC) to identify possible PCR inhibitions from sample processing.</i>
<i>Date of Provisional Authorisation</i>	14 October 2020

---