
Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Camtech Diagnostics Pte Ltd
<i>Name of test</i>	Camtech COVID-19 Antigen Rapid Test
<i>Intended purpose (As per manufacturer's information for use)</i>	<p><i>The Camtech COVID-19 Antigen Rapid Test Kit is used for the qualitative detection of novel coronavirus (SARS-CoV-2) antigen in human throat swab samples in vitro. The kit should only be used as a supplementary indicator or used in conjunction with nucleic acid detection in the diagnosis of suspected COVID-19 cases. It cannot be used as a sole basis for the diagnosis and exclusion of pneumonitis patients infected by the novel coronavirus and is not suitable for screening of the general population. Negative test results do not exclude the possibility of virus infection. If test result is negative and the patient has clinical symptoms, it is recommended to use virus isolation and culture for confirmation, and a comprehensive diagnosis by the attending physician.</i></p> <p><i>The kit is to be used in hospitals or by healthcare professionals.</i></p>
<i>Date of Provisional Authorisation</i>	28 August 2020
