

---

## Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Cellbae Pte Ltd
<i>Name of test</i>	TEPAT 1.0 SARS-CoV-2 Multiplex RT-PCR Kit
<i>Intended purpose (As per manufacturer's information for use)</i>	<p><i>The TEPAT 1.0 SARS-CoV-2 Multiplex RT-PCR Kit is a real-time reverse transcriptase polymerase chain reaction (RT-PCR) technology for the qualitative and simultaneous detection of 2 genetic targets of SARS-CoV-2-specific RNA and 1 genetic target of human-specific mRNA from nasopharyngeal swabs.</i></p> <p><i>A positive result from the test may indicate the presence of SARS-CoV-2-specific RNA in the test sample. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of SARS-COV-2 infection.</i></p>
<i>Date of Provisional Authorisation</i>	14 September 2020

---