
Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Credo Diagnostics Biomedical Pte. Ltd
<i>Name of test</i>	VitaPCR™ SARS-CoV-2 Gen 2 Assay
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>The VitaPCR™ SARS-CoV-2 Gen 2 Assay performed on the VitaPCR™ Instrument is a rapid molecular in vitro diagnostic test utilizing a real-time reverse transcription polymerase chain reaction (RT-PCR) amplification technology for the qualitative detection of Coronavirus Disease 2019 (COVID-19) viral RNA in nasopharyngeal (NP) or oropharyngeal (OP) swabs from patients with signs and symptoms of respiratory infection.</i>
<i>Date of Provisional Authorisation</i>	16 Jun 2020
