
Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	DSO National Laboratories
<i>Name of test</i>	RESOLUTE 1.0 LITE
<i>Intended purpose (As per manufacturer's information for use)</i>	<p><i>RESOLUTE 1.0 LITE is a real-time RT-PCR test intended for the qualitative detection of nucleic acid from the SARS-CoV-2 in nasopharyngeal swabs specimen from individuals with signs and symptoms of infection who are suspected of COVID-19.</i></p> <p><i>RESOLUTE 1.0 LITE is intended for the direct amplification of Coronavirus SARS-CoV-2 RNA from nasopharyngeal swabs (NPS) in Universal Transport Medium (UTM). RNA extraction is not required for the use of this kit.</i></p>
<i>Date of Provisional Authorisation</i>	25 Jun 2020
