
Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	I-Test Pte Ltd
<i>Name of test</i>	i-Test COVID -19 IgM/IgG Antibody Rapid Test (Colloidal Gold)
<i>Intended purpose (As per manufacturer's information for use)</i>	<p><i>This kit is used for in-vitro qualitative determination of Novel coronavirus (COVID-2019) Antibody in human serum or plasma or whole blood.</i></p> <p><i>As there may be a few days' delay between the viral infection and development of antibodies in an individual, any negative result from this test should not be used as a sole basis to rule out COVID-19 infection. The negative results have to be carefully considered together with clinical presentations. Any positive result from this test should be confirmed with supplementary laboratory testing (e.g. RT-PCR).</i></p>
<i>Date of Provisional Authorisation</i>	20 May 2020
