

---

## Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	IceSing Bioscience Pte Ltd
<i>Name of test</i>	DirectDetect SARS-CoV-2 Detection Kit
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>The DirectDetect SARS-CoV-2 Detection Kit is a real-time RT-PCR test intended for the qualitative detection of ORF1ab and N gene of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) in oropharyngeal swabs from individuals that are suspected Corona Virus Disease 2019 cases, suspected clustered cases, or that needs to be infectious diagnosis or differential diagnosis.</i>
<i>Date of Provisional Authorisation</i>	12 January 2021

---