

---

## Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Illumina Singapore Pte Ltd
<i>Name of test</i>	Illumina COVIDSeq Test
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>The Illumina COVIDSeq Test is a Next-Generation Sequencing (NGS) in vitro diagnostic test on the Illumina NovaSeq 6000 Sequencing System, NextSeq 2000 Sequencing System, NextSeq 500 Sequencing System, NextSeq 550 Sequencing System, or NextSeq 550Dx Instrument intended for the qualitative detection of SARS-CoV-2 RNA from nasopharyngeal (NP) swabs, oropharyngeal (OP) swabs, anterior nasal swabs, mid-turbinate nasal swabs, nasopharyngeal wash/aspirates, and nasal aspirates from individuals suspected of COVID-19 by their healthcare provider.</i>
<i>Date of Provisional Authorisation</i>	3 September 2020

---