
Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Link Trust Laboratories Pte Ltd
<i>Name of test</i>	DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection kit
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>The DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit is a real-time reverse transcriptase (RT)-PCR test intended for the qualitative detection of nucleic acid from SARS-CoV-2 in nasopharyngeal swabs, oropharyngeal (throat) swabs, anterior nasal swabs, mid-turbinate nasal swabs, nasal aspirates, nasal washes, bronchoalveolar lavage (BAL) fluid and sputum from individuals suspected of COVID-19 by their healthcare provider.</i>
<i>Date of Provisional Authorisation</i>	14 August 2020
