
Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Roche Diagnostics Asia Pacific Pte Ltd
<i>Name of test</i>	Cobas® SARS-CoV-2
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>cobas® SARS-CoV-2 for use on the cobas® 6800/8800 Systems is a real-time RT-PCR test intended for the qualitative detection of nucleic acids from SARS-CoV-2 in clinician-instructed self-collected nasal swab specimens (collected on site), and clinician-collected nasal, nasopharyngeal and oropharyngeal swab samples from patients with signs and symptoms suggestive of COVID-19 (e.g., fever and/or symptoms of acute respiratory illness).</i>
<i>Date of Provisional Authorisation</i>	19 March 2020
