
Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Tigermed Asia Pacific Private Limited
<i>Name of test</i>	GenBody COVID-19 Ag
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>GenBody COVID-19 Ag is an immunochromatographic assay kit for the qualitative detection of SARS-CoV-2 antigen in nasopharyngeal and oropharyngeal swab from human. The negative results from this test do not preclude the possibility of COVID-19 infection. The negative results will have to be confirmed with molecular assays.</i>
<i>Date of Provisional Authorisation</i>	17 September 2020
