

---

## Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Unison Collaborative Pte Ltd
<i>Name of test</i>	Panbio™ COVID-19 Antigen Self-test
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>The Panbio™ COVID-19 Antigen Self-Test is a single-use, in vitro (outside the body) visually read rapid immunoassay that uses a human nasal swab specimen for the qualitative detection of nucleocapsid protein SARS-CoV-2 antigen (Ag). The Panbio™ COVID-19 Antigen Self-Test is intended to be used manually by untrained lay users (self testing) in a private setting to aid in the diagnosis of an active SARS-CoV-2 infection. Children under 14 years should be supported by an adult.</i>
<i>Date of Provisional Authorisation</i>	31 May 2021

---