

---

## Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Unison Collaborative Pte Ltd
<i>Name of test</i>	RightSign COVID-19 IgG/IgM Rapid Test Cassette
<i>Intended purpose (As per manufacturer's information for use)</i>	<p><i>The COVID-19 IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of IgG and IgM antibodies to SARS-CoV-2 in human whole blood, serum, or plasma.</i></p> <p><i>The results from this test is not to be used for confirmatory testing or as sole basis for diagnosis. The results will have to be interpreted together with clinical presentation and are to be confirmed with supplemental testing (e.g. RT-PCR).</i></p>
<i>Date of Provisional Authorisation</i>	11 June 2020

---