
Provisional Authorisation for COVID-19 Tests

<i>Applicant</i>	Veredus Laboratories Pte Ltd
<i>Name of test</i>	VereRT™ ZeroPrep™ COVID-19 PCR Kit
<i>Intended purpose (As per manufacturer's information for use)</i>	<i>VereRT™ ZeroPrep™ COVID-19 PCR Kit is a Reverse Transcription Polymerase Chain Reaction (RT-PCR)-based In Vitro Diagnostic (IVD) test intended for the qualitative detection of nucleic acid belonging to SARS-CoV-2 also known as 2019-nCoV (2019 novel Coronavirus) in nasopharyngeal swab specimens. Testing with VereRT™ ZeroPrep™ COVID-19 PCR Kit is intended for use directly with Viral Transport Medium (VTM) or Universal Transport Medium (UTM) containing the swab specimen without the need for viral RNA extraction.</i>
<i>Date of Provisional Authorisation</i>	14 August 2020
