


Quality CONTROL in Clinical Chemistry



Certified Reference Materials from HSA

27 Jul 2017

The HSA CRMs serve as independent quality controls for clinical laboratories. The CRMs are produced in accordance with the requirements of ISO 17034:2016 and ISO Guide 35:2006.

Every CRM is accompanied by a Certificate of Analysis, containing traceable certified values for each analyte.

To order, email us at:

HSA_GML@HSA.gov.sg

Chemical Metrology Laboratory
Applied Sciences Group
Health Sciences Authority
1 Science Park Road
#01-05/06, The Capricorn
Singapore Science Park II
Singapore 117528

*The CRMs are stable for a long term when stored below -60°C and can be used for up to a month when stored at -20°C.

^The CRM is stable for a long term when stored at -20°C.

CRMs Available

HRM-3002B*

Creatinine • Glucose • Total Cholesterol • HDL-Cholesterol • LDL-Cholesterol • Total Glycerides • Urea • Uric Acid

(Matrix: Human Sera)

HRM-2005A*

Calcium • Potassium • Sodium • Magnesium • Iron • Chloride

(Matrix: Human Sera)

HRM-3003A^

HRM-3003B^

Glycated haemoglobin (HbA1c)

(Matrix: Human Blood)

HRM-2011A*

Sodium • Chloride • Copper • Selenium • Phosphorous

(Matrix: Human Sera)

